SCOOTER ELETTRICO
MOBILITY 250 E MOBILITY 260

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INDICE

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Via _____ Località _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com**MADE IN P.R.C.**

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Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com

ELECTRIC SCOOTER
MOBILITY 250 AND MOBILITY 260

INSTRUCTION MANUAL

INDEX

1. CODES.....	PAG.3
2. INTRODUCTION.....	PAG.3
3. INTENDED USE.....	PAG.3
4. DECLARATION OF CONFORMITY.....	PAG.3
4.1 Applicable regulations and directives.....	pag. 4
5. GENERAL WARNINGS.....	PAG.4
6. SYMBOLS.....	PAG.4
7. GENERAL DESCRIPTION.....	PAG.5
8. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY.....	PAG.7
9. SAFETY.....	PAG.11
10. WARNINGS ON USE.....	PAG.12
10.1 General precautions.....	pag. 12
10.2 Changes.....	pag. 12
10.3 Before every use.....	pag. 13
10.4 Weight limits.....	pag. 13
10.5 Temperature.....	pag. 13
10.6 Tyre Inflation.....	pag. 13
11. FUNCTIONING OF THE SCOOTER.....	PAG.13
11.1 Control panel.....	pag. 13
11.2 Steering adjustment.....	pag. 15
11.3 Adjustment of the handlebar.....	pag. 16
11.4 Armrest adjustment.....	pag. 16
11.5 Seat adjustment.....	pag. 17
11.6 Wheel locking/unlocking function.....	pag. 18
12. STARTING THE SCOOTER FOR THE FIRST TIME.....	PAG.18
13. ASSEMBLY AND DISASSEMBLY.....	PAG.19
14. DRIVING.....	PAG.19
14.1 Driving up and down.....	pag. 19
14.2 Basic driving.....	pag. 20
14.3 Steering.....	pag. 20
14.4 Control through narrow passages.....	pag. 20
14.5 Driving up and down slopes.....	pag. 21
14.6 Engine braking system.....	pag. 22
14.7 Stopping brake.....	pag. 23
14.8 Thermal protection.....	pag. 23
15. TRANSPORT.....	PAG.23
16. CLEANING.....	PAG.23
16.1 Cleaning the tyres.....	pag. 23
16.2 Cleaning the body.....	pag. 23
16.3 Cleaning the seat.....	pag. 23
17. GENERAL MAINTENANCE INSTRUCTIONS.....	PAG.24
17.1 Controls and reminders.....	pag. 24
17.2 Wheel replacement.....	pag. 25
17.3 Control panels, battery charger and rear electronics.....	pag. 25
17.4 Storing the scooter.....	pag. 25
17.5 Battery and charging.....	pag. 25
17.6 Charging the battery.....	pag. 26
17.7 If the scooter does not work.....	pag. 27
18. PROBLEMS AND SOLUTIONS.....	PAG.28
19. CONDITIONS OF DISPOSAL.....	PAG.29
19.1 General terms of disposal.....	pag. 29
19.2 Instructions for correct disposal in accordance with European directive 2012/19/UE: ..	pag. 29
19.3 Disposal of exhausted batteries - (Directive 2006/66/EC):.....	pag. 29
20. TECHNICAL SPECIFICATIONS.....	PAG.29
20.1 CN250 technical specifications.....	pag. 29
20.2 CN260 technical specifications.....	pag. 30
21. WARRANTY.....	PAG.32
21.1 Warranty on rechargeable batteries (if provided).....	pag. 32
22. REPAIRS.....	PAG.32
23. REPLACEMENT PARTS.....	PAG.32
24. NON-LIABILITY CLAUSE.....	PAG.32

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices

1. CODES

CN250 Mobility Scooter 250

CN260 Mobility Scooter 260

2. INTRODUCTION

Thank you for choosing an ARDEA MOBILITY scooter by Moretti. Moretti electric scooters have been designed and constructed to satisfy all your needs for practical, correct and safe use. This manual contains useful suggestions for using your scooter properly and safely. Please read this manual in full before you use the scooter. Should you have any queries, please contact your dealer for suitable advice and assistance.

NB Check to ensure that no parts of the product have been damaged during shipment. Do not use the product if it is damaged and contact the retailer for further instructions.

3. INTENDED USE

The electric scooter is intended for people with mobility impairments.

WARNING!

- Do not use the product for any purpose other than that specified in this manual.
- The electric scooter must always be used by people of sound body and mind and not under the influence of drugs or alcohol.
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than that specified in this manual.
- The manufacturer reserves the right to modify the device and its relevant manual without prior notice, in order to improve the product's features.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of ELECTRIC SCOOTER - ARDEA MOBILITY complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

4.1 Applicable regulations and directives

In order to satisfy safety standards for users, Moretti S.p.A. complies with the following standards:

- EN 12184:2014;
- EN 60601-1:2007;
- EN 60601-1-2:2006.

5. GENERAL WARNINGS

WARNING!

After using the scooter even for a short distance, do not touch the engine: risk of burns from overheating.

WARNING!

Take extreme care in the presence of moving parts that could trap limbs and cause injury.

- Please consult this manual carefully for correct use of the device.
- Always consult your doctor or therapist for correct use of the device.
- Keep the packaged product away from sources of heat, as the packaging is cardboard.
- The lifetime of the device depends on wear and tear of non-repairable and/or non-replaceable parts.
- Always take care when children are present.
- Maximum user weight: 203 kg.
- Classification: class B (EN 12184).
- The product is not suitable for the visually impaired.
- The device must not be used by children under 12 years of age.
- The scooter is not designed for use as a seat in a motor vehicle.
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

6. SYMBOLS

Product code

Serial number

Unique Device Identification

CE mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Attention, should read the instructions

Production dates

Waste Electrical & Electronic Equipment (WEEE)

B Type

IPX4

Protection class

Max loading weight

7. GENERAL DESCRIPTION
CN250

Illustrative images

CN260

ENGLISH

Illustrative images

8. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY

The CN250 and CN260 Scooters are designed for use in the electromagnetic conditions specified in the following tables. The user of the CN250 and CN260 Scooters must make sure that it is used in the conditions specified.

Warnings related to risks from electromagnetic interference:

Interference can cause accidental movement and/or erratic vehicle control.

Electromagnetic energy emitted by sources such as:

- radio transmitting stations;
- TV broadcasting stations;
- amateur radio stations;
- shop alarms;
- mobile phones and cordless telephones can interfere with electric wheelchairs and electric scooters.

Interference can cause the electromagnetic brake to release, the device to move in an unintended direction. They can also permanently damage the permanently damage the device's electronic control unit. In the everyday environment there are a number number of relatively intense electromagnetic sources. The intensity of EM energy can be measured in volts per metre (V/m). Your scooter is equipped with protection against interference electromagnetic interference tested and certified in accordance with current international requirements, up to a certain intensity. This is called the 'immunity level'.

Table 1

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC EMISSIONS		
The scooters are intended for use in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.		
EMISSION TESTS	CONFORMITY	ELECTROMAGNETIC ENVIRONMENT - GUIDE
RF emissions CISPR 11	Group 1	The scooters use RF energy only for internal functions and for charging the battery. Their RF emissions are therefore very low and in all likelihood do not interfere with nearby electronic equipment.
RF emissions CISPR 11	Class B	The scooters are suitable for use in all buildings, including the home and buildings connected directly to the low voltage public power grid for residential applications.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage flutter/flicker IEC 61000-3-3	Conforming	

Table 2

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC EMISSIONS			
The scooters are intended for use in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.			
IMMUNITY TESTS	TEST LEVEL IEC 60601	CONFORMITY LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDE
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Information for using the device and charging the battery. The floor must be wood, concrete or ceramic. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Electrical Fast Transient/ Burst Immunity IEC 61000-4-4	±2 kV for power lines	±2 kV for power lines	The quality of the mains supply must be that of a typical commercial or hospital environment.
Over voltages IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	The quality of the mains supply must be that of a typical commercial or hospital environment.
Voltage drop outs, brief interruptions and variations of voltage on the power supply lines. IEC 61000-4-11	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	The quality of the mains supply must be that of a typical commercial hospital when charging the battery.
Mains frequency magnetic field (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should have the typical characteristics of a commercial or hospital environment. (Applies both to using the device and charging the battery).
Note: UT is the AC mains voltage prior to application of the test level.			

Table 3

MANUFACTURER'S GUIDE AND DECLARATION - ELECTROMAGNETIC IMMUNITY			
The scooters are intended for use in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.			
IMMUNITY TESTS	TEST LEVEL IEC 60601	CONFORMITY LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDE
conducted RF IEC 61000-4-6	3Vrms From 150 kHz to 80 MHz	1 Vrms	Portable and mobile RF communications devices should not be used near any part of the scooters, including the cables, unless the recommended separation distances are observed, calculated by the equation applicable to the transmitter frequency. Recommended distances: $d = [3,5/V1] \sqrt{P}$ from 0 Hz to 80 MHz $d = [3,5/E1] \sqrt{P}$ from 80 MHz to 800MHz, $d = [7/E1] \sqrt{P}$ from 800 MHz to 2.5 GHz, where "P" indicates the maximum rated output of the transmitter in Watt (W) according to the transmitter manufacturer and "d" is the recommended separation distance in metres (m). The intensity of the field of fixed RF transmitters, as determined by an electromagnetic survey of the site, a could be lower than the level of conformity of each frequency range. b There may be interference near devices with the following symbol: 
radiated RF IEC 61000-4-3	20 V/m From 80 MHz to 2.5 GHz	1 V/m	
NOTE 1: at 80 MHz and 800 MHz, apply the separation distance for the higher frequency band. NOTE 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a The intensity of the field for fixed transmitters, such as base stations for radio-telephones (cellulars and cordless) and terrestrial mobile radios, amateur radio equipment, AM and FM radio transmitters and TV transmitters may not be known precisely in theory. To establish an electromagnetic environment caused by fixed RF transmitters, one must run an on-site electromagnetic inspection. If the field intensity measured in the place where Nefti is used exceeds the applicable RF level of conformity indicated above, normal functioning of the scooters must be monitored. If abnormalities are observed, supplementary measures may be necessary, such as re-orientation or repositioning of the scooters. b The field intensity over the frequency range from 150 kHz to 80 MHz should be less than 3 V/m.			

ENGLISH

Table 4

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RADIO COMMUNICATIONS EQUIPMENT AND THE SCOOTERS						
<p>The scooters are intended to operate in an electromagnetic environment in which RF disturbances are under control. The client or user of the scooters can contribute to preventing EM interference by ensuring a minimum distance between the mobile and portable RF communications devices (transmitters) and the scooter battery charger during charging and the internal scooter devices during their operation as indicated below, in relation to the maximum output power of the radio communication equipment.</p>						
TESTS	SEPARATION DISTANCE AS A FUNCTION OF THE FREQUENCY OF THE TRANSMITTER (M)					
Maximum rated output power of transmitter (W)	During battery charging	During use of the CN250-CN260 Scooters	During battery charging	During use of the CN250-CN260 Scooters	During battery charging	During use of the CN250-CN260 Scooters
	from 150 kHz to 80 MHz $d=[3,5/V1]\sqrt{P}$	from 150 kHz to 80 MHz $d=[3,5/V1]\sqrt{P}$	from 80 MHz to 800 MHz $d=[3,5/E1]\sqrt{P}$	from 80 MHz to 800 MHz $d=[3,5/E1]\sqrt{P}$	from 800 MHz to 2.5 GHz $d=[7/E1]\sqrt{P}$	from 800 MHz to 2.5 GHz $d=[7/E1]\sqrt{P}$
0.01	0.35	0.35	0.35	0.35	0.7	0.7
0.1	1.1	1.1	1.1	1.1	2.2	2.2
1	3.5	3.5	3.5	3.5	7	7
10	11	11	11	11	7.3	7.3
100	35	35	35	35	70	70
<p>For transmitters whose maximum rated output power is not listed above, the recommended separation distance d, in metres (m), can be calculated using the equation applicable to the transmitter frequency, where P is the maximum rated output power of the transmitter, in Watt (W), as given by the transmitter's manufacturer.</p> <p>NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>						

9. SAFETY

 <p>Safety Belt</p> <p>Always use the safety belt and keep your feet resting on the scooter at all times</p>	 <p>Do not use the scooter under the influence of alcohol</p>
 <p>Do not use electronic radio transmitters such as walkie-talkies or mobile telephones</p>	 <p>Make sure there are no obstacles nearby when you are driving the scooter</p>
 <p>Do not make tight curves or brake suddenly when you are driving the scooter</p>	 <p>Do not drive the scooter in traffic</p>
 <p>Do not attempt to mount curbs higher than those indicated in the technical specifications</p>	 <p>Do not put your hands or legs out when driving the scooter</p>
 <p>Do not drive the scooter on slippery roads or when it is snowing.</p>	 <p>Do not allow unsupervised children to play near the scooter when the batteries are charging.</p>

WARNING!
Do not use the scooter without first reading and understanding this manual completely.

ENGLISH

10. WARNINGS ON USE

10.1 General precautions

1. Do not use the scooter on public highways. Please be aware that it may be difficult for drivers to see you when you are sitting on the scooter. Always use pedestrian routes. Comply with the rules of pedestrian traffic. Wait for the way to be clear, then proceed with extreme caution.
2. Before getting on and off your scooter, always check that it is switched off and that the battery charger is disconnected, to avoid injury to yourself and to others.
3. Always check that the drive wheels are in drive mode before driving. Do not switch off the scooter when it is still moving, as it would stop very suddenly.
4. Do not use this product or any other optional accessory available without first reading and understanding this manual completely. If you do not understand the warnings, the precautions or the instructions, contact the retailer or a healthcare operator before attempting to use the device, to avoid injury or damage.
5. In certain cases, such as several medical conditions, the scooter user will need practice in use of the scooter in the presence of a qualified assistant, who may be a relative or a professional qualified to assist the user in use of the scooter during day-to-day activities.
6. Avoid lifting or moving the scooter by gripping one of its removable parts, including the armrests, seat or fairing. This could cause personal injury and/or damage the scooter.
7. Avoid using the scooter beyond its limits, as described in this manual.
8. Do not sit on the scooter when it is inside a moving vehicle.
9. Keep your hands away from the wheels (tyres) while driving the scooter. Please remember that loose clothing could get caught in the wheels.
10. Consult your doctor if you are taking drugs or if you have specific physical limitations. Certain drugs and certain limitations could prevent safe use of the scooter.
11. Always check that the wheel locking/unlocking lever is locked or unlocked.
12. Do not remove the roll-over protection wheels if they are present on the scooter.
13. Contact with tools could cause an electric shock; do not connect an electrical extension cord to the AC / DC converter or to the battery charger.
14. The scooter can be driven on a maximum estimated slope of 12°/21%. However, take great care when using the scooter on slopes over 10°/17%.
15. Do not drive on a slope over the limit of the scooter.
16. Avoid reversing off steps, pavements or other obstacles. This could cause the scooter to fall/tip over.
17. Always reduce speed and maintain a stable centre of gravity on tight curves. Do not use the maximum power on tight curves.
18. Using the scooter in the rain, snow or foggy conditions and on icy or slippery surfaces could damage the electrical wiring.
19. Never sit on the scooter when it is being used with any type of lift or lifting device. The scooter is not designed for this use. The manufacturer is not liable for any damage or injury caused by improper use.

10.2 Changes

The scooters in the Ardea Mobility range have been designed and constructed to satisfy all your needs for practical, correct and safe use. Do not alter, add, remove or disable any part or function of the scooter under any circumstances. This could cause damage to the scooter and/or personal injury.

1. The control parameter must only be altered by authorised technicians, for safety reasons.
2. It is recommended to carry out a safety check before every use of the scooter to ensure that it is functioning safely.

10.3 Before every use

1. Check the tyres are correctly inflated, if present.
2. Check all the electrical connections and make sure they are properly connected and not corroded.
3. Check all the cables and make sure they are connected properly.
4. Check the brakes.

10.4 Weight limits

1. Check the information on maximum capacity in the table of the technical specifications in this manual or on the label on the device. The power of the scooter is assessed according to the capacity.
2. Comply with the weight limits indicated for your scooter. Exceeding the maximum capacity invalidates the warranty. The manufacturer cannot be held liable for injury or damage to property caused by failure to comply with the weight limit.
3. Do not carry passengers on the scooter, as this could alter the centre of gravity and cause it to tip over or fall.

10.5 Temperature

1. Certain parts of the scooter are affected by changes in temperature. The controller only functions at a temperature between -25 C° and 50 C°.
2. The batteries could freeze at extremely low temperatures and prevent your scooter from operating. In extremely hot climates, the scooter might only function at lower speeds, due to a safety feature of the controller that prevents damage to the engine and other electrical components.

10.6 Tyre Inflation

1. If the scooter is equipped with tyres, the air pressure must be checked air pressure at least once a week.
2. Correct inflation pressure prolongs the life of the tyres and ensures better performance while riding.
3. Do not underinflate or overinflate your tyres. It is extremely important that tyre pressure is maintained between 30-25 psi (2-2.4 bar) at all times.
4. Inflating tyres with an unregulated air source can cause them to burst.

11. FUNCTIONING OF THE SCOOTER

Your scooter is simple to use. However, it is recommended to read the following instructions carefully, to familiarise yourself with your new vehicle.

11.1 Control panel


WARNING!


Before switching on the scooter, always be aware of your surroundings when selecting the speed. It is recommended to select the slowest speed for use indoors. It is recommended to select a speed that allows you to control the scooter safely for use outdoors. The steps listed below are necessary to operate your scooter safely.

A) Speed regulator

- Press the button with the rabbit figure to increase the speed.
- Press on the button with the turtle figure to decrease the speed.

Illustrative images

Illustrative images

B) Key

- The key serves to switch the scooter on and off.
- Insert the key.
- Turn the key through 90° in a clockwise direction to switch on the scooter. It will beep and the needle on the battery indicator will move to indicate that the scooter is switched on.
- Return the key to the vertical position to switch off the scooter.

NOTE: Always remove the key before getting off the scooter and before lowering the handlebar.

C) Horn

- Press one of the two keys with the trumpet figure to sound the horn.

D) Control lever

- This allows the scooter to be used with one hand only.
- It controls the forward and reverse speed.

DRIVING

- to drive forwards:
 - use the fingers of your right hand to pull back the accelerator lever on the right;
 - use your left thumb to push the accelerator lever on the left.

NOTE: Always bring the scooter to a complete stop before changing direction from forward to reverse or from reverse to forward.

- to drive in reverse:
 - use the fingers of your left hand to pull back the accelerator lever on the left;
 - use your right thumb to push the accelerator lever on the right.

NOTE: When the accelerator control lever is completely released, it returns automatically to the central stop position and activates the scooter brakes, bringing the scooter to a complete stop. You will hear a "click" when the brake is on.

E) Battery indicator

The battery charge level will be displayed on the LCD screen once the scooter is switched on.

F) Four direction indicators

Press the button once to activate simultaneous flashing of both direction indicators; flashing is accompanied by an acoustic signal. Press again to deactivate them.

G) Switching on the lights

Press the button once to switch on the lights.
Press the button again to switch off the lights.

H) Direction indicators

Press the buttons to activate the direction indicator.

1. Press the left-hand button to activate the left direction indicator.
2. Press the right-hand button to activate the right direction indicator.
3. The direction indicator will switch off automatically after 15 sec. Press the button again to deactivate it sooner.

I) Battery charger socket

The battery charger socket is positioned on the steering to avoid having to lean down to charge the batteries.

WARNING!

- The battery chargers are selected precisely for specific applications and are particularly suited to the type, dimensions and chemical formulation of the specific batteries. For the safest and most efficient charging of your scooter's batteries, it is recommended to use the battery charger supplied with the scooter. Any charging method that charges the batteries individually is prohibited.
- Do not attempt to disassemble or remount the batteries on the scooter.
- If the battery charger has not been tested and approved for external use, do not expose it to adverse weather conditions. If the battery charger is exposed to adverse or extreme weather conditions, adaptation to the difference in environmental conditions must be allowed before using it in closed environments.

L) Driving mode selector

Press the "H/L" button to switch to alternate between the two functions:

- "H" is the fast driving mode
- "L" is the slow driving mode

M) USB port

11.2 Steering adjustment

The steering adjustment lever is positioned in front of the steering handle on the opposite side to the handbrake.

1. Pull the steering adjustment lever towards the steering handle to disengage it; **(Fig. 1)**
2. With the steering adjustment lever disengaged, continue to hold it and simultaneously push or pull the steering into the desired position; **(Fig. 1)**
3. Release the steering adjustment lever to lock the steering in position.

WARNING!

- The following situations may influence steering and stability during scooter functioning:
- attaching or tying a pet's lead to the steering
 - carrying passengers (including pets)
 - attaching anything to the steering
 - towing or being pushed by another vehicle

IMPORTANT: keep both hands on the steering and your feet well-positioned on the footrest at all times when operating your scooter. This driving position allows you to have maximum control of the vehicle

11.3 Adjustment of the handlebar

The handlebar adjustment lever is positioned in front of the steering handle on the same side as the handbrake.

1. Pull the handlebar adjustment lever towards the handlebar handle to disengage it; **(Fig. 2)**
2. With the handlebar adjustment lever disengaged, continue to hold it and simultaneously push or pull the handlebar into the desired position; **(Fig. 2)**
3. Release the handlebar adjustment lever to lock the handlebar in position.

(Fig. 1) Illustrative image

(Fig. 2) Illustrative image

11.4 Armrest adjustment

ADJUSTING THE SIZE OF THE ARMRESTS

(Fig. 3) Illustrative image

1. Loosen the adjustment knob on the back of the seat (Fig. 3);
2. Slide the armrest in or out to obtain the desired width;
3. Tighten the adjustment knob again.

WARNING!

Ensure that the adjustment knobs are firmly screwed onto the tube inserted in the armrest. Do not extract the armrest too far.

(Fig. 4) Illustrative image

TIPPING THE ARMRESTS

Raise the armrests to assist getting on and/or off the scooter.


SLOPING THE ARMRESTS

Turn the adjustment knob in an anti-clockwise direction to increase the slope of the armrests. Turn the adjustment knob in a clockwise direction to reduce the slope of the armrest (Fig. 4 - point 1).

ADJUSTING THE ARMREST HEIGHT

Turn the knob in an anti-clockwise direction (Fig. 4 - point 2), pull the knob outwards and move the armrest vertically to increase or reduce the armrest height. Release the knob and turn it in a clockwise direction once you have chosen the desired height.

11.5 Seat adjustment

(Fig. 5) Illustrative image

SEAT DEPTH LEVER

The depth of the seat can be adjusted. Pull the lever (Fig. 5 - point 1) and position the seat in the desired position. Release the lever.


SEAT ROTATION LEVER

The seat can be rotated through 360 ° and locked at every 45 °. Press the lever down and rotate the seat to the desired position (Fig. 5 - point 2). Release the lever to lock the seat in position.

ADJUSTMENT OF THE SEAT HEIGHT

(Fig. 6) Illustrative images

(Fig. 7)

(Fig. 8)

The seat can be adjusted using the seat post height adjustment.

1. Fold down the backrest by pulling up the lever on the left-hand side of the seat, then stand behind the scooter and, after pushing down the lever on the underside of the seat (right-hand side), use both hands to lift the seat off the seat post (Fig. 6);
2. Remove the bolt from the seat post rod, then adjust the seat post to the desired height, reinsert and tighten the bolt completely (Fig. 7);
3. Reinsert the seat in the seat post (Fig. 8).

WARNING!

Do not sit on the scooter or try to move it without first checking that the frame and the rear trolley are firmly inserted. The frame could accidentally detach from the trolley and cause injury or damage.

WARNING!

- Never sit on the scooter when the “free wheel” lever is inserted.
- Never disengage the stopping brake when you are using the scooter.

WARNING!

Always check that the stopping brake of the scooter is inserted before inserting the key in the ignition and turning it to the ON position.

11.6 Wheel locking/unlocking function

The scooter has a lever to lock/unlock the wheels.

(Fig. 9) Illustrative images


- To disengage the stop brake and unlock the scooter wheels, push the lever backwards towards the letter N (Fig. 9).
- To engage the brake and lock the scooter wheels, push the lever forwards towards the letter D (Fig. 9).

WARNING!

- If the scooter is pushed too fast with the wheels unlocked, the engine will act like a generator and make the scooter difficult to push;
- Do not place the scooter in unlocked wheels mode when you are on a slope;
- Never disengage the stopping brake when you are using the scooter;
- Always check that the stopping brake of the scooter is inserted before inserting the key in the ignition and turning it to the ON position.

12. STARTING THE SCOOTER FOR THE FIRST TIME

- Making sure to connect the battery cables to the motor cables as shown in the pictures;
- Ensure that the parking brake is in position D, otherwise as soon as the scooter is started it will produce an intermittent noise to warn of the problem;
- Adjust the seat, backrest and armrests according to your postural needs, as previously explained;
- Adjust the handlebars, as previously explained;
- Ensure that the rear view mirrors are properly aligned;
- Set the driving speed to minimum, for the first ride.

13. ASSEMBLY AND DISASSEMBLY

WARNING!

- Take care over any pinching points when attaching the rear frame to the front frame. During assembly, check the warning labels at the pinching points;
- Pinching and crushing risk! There is an increased risk of pinching or crushing, due to the weight of the components (such as the batteries), when preparing the scooter for transport or during maintenance activities. Always take extreme care over every operation. Use the help of another person if possible, particularly when storing the parts to be transported;
- Check that all the scooter components have been correctly assembled. After assembly, check that all locking devices are correctly inserted by raising the rear frame, to check that the transmission is locked in position;
- Incorrect assembly can cause injury! Check that all the scooter components have been correctly assembled. After assembly, check that all locking devices are correctly inserted;
- Check that all basic components are present.

14. DRIVING

14.1 Driving up and down

NOTE: For your first driving session, check that the scooter is being driven on a flat surface and will continue to be driven on a flat surface.

BEFORE YOU GET ON THE SCOOTER

- Check that it is switched off. This will avoid accidentally activating the paddle controls and causing injury to yourself or to others;
- Check that the scooter brake handle is in the engaged position;
- Check that the wheel unlocking lever is in the engaged position.

WARNING!

Keep your weight towards the centre of the platform when getting on/off the scooter. Shifting your weight towards the edge of the platform could cause instability!

GETTING ON THE SCOOTER

- Position the seat so that you can get on safely and easily;
- Return the handlebar to the vertical position;
- Carefully place one foot at the approximate centre of the platform and sit on the seat comfortably and safely;
- Fasten the safety belt, if the scooter has one;
- Lower or reposition the armrests in place;
- Insert the key in the ignition;
- Turn the key to the ignition position, as shown in **Fig. 10**. **Fig. 11** shows the key in the switched-off position.

(Fig. 10) Illustrative image

(Fig. 11) Illustrative image

GETTING OFF THE SCOOTER

- Check that the scooter is switched off and the key has been removed from the ignition;
- Return the handlebar to the vertical position;
- Lift or remove the armrests;
- Undo the safety belt;
- Carefully place one foot on the ground, transfer your weight onto that leg and slowly get up;
- Step away from the scooter.

14.2 Basic driving

- Check that you are seated safely and correctly on the scooter;
- Rotate the speed selector completely in an anti-clockwise direction to set the minimum speed;
- Insert the key in the ignition;
- Turn the key in a clockwise direction to the “On” position;
- Position your hand on the handles;
- To drive forwards, pull the right side of the accelerator control lever backwards (or push the left part of the accelerator control lever forwards);
- To drive in reverse, pull the left side of the accelerator control lever backwards (or push the right part of the accelerator control lever forwards);
- Pull the accelerator control lever to accelerate slightly;
- Release the accelerator control lever to allow the scooter to slow down gently and come to a complete stop;
- Practice these two basic functions until you feel you have control of the scooter.

14.3 Steering

- Position both hands on the handlebar handles and turn the handlebar right to move to the right;
- Turn the handlebar left to move to the left;
- Check that the area is sufficiently free of obstacles when turning the scooter, so that the rear wheels can pass any obstacle.

WARNING!

Turning the scooter sharply and at excessive speed could cause one of the rear wheels to leave the ground and the scooter to tip over. Avoid this risk by decelerating and steering in a wide arc around corners and obstacles.

STEERING NEAR A NARROW PASSAGE

To steer near a narrow passage, such as when entering or exiting a door or changing direction:

1. Bring the scooter to a complete stop;
2. Set the minimum speed;
3. Turn the handlebar in the direction in which you wish to go.

STEERING IN REVERSE

Be extremely careful when reversing.

1. Use your right hand to push the accelerator control lever forwards or use your left hand to pull the accelerator control lever backwards;
2. Turn the handlebar left to move to the left in reverse;
3. Turn the handlebar right to move to the right in reverse.

NOTE: The scooter speed in reverse is equal to 50% of the speed set on the speed selector.

14.4 Control through narrow passages

When you use a scooter to increase your mobility, you are bound to come across obstacles that require a certain experience in manoeuvring the scooter safely and easily. A list of common obstacles you may find in day-to-day use of your scooter is provided below. Several suggestions

on getting around them are listed with the obstacles. Learn them and follow the suggestions and you will find it surprisingly easy to control your scooter, manoeuvring it through doorways, up and down ramps, on and off pavements, onto grass and gravel and downhill.

RAMPS

When you go up a ramp, onto a pavement or a sloping surface:

- Lean forwards on the seat to move your own centre of gravity forwards and guarantee maximum stability and safety.

If there is a bump on the ramp, it is necessary to have a good grip on the road surface:

- Manoeuvre the scooter to allow the front wheels to make wide curves around the corners of the ramp.
- This will allow the rear wheels of the scooter to make a wide arc around the corner and steer clear of the obstacles.

If necessary, stop the scooter before driving onto the ramp:

- Before restarting it, press gently on the accelerator control lever.
- Accelerate gently after stopping on any sloping surface.

DRIVING OFF A RAMP

- Keep the scooter speed selector completely turned in a clockwise direction and set at the minimum speed.
- If it is necessary to stop, release the accelerator control lever slowly and constantly.

PAVEMENTS

- Never drive on/off a pavement that is higher than indicated in the technical specifications;
- Always drive carefully off a pavement;
- Approach the pavement in a manner that both rear wheels of the scooter go over the curb at the same moment;
- Do not drive off pavements slantwise, otherwise the scooter will tip over;
- Drive slowly down the curb to avoid jolts. Use as little power as possible.

GRASS AND GRAVEL

The scooter functions well on grass and gravel, and also on hills, but it is necessary to comply with the operating parameters indicated in this manual, in the technical specifications chapter. If any situation causes doubt, avoid it.

- Feel free to use the scooter on lawns and in parks;
- Avoid long or high grass, as it could wrap around the scooter axles;
- Avoid loose gravel.

14.5 Driving up and down slopes

DRIVING UP A SLOPE

- To guarantee maximum stability, lean forwards on the scooter seat when you go up ramps or hills, onto pavements or any low rise;
- Drive carefully when attempting to drive up any slope, including ramps for the disabled;
- Always remain perpendicular to the slope when getting onto/off a sloping surface;
- Never cut sideways across a sloping surface in any direction;
- Do not attempt to drive down a surface covered in snow, ice, freshly cut grass, leaves or other potentially hazardous materials;
- Do not reverse down a slope;
- Never drive the scooter up or down a slope that is higher than the recommended one. See the chapter "Technical specifications". Failure to do so could cause serious injury or death;
- Try to keep your scooter moving when driving up or down a slope. If you have to stop, start up again and accelerate slowly and carefully.

**WARNING!**

If the scooter starts moving more quickly than you consider safe on a downward slope, release the accelerator control lever and bring the scooter to a stop. As soon as you feel you have regained control of the scooter, push the accelerator control lever and continue driving down the rest of the slope with care.

GOING DOWN A SLOPE

- Set the minimum speed;
- If it is possible to do so safely, drive forwards down any ramp, low rise or sloping surface. The manufacturer advises against reversing on slopes, ramps, pavements and low rises. Reversing on a sloping surface can be extremely dangerous. However, if you do need to reverse, follow one of two procedures.

NOTE: When one of the following procedures is used to reverse down a sloping surface, the manufacturer strongly advises using the assistance of another person.

Procedure 1: In operation

- Turn the ignition to the "Off" position;
- Get off the scooter;
- Turn the ignition to the "On" position;
- Remain standing beside the scooter and carefully accelerate the controls at the minimum speed;
- Slowly and carefully accompany the scooter along the sloping surface;
- Carefully get back on the scooter and return to normal functioning.

Procedure 2: Emergency (not in operation)

- Turn the ignition to the "Off" position;
- Get off the scooter;
- Set the wheel unlocking lever to "Unlocked";
- Remain standing alongside the scooter and accompany it by hand along the entire sloping surface;
- As soon as you reach the flat surface at the bottom of the slope, position the wheel unlocking lever in the locked position.

**WARNING!**

When the scooter is in unlocked wheels mode, the stopping brake is released. The weight of the scooter could cause you to lose control on a slope. If you do not feel able to drive the scooter on a slope, ask for assistance or do not proceed.

14.6 Engine braking system

The scooter has a system that uses the engine to assist braking. The engine braking system is designed to function when the key is in both the "On" position and the "Off" position. When the key is in the "On" position, the wheel lock is in the locked position and the scooter is operating, the engine will contribute to slowing down the scooter as soon as you remove your hand from the accelerator lever.

When the key is in the "Off" position and the wheel lock is in the unlocked position, the engine braking system will prevent you from pushing the scooter too quickly (such as on a sloping surface) and you will notice this as you push it. The scooter will move freely until you reach a certain speed, then you will feel resistance, as the engine braking system has been activated.

14.7 Stopping brake

Your scooter has an automatic parking brake included in the electromechanical brake. The scooter stops when the engine is in drive and the ignition is switched off or when the ignition is switched on and the hand lever is in the neutral position. If the scooter is in free wheel mode (the engine is disengaged), it is possible to use the manual stopping brake function by placing the lever in the drive position.

14.8 Thermal protection

The controller on your scooter has a thermal restore safety system. An incorporated circuit controls the controller and engine temperature. If the controller and the engine overheat, the controller will remove power to allow the electrical components to cool down. Although your scooter returns to normal speed when the temperature returns to safe levels, it is advisable to wait five minutes before restarting the scooter, to allow complete cooling of all components.

15. TRANSPORT

- Disassemble or fold the seat and the handlebar as far as possible in the loading area of the vehicle used for transport;
- The choice to disassemble or lower the seat and the handlebar depends on the size and shape of the trunk of the vehicle used for transport.
- Do not lift the scooter using the plastic parts of the body or the handlebar. Breakage of these parts is not covered by the warranty.
- Do not lift the rear compartment using the tyres or the wheels. The compartment could rotate and cause injury or damage.
- IT IS advisable to use removable sheets or other types of covering to protect the scooter during transport.

WARNING!

- If the scooter and its components are not correctly and safely stored, they could move or be carried by the air and cause injury or damage;
- Do not sit on the scooter during transport. When transporting the scooter, ensure it is firmly secured with an approved anchoring/attachment system.

16. CLEANING

16.1 Cleaning the tyres

Clean the tyres with normal kitchen detergents and a damp cloth. Do not use solvents on the tyres. Solvents could damage or soften the tyre material.

16.2 Cleaning the body

- Clean the scooter body with a damp cloth. Do not wash the scooter with a hose. Dry with a soft and clean cloth;
- Use cold water mixed with a specific soap to remove dirt;
- Clean by hand with a soft cloth.

16.3 Cleaning the seat

Clean with a delicate soap or detergent and a damp cloth. A detergent for vinyl can also be used.

17. GENERAL MAINTENANCE INSTRUCTIONS


The scooter requires periodic maintenance. Incorrect maintenance of the scooter could lead to a higher number of technical problems, make it less flexible and invalidate the terms of the warranty. Preventive maintenance is important. The main operations to perform are listed below.

INSPECTION	DAILY	MONTHLY	EVERY SIX MONTHS	ANNUAL	PERFORMED BY
Brake functioning	X				User
Tyre conditions	X				User
Battery level check	X				User
Check on reversibility of the front wheels	X				User
Cleaning		X			User
Cable check		X			User
Battery terminals check			X		User
Engine brushes			X		Retailers
Frame stability check				X	Retailers
Wheel bearing lubrication				X	Retailers


Failure to perform maintenance on the brushes could invalidate the warranty on your scooter. To check and replace the engine brushes:

- Unscrew the engine brush cap;
- Remove the brushes;
- Check wear of the brushes;
- Replace the brushes if necessary.

Brush cap

New brush

Worn brush

17.1 Controls and reminders

- Keep the controller clean and protect it from water and rain. Never expose the device to direct contact with water.
- Keep the wheel clean of residue, hair, sand and carpet fibres.

- Inspect the tyre tread. If it is less than 1 mm (1/32”), have the tyres replaced by your local retailer.
- All the upholstery can be washed with warm water and neutral soap. Check the seat and backrest for rips and/or tears every now and then. Replace them if necessary. Do not store the product in damp environments, as this can cause mould to form and quickly damage the upholstery.
- Inspect and lubricate all moving parts. Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains. Always check that all nuts and bolts are tight.
- Check both the length of braking and the retention of the brake lock from a standstill.

17.2 Wheel replacement

- If the scooter has a solid tyre, replace the entire wheel.

WARNING!

The wheels must be replaced by an authorised retailer or a qualified technician in a workshop

17.3 Control panels, battery charger and rear electronics

- Protect these components from damp.
- If exposed to damp, dry them off completely before you use the scooter again.

17.4 Storing the scooter

If you do not intend to use the scooter for a lengthy period of time, it is advisable to:

- Completely charge the batteries before you store it;
- Disconnect the batteries from the scooter;
- Keep the scooter in a warm, dry place;
- Do not keep the scooter in places exposed to extreme temperatures;
- Operating conditions from -25°C to +50°C;
- Storage conditions from -40°C to +65°C;
- Batteries that are completely discharged, charged rarely, stored at extreme temperatures or stored without being completely charged may suffer permanent damage, resulting in their unreliability and limited durability. It is advisable to charge your scooter batteries periodically in the case of lengthy storage, to guarantee proper performance.

17.5 Battery and charging

Battery maintenance is the most important part of scooter maintenance. Keeping the batteries completely charged contributes to prolonging the battery duration. Follow the instructions below to help store the batteries in optimal conditions.

WARNING!

New batteries must be completely charged before using your scooter for the first time. Charge new batteries for 12 hours, even if the battery indicator already shows a full charge. Basic condition for maximising battery performance.

- Keep the batteries completely charged for daily use. It is advisable to connect the external battery charger after every use and charge for 6-8 hours;
- If the scooter will never be used for more than one week, charge the batteries completely and then disconnect them from the scooter.


CHARGING THE BATTERIES

List of charging guidelines to maximise battery duration:

1. Only use the battery charger supplied with your scooter;
2. NEVER use a car battery charger or a wet one;

3. Avoid deep discharges and never allow the batteries to discharge completely;
4. Do not leave the batteries at low charge for lengthy periods. Charge a discharged battery as soon as possible;
5. Charge the batteries completely on a regular basis;
6. Store the batteries completely charged;
7. Check the batteries once a month and charge them if necessary.

The battery charger is external. Follow the steps below to charge the batteries.

Illustrative image

- Position the scooter near a normal electrical socket on the wall.
- Remove the key to switch off the scooter.
- Rotate the battery charger compartment cover.
- Insert the XLR connector of the battery charger in the battery charger port.
- Insert the other end of the power cable in the wall socket.
- The LED on the battery charger will be YELLOW during charging.
- The LED on the battery charger will be GREEN when charging is complete.
- The battery capacity is displayed when charging is complete.
- Only disconnect the battery charger cable from the wall socket when the batteries are completely charged.

Only charge the batteries with the key in the OFF position (switched off).

BATTERY MAINTENANCE

Continuous cycle GEL or AGM and SLA batteries are used for the scooters.

- These batteries do not require maintenance;
- There is no risk of leakage or spills, so these batteries can be transported safely on aircraft, buses, trains, etc.;
- The battery lifetime should be prolonged if you follow the procedures indicated in this manual.

Note: Do not use car batteries on this scooter, as they are not designed to be totally discharged and are not safe for use on a scooter. A battery's lifetime often depends on its care.

IF THE BATTERY TERMINALS ARE CORRODED

- Corrosion can cause a poor electrical connection or operating problems;
- Clean corroded batteries with a stiff brush and a mixture of sodium bicarbonate and water.

17.6 Charging the battery

The battery charger functions at the normal voltage of a socket (alternating current) and converts it into V DC (direct current). The batteries use direct current to operate the scooter. When the batteries are completely charged, the amperage of the battery charger is almost at zero. This allows the battery charger to maintain the charge but not overcharge the battery.

Notes:

- The batteries cannot be charged if they have discharged almost to zero voltage.
- Always charge the batteries in a well-ventilated place.
- The battery charger is only for internal use. Protect it against humidity.
- To obtain maximum performance, it is recommended to replace both batteries simultaneously when they are discharged.
- All batteries slowly lose their energy charge if they are not used for long periods. After 3 months of non-use, a loss of charge of approximately 10% can occur. Therefore, in the case of an expected long period of non-use, it is recommended to fully recharge the batteries for at least 10 hours once a month and immediately afterwards disconnect them, so that they are not connected to the motor/control unit.

17.7 If the scooter does not work

- Check that the wheel unlocking lever is in the locked position;
- Check the automatic main switch. Reset the automatic switch if necessary;
- Check that the speed selector is set to the desired mode;
- Check that the ignition is in the "On" position.

If none of the above procedures solves the problem, contact an authorised retailer.

AUTOMATIC MAIN SWITCH

WARNING!

Do not attempt to carry out do-it-yourself electrical repairs. Consult your local retailer.

If the scooter should stop working, for no apparent reason, it could be that the automatic main switch has tripped.

POSSIBLE REASONS FOR THE AUTOMATIC MAIN SWITCH TRIPPING

- Driving on a steep incline;
- Overtaking on a pavement;
- Depleted batteries;
- Exceeding the maximum weight capacity.

When you switch on the scooter, the battery voltage decreases and the battery current increases to satisfy the needs of the engine or the other electrical devices on the scooter. This could result in heavy absorption of electric current that causes the automatic main switch to trip.

Solutions:

- Charge the scooter batteries. Consult chapter 16.5 (Battery and charging) of this manual;
- If the problem persists, contact your authorised retailer to carry out a battery charging test;
- If the batteries are working properly, the problem could be the battery charger. Contact your authorised retailer.

RESTARTING THE AUTOMATIC MAIN SWITCH.

- If the automatic main switch trips due to the batteries being depleted or due to a temporary overload, restart the automatic main switch;
- Wait 10 minutes for the engine control panel to start operating again;
- Check that the ignition is in the "Off" position;
- Press the restart button on the automatic main switch.

NOTE: If the automatic main switch continues to trip, there is probably a basic electrical fault that requires repair by qualified personnel.

18. PROBLEMS AND SOLUTIONS

The scooters have a controller that continually controls the operating conditions of your scooter. If it detects a problem, the ON / OFF button flashes. You must count the number of flashes and consult the list to check which type of error has been detected.

18.1 CN250 - CN260 problems and solutions

Number of flashes	Error	Notes
1 Flash	Low battery voltage	The battery needs charging or there is a battery connection fault. Check the battery connections. If the connections are correct, try charging the battery.
2 Flashes	Engine disconnected	There is an engine connection fault. Check all the connections between the engine and the controller.
3 Flashes	Problems on the engine circuits	The engine has a short circuit to a battery connection. Contact the assistance centre.
4 Flashes	Free wheel switch problem	The free wheel lever is activated or the manual brake disengage mechanism is enabled. Check the lever position.
5 Flashes	-	Not used.
6 Flashes	Battery charger connected	The S-Drive controller is disabled. This may occur because the battery charger is connected or the lever is not in the drive position.
7 Flashes	Accelerator problem	Indicates a problem on the command lever. Check that the lever is in the idle position before switching on the scooter.
8 Flashes	Problem in the controller	Indicates an error in the controller. Check that all the connections are correct.
9 Flashes	Stopping brake error	The stopping brake is not properly connected. Check the brake and the engine connections. Check that the controller connections are secure.
10 Flashes	High battery voltage	Excessive voltage has been applied to the controller, usually as a result of poor battery connections. Check the battery connections.

Note: If technical problems occur, it is recommended to check the device with a local retailer before trying to solve the problems yourself.

The following symptoms may indicate that your scooter has a serious problem. Contact your local retailer if any one of the following cases occurs:

1. Engine noise;
2. Worn cables and connections;
3. Bent or broken connectors;
4. Uneven wear on any one of the tyres;
5. Sudden movements;
6. The scooter pulls more on one side;
7. Bent or broken wheel units;
8. The scooter does not switch on;
9. The scooter switches on, but does not move.

19. CONDITIONS OF DISPOSAL 

19.1 General terms of disposal

Never dispose of the product as normal domestic waste. Dispose of the product at a sorted waste collection centre for recycling.

19.2 Instructions for correct disposal in accordance with European directive 2012/19/UE:

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service.

Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment.

19.3 Disposal of exhausted batteries - (Directive 2006/66/EC): 

This symbol on the product indicates that the batteries must not be considered as normal domestic waste. Ensuring that the batteries are discarded correctly contributes to preventing potentially negative consequences for health and the environment that would otherwise be caused by their inappropriate disposal. Recycling materials helps to conserve natural resources. Take depleted batteries to the collection points indicated for recycling. You can contact your local council, the local waste disposal service or the shop where you purchased the device for more detailed information on disposal of depleted batteries or the product.

20. TECHNICAL SPECIFICATIONS

20.1 CN250 technical specifications

CODE - MODEL	CN250
MAXIMUM CAPACITY	203 Kg
REAR WHEELS	350x100 mm
FRONT WHEELS	330x90 mm
ROLL-OVER PROTECTION WHEELS	Optional features
MAXIMUM SPEED	15 Km/h
BATTERY SPECIFICATIONS	12V 100Ah*2
AUTONOMY*	>45 Km
BATTERY CHARGER TYPE	8 Amp off-board, 220V 50Hz
CONTROLLER	PG S-200A
ENGINE TYPE	1650W
WEIGHT WITH BATTERIES	143.6 Kg
WEIGHT WITHOUT BATTERIES	100.6 Kg
STEERING RADIUS	1750 mm
SUSPENSION	Yes
LENGTH	1460 mm

WIDTH	700 mm
HEIGHT	1350 - 1410 mm
SEAT WIDTH	530 mm
SEAT HEIGHT FROM PLATFORM	485 mm
SEAT HEIGHT FROM THE GROUND	680
SEAT DEPTH	485 mm
BACKREST HEIGHT	490 mm
WHEELBASE	1055 mm
HEIGHT FROM THE GROUND	135 mm
MAXIMUM SLOPE ESTIMATED SURMOUNTABLE**	12° - 21%
SURMOUNTING OF AN OBSTACLE	120 mm

20.2 CN260 technical specifications

CODE - MODEL	CN260
MAXIMUM CAPACITY	203 Kg
REAR WHEELS	350x100 mm
FRONT WHEELS	350x100 mm
ROLL-OVER PROTECTION WHEELS	Optional features
MAXIMUM SPEED	15 Km/h
BATTERY SPECIFICATIONS	12V 120Ah*2
AUTONOMY*	>54 Km
BATTERY CHARGER TYPE	8 Amp off-board, 220V 50Hz
CONTROLLER	PG S-200A
ENGINE TYPE	1650W
WEIGHT WITH BATTERIES	156.3 Kg
WEIGHT WITHOUT BATTERIES	96.3 Kg
STEERING RADIUS	1660 mm
SUSPENSION	Yes
LENGTH	1600 mm
WIDTH	820 mm
HEIGHT WITH CANOPY	1700 mm
MIN. HEIGHT WITHOUT CANOPY	715 mm
MAX HEIGHT WITHOUT CANOPY	770 mm
SEAT WIDTH	530 mm

SEAT HEIGHT TO THE PLATFORM	475 mm
SEAT HEIGHT TO THE GROUND	720 mm
SEAT DEPTH	485 mm
BACKREST HEIGHT	500 mm
WHEELBASE	1130 mm
HEIGHT FROM THE GROUND	150 mm
MAXIMUM SLOPE ESTIMATED SURMOUNTABLE**	12° - 21%
SURMOUNTING OF AN OBSTACLE	120 mm

The scooter is tested according to standard ISO 7176-16:2012 for combustion resistance, but it is recommended to avoid naked flames close to the scooter and smoking when sitting on the scooter seat. The electrical system of this scooter satisfies standard ISO 7176-14:2008.

*** Scooter autonomy may vary in relation to:**

- Weight of the user;
- Slopes along the route;
- Battery wear;
- Driving style;
- Charging modes.

**** The term rated slope is defined in standard EN 12184:2014**

21. WARRANTY

Moretti products are guaranteed for 2 (two) years from the date of sale against material and manufacturing defects, subject to the following limitations. The warranty is voided by improper use, abuse, modifications to the product and failure to follow the instructions. The intended use of the product is given in the user manual. Moretti is not liable for damage, injury or any other consequences resulting from installation or use which are not scrupulously conforming with the instructions given in the installation, assembly and user manual. Moretti does not guarantee its products against damage or defects in the following circumstances: natural disasters, unauthorised repair or maintenance, improper electric power supply (as applicable), use of parts or components not supplied by Moretti, failure to follow the guidelines and instructions for use, tampering, shipping damage (other than the original shipping by Moretti), or failure to perform maintenance as indicated in the manual. Components subject to wear and tear are not covered by this warranty if the damage is caused by normal use of the product.

21.1 Warranty on rechargeable batteries (if provided)

The original and replacement batteries are covered by a ninety-day warranty on performance and a six-month warranty on manufacturing defects or as required by law. If completely charged batteries are left unused for more than three consecutive months, the warranty is automatically invalidated. If completely discharged batteries are left unused for more than three consecutive days, the warranty is automatically invalidated.

22. REPAIRS

Repairs under warranty If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti reseller or its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

Repairing a product not covered by the warranty Product not covered by warranty may be returned for repair only if authorised in advance by Moretti customer service. The costs of labour and shipping incurred by repairs not covered by the warranty are borne by the client or reseller in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for 6 (six) months from the day of reception of the repaired product.

Non-defective products The client will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the client at his expense.

23. REPLACEMENT PARTS

Use only spare parts and accessories listed in the Moretti general catalogue. Moretti original spare parts are guaranteed for 6 (six) months from the day of delivery.

24. NON-LIABILITY CLAUSE

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti makes no declaration, guarantee or condition, express or implicit, including any future declaration, guarantee or condition of sale, suitability for a given purpose, non-violation and non-interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit guarantee under the law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit guarantee or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice.

WARRANTY CERTIFICATE

ENGLISH

Product _____

Date of purchase _____

Reseller _____

Street _____ **Town** _____

Sold to _____

Street _____ **Town** _____

Via Bruxelles, 3 - Melegnano 20122 Cavriago (Arezzo) - ITALY - Tel. +39 055 96 21 11
www.morettispa.com email: info@morettispa.com
MADE IN P.R.C.

** Please consult our website for the latest available version of the user manual*

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