ITALIANO

MA CMB40-CM850-CMB60 01 A*_ITA_10-2023

**CARROZZINE ELETTRICHE
PIEGHEVOLI MOBILITY 840
MOBILITY 850 E MOBILITY 860**

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INDICE

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ENGLISH

MA_CM840-CM850-CM860_01 A*_ENG_10-2023

**MOBILITY 840. MOBILITY 850
AND MOBILITY 860 FOLDING
ELECTRIC WHEELCHAIRS**

INSTRUCTION MANUAL

CONTENTS

1. CODES	PAG.3
2. INTRODUCTION	PAG.3
3. INTENDED USE	PAG.3
4. EU DECLARATION OF CONFORMITY	PAG.3
5. GENERAL PRECAUTIONS	PAG.4
6. SYMBOLS USED	PAG.4
8. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY	PAG.4
9. SAFETY	PAG.8
10. WARNINGS ON USE	PAG.9
10.1 General warnings.....	pag. 9
10.2 Changes	pag. 9
10.3 Before every use	pag. 10
10.4 Weight limits	pag. 10
11. GENERAL VIEW OF THE PARTS	PAG.10
12. PLACING THE WHEELCHAIR IN OPERATION	PAG.11
12.1 Installing the controller and the connection cables	pag. 11
12.3 Switching on and using for the first time	pag. 12
13. OPERATION	PAG.13
13.4 Recommendations for safe use	pag. 18
13.5 Other information for use	pag. 18
13.6 Free wheel - neutral.....	pag. 19
14. BATTERY & CHARGING	PAG.19
14.1 Attention.....	pag. 20
14.2 Indications for battery charging and health	pag. 20
15. REPLACING THE WHEELS	PAG.21
16. MAINTENANCE AND REPAIR OF THE ELECTRIC WHEELCHAIR	PAG.21
16.1 Daily checks.....	pag. 21
16.2 Weekly checks.....	pag. 22
17. CARE AND MAINTENANCE	PAG.22
17.1 Cleaning and disinfection.....	pag. 22
17.2 Battery terminal connections.....	pag. 23
17.3 Cables	pag. 23
17.4 Axle supports and motor/differential	pag. 23
17.5 Control panels, battery charger and rear electronics.....	pag. 23
17.6 Storing the electric wheelchair	pag. 23
18. TROUBLESHOOTING	PAG.23
18.1 Self-diagnosis.....	pag. 24
19. DISPOSAL CONDITIONS	PAG.24
19.1 General conditions of disposal	pag. 24
19.2 Instructions for correct disposal in accordance with European Directive 2012/19/EU:.....	pag. 24
19.3 Disposal of depleted batteries - (Directive 2006/66/EC):	pag. 25
20. REPLACEMENT PARTS AND ACCESSORIES	PAG.25
21. TECHNICAL SPECIFICATIONS	PAG.25
22. WARRANTY	PAG.27
22.1 Warranty on rechargeable batteries (if provided).....	pag. 27
22.2 Repairs	pag. 27
23. SPARE PARTS	PAG.28
24. NON-LIABILITY CLAUSES	PAG.28

CE Class I medical device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 on medical devices

1. CODES

CM840 Mobility 840 folding electric wheelchair

CM850 Mobility 850 folding electric wheelchair

CM860 Mobility 860 folding electric wheelchair

2. INTRODUCTION

Thank you for choosing an ARDEA MOBILITY electric wheelchair by Moretti. Moretti electric wheelchairs have been designed and manufactured to satisfy all your needs for practical, correct and safe use. This manual contains useful suggestions for using your device properly and important recommendations for your safety. Please read this manual carefully before you use the electric wheelchair. Should you have any queries, please contact your retailer for suitable advice and assistance.

NB Check to ensure that no parts of the wheelchair have been damaged during shipment. Do not use the product if it is damaged and contact the retailer for further instructions.

3. INTENDED USE

The electric wheelchair is intended for people with mobility impairments

WARNING!

- Do not use the product for any purpose other than the one specified in this manual.
- The electric wheelchair must always be used by people of sound body and mind and not under the influence of drugs or alcohol.
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than the one specified in this manual.
- The manufacturer reserves the right to modify the device and its relevant manual without prior notice, in order to improve the product's features.

4. EU DECLARATION OF CONFORMITY

MORETTI SpA declares, under its own exclusive responsibility, that the products manufactured and sold by MORETTI SpA in the ELECTRIC WHEELCHAIR product family conform with the applicable provisions of Regulation (EU) 2017/745 of 5 April 2017 on MEDICAL DEVICES.

To this end, MORETTI SpA guarantees and declares as follows, under its own exclusive responsibility:

1. The devices in question satisfy the general safety and performance requirements set out in Annex I to Regulation (EU) 2017/745, as prescribed by Annex IV to said regulation.
2. The devices in question ARE NOT MEASUREMENT TOOLS.
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATIONS.
4. The devices in question are sold in NON-STERILE packing.
5. The devices in question should be considered as belonging to class I, in accordance with the rules set out in Annex VIII to said regulation.
6. MORETTI SpA maintains and provides to the competent authorities, for at least ten years from the date of manufacture of the last production lot, the technical documentation proving conformity with Regulation (EU) 2017/745.

Note: The complete product codes, the manufacturer's single registration number (SRN), the basic UDI-DI code and any references to standards used are indicated in the EU Declaration of Conformity that MORETTI SPA issues and provides through its own channels.

5. GENERAL PRECAUTIONS

- Please consult this manual carefully for correct use of the device
- The wheelchair is suitable for use both outdoors and indoors, in hospitals, centres for the elderly, at home or in similar circumstances.
- Keep the packaged product away from sources of heat, as the packaging is cardboard
- Avoid exposure of the device to direct sunlight or other sources of heat
- The lifetime of the device depends on wear and tear of non-repairable and/or non-replaceable parts and is, in any case, no longer than 10 years
- Always take care when children are present
- Always consult your doctor or therapist for correct use of the device.
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.
- Use the wheelchair after opening it and always have it used by one person only.

6. SYMBOLS USED

	Product code
	Device unique identifier
	Serial number
	CE marking
	Manufacturer
	Production lot
	Read the instruction manual
	Medical Device
	Maximum supported weight
	Disposal conditions
	Warning
	Production date
	Disposal of the product according to Directive EC/19/2012
IPX4	Protection against dust and water spray
	Type BF applied part

8. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY

The rapid growth of electronics, especially in communications, has saturated the environment with electromagnetic radiation (EM) emitted by television, radio and communication signals. These EM waves are invisible and are stronger the closer you are to their source. All electrical conductors act as antennas for EM signals and, to varying extents, all electric wheelchairs and scooters may be affected by electromagnetic interference (EMI). EMI may cause unintentional movements and/or irregular control of the vehicle.

The electromagnetic energy emitted by sources such as:

- radio stations;
- TV stations;
- amateur radio stations;
- shop alarm systems;
- mobile telephones and cordless telephones can interfere with electric wheelchairs and electric scooters.

The interference can cause the electromagnetic brake to release, unintentional movement of the device and movement in the wrong direction. It can also permanently damage the device's electronic control unit. In our everyday environment, there are a certain number of intense sources of EMI.

The intensity of the EM energy can be measured in Volts per metre (V/m). Your electric wheelchair is equipped with a protection against EMI which has been tested and certified as compliant with current international regulations, up to a certain intensity. This is called the "level of immunity" The greater the level of immunity, the better the device is protected. At this time, current technology can provide a level of immunity of at least 20 V/m; this is sufficient protection against the most common sources of EMI. However, we consider that the following precautions can further minimise the risk.

- A) Do not use portable radio transmitters, CB radios and avoid switching on mobile telephones or other personal communications equipment when the electric wheelchair is switched on.
- B) Do not approach powerful radio and TV transmitters too closely.
- C) If an unintended movement or unintended release of the brake occurs, turn the key to the OFF position as soon as you can do so safely.
- D) Installing additional electrical equipment or accessories and modifying the motor of your electric wheelchair may make it more susceptible to EMI and hence make its use hazardous. Note that doing so also voids the warranty.
- E) In any case, report any unintentional movement of the vehicle to your retailer or to MORETTI SPA, note it down in this manual and note down the fact that there is a source of electromagnetic radiation in your vicinity, if any.

8.1 Electromagnetic compatibility

Table 1

Manufacturer's guide and declaration - electromagnetic emissions		
The device is intended to be used in the following electromagnetic conditions. The client or user of the device must ensure that it is used in these conditions.		
Emission tests	Conformity	Electromagnetic environment - guide
RF Emissions IEC60601	Group 1	The device uses RF energy only for internal functioning. Its RF emissions are therefore very low and in all likelihood do not interfere with nearby electronic equipment.
RF Emissions IEC60061	Class B	The device is intended to be used in all buildings, including the home and buildings connected directly to the low voltage public power grid for residential applications.
Harmonic emissions IEC 60601	Applicable	
Voltage flutter/flicker emissions IEC 60601	Conforming	

Table 2

Manufacturer's guide and declaration - electromagnetic emissions			
The device is intended to be used in the following electromagnetic conditions. The client or user of the device must ensure that it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
Electrostatic discharge (ESD) ISO 7176 IE C60601	+ 6 kV contact + 8 kV air	+ 6 kV contact + 8 kV air	Information for using the device and charging the battery. The floor must be wood, concrete or ceramic. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst immunity ISO 7176 IE C60601	+2KV for power lines +1 kV for I/O lines	2KV for power lines +1 kV for I/O lines	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery
Over voltages ISO7176 IEC60601	1 kV differential mode +2 KV common mode	1 kV differential mode +2 KV common mode	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery
Voltage drop outs, brief interruptions and variations of voltage on the power supply lines. ISO 7176 IE C60601	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery.
Mains frequency magnetic field (50/60 Hz) ISO 7176 IE C60601	30 A/m	30 A/m	Magnetic fields at the mains frequency should have the typical characteristics of a commercial or hospital environment. (Applies both to using the device and charging the battery).
Note: UT is the AC mains voltage prior to application of the test level.			

Table 3

Manufacturer's guide and declaration - electromagnetic immunity			
The device is intended to be used in the following electromagnetic conditions. The client or user of the device must ensure that it is used in these conditions. Make sure it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
conducted RF ISO7176 IEC60601	3 V eff. From 150 KHz to 80 MHz	3 V eff.	Portable and mobile RF communications devices should not be used near any part, including the cables, unless the recommended separation distances are observed, calculated from the equation applicable to the transmitter frequency. Recommended distances $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ From 80 MHz to 800 Mhz $d = 2.3\sqrt{P}$ from 800 MHz to 1.0 GHz $d = 0.2 \sqrt{P}$ 26 MHz to 800 MHz $d = 0.4 \sqrt{P}$ 800 MHz to 2.5 GHz during charging operations where "P" indicates the maximum rated output of the transmitter in Watt (W) according to the transmitter manufacturer and "d" is the recommended separation distance in metres (m). The intensity of the field generated by fixed RF transmitters, as determined in an on-site electromagnetic inspection, might be below the level of conformity in each frequency band. Interference may be experienced in the proximity of equipment marked with this symbol: 
radiated RF ISO7176 IEC60601	3 V/m From 80 MHz to 1.0 GHz 20 V/m From 26 MHz to 2.5 GHz	3 V/m 20 V/m	

ENGLISH

NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a The intensity of the field for fixed transmitters, such as base stations for radio-telephones (mobile and cordless) and terrestrial mobile radios, amateur radio equipment, AM and FM radio transmitters and TV transmitters cannot be known precisely in theory. An on-site electromagnetic inspection should be performed to establish an electromagnetic environment caused by fixed RF transmitters. If the field intensity measured in the place where Nefti is used exceeds the applicable RF level of conformity indicated above, normal functioning of the device must be monitored.

If abnormalities are observed, supplementary measures may be necessary, such as re-orientation or repositioning.

^b The field intensity over the frequency range from 150 kHz to 80 MHz should be less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile radio communications equipment					
The device is intended to operate in an electromagnetic environment in which RF disturbances are under control. The client or user can contribute to preventing EM interference by ensuring a minimum distance between the mobile and portable RF communications devices (transmitters) and the battery charger during charging and the entire device during its operation as recommended below, in relation to the maximum output power of the radio communication equipment					
Tests	Separation distance as a function of the frequency of the transmitter m				
Maximum rated output power of transmitter W		During battery charging	During battery charging	During use	During use
	from 150 kHz to 80 MHz d = 1.2√P	from 80MHz to 800MH d = 1.2√P	from 800 MHz to 2.5 GHz d = 2.3 √P	from 26 MHz to 800MHz d= 0.2 √P	from 800 kHz to 2.5 GHz d = 0.4√P
0.01	0.12	0.12	0.23	0.2	0.04
0.1	0.38	0.38	0.73	0.06	0.13
1	1.2	1.2	2.3	0.2	0.4
10	3.8	3.8	7.3	0.63	1.26
100	12	12	23	2	4
For transmitters for which the maximum rated output power is not listed above, the recommended separation distance d, in metres (m), can be calculated using the equation applicable to the transmitter frequency, where P is the maximum rated output power of the transmitter, in Watt (W), as given by the transmitter's manufacturer.					
NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.					
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection from structures, objects and people.					

9. SAFETY

- Always use the safety belt and keep your feet resting on the footplate at all times.
- Do not use the electric wheelchair under the influence of alcohol.
- Do not use electronic radio transmitters such as walkie-talkies or mobile telephones.
- Make sure there are no obstacles nearby when you are driving the electric wheelchair.
- Do not make a sharp turn or brake suddenly when you are driving the electric wheelchair.
- Do not drive the electric wheelchair in traffic.
- Do not attempt to mount curbs higher than those indicated in the technical specifications.
- Do not put your hands or legs outside the electric wheelchair while driving.

- Do not drive the electric wheelchair on slippery roads or when it is snowing.
- The device is intended for use by one person at a time
- Do not overload the wheelchair above the maximum limit (see technical characteristics)

WARNING!

The wheelchair must only be transported in a car when in the folded position and placed in the trunk. The wheelchair is not designed for use inside moving vehicles.

10. WARNINGS ON USE

10.1 General warnings

WARNING!

Do not use the electric wheelchair on public highways. Please be aware that it may be difficult for drivers to see you when you are sitting on the electric wheelchair. Always use pedestrian routes. Comply with the rules of pedestrian traffic. Wait for the way to be clear, then proceed with extreme caution.

1. Before getting into or out of your electric wheelchair, always check that the device is switched off and that the battery charger is disconnected, to avoid injury to yourself or to others.
2. Always check that the drive wheels are in drive mode before driving. Do not switch off the electric wheelchair when it is still moving, as it would stop very suddenly.
3. Do not use this product or any other optional accessory available without first reading and understanding this manual completely. If you do not understand the warnings, the precautions or the instructions, contact the retailer or a healthcare operator before attempting to use the device, to avoid injury or damage.
4. In certain cases, such as certain medical conditions, the user of the electric wheelchair will need to practice using it in the presence of a qualified assistant, who may be a relative or a professional qualified to assist the user in use of the electric wheelchair during various everyday activities.
5. Avoid lifting or moving the electric wheelchair by gripping any one of its removable parts, including the armrests, seat or fairing. This could cause personal injury and/or damage the wheelchair.
6. Avoid using the electric wheelchair beyond its limits, as described in this manual.
7. Do not sit on the electric wheelchair when it is inside a moving vehicle.
8. Keep your hands away from the wheels while driving the electric wheelchair. Please remember that loose clothing could get caught in the wheels.
9. Consult your doctor if you are taking drugs or if you have specific physical limitations. Certain drugs and certain limitations could affect the ability to use the electric wheelchair safely.
10. Always check that the wheel locking/unlocking lever is locked or unlocked.
11. Do not remove the roll-over protection wheels if they are present on the electric wheelchair.
12. Contact with tools could cause an electric shock; do not connect an electrical extension cord to the AC / DC converter or to the battery charger.
13. The maximum slope at which the wheelchair can operate is 8°/14° and great care should be taken when turning on a slope.
14. Do not drive on a slope over the limit of the electric wheelchair.
15. Avoid reversing off steps, pavements or other obstacles. This could cause the electric wheelchair to fall or tip over.
16. Always reduce speed and maintain a stable centre of gravity on sharp bends. Do not use the maximum power on sharp bends.
17. Using the electric wheelchair in the rain, snow or foggy conditions and on icy or slippery surfaces could damage the electrical wiring.
18. Never sit in the electric wheelchair when it is being used with any type of lift or lifting device.

The electric wheelchair is not designed for this use.

The manufacturer is not liable for any damage or injury caused by improper use.

19. The environmental conditions suitable for use of the electric wheelchair are a temperature between 5-40°C and humidity between 10%-80%.

10.2 Changes

The electric wheelchairs in the ARDEA MOBILITY range by Moretti have been designed and manufactured to satisfy all your needs for practical, correct and safe use. Do not alter, add, remove or disable any part or function of the electric wheelchair under any circumstances. This could cause damage to the electric wheelchair and/or personal injury.

1. The control parameter must only be altered by authorised technicians, for safety reasons.
2. It is recommended to carry out a safety check before every use of the electric wheelchair, in order to ensure that it is functioning safely.

10.3 Before every use

1. Check all the electrical connections and make sure they are properly connected and not corroded.
2. Check all the cables and make sure they are connected properly.
3. Check the brakes.

10.4 Weight limits

1. Check the information on maximum capacity in the table of the technical specifications in this manual or on the label on the device. The power of the electric wheelchair is assessed according to the capacity.
2. Comply with the weight limits indicated for your electric wheelchair. Exceeding the maximum capacity invalidates the warranty. The manufacturer cannot be held liable for injury or damage to property caused by failure to comply with the weight limit.
3. Do not carry passengers on the electric wheelchair, as this could alter the centre of gravity and cause tipping over or falls.

11. GENERAL VIEW OF THE PARTS

CM840

1. Main structure with motor
2. Controller
3. Controller attachment bar
4. Footplate
5. Motor
6. Padded seat
7. Backrest
8. Front wheel
9. Rear wheel
10. CPU
11. Batteries

CM850

1. Main structure with motor
2. Controller
3. Controller attachment bar
4. Footplate
5. Motor
6. Mudguard
7. Padded seat
8. Backrest
9. Front wheel
10. Rear wheel
11. Shock absorbers
12. CPU
13. Batteries

CM860

1. Main structure with motor
2. Controller
3. Controller attachment bar
4. Footplate
5. Motor
6. Mudguard
7. Padded seat
8. Backrest
9. Front wheel
10. Rear wheel
11. Shock absorbers
12. CPU
13. Batteries

12. PLACING THE WHEELCHAIR IN OPERATION

Before starting to prepare the wheelchair, carefully unpack it and position it on the ground in the closed position, in a manner that all four wheels are on the ground. Then open it by gently pulling the backrest upwards.

12.1 Installing the controller and the connection cables

The controller of your electric wheelchair can be positioned on either the left or the right, so **decide which hand the user will use** for the controller and assemble it on that side. Insert the support bar into the specific guide positioned underneath the armrest on the chosen side. Check that the distance of the drive joystick is the best one for use.

- Loosen the knob

- Slide the controller bar into the desired position.
- Tighten the knob

Connect the controller to the power cable and check that the cables of the motor and of the control unit underneath the structure are connected correctly (Fig. 1)


(Fig. 1)


12.2 Installation and operation of the batteries

CM840

- To remove the battery, grip the handle and pull firmly upwards. In order to reposition it, first insert the rear part in the compartment and then press firmly on the front part to reinsert it (Fig. 2-3)


(Fig. 2)

(Fig. 3)

CM850

- To remove the battery, rotate the lever to the left to unlock the locking system and then lift the top part of the battery (Fig. 4);
- Press the red button  on the battery to check the charge status (Fig. 5);
- The batter has a general on/off switch, which must be set to the on position during periods of use. To reinsert the battery, first position the rear part and then the top part, pressing down until it clicks into place. (Fig. 6)

(Fig. 4)

(Fig. 5)

(Fig. 6)

CM860

- To remove the battery, lift the lock button, grip the black plastic head of the battery and carefully extract it from the compartment (Fig. 7); to reinsert it, lift the lock button again and gently push the battery completely into its compartment (Fig. 8);

illustrative smage


- Check the battery charge status by pressing the blue button on the black head (Fig. 9).

(Fig. 7)

(Fig. 8)

(Fig. 9)

12.3 Switching on and using for the first time

Press the On button to switch on the wheelchair and set the driving speed to the minimum, then try moving the wheelchair in a large and empty space. The speed can be increased once the user is familiar with the wheelchair

WARNING!

- Before switching on the power supply, always be aware of your surroundings when selecting the desired speed.
- It is recommended to select the slowest speed for use indoors.
- It is recommended to select a speed that allows you to control the wheelchair safely for use outdoors.

13. OPERATION

The electric wheelchair is simple to use. However, it is recommended to read the following instructions carefully, to familiarise yourself with your new vehicle. The following steps are necessary to ensure that the wheelchair operates safely. To guarantee safe and easy use, both armrests on the electronic wheelchairs are liftable:

CM840-CM850

Rotate the lever upwards to unlock the mechanism and lift the armrest (Fig. 10); lower the armrest completely and rotate the lever downwards to reposition it (Fig. 11).

(Fig. 10)

(Fig. 11)


CM860

Lift the lock button upwards to unlock the armrest, then lift it upwards; lower the armrest completely and press lightly until you hear a click to reposition the armrest and lock it in place (Fig. 12).


illustrative smage

(Fig. 12)


(Fig. 13)

(Fig. 14)

On the CM860 model, it is possible to adjust the depth of the footrest by up to 10cm, for a total of three steps; to do so, use both hands to lift the lock buttons on the side of the footrest and pull forwards or backwards, as needed (Fig. 13). It is also possible to adjust the backrest inclination by up to around 30° from the erect position; to adjust it, press the lever behind the backrest and pull it backwards or upwards, as needed.

The CM860 electric wheelchair comes equipped with an armrest widening kit, to increase the space between the armrests; it can be applied on only one or on both sides, as needed.


(Fig. 15)

Unscrew the two side screws of the armrest shown in Fig. 16 and remove the armrest, leaving the black plastic shims in position.

(Fig. 16)

(Ref. Fig. 15)

(Fig. 17)

(Fig. 18)

- Apply shim A for the top part and shim B for the bottom part. (Fig. 17)
- Insert the armrest from the top hole. (Fig. 18)

(Ref. Fig. 15)

(Fig. 19)

(Fig. 20)

- Tighten the pieces again, using screws A for the top part (Fig. 19) and screws B for the bottom part (Fig. 20)

13.1 Controller a Joystick

The controller has a multi-purpose visual display that can provide four types of information:

Controller

1. On button
2. Off button
3. Battery charge level bar
4. Horn
5. Increase speed button
6. Reduce speed button
7. Speed level bar

The joystick controls the direction and speed of the wheelchair. Operate the wheelchair by pressing gently on the joystick in the direction you intend to take. Returning the joystick to the neutral position (central) reduces the speed and stops the wheelchair by automatically applying the electromechanical brakes.

1. ON / OFF status

When the wheelchair is switched on, the controller LED is lit up. When the LED is not lit up, the controller is switched off.

Notes:

- Press the ON / OFF button, then wait two seconds before using the joystick. This is a safety feature that prevents a sudden start-up.
- Gentle use of the joystick results in smoother movements in speed and direction, whereas brusque use of the joystick results in sudden movements in speed and direction.
- In an emergency, release the joystick. The wheelchair will stop immediately.

2. Battery status indicator

The battery status indicator is formed of seven segments (one red, four yellow and three green). It allows you to monitor the battery charge. The battery indicator indicates the approximate quantity of the remaining charge.

- The red, yellow and green LEDs lit up together indicate that the battery is completely charged.
- The red and yellow LEDs lit up together indicate that the battery needs charging as soon as possible.
- The red LEDs lit up on their own indicate that the battery needs charging as soon as possible, as it is almost flat. This could cause a loss of power in the wheelchair.

13.2 Closing and moving the electric wheelchair

The CM840-850-860 electric wheelchairs can be folded for easier transport and storage:

WARNING!

- Never transport the wheelchair with a fastening system, as it is not compatible with these systems.
- Never sit on the wheelchair when it is being transported in a moving vehicle.
- Do not place the folded wheelchair on the front seat with the driver.

CM840

- To close the wheelchair, press the red lever (Fig. 21) underneath the backrest downwards and gently push the backrest forwards until it is completely folded in the forward position. The wheelchair has wheels for transport, but it is not stable in the erect position, so it is not advisable to leave it in this position without supports or locks (Fig. 22).
- It is recommended to use the backrest handle for transporting the wheelchair when closed, pivoting on the rear wheels.
- To reopen the wheelchair, rest it in the forward position with all four wheels on the ground and lift the backrest until it is fully extended.

(Fig. 21)

(Fig. 22)

CM850

- To close the wheelchair, press the black lever (Fig. 23) downwards with respect to the backrest and gently push the backrest forwards until it is completely folded in the forward position, then grip the footplate and lift the wheelchair until the special shims are resting on the ground (Fig. 24). The wheelchair must be in D (drive) mode to do this;
- However, it is possible to transport the folded wheelchair in M (manual) mode, using the footplate as a handle, although care should be taken when lifting from the ground and repositioning;
- To reopen the wheelchair, set D (drive) mode again, rest it in the forward position with all four wheels on the ground and lift the backrest until it is fully extended.


(Fig. 23)

(Fig. 24)

CM860

- To close the wheelchair, press the red lever (Fig. 25) underneath the backrest downwards and gently push the backrest forwards until it is completely folded in the forward position. The wheelchair has wheels for transport, but it is not stable in the erect position, so it is not advisable to leave it in this position without supports or locks (Fig. 26).
- It is recommended to use the backrest handle for transporting the wheelchair when closed, pivoting on the rear wheels.
- To reopen the wheelchair, rest it in the forward position with all four wheels on the ground and lift the backrest until it is fully extended.

(Fig. 25)

(Fig. 26)

IT is also possible to remove both wheel-motor blocks in order to further reduce the weight and dimensions:

(Fig. 27)

(Fig. 28)

(Fig. 29)

- Disconnect the motors from the batteries. (Fig. 27)
- Press the metal head to remove the safety pin. (Fig. 28)
- Slide the block outwards to remove it. (Fig. 29)
- Perform the procedure in reverse to reassemble it.

13.3 Using the electric wheelchair

Before you get into your electric wheelchair, make sure you are on a flat surface with no obstacles. Although your wheelchair is designed for use on slopes, it is safer to practise using it on a flat surface.

Several checks to make before every use are indicated below:

- Check functioning of the horn
- Check that the wheels are in good condition
- Check the battery charge status of your wheelchair
- The motor unlocking levers must be in position D (DRIVE).
- It is recommended to select a low speed for the first few times you use it.

Follow these steps for use of the electric wheelchair:

1. Sit on your wheelchair and press the ON/OFF button. All the LED indicator lights of the battery charger must be lit up with a fixed light.
2. Rest your forearm firmly on the armrest and place your hand on the controller. After switching on the wheelchair, wait two seconds before activating the joystick. This is a safety feature to prevent sudden start-ups.
3. The joystick central lever controls all movements of the wheelchair. If you move the joystick forwards, the wheelchair moves forwards. If you move the joystick backwards, the wheelchair moves backwards. If you move the joystick to the left, the wheelchair moves to the left. If you move the joystick to the right, the wheelchair moves to the right. When you release the joystick in the vertical or central position, the automatic brake is activated and the wheelchair stops automatically. The joystick also operates as an accelerator. The further you push the lever, the faster your wheelchair will move. The speed depends on the position of the speed control button.
4. Move the joystick central lever gently, applying stable and controlled pressure.
5. In an emergency, release the joystick central lever and the wheelchair will stop automatically.
6. Never drive at a speed beyond your ability to control the wheelchair.
7. Practise driving in a place where there are no obstacles. Start at the lowest speed and move forwards and backwards. Do this several times. When you feel secure, increase the speed gradually by pressing the speed button.
8. The number of LEDs lit up on the speed control indicator identifies the level of speed set.
9. If the battery charge indicator only has two LEDs lit up, the batteries must be charged as soon as possible.
10. When you have finished travelling or before putting the wheelchair away, always switch off the controller by pressing the ON/OFF button.

13.4 Recommendations for safe use

- Before changing direction from forwards to backwards or backwards to forwards, bring the control lever into the central position and wait for the wheelchair to come to a complete stop.
- Only use the wheelchair where it would also be safe to walk.
- If you need to go up or down a slope, check that the gradient does not exceed 8°/14% and that the surface has a good grip and is not slippery.
- Use a slow speed to reverse, go up a slope, down a slope, turn a corner or on uneven surfaces.
- NEVER use the leg rests to lift the wheelchair. The leg rests are movable and therefore cannot be used for lifting.

13.5 Other information for use

DRIVING UP A SLOPE

If you need to use the wheelchair on the maximum gradient, you may need to travel at a higher speed. Press the speed button to increase the speed. Be extremely careful when travelling on steep slopes.

DRIVING DOWN A SLOPE


Proceed slowly down steep descents. Press the speed button to reduce the speed. This allows better control of driving. However, your wheelchair will not accelerate on its own down a slope, due to the automatic braking system, which starts functioning if you attempt to drive too fast. In general, you should drive in a linear manner on slopes, without sudden changes in direction, and never stop the wheelchair on a slope that exceeds the safety limit.

GOING OVER OBSTACLES

Go over obstacles slowly and in a forwards direction. Do not attempt to go over obstacles from the side. This will allow you to go over obstacles up to 5 cm. Do not go over obstacles higher than 5 cm. Your wheelchair has a roll-over protection system. These devices may interfere when you are ascending onto or descending from a sidewalk. Take maximum care when going over these obstacles.

13.6 Free wheel - neutral

The electromechanical brakes are activated when the wheelchair is not in use or is switched off. The wheelchair has a "free wheel" (Manual) function that is activated by turning the levers from the wheel in drive position (Drive) to the free wheel position. Before switching on for use, check that both red levers are in the wheel in drive position. (Fig. 20)

(Fig. 20)

**WARNING!**

- Never use the free wheel position on a slope. Never use the free wheel position when driving the wheelchair.
- Always remember to place the lever in “wheel in drive” position before switching on the wheelchair

13.7 Electromechanical brakes

The electric wheelchair has electromechanical brakes. This is an automatic magnetic disc, also known as a safety brake. The electromechanical brakes are automatic and are engaged when the electric wheelchair is switched on but stationary, even when the wheelchair is on a slope. The electromechanical brakes are also engaged when the electric wheelchair is switched off, but the motor levers are in the drive position (vertical).

Note: Please consult the section on maintenance and repair to ensure that the brakes are in good condition.

14. BATTERY & CHARGING

The electric wheelchair comes supplied with two 24 V 6AH lithium batteries and they can be charged by connecting the battery charger directly on the controller, or charged singly by connecting the power cable directly onto the battery, in the specific port. In the first case, the electric wheelchair must be switched off. The charging time is three hours for each individual battery and six hours for charging both batteries from the controller. Once the battery is fully charged, the red light of the battery will turn green.

Note: The battery lifetime is also influenced by the weight that the wheelchair must carry

14.1 Attention

1. Read all the instructions and warnings before using the battery charger.
2. Use the battery charger in a well-ventilated area.
3. To avoid injury, only charge lead-acid or gel batteries.
4. Disconnect the power after charging.

NB: only use the battery charger provided by the supplier of the electric wheelchair. Use of any other type of battery charger can be dangerous and requires the manufacturer's approval.

- Always charge the batteries in a well-ventilated place.
- The battery charger is only for internal use. Protect against the damp.
- For better performance, it is recommended to replace both batteries at the same time if they are flat
- If the wheelchair is not used for a lengthy period of time, it is advisable to charge the batteries at least once a month to avoid their deterioration.
- Any method of individual charging of the batteries is prohibited.

14.2 Indications for battery charging and health

- Never charge each battery for more than three hours.
- Switch off the main unit during charging.
- Do not leave unsupervised children near the wheelchair when it is charging.
- Do not attempt to open the battery compartment.
- The red lights indicates that the unit is switched on and charging; green indicates that the unit is fully charged.
- When the battery is fully charged, remove the cable, then extract the charging connector and the power socket.
- If the wheelchair is not used for more than two months, it is recommended to charge the battery completely before use.

- The wheelchair has an overcurrent protection system that cuts off the current from the batteries if there is a sudden increase in voltage, switching off the system and guaranteeing that the motor and the battery are not damaged.
- The life cycle of a lithium battery is measured in charging cycles, with an average of 500 charging cycles.
- Check the battery light every day if the wheelchair is used frequently.
- Do not use battery chargers other than the one provided
- If the battery is damaged, do not handle it without gloves and protective goggles.
- If you need to dispose of the battery, find the nearest hazardous material disposal centre and call to obtain instructions.

Perform the following battery maintenance procedures before storing the wheelchair for lengthy periods of time:

- Completely charge the battery before storing it.
- Remove the battery from the wheelchair.
- Store the battery in dry place at room temperature.
- Avoid extreme temperatures and do not allow the battery to freeze.
- During lengthy periods of non-use, charge the batteries once a month to stop them from being damaged

WARNING!
 Do not charge the individual batteries directly from the charger for more than three hours

WARNING!
 Do not attempt to charge the battery if it is frozen due to extreme weather conditions. Leave it to defrost for several hours.

WARNING!
 Disconnect the battery charger power cable immediately when the light turns green.

WARNING!
 Do not charge both batteries for more than six consecutive hours from the controller and for longer than three hours if only one is inserted.

15. REPLACING THE WHEELS

WARNING!
 The wheels must be replaced by an authorised retailer or a qualified technician in a workshop.

16. MAINTENANCE AND REPAIR OF THE ELECTRIC WHEELCHAIR

The ARDEA MOBILITY range of devices by Moretti are CE marked and carefully checked before being placed on the market. For the safety of the patient, your device should be checked at least once every year by the manufacturer or an authorised laboratory to verify its suitability for use. Use only original replacement parts and accessories for repairs.

MAINTENANCE	Daily	Weekly	Monthly	Annual
ELECTRICAL SYSTEM				
Battery indicator - Check the battery indicator to see whether the battery needs charging	•			
Controller / Display - Check for damage and exposed cables			•	
Check that the sockets and connection cables are not loose			•	

Check that the batteries are completely charged before daily use	•			
Check that all bolts are tight	•			
Check that the lighting system is working properly (if present)	•			
WHEELS				
The front and rear wheels must rotate smoothly without any obstruction		•		
The front and rear wheels must rotate without wobbling			•	
Visually inspect the tread of the tyre. If it is less than 1 mm (1/32"), replace it with a new tyre (available from your local retailer)			•	

16.1 Daily checks

1. With the controller switched off, check the joystick. Check to ensure that it is not bent or damaged and that it returns to the centre on being released. Visually inspect the rubber sleeve at the base of the joystick to check for any damage. Do not touch it or try to repair it.
2. Visually inspect the controller cables. Check to ensure that there are no damaged or exposed wires.
3. Check the battery status indicator on the controller to see whether the batteries need to be charged.

16.2 Weekly checks

1. Disconnect the controller and inspect the connections to check for any signs of corrosion
2. Ensure that all parts of the controller are securely attached to the electric wheelchair. Do not over-tighten the screws.
3. Check the brakes. This test must be carried out on a flat surface with at least one metre of space around your electric wheelchair

To check the brakes (the electric wheelchair should move slightly during performance of this test):

- Switch on the controller and lower the speed and the regulation knob
- After one second, check that the battery level indicator remains on
- Slowly push the joystick forwards until you hear the parking brakes go click, then release the joystick immediately. You should hear each parking brake operate within a few seconds of releasing the joystick.
- Repeat the check on the brake for the backwards, left and right positions of the joystick

CHECKS & REMINDERS:

- Keep the controller clean and protect it from rain or water. Never expose the electric wheelchair to direct contact with water.
- Keep the wheel clean of residue, hair, sand and carpet fibres.
- All the upholstery can be washed with warm water and neutral soap. Check the seat and backrest for rips and/or tears every now and then. Replace them if necessary.
Do not store your electric wheelchair in damp environments, as this can cause mould to form and quickly damage the upholstery.
- Inspect and lubricate all moving parts. Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains. Always check that all nuts and bolts are tight.

PERIODIC CHECKS

1. Keep the controller clean and protect it from rain or water. Never expose the wheelchair to direct contact with water.
2. Keep the wheels free of fluff, hair, sand and textile fibres.

3. Visually inspect the tread of the wheels. If it is less than 1/32", replace them.
4. All the upholstery can be washed with warm water and neutral soap. Check the seat and backrest for rips and/or tears every now and then, and replace them if necessary.
Do not store your wheelchair in damp environments, as this can cause mould to form and quickly damage the upholstery.
5. Inspect and lubricate all moving parts. Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains. Always check that all nuts and bolts are tight.

17. CARE AND MAINTENANCE

Your electric wheelchair requires a minimum amount of care and maintenance. The following parts require checking and/or care and maintenance.

17.1 Cleaning and disinfection

- Use a damp cloth and non-abrasive neutral detergent to clean the plastic and metal parts of the electric wheelchair. Avoid using products that could scratch the surface of your electric wheelchair.
- DO NOT immerse the device in water for any reason!
- If necessary, clean the product with an approved disinfectant. Make sure the disinfectant is safe for the product before you use it.
- Read all the safety instructions for correct use of the disinfectant and/or detergent before using it on your product. Failure to do so may result in irritation of the skin or premature damage to the upholstery and/or finishes of the electric wheelchair.
- Never use acid or alkali based products or solvents (such as acetone or diluent) to clean it.

17.2 Battery terminal connections

- Make sure the connectors are tightly connected and not corroded.
- The batteries must be correctly housed in their compartments

17.3 Cables

- Regularly check the electrical connections.
- Regularly check the cable insulation, including the battery charger cable, for wear or damage.
- Have any damaged connector, connection or insulation repaired or replaced by an authorised retailer before using the electric wheelchair.

17.4 Axle supports and motor/differential

These components are pre-lubricated and sealed and do not require further lubrication.

REPLACING A WHEEL

- If the electric wheelchair has tyred wheels and one of the tyres is flat, replace the inner tube.
- If the electric wheelchair has a solid rubber tyre, replace the entire wheel. Contact an authorised retailer for details on spare wheels for your electric wheelchair.

17.5 Control panels, battery charger and rear electronics

- Protect these components from damp.
- If exposed to damp, leave them to dry completely before you use the electric wheelchair again.

17.6 Storing the electric wheelchair

If you do not plan to use the electric wheelchair for a lengthy period of time, it is advisable to:

- Completely charge the batteries before you store it.
- Disconnect the batteries from the electric wheelchair.
- Keep the electric wheelchair in a warm and dry place.
- Do not store the electric wheelchair in places where it could be exposed to extreme temperatures.

- Operating conditions 5-40°
- Completely charge the batteries around every two months of non-use

Batteries that are completely flat, charged rarely, stored at extreme temperatures or stored without being completely charged may suffer permanent damage, resulting in their unreliability and limited durability. It is advisable to charge the batteries of your electric wheelchair periodically in the case of lengthy storage, to guarantee proper performance.

18. TROUBLESHOOTING

Note: If technical problems occur, it is recommended to check the device with a local retailer before trying to solve the problems yourself.






The following symptoms may indicate that your electric wheelchair has a serious problem. Contact your local retailer if any one of the following cases occurs:

1. Motor noise
2. Damaged safety belts
3. Bent or broken connectors
4. Uneven wear on any one of the wheels
5. Sudden movements
6. The electric wheelchair pulls more to one side
7. Bent or broken wheel units
8. The wheelchair does not switch on
9. The wheelchair switches on, but does not move

Problem	Possible causes	Possible solutions
The wheelchair does not switch on	The batteries are not correctly connected to the CPU. The CPU is not correctly connected to the joystick. Completely flat batteries. Batteries in Off mode (CM850)	Position and connect the batteries correctly. Check all the connections between the joystick, the CPU and the batteries. Charge the batteries. Check that the batteries are in On mode (CM850)
The wheelchair is making strange noises or vibrating during steering	The speed is too low and is causing excess friction. A motor might be damaged.	Increase the speed Replace the motors
The batteries do not charge	Damaged battery charger Damaged battery The charging time ends before the battery is completely charged.	Replace the battery charger Replace the batteries Always make sure to charge the batteries completely in a single charging session
The controller cables are not connected	The connector pins (male) on the controller are bent or misaligned with the connector holes (female) on the CPU connector.	Request assistance for repair or replacement of the damaged part

18.1 Self-diagnosis

- The controller lights also act as a diagnostic tool and flash in sequences that indicate specific problems.
- If the warning light indicates a fault, switch off the wheelchair and switch it back on again. If the fault has not been eliminated after switching it back on again, the warning light repeats the error message.
- Consult the following table to identify potential problems and the relative solutions
- Please contact the retailer if the method described above does not help to solve the problem

 <p>1</p>	<p>1st light flashing: Insufficient battery charge</p>
 <p>2</p>	<p>2nd light flashing: no signal from the joystick</p>
 <p>3</p>	<p>3rd light flashing: batteries charging</p>
 <p>4</p>	<p>4th light flashing: the joystick is moving during switching on, or the motor cables have not been correctly connected</p>
 <p>5</p>	<p>5th light flashing: the red levers of the motor are not both in Drive mode, or are damaged</p>

19. DISPOSAL CONDITIONS

19.1 General conditions of disposal 

Never dispose of the device as normal domestic waste. Dispose of the device at a municipal sort-ed waste collection centre for recycling.

19.2 Instructions for correct disposal in accordance with European Directive 2012/19/EU:

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling of the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment separately.

19.3 Disposal of depleted batteries - (Directive 2006/66/EC): 

This symbol on the product indicates that the batteries must not be considered as normal domestic waste. Ensuring that the batteries are discarded correctly contributes to preventing potentially negative consequences for health and the environment that would otherwise be caused by their inappropriate disposal. Recycling materials helps to conserve natural resources. Take depleted batteries to the collection points indicated for recycling. You can contact your local council, the local waste disposal service or the shop where you purchased the device for more detailed information on disposal of depleted batteries or the product.

20. REPLACEMENT PARTS AND ACCESSORIES

Use only replacement parts and accessories listed in the Moretti general catalogue.

21. TECHNICAL SPECIFICATIONS

Code	CM840
Material	Aluminium
PERFORMANCE	
Maximum speed	6 km/h
Steering radius	800 mm
Maximum safe slope	8°/14%
Max. step	40 mm
Autonomy	20 km
Charging time	6h
TECHNICAL SPECIFICATIONS	
Battery specifications	Two 24V 6AH lithium batteries
Battery charger type	5A
Controller electronics	Brushless-24V/30A
Motor	Two 180W brushless motors
Brake	Magnetic
Drive	Two motors on rear wheels
Type A	EN12184:2022

DIMENSIONS* AND WEIGHTS

Closed: Length: 33 cm Width 60 cm Height: 91 cm

91	97	44	97	72	56 front 52 back	44	45	48	60	27	29.5	20	18	-	120 kg

* Measurements in cm

Code	CM850
Material	Aluminium
PERFORMANCE	
Maximum speed	6 km/h
Steering radius	800 mm
Maximum safe slope	8°/14%
Max. step	40 mm
Autonomy	20 km
Charging time	6h
TECHNICAL SPECIFICATIONS	
Battery specifications	Two 24V 6AH lithium batteries
Battery charger type	5A
Controller electronics	Brushless-24V/30A
Motor	Two 250W brushless motors
Brake	Magnetic
Drive	Two motors on rear wheels
Type A	EN12184:2022

DIMENSIONS* AND WEIGHTS

Closed: Length: 41 cm Width: 65 cm Height: 91 cm

100	80	44	97	72	54 front 50 back	46	45	57	65	28.5 kg	31 kg	32	18	-	150 kg

* Measurements in cm

22.2 Repairs

- Repairs under warranty

If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client on whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti retailer or on its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

- Repairing a product not covered by the warranty

A product not covered by warranty may be returned for repair only if authorised in advance by the Moretti customer service. The costs of labour and shipping incurred for repairs not covered by the warranty are borne by the customer or retailer in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for six months from the day on which the repaired product is received.

- Non-defective products

The customer will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the customer at their own expense.

23. SPARE PARTS

Moretti original replacement parts are guaranteed for six months from the day of their delivery.

24. NON-LIABILITY CLAUSES

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti provides no other declaration, guarantee or condition, express or implicit, including any future declarations, guarantees or conditions of sale, suitability for a given purpose, non violation and non interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit warranty provided by law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit warranty or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice

WARRANTY CERTIFICATE

Product _____

Date of purchase _____

Reseller _____

Street _____ **Town** _____

Sold to _____

Street _____ **Town** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) - ITALY - Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

** Consult our website for the most recent version of the user manual available*

ENGLISH

WARRANTY CERTIFICATE

Product _____

Date of purchase _____

Reseller _____

Street _____ **Town** _____

Sold to _____

Street _____ **Town** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) - ITALY - Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com

MADE IN P.R.C.

** Consult our website for the most recent version of the user manual available*

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