ITALIANO

MA LAD47x-LAD48x-LAD49x Piuma Up 02 A\_05 2021

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CODICE PRODOTTO	MATERIALE SINGOLO ELEMENTO	MATERIALE BASE MATERASSO	MATERIALE COPERTA
LAD470	NYLON/PVC	NYLON/PVC	NYLON/PU imbottita
LAD475	NYLON/PVC	NYLON/PVC	NYLON/PU imbottita
LAD480	NYLON/PVC	NYLON/PVC	NYLON/PU imbottita
LAD471	NYLON/TPU	NYLON/PVC	NYLON/PU imbottita
LAD476	NYLON/TPU	NYLON/PVC	NYLON/PU imbottita
LAD481	NYLON/TPU	NYLON/PVC	NYLON/PU imbottita
LAD491	NYLON/TPU	NYLON/PVC	NYLON/PU imbottita

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**LEVITAS®**

MA LAD47x-LAD48x-LAD49x Piuma Up 02 A\_05 2021

**PIUMA UP SERIES OF ALTERNATING PRESSURE  
MATTRESS AND COMPRESSOR SYSTEMS**

---

**INSTRUCTION MANUAL**

ENGLISH

## INDEX

<b>1. CODES</b> .....	<b>PAG.3</b>
<b>2. INTRODUCTION</b> .....	<b>PAG.3</b>
<b>3. INTENDED USE</b> .....	<b>PAG.3</b>
<b>4. EU DECLARATION OF CONFORMITY</b> .....	<b>PAG.4</b>
<b>5. APPLICABLE REGULATIONS AND DIRECTIVES</b> .....	<b>PAG.4</b>
<b>6. GENERAL PRECAUTIONS</b> .....	<b>PAG.4</b>
6.1 Precautions concerning electromagnetic compatibility.....	pag. 4
<b>7. SYMBOLS USED</b> .....	<b>PAG.5</b>
<b>8. CONTENTS OF THE KIT</b> .....	<b>PAG.5</b>
<b>9 GENERAL DESCRIPTION</b> .....	<b>PAG.6</b>
9.1 Mattress and compressor system LAD470, LAD471, LAD475, LAD476.....	pag. 6
9.1.1 Mattress.....	pag. 6
9.2 Compressor operations.....	pag. 7
9.2.1 Pressure regulation knob.....	pag. 7
9.2.2 Normal pressure indicator.....	pag. 7
9.2.3 Low pressure indicator.....	pag. 7
9.2.4 Alternate/static switch.....	pag. 7
9.2.5 Alarm deactivation.....	pag. 7
9.2.6 CPR valve.....	pag. 8
9.3 Mattress and compressor system LAD480, LAD481, LAD491.....	pag. 8
9.3.1 Mattress.....	pag. 8
9.3.2 Compressor and control panel.....	pag. 8
9.4 Compressor operations.....	pag. 8
9.4.1 Pressure regulation levels.....	pag. 8
9.4.2 Treatment.....	pag. 9
9.4.3 Alarm deactivation.....	pag. 9
<b>10. INSTALLATION</b> .....	<b>PAG.10</b>
10.1 Installation of the compressor and the mattress.....	pag. 10
10.2. Only for the digital compressor model.....	pag. 11
<b>11. BEFORE EVERY USE</b> .....	<b>PAG.11</b>
<b>12. WARNINGS ON USE</b> .....	<b>PAG.11</b>
<b>13. CLEANING</b> .....	<b>PAG.12</b>
<b>14. STORAGE</b> .....	<b>PAG.13</b>
<b>15. MAINTENANCE</b> .....	<b>PAG.13</b>
15.1 General.....	pag. 13
15.2 Replacement of the air filter.....	pag. 14
15.3 Element replacement.....	pag. 14
15.4 Fuse replacement.....	pag. 14
<b>16. CONDITIONS OF DISPOSAL</b> .....	<b>PAG.14</b>
16.1 General disposal conditions.....	pag. 14
16.2 Instructions for correct disposal in accordance with European directive 2012/19/UE:.....	pag. 14
<b>17. REPLACEMENT PARTS AND ACCESSORIES</b> .....	<b>PAG.14</b>
<b>18. TROUBLESHOOTING</b> .....	<b>PAG.15</b>
<b>19. DECLARATION ELECTROMAGNETIC COMPATIBILITY</b> .....	<b>PAG.15</b>
19.1 Manufacturer's guide and declaration - electromagnetic emissions.....	pag. 15
19.2 Manufacturer's guide and declaration - electromagnetic immunity.....	pag. 16
<b>20. TECHNICAL SPECIFICATIONS</b> .....	<b>PAG.18</b>
20.1 Compressor.....	pag. 18
20.2 Mattress.....	pag. 19
<b>21. WARRANTY</b> .....	<b>PAG.19</b>
22.1 Repairs under warranty.....	pag. 19
22.2 Repairing a product not covered by the warranty.....	pag. 19
22.3 Non-defective products.....	pag. 19
<b>23. SPARE PARTS</b> .....	<b>PAG.20</b>
<b>24. NON-LIABILITY CLAUSE</b> .....	<b>PAG.20</b>

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL  
of 5 April 2017 on medical devices

### 1. CODES

**LAD470** Piuma Up3 mattress and compressor system - PVC mattress - height 12.7 cm

**LAD471** Piuma Up3 mattress and compressor system - TPU mattress - height 12.7 cm

**LAD475** Piuma Up3 mattress and compressor system - PVC mattress - height 20.3 cm

**LAD476** Piuma Up3 mattress and compressor system - TPU mattress - height 20.3 cm

**LAD480** Piuma Up4 mattress and compressor system - PVC mattress - height 20.3 cm

**LAD481** Piuma Up4 mattress and digital compressor system - TPU mattress - height 20.3 cm

**LAD491** Piuma Up4 mattress and digital pressure compressor system - TPU mattress - height 20.3 cm - length 208 cm

### 2. INTRODUCTION

Thank you for choosing the Piuma UP mattress and compressor system in the LEVITAS range by Moretti. Moretti mattress and compressor kits have been designed and manufactured to satisfy all your needs for practical, correct and safe use. This manual contains useful suggestions for using your device properly and safely. Please read this entire manual carefully before you use the mattress and compressor kit. Should you have any queries, please contact your dealer for suitable advice and assistance.

#### NOTE, CAUTION AND WARNING!

 **NOTE:** Indicate any recommendations.

**CAUTION:** indicate the correct operating or maintenance procedures, in order to avoid damage to the device, the equipment or other items.

**WARNING:** draw attention to a potential hazard that requires correct procedures or practices to avoid injury.

### 3. INTENDED USE

The mattress and compressor kit is designed to prevent and treat the bedsores that can occur in any situation involving lengthy hospitalisation and/or being laid up at home.

#### WARNING!

- Do not use the product for any purpose other than that specified in this manual
- The device must be installed by personnel qualified to perform general nursing procedures and who have received appropriate training in the prevention and treatment of bedsores
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than that specified in this manual
- The manufacturer reserves the right to modify the device and this manual without prior notice, in order to improve its features

#### 4. EU DECLARATION OF CONFORMITY

MORETTI SpA declares, under its own exclusive responsibility, that the products made and sold by MORETTI SpA in the PIUMA UP ANTI-DECUBITUS product family conform with the applicable provisions of Regulation (EU) 2017/745 of 5 April 2017 on MEDICAL DEVICES.

To this end, MORETTI SpA guarantees and declares as follows, under its own exclusive responsibility:

1. The devices in question satisfy the general safety and performance requirements set out in Annex I to Regulation (EU) 2017/745, as prescribed by Annex IV to said regulation.
2. The devices in question ARE NOT MEASUREMENT INSTRUMENTS.
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATIONS.
4. The devices in question are sold in NON STERILE packing.
5. The devices in question should be considered as belonging to class I, in accordance with the rules set out in Annex VIII to said regulation.
6. MORETTI SpA maintains and provides to the competent authorities, for at least ten years from the date of manufacture of the last production lot, the technical documentation proving conformity with Regulation (EU) 2017/745.

**⚠ NOTE:** The complete product codes, the manufacturer's single registration number (SRN), the basic UDI-DI code and any references to standards used are indicated in the EU Declaration of Conformity that MORETTI SPA issues and provides through its own channels.

#### 5. APPLICABLE REGULATIONS AND DIRECTIVES

This device has been tested and approved according to the following directives and standards:

EN 60601-1:2005+A1:2012  
EN 60601-1-2:2007

#### 6. GENERAL PRECAUTIONS

- Please consult this manual carefully for correct use of the device
- Always consult your doctor or therapist for correct use of the device
- Keep the packaged product away from sources of heat, as the packaging is cardboard
- The lifetime of the device depends on wear and tear of non-repairable and/or non-replaceable parts
- Always take care when children are present
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

##### 6.1 Precautions concerning electromagnetic compatibility

This device has been tested and found to comply with the limits for medical devices according to standard EN 60601-1-2:2007.

These limits are established to provide reasonable protection against harmful interferences in typical medical applications. This device generates, uses and may irradiate energy in radio frequency and, if not installed and used in accordance with the instructions, may cause harmful interferences for other nearby devices. However, there are no total guarantees against electromagnetic interferences occurring in certain circumstances. If this device should cause harmful interferences to other devices, which could be caused by switching it on or off, it is recommended that the user take the following steps:

- Redirect and reposition the receiving device
- Increase the distance separating the devices
- Connect the device to a different socket or circuit to the other connected devices
- Consult the manufacturer or an authorised assistance centre

**7. SYMBOLS USED**

-  **REF** Product code
-  **UDI** Unique device identification
-  **SN** Serial number
-  Date of manufacture
-  **MD** Medical Device
-  **CE** CE Mark (EU conformity)
-  Read the instruction manual
-  **LF** LF type
- IP21** Protection against solid foreign bodies of 12.5 mm or higher  
Protection against vertically falling water drops
-  Disposal of the product according to directive EC/19/2012
-  Warning
-  Insulation Class II
-  Manufacturer
-  Conditions of disposal
-  Maximum washing temperature 60°C; normal process
-  Do not bleach
-  Do not iron
-  Do not dry clean
-  Tumble drying possible; low temperature

**8. CONTENTS OF THE KIT**

This product is designed for:

- Assistance in reducing the occurrence of bedsores, while simultaneously optimising the patient's comfort
- Long-term home care of patients suffering from bedsores
- Pain management based on the doctor's prescriptions

**⚠ NOTE:** This device is not suitable for use in the presence of flammable anaesthetic mixtures with air or with pure oxygen or nitrous oxide.

**LAD470** Piuma Up3 mattress and compressor system - PVC mattress - height 12.7 cm

**LAD471** Piuma Up3 mattress and compressor system - TPU mattress - height 12.7 cm

**LAD475** Piuma Up3 mattress and compressor system - PVC mattress - height 20.3 cm

**LAD476** Piuma Up3 mattress and compressor system - TPU mattress - height 20.3 cm

**LAD480** Piuma Up4 mattress and compressor system - PVC mattress - height 20.3 cm


**LAD481** Piuma Up4 mattress and digital pressure compressor system - TPU mattress - height 20.3 cm

**LAD491** Piuma Up4 mattress and digital pressure compressor system - TPU mattress - height 20.3 cm - length 208 cm


## 9 GENERAL DESCRIPTION

### 9.1 Mattress and compressor system LAD470, LAD471, LAD475, LAD476

#### 9.1.1 Mattress


#### 9.1.2 Compressor

Control panel:

1. pressure regulation knob
2. on/off button
3. treatment button
4. clasps

Control panel:

1. quick connection
2. clasps
3. fuses
4. air filter
5. power cable

## 9.2 Compressor operations

### 9.2.1 Pressure regulation knob

This controls the air outlet pressure. Turn clockwise to increase the pressure. Turn anticlockwise to reduce the pressure. Higher pressures will support heavier patients. When the mattress has reached the desired pressure level, the normal pressure LED lights up.

 **NOTE:** It is possible to check whether the pressure is suited to the patient by running your hand underneath the air cells at the level of the patient's buttocks. The air cells will alternatively inflate and deflate. You should feel a slight contact with the buttocks when the cells underneath the buttocks deflate.

 **NOTE:** Each time the mattress is inflated, it is recommended to set the pressure knob to "Max" to speed up inflation. The air mattress can be inflated to the desired degree of firmness to ensure the level of comfort of the mattress.

### 9.2.2 Normal pressure indicator

When the green LED lights up, the pressure inside the air mattress has reached the desired pressure setting (pre-set on the pressure regulation knob).

### 9.2.3 Low pressure indicator

When the mattress is in the start-up phase, the low pressure LED remains lit until the appropriate pressure is reached (as a function of the pressure selected on the knob). In any other case, the low pressure LED is a warning, indicating that the pressure in the mattress is unusually low. In this case, check that all the connections are correctly installed, according to the installation instructions.

 **NOTE:** if the pressure level is always low, check for possible leaks (components or air tubes). If necessary, replace any damaged tubes or contact a local retailer qualified to carry out repairs.

### 9.2.4 Alternate/static switch

The control for the alternate/static function allows you to choose between the alternating pressure and the static pressure mode. With the alternating pressure mode, the air cells that alternate are partially deflated and inflated, avoiding prolonged pressure on every single point underneath the patient and thus preventing bedsores. With the static pressure mode, all the air cells are equally inflated. Mode used for treating the patient and for positioning or repositioning on the bed.

### 9.2.5 Alarm deactivation

When low pressure occurs, the LED lights up and the alarm buzzer sounds. Press the button to deactivate the acoustic alarm and check for any air leaks.


 **NOTE:** Press the  "button" again to reactivate the acoustic alarm.

### 9.2.6 CPR valve

Pull the CPR valve strap to rapidly deflate the mattress.


Ensure that the CPR valve connector is correctly connected after use.

Remove the tube connections on the pump side to deflate the mattress more quickly


## 9.3 Mattress and compressor system LAD480, LAD481, LAD491

### 9.3.1 Mattress

### 9.3.2 Compressor and control panel

## 9.4 Compressor operations

### 9.4.1 Pressure regulation levels

The pressure regulation levels control the air pressure output. When they are raised, the output pressure increases. Lower them to reduce the air pressure. Heavier patients will need a higher output pressure. Identify the correct pressure setting as follows: Start by inflating to the maximum pressure level, have the patient lie down on the mattress, then gradually decrease the pressure until the patient feels comfortable, without touching the bottom. At this point, the carer inserts their hand palm upwards between the area under the patient's buttocks and the

cell. If the palm can touch the patient's buttocks, increase the pressure until the skin is around 1.5 cm from the palm. Always leave at least three centimetres of space between the area of the patient's buttocks and the underlying air cells, in order to avoid them touching the bottom.

#### 9.4.2 Treatment

##### A. Max Inflation

Press the switch to allow the compressor to reach its maximum operating pressure. Once the maximum pressure is reached, the compressor will automatically pass into alternate mode. This function can also be used to inflate the mattress completely when getting the patient onto/off the bed, to obtain better support. The alternate function is automatically restored after twenty minutes.

##### B. Alternate

In Alternate Treatment Mode, the mattress system will alternate inflation of the different elements. The four treatment cycles to choose from are 10-15-20-25 minutes.

##### C. Static

This interrupts the alternate mode. In static mode, the cell pressure level will be reduced with respect to the same pressure level in alternate mode. This special feature serves to protect patients at high risk of bedsores from additional pressure, which could cause the bedsores.

##### D. Seat Inflate

This function offers additional support to a seated patient. This additional feature can be selected in static or alternate mode. If the Seat Inflate function is selected with the Static mode, this function will last for twenty minutes, after which it will automatically return to the Alternate and Seat Inflate function.

##### E. Control panel lock/unlock

Press the lock/unlock button to lock the control panel, so that users cannot choose any function on the panel. Press the button again for around three seconds to unlock the panel.

#### 9.4.3 Alarm deactivation

##### A. PFA (power failure alarm)

In the case of a power failure, the Power Failure LED lights up and a buzzer sounds. Press the deactivate button to deactivate both the buzzer and the LED.

##### B. Low pressure indicator

When the low pressure LED lights up, the pressure inside the air mattress is lower than normal. Refer to the problem solving section.

##### C. Maintenance

This element switches on when there is a mechanical failure. Contact technical assistance for repairs.


## 10. INSTALLATION

Unpack the product and check the device for any damage that could have occurred during shipment. Contact the retailer immediately in the case of damage.

Contents of the box:

- 1 mattress (might not be included if only the compressor unit is purchased)
- 1 Compressor
- 1 User manual

### 10.1 Installation of the compressor and the mattress


1. Position the mattress or padding on the bed structure. Check for the indication of the feet end.


**WARNING:** The mattress cover must be positioned on the mattress underneath

2. Hang the compressor on the bed frame (feet end) and adjust the brackets to obtain the best erect position of the compressor, or position the compressor on a flat surface.


3. Connect the air tube connectors from the mattress to the air unit of the compressor. When you hear or note a "click", the connection is complete and secure.

**NOTE:** Check and ensure that the air tubes are not bent or trapped underneath the mattress.

4. Connect the plug to an electric socket.

- ⚠ **NOTE:** 1. Ensure that the compressor unit is suitable for the local voltage. Do not position the device in a manner that makes it difficult to use the circuit breaker.  
2. The plug is also used to disconnect the device.

**WARNING:** The pump provided with the mattress can only be used for mattresses recommended by the manufacturer. Do not use it for any other purpose. (part applied: air mattress)

#### 5. Switch on

- ⚠ **NOTE:** The unit can be disconnected to switch off the device.

⚠ **NOTE:** (For models WITHOUT the low air leak function) during power failures, the quick connector can be covered with the transport cap to maintain the air pressure inside the cells.

⚠ **NOTE:** After installation, ensure that the cable or any excess tubes are removed from the transit area, to avoid tripping. Position all DEVICES in a manner to ensure doctors and carers have easy access to the patient.

#### 10.2. Only for the digital compressor model

Based on the patient's weight, set the required pressure using the weight settings (40~200 kg).

⚠ **NOTE:** When the headboard is raised, the weight setting must be increased by at least one step.

#### *Cycle time*

There are four cycle times to choose from: 10 minutes, 15 minutes, 20 minutes, 25 minutes. Press the button to select.

The pump starts up, always in Auto Firm mode. No other functions on the panel can be used during initial inflation. When the mattress is completely inflated, the pump automatically passes into alternate mode.

Press the Mode button to change the mode. The sequence is: Static, Seat Inflate + Alternate, Seat inflate + Static, Max Firm, Alternate.

#### 11. BEFORE EVERY USE

Check correct assembly of the device, particularly connection of the tubes between the mattress and the compressor.

#### 12. WARNINGS ON USE

DANGER - To reduce the risk of electrocution

1. Disconnect the device after every use
2. Do not use the device while washing the patient
3. Do not position the device in a place from which it could fall
4. Do not position the device near or allow it to fall into water or any other liquid
5. Add disconnect immediately at the end

WARNINGS - To reduce the risk of electrocution, fire or injury

1. Never leave the device unsupervised when it is connected to the mains current
2. Assess the risk of trapping the patient's limbs according to current protocols at the structure and monitor the patient appropriately
3. Do not use this system with patients who have spinal injuries

4. Carefully supervise use of this device on children or with children nearby, as there is an electrocution risk or suffocation risk due to swallowing of detached parts of the device
5. Only use this device in the manner indicated and described in this manual  
Do not use accessories not recommended by the manufacturer
6. Never use this device if the power cable or plug is damaged, if it is not working properly, has been damaged, has been dropped or has fallen into water. Take the device to a service centre for checking or repair
7. Keep the cable away from heated surfaces
8. Never block the air vents on the device. Never position the device on a soft surface, such as a bed or sofa, as the air vents could be obstructed. Keep the air intakes free of threads, hair and other similar particles
9. Do not insert or allow any object to fall into any opening or tube on the device
10. Do not alter the device without the manufacturer's authorisation
11. The mattress covers have passed the skin sensitivity and irritation tests  
However, if you suspect that there could have been or there is an allergic reaction, immediately seek medical advice
12. Do not leave overly long tubes in the areas above the bed, to avoid the strangulation risk

### 13. CLEANING

Before using the device, it is recommended to carry out the cleaning procedure described below. Pass a damp cloth with mild detergent over the compressor and keep it away from dust. Always choose a detergent that contains no chemical agents that could damage the plastic surface of the compressor casing.

 **WARNING:** do not get the compressor wet or immerse it in any type of liquid. Wipe the mattress with a cloth dampened in warm water (maximum 65°C) and a mild detergent. The covering can be washed in sodium hypochlorite diluted in water. Leave all parts to dry in the air before use.

 **WARNING:** do not use products containing alcohol/phenolic components

 **WARNING:** leave the mattress to dry in the air after cleaning, but do not expose it to direct sunlight for a prolonged period

**WARNING: All operations must be performed after switching off the compressor and unplugging the device!**

PRODUCT CODE	SINGLE ELEMENT MATERIAL	MATTRESS BASE MATERIAL	COVER MATERIAL
LAD470	NYLON/PVC	NYLON/PVC	NYLON/PU padded
LAD475	NYLON/PVC	NYLON/PVC	NYLON/PU padded
LAD480	NYLON/PVC	NYLON/PVC	NYLON/PU padded
LAD471	NYLON/TPU	NYLON/PVC	NYLON/PU padded
LAD476	NYLON/TPU	NYLON/PVC	NYLON/PU padded
LAD481	NYLON/TPU	NYLON/PVC	NYLON/PU padded
LAD491	NYLON/TPU	NYLON/PVC	NYLON/PU padded

Maximum washing temperature 60°C; normal process

Do not bleach

Do not iron

Do not dry clean

Tumble drying possible; low temperature

#### 14. STORAGE

1. Position the mattress on a flat surface, turn it over, roll it up and store it.
2. Roll the mattress from top to bottom

 **NOTE:** Do not fold, crease or stack the mattress.

#### 15. MAINTENANCE

The devices in the LEVITAS range by Moretti have been carefully checked and CE marked before being placed on sale. The devices are intended to offer safe and reliable functioning when used or installed according to the instructions of the manufacturer, who recommends having the system inspected and repaired by authorised technicians if there are signs of wear and tear or malfunctioning of the device. If not, the devices should be serviced and checked every two years.

##### 15.1 General

1. Check the power cable and plug for any abrasions or excessive wear and tear.
2. Check the mattress cover for signs of wear and tear or damage.
3. Disconnect the air tube from the mattress. Then check the air flow from the two air outlets on the compressor. These outlets should blow out air alternately when the compressor is set to "alternate" mode.
4. Check the air tubes to identify any faults or breakages. Contact the local distributor for replacement.

### 15.2 Replacement of the air filter


1. Disconnect the compressor from the mains electricity.
2. Remove the air filter cover, which is located on the back of the compressor.
3. Replace the air filter and put the cover back on.

### 15.3 Element replacement

1. Disconnect the compressor<sup>2</sup> from the mains electricity<sup>3</sup>
2. Remove the ring on the air connector and disconnect the tube
3. Remove the damaged element and install the new one

### 15.4 Fuse replacement

1. Disconnect the compressor from the mains electricity
2. Use a screwdriver to remove the fuse holder
3. Insert a new fuse and put the protective cover back on

## 16. CONDITIONS OF DISPOSAL

### 16.1 General disposal conditions

Never dispose of the product as normal domestic waste. Dispose of the product at a sorted waste collection centre for recycling.

### 16.2 Instructions for correct disposal in accordance with European directive 2012/19/UE:

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment.

## 17. REPLACEMENT PARTS AND ACCESSORIES

Use only replacement parts and accessories listed in the Moretti general catalogue.

## 18. TROUBLESHOOTING

If the following solutions to your problems are not sufficient, please contact your retailer for qualified technical assistance.

PROBLEM	SOLUTION
There is no power supply	<ul style="list-style-type: none"> <li>• Check that the plug is connected to the mains.</li> </ul>
The alarm is on	<ul style="list-style-type: none"> <li>• Check whether the power supply switches off suddenly</li> <li>• Check that the CPR is insulated</li> <li>• Check that there is a firm connection between the air tube connector and the compressor unit</li> <li>• Check that all tube connections along the mattress are firm</li> </ul>
The patient touches the bottom	<ul style="list-style-type: none"> <li>• The pressure settings might be inadequate for the patient. Adjust the comfort range to one or two levels higher and wait several minutes until the greatest level of comfort is reached</li> </ul>
The mattress is floppy	<ul style="list-style-type: none"> <li>• Check that all the snap fasteners and straps on the mattress are properly attached.</li> <li>• Check that the mattress is attached to the bed structure with the straps</li> </ul>
No air is coming out of several openings on the air tube connector	<ul style="list-style-type: none"> <li>• This is normal when the device is being used in alternate mode. The air outlets are activated alternately to produce air during the cycle time</li> </ul>

**⚠ NOTE:** If the pressure level always remains low, check for possible leaks. If necessary, replace any damaged tubes or contact a local retailer qualified to carry out repairs.

## 19. DECLARATION ELECTROMAGNETIC COMPATIBILITY

### 19.1 Manufacturer's guide and declaration - electromagnetic emissions

The Piuma UP system is intended to be used in the following electromagnetic conditions.

The client or user of the Piuma UP system must ensure that it is used in these conditions.

EMISSION TEST	CONFORMITY	ELECTROMAGNETIC ENVIRONMENT - GUIDE
Emissions RF CISPR 15	Conforming	The alternating pressure mattress is not suitable for interconnection with other equipment
Harmonic emissions IEC 1000-3-2	Class C	Class C in reference to IEC 61000-3-2. Conforms with IEC 61000-3-3. The alternating pressure mattress is suitable for use in all buildings, including homes and buildings connected directly to the low-voltage mains electricity network for residential buildings
Voltage flutter/flicker IEC 61000-3-3	Conforming	

### 19.2 Manufacturer's guide and declaration - electromagnetic immunity

The Piuma UP system is intended to be used in the following electromagnetic conditions. The client or user of the Piuma UP system must ensure that it is used in these conditions.

Immunity test	Test level IEC60601	CONFORMITY level	ELECTROMAGNETIC environment - guide
Electrostatic discharge (ESD) IEC 61000-4-2	6kV contact ± 8kV air	± 6kV contact ± 8kV air	The floor must be wood, concrete or tile. If the floor is covered with a synthetic material, the relative humidity must be at least 30%
Electrical fast transient/ burst immunity IEC 61000-4-4	± 2kV for power lines ± 1kV for I/O lines	± 2kV for power lines ± 1kV for I/O lines	The quality of the mains voltage must be that of a typical commercial or hospital environment.
Surge voltage IEC61000-4-5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line	The mains power must be that of a typical commercial or hospital environment.
Immunity test	Test level IEC60601	CONFORMITY level	ELECTROMAGNETIC environment - guide
Voltage drop outs, brief interruptions and variations of voltage on the power supply lines IEC61000-4-11	<5% UT (>95% UT drop out) for 0.5 cycles  40% UT (60% UT drop out) for 5 cycles  70 % UT (30 % UT drop out) for 25 cycles  <5% UT (>95% UT drop out) for 5 seconds	<5% UT (>95% UT drop out) for 0.5 cycles  40% UT (60% UT drop out) for 5 cycles  70 % UT (30 % UT drop out) for 25 cycles  <5% UT (>95% UT drop out) for 5 seconds	The quality of the mains supply must be that of a typical commercial or hospital environment. If the user of the Piuma UP system requires continuous functioning, even during a power failure, it is recommended to supply power to the Piuma UP system with a UPS or with batteries.
Mains frequency magnetic field (50/60Hz) IEC61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should have the typical characteristics of a commercial or hospital environment.
 NOTE: UT is the AC mains voltage prior to application of the test level.			

The Piuma UP system is intended to be used in the following electromagnetic conditions. The client or user of the Piuma UP system must ensure that it is used in these conditions.

Immunity test	Test level IEC60601	CONFORMITY level	ELECTROMAGNETIC environment - guide
conducted RF IEC 61000-4-6	3Vrms/Veff From 150 KHz to 80 MHz	3 Vrms	<p>The Piuma UP system must be used in shielded places with a minimum RF shielding efficiency and, for each cable that enters the shielded room, with a minimum RF filter attenuation</p> <p>Recommended separation distances</p> <p><math>d = 1.167 \sqrt{P}</math></p> <p><math>d = 1.167 \sqrt{P}</math></p> <p>80MHz - 800 MHz</p> <p><math>d = 2.333 \sqrt{P}</math></p> <p>800MHz - 2.5 GHz</p> <p>where P is the maximum output of the transmitter in Watt (W) according to the transmitter manufacturer and "d" is the recommended separation distance in metres (m).</p> <p>The resistances of the RF fixed transmitter fields, as determined from a detection site (NOTE 3) must be lower than the level of conformity of each frequency range. (NOTE 4)</p> <p>Interference may occur in the proximity of equipment marked with this symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m From 80 MHz to 2.5 GHz	3 V/m	

 **NOTE:** The maximum frequency range is applied at 80MHz and 800MHz

 **NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection from structures, objects and people.

 **NOTE:** The intensity of the field for fixed transmitters, such as base stations for radio-telephones (mobile and cordless telephones) and terrestrial mobile radios, amateur radio equipment, AM and FM radio transmitters and TV transmitters cannot be known precisely in theory. To assess an electromagnetic environment caused by fixed RF transmitters, consider performing an on-site electromagnetic inspection. If the field intensity measured outside the shielded room, in the place where a Piuma UP system is being used, exceeds the applicable level of conformity indicated above, normal functioning of the Piuma UP system must be monitored. If abnormal functioning is observed, additional measures may be necessary, such as changing the position of the Piuma UP system or using a shielded room with more efficient RF shielding and greater filter attenuation

 **NOTE:** Above the frequency range of 150KHz and up to 80MHz, the resistance field should be under 3 V/m

Recommended separation distances between portable and mobile radio communications equipment and the Piuma UP system

The Piuma UP system is intended to operate in an electromagnetic environment in which radiated RF disturbances are under control. The client or operator of the Piuma UP system can contribute to preventing electromagnetic interference by ensuring a minimum distance between the mobile and portable RF communications devices (transmitters) and the Piuma UP system as recommended below, in relation to the maximum output power of the radio communication equipment.

For transmitters specified for a maximum output power not listed above, the recommended separation distance "d" in metres (m) can be calculated using the equation applicable to the transmitter frequency, where "P" is the maximum rated output power of the transmitter in Watt (W), as indicated by the transmitter manufacturer.

Maximum output power of the transmitter (specified W)	Separation distance at the transmitter frequency (m)		
	From 150 KHz up to 80MHz $d=1.2\sqrt{P}$	From 80 MHz to 800 MHz $d=1.2\sqrt{P}$	From 800MHz up to 2,5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

**⚠ NOTE:** At 80MHz and 800MHz, the higher frequency band is applied

**⚠ NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection from structures, objects and people

## 20. TECHNICAL SPECIFICATIONS

### 20.1 Compressor

CHARACTERISTICS	LAD470_LAD471_LAD475_LAD476	LAD480_LAD481_LAD491
Power supply	~ 220V/50Hz	
T1A fuse	250V	
Power cable	5 m	
Pressure supplied	55~120mmHg	50~120mmHg
Air outlets	7.5/8.5LPM	8/9LPM
Cycle time	12 min	adjustable 10-15-20-25 min
Noise level	30-35 dB	20-25 dB
Compressor dimensions	30x13.5x11 cm	30x13.5x11 cm
Ambient temperature	- 20°C + 70°C	
Ambient humidity	10%/100%	
Atmospheric pressure	500/1060 hPa	
Classification	Class II, Type LF, IP21 (keep dry) Parts applied: air mattress Cannot be used in the presence of flammable anaesthetic mixtures (AP/APG protection not present)	Class I, Type LF, IP21 (keep dry) Parts applied: air mattress Cannot be used in the presence of flammable anaesthetic mixtures (AP/APG protection not present)

## 20.2 Mattress

CHARACTERISTICS	MATTRESS A ELEMENTS H 12.7	MATTRESS A ELEMENTS H 20.3
Modell	A elements	
Material	NYLON/PVC	NYLON/TPU
Mattress dimensions	196x86x12.8 cm	200x86x21 cm
Single cell height	~ 12.7 cm	~ 20.3 cm
Micro-perforated elements	8	
Pillow function	first three elements	
CPR valve	DEFLATION TIME AROUND 15 SEC.	
Maximum patient weight	160 Kg	200 Kg

## 21. WARRANTY

Moretti products are guaranteed for 2 (two) years from the date of sale against material and manufacturing defects, subject the following limitations. The warranty is voided by improper use, abuse, modifications to the product and failure to follow the instructions. The intended use of the product is given in the user manual. Moretti is not liable for damage, injury or any other consequences resulting from installation or use which are not scrupulously conforming with the instructions given in the installation, assembly and user manual. Moretti does not guarantee its products against damage or defects in the following circumstances: natural disasters, unauthorised repair or maintenance, improper electric power supply (as applicable), use of parts or components not supplied by Moretti, failure to follow the guidelines and instructions for use, tampering, shipping damage (other than the original shipping by Moretti), or failure to perform maintenance as indicated in the manual.

Components subject to wear and tear are not covered by this warranty if the damage is caused by normal use of the product.

## 22. REPAIRS

### 22.1 Repairs under warranty

If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti reseller or its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

### 22.2 Repairing a product not covered by the warranty

A product not covered by warranty may be returned for repair only if authorised in advance by Moretti customer service. The costs of labour and shipping incurred by repairs not covered by the warranty are borne by the client or reseller in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for 6 (six) months from the day of reception of the repaired product.

### 22.3 Non-defective products

The client will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the client at his expense.

**23. SPARE PARTS**

Moretti original replacement parts are guaranteed for 6 (six) months from the day of delivery.

**24. NON-LIABILITY CLAUSE**

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti makes no declaration, guarantee or condition, express or implicit, including any future declaration, guarantee or condition of sale, suitability for a given purpose, non violation and non interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit guarantee under the law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit guarantee or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice.

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