# ULTRASOUND THERAPY UNIT

MANUALE D'USO

IT

## SOMMARIO

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EN

# SUMMARY

IMPORTANT.....	1
WARNINGS AND PRECAUTIONS.....	2
DESCRIPTION.....	4
INTENDED USE .....	5
TARGET USER .....	5
BLOCK DIAGRAM.....	6
STANDARD AND OPTIONAL COMPOSITION .....	6
PHYSICAL ELEMENTS.....	7
<i>Speed</i> .....	8
<i>Wavelength</i> .....	8
<i>Frequency</i> .....	8
<i>Amplitude</i> .....	8
<i>Power</i> .....	8
<i>Intensity</i> .....	9
<i>Modulation</i> .....	9
TRANSDUCERS .....	9
INTERACTIONS BETWEEN TISSUE AND ULTRASOUNDS.....	10
BIOPHYSICAL EFFECTS.....	10
<i>Mechanical effects</i> .....	10
<i>Chemical effects</i> .....	10
<i>Cavitation effects</i> .....	10
<i>Heat effects</i> .....	10
METHODS OF APPLICATION .....	11
<i>By immersion</i> .....	11
<i>Direct contact</i> .....	11
<i>Fixed head</i> .....	11
<i>Mobile head</i> .....	12

CONTROINDICATIONS .....	13
TECHNICAL CHARACTERISTICS .....	14
HARDWARE REQUIREMENTS .....	14
MEANING OF GRAPHICS SYMBOLS .....	15
BOX LABEL .....	16
DESCRIPTION OF THE UNIT .....	16
DESCRIPTION OF THE CONTROLS AND PARTS .....	17
<i>Unit power supply socket</i> .....	17
<i>Power switch</i> .....	17
<i>Alphanumerical display</i> .....	17
<i>Control keyboard</i> .....	17
<i>Warning/indication lights</i> .....	18
<i>Output connectors</i> .....	18
PUTTING INTO SERVICE .....	18
ACCESSORIES INSTALLATION .....	19
APPLICATIONS .....	20
APPLICATIONS TABLE .....	20
POINTS OF APPLICATION .....	22
TRATMENTS .....	24
FREE PROGRAM .....	24
Pre-set program .....	27
PREVENTIVE MAINTENANCE .....	28
<i>Cleaning the Appliance</i> .....	29
<i>Cleaning and Disinfection of the transducer handpiece</i> .....	29
CORRECTIVE MAINTENANCE .....	29
TROUBLESHOOTING .....	30
WARRANTY .....	30

## IMPORTANT

The information contained herein is subject to modifications without forewarning. No part of this document can be photocopied, reproduced, or translated into another language.

It is very important that this instruction manual is kept with the appliance for future consultation. If the appliance should be sold or transferred to another user, ensure that the manual accompanies the appliance, in a way that the new user can be informed of the functioning regarding the appliance and all relative warnings.

If after-sales assistance is required, contact LED SpA.

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[www.led.it](http://www.led.it)

## WARNINGS AND PRECAUTIONS

It is very important that this manual of instructions is preserved together the unit for any future consultation. If the unit will be sold or transferred to other consumer, make sure that the manual comes together, so that the new consumer could be aware of the operation of the unit and of the relative instructions.

**THESE WARNINGS HAVE BEEN WRITTEN FOR YOUR SAFETY AND FOR THE OTHERS ONE, SO PLEASE PAY ATTENTION BEFORE YOU INSTALL AND USE THE UNIT.**

After removing the packing, make sure of the equipment's integrity. If there is a doubt do not use it and turn to qualified staff.

- The electric security of the equipment is assured only when it is correctly connected to an electric installation conformable to the standard norms. It is necessary to verify this fundamental requisite of security. If there is a doubt, ask for a control by qualified staff.
- Before connecting the equipment be sure that the rating (on the back panel) is in conformity to the electric system.
- If there is incompatibility between the wall socket and the supply cable of the equipment substitute it with another one. In general, it is not advised the use of adapters, multiple wall sockets and/or extensions. If their use is essential, it is necessary to use only simple adapters and extensions in conformity to the norms of safe.
- The use of any type of electric equipment requires the observance of some fundamental rules. In particular:
  - Do not touch the equipment with wet hands or feet.
  - Do not use the equipment with barefoot.
- Do not expose the equipment to atmospheric agents (rain, sun, etc.)
- Do not keep the equipment uselessly connected. Switch off the general interrupter when the equipment is not used.
- This equipment is destined only to the use for which it has been designed. Every other use is improper and so dangerous. The manufacturer cannot be considered responsible for possible damages caused by improper or unreasonable uses.
- It is dangerous to modify the characteristics of this equipment.

- Before carrying out any type of cleaning or maintenance, disconnect the equipment through the plug or switching off the mains switch.
- If there is a damage or the equipment doesn't work well, switch it off. For the repairing turn only to a qualified Technical Assistance centre and ask for the original spare parts. The missed respect of these things can be dangerous for the safety of the equipment and for the user.
- It is recommended to carry out the treatments after a medical consultation.
- The device emits electromagnetic radiation and during its use could interfere with specific investigations or treatments.
- This equipment has been manufactured to satisfy the present electromagnetic compatibility. If the functioning of the unit is disturbed or interferes with the normal operation of other electrical and electronic equipment, it is recommended to power the equipment from a different electrical plug and/or try to dispose the unit in a different way until the interference ends.
- Do not use mobile phone near to the equipment to avoid possible interference with it.
- It is recommended to not carry out this therapy at the same time of TENS and/or Holter.
- Many components used inside this equipment are sensible to the electrostatic charges. If you need to handle any electronic components, to avoid eventual electrostatic energy that you could be carrying, touch a grounding surface to discharge the energy. If it is possible, wear a bracelet connected to the ground. If you do not observe these precautions, you could cause permanent damages to the electronic components of the unit.
- When the unit is working do not put the magnets near any type of screens (monitor, tv, etc.), this could cause the images and colours distortion.
- Before using the unit, verify the conditions of cables and equipment to individualize electric dangers.
- The unit does not fit explosive or gas environments.

- Any serious accident occurring in relation to the device must be reported to LED SpA (via Selciatella n.40, 04011 Aprilia (LT) ITALY) and to the competent authority:  
Ministry of Health – General Directorate of Medical Devices and Pharmaceutical Service  
Viale Giorgio Ribotta, 5 – Rome  
E-mail: segr.dgfdm@sanita.it  
Tel.: +39 06 5994 3199 / +39 06 5994 3207

## DESCRIPTION

**unisonic** is an appliance able to pilot an ultrasound transducer hand piece with work frequency of 1 to 3 MHz. Output mode is both continuous and pulsed with modulation that can be set at 25, 50, 75 and 100%.

The **unisonic** ultrasound therapy appliance is destined for temporary use, for surgery and/or home therapeutic applications.

The main fields of application are the pathologies of the locomotory system such as:

- Chronic degenerative tendonitis, tenosynovitis, epicondylitis, enthesitis;
- Fibrous sequelae of distorsions;
- Capsulitis, joint stiffness;
- Scars, keloid adhesions;
- Dupuytren's contracture;
- Ischialgia and neuritis in general;
- Periarthritis;
- Non-surgical treatment of lipomas and localised excess fat (ULTRASONIC HYDROLIPOCLASIA).

The **Unisonic** appliance is completely managed and controlled by a microprocessor.

The appliance envisions functioning with selection of parameters, or the use of programs memorised internally.

The appliance is supplied with a multi-frequency transducer hand piece able to adapt itself to the selected work frequency (1 to 3 MHz). The supplied transducer hand piece is unmistakably associated to the appliance.

The supplied transducer hand piece is compatible with the use of immersion treatments (IP67).

## LEM MEANING

LEM is the brand owned by MORETTI SpA dedicated to the line of electro-medical devices for professional and home use in the healthcare environment.

## INTENDED USE

Active medical device and accessories for temporary use, for professional and home use in the healthcare environment for ultrasound therapy used for the treatment of inflammatory and painful conditions in osteo-articular pathologies.

## TARGET USER


Device for home use: lay person.

The minimum requirements required for the lay person are:

- must be older than 18 years of age.
- must have good knowledge of the language in which the manual is written and used for the medical device.
- must not have any serious deficits and must not have any physical impairments that prevent the activation of the controls on the medical device.
- must have a sound mental constitution (e.g., no use of alcohol or drugs, no mental illnesses leading to hallucinations, loss of balance and similar).
- must be able to assess dangerous situations and act calmly and thoughtfully.
- must know the mathematical symbols.
- must have at least an eighth-grade education.
- must have good eyesight to read and good hearing.

## BLOCK DIAGRAM

MEDICAL DEVICE

## STANDARD AND OPTIONAL COMPOSITION

CODE	DESCRIPTION	QUANTITY
00100.02	Power Cord	■/1
LEM269US138	Ultrasound Handpiece	■/1
LEM269PMDX1	Handpiece support	■/1

00100.02

LEM269US138

LEM269PMDX1

## PHYSICAL ELEMENTS

The ultrasonic waves are sound waves with a frequency higher than can be perceived by the human ear. The frequencies that can be perceived by man are those within 16 and 20.000 Hz.<sup>1</sup> The frequencies below 16Hz are called infrasound, while those above 20kHz are called ultrasounds. Like all sound waves, the ultrasounds need a solid, liquid, or gaseous means for propagation. This takes place when some of the particles of the means start to oscillate when they undergo disturbance exciting the particles nearby due to molecular interactions. This causes compressions and decompressions that propagate through the means.

---

<sup>1</sup> 1 Hertz = 1 oscillation per second

The sound waves are longitudinal waves in so much as the direction of the oscillations of the molecules of the means is parallel to the direction of propagation of the wave.

### *Speed*

The propagation speed of a sound wave depends on the elastic features of the means of transmission. For this reason, it is greater in solids than in liquids and is greater in liquids than in gases.

<i>Means</i>	<i>Speed (m/sec)</i>
Water	1430
Aluminium	6400
Air	330
Grease	1450
Bone	4080
Blood	1570
Soft tissue	1540

### *Wavelength*

The wavelength is the shortest distance between any two points in which the particles oscillate in phase agreement.

### *Frequency*

The movement of a particle from a position and its return to the start position is defined as a cycle. The frequency of a wave means the number of cycles in the unit of time. The unit of measurement of frequency is Hertz (Hz).

### *Amplitude*

The amplitude of a wave is the maximum value of the oscillation.

### *Power*

The power is the energy transmitted by one wave in the unit of time.

### *Intensity*

The intensity of a sound wave is the energy that passes the unit of surface in the unit of time. The unit of measurement used is  $W/cm^2$ .

### *Modulation*

This parameter is important in pulsed emissions, in other words when disturbance of duration  $m$  repeats itself after time  $s$ . A relationship is defined between the duration  $m$  of the impulse and the interval of time given by the sum of the duration of the impulse plus the duration of the interval between two successive impulses ( $m+s$ ). It is expressed as a percentage (%).

## TRANSDUCERS

Devices that can convert a physical quantity of a given nature into a different nature are called transducers. In our case the transducer can change electrical energy into sound energy.

Very often the devices (transducers) can function as Ultrasound generators (US) and as receivers. The latter devices are used in diagnostics.

Generally piezo-electric transducers are used which, excited by an alternating current, create vibrations transmitted to the means of contact. The vibration frequency of the piezo-electric crystal determines the vibration frequency of the ultrasounds. In therapeutic use it is necessary to remember that the depth of penetration is inversely proportional to their frequency (1 MHz about 4/5 cm – 3MHz circa 1/1.5 cm). Generally, we can say that the frequencies of 1 MHz have a field of application in physiotherapy treatments while those at 3 MHz in beauty treatments.

The transducer treatment head is, with the scope of preventing reflection of the band, generally put into contact with the skin using a GEL that is specific for ultrasounds. This gel must have the following requisites:

- it must not be easily absorbed by the skin;
- it must not irritate the dermis;

- it must be chemically inert;
- it must not have gas bubbles.

## INTERACTIONS BETWEEN TISSUE AND ULTRASOUNDS

The intensity of an ultrasound wave, which travels through a tissue, decreases with the increase of the distance travelled. This attenuation of the beam takes place because the energy is partly absorbed and partly dispersed due to diffusion and reflection phenomena. The latter is verified when the wave meets a separation surface between two means in which due to the different elastic features, the sound propagates itself with different speeds. In correspondence of this separation a part of the wave is reflected, and a part transmitted according to the means, however changing direction (refraction). Clearly, the amount of reflection depends on the acoustic impedance of the two means. For example, there is a greater reflection in the passage from a soft tissue to bone and vice versa.

## BIOPHYSICAL EFFECTS

The interaction between ultrasounds and biological tissues produces heat and non-heat effects (mechanical, chemical, cavitation).

### *Mechanical effects*

This effect is due to the movement of the particles of the tissue through which the ultrasound wave travels. In presence of non-homogeneity pressure variations are obtained that determine a movement of the liquids, an increase in permeability of the membrane and the breaking up of the tissues due to separation of the collagen fibres. This movement, apart from rapidity, is more intense than a normal movement created by a manual or instrumental massage, because it has different concentrations of different strengths, or due to the frequency. The pressure and traction phenomenon, with consequent contraction and dilation phenomenon take place in the intimate fibres of the tissues and are separated by very brief intervals, while in a normal massage the mass of tissue is moved.

### *Chemical effects*

The noteworthy acceleration force of the tissue particles through which the ultrasound waves pass, cause the modification of the local PH and the permeability of the cell membranes with molecular changes.

### *Cavitation effects*

On passage of the ultrasound wave through liquids, minute air bubbles may form due to compression-decompression phenomena. On reaching a determined volume these bubbles, following the pressure variations of the field, may become concentrated and oscillate in a stable beam (stable cavitation), or in a turbulent beam where they explode due to great volume variations (unstable cavitation).

The cavitation causes an increase in the flow of liquids in proximity of the oscillations of the bubbles and can be responsible for some therapeutic effects such as cell damage. At therapeutic doses the destructive reactions would only occur in the presence of low cell concentration with low viscosity of the means, e.g., the level of the eye or uterus.

### *Heat effects*

When an ultrasonic wave passes through biological tissue, part of the mechanical energy is converted into heat. The increase of the temperature depends on the energy absorbed and is proportional to the intensity of the ultrasound beam and the absorption co-efficient of the tissue. The increase in temperature varies in relation to the possibility of cooling of the tissue due to blood circulation or movement of the transducers.

The generation of a local diathermy in ultrasound treatment is not immediate. The appearance of the heat phenomenon and most of all its perception by the patient assumes great importance for this non-immediacy detail regarding the "dose" concept during treatment, generally by means of a circular massage on successive delimited tracts. The mentioned correct "dose" phenomenon should, in fact, appear as a pleasant sensation of heat after a few minutes of treatment. The lack of or excessive delay of the manifestation of the phenomenon, denotes a dose-base insufficiency or incorrect location (on healthy tissue, for example, rather than

sick). The immediate appearance or excessive or irritating sense of warmth indicates excessive dose-base and the pain in this case must be interpreted as a synonym of excessive heat, as also the probability of an error in diagnosis in the evaluation of the state of acuity of a pathological phenomenon, in the sense that it is much more acute than envisioned.

## METHODS OF APPLICATION

The therapy with ultrasounds can be applied in two different ways: direct contact or in water (by immersion).

### *By immersion*

This method is used when the surfaces to treat are irregular or small or when the area is painful to prevent direct contact. The part to be treated is placed in a container with water at a temperature of 37°C, together with the treatment head at a max. distance of 2 –3 cm. The average intensity for this type of treatment is 2 W/cm<sup>2</sup> for 10-15 minutes. The frequency of the treatments is generally daily for a total of 10 sittings.

### *Direct contact*

In treatments using direct contact. A conducting GEL is used to favour transmission, adherence, slipping and elimination of air between the skin and the transducer.

Two techniques are used for direct contact treatments: fixed or mobile head.

### *Fixed head*

The treatment head is positioned on the area to be treated held by a support that maintains the treatment head in contact with the skin for the entire duration of the therapy. The lack of movement and the small treatment area cause rapid temperature increases, therefore a lower power intensity must be used (0.5-1 W/cm<sup>3</sup>) with pulsed supply (25%). Average treatment time is 5-10 minute daily for a total of ten sittings.

### *Mobile head*

The energy is distributed by direct massage, practically by moving the head in a way to follow the principle of the same distribution of the dose of energy in the space. With the head well adherent in all of its points, possibly without air bubbles, which would reflect some of the waves with a localised heating effect of the skin in those small points, the massage must be performed slowly starting from one end of the area moving towards the other end after having touched the points of interest. To do this, the best technique is that of tracing concentric circles, almost overlapping each other, by performing a very restricted action it is certain that the area will be saturated with energy, abandoning it to pass to the surrounding area. This method is called "micro massage" and must be performed slowly, at the speed of one rotation per minute for a circle not much bigger than the diameter of the head. Pressure must be homogeneous and not exaggerated. After a few minutes of circular movements, with the correct dose, the patient detects a warm feeling in the treated area, which represents the certainty to have practiced the correct intensity-time ratio in that area, which is a very important factor. The time passing before feeling the heat is naturally variable and depends on the power, the speed of the massage, local anatomic conditions, the pathological state and the various neurovegetative reflexes of the patient. The sensation of immediate pain must be avoided, it being an index that the vibrations pass from a less dense means to a denser means (bone), causing excessive refraction, which determines irradiation of the periosteal. In these cases, vary the beam incidence and reduce the power. The average intensity for use for mobile head treatment is  $1 \text{ W/cm}^2$  for 10-15 minutes. The frequency of treatments is generally daily for a total of 10 sittings.

## CONTRAINDICATIONS

The application of ultrasound therapy is contraindicated in subjects:

- carriers of pacemaker or of other electronic installed device, or in proximity of such carrier patient.
- with cardiac illness.
- with serious unbalances of the arterial pressure.
- with serious illness of the nervous system.
- with serious renal insufficiencies.
- in state of pregnancy.
- haemorrhagic.
- epileptic.
- carriers of metallic prosthesis.
- with troubles of the thermic sensibility.

markedly asthenic.

- with troubles of the thermic sensibility.
- markedly asthenic.

It is besides contraindicated its application:

- in the tumoral focuses.
- on noble organs, like heart, vessels, gonads.
- on neoplastic lesions.
- in progress of contagious infectious illnesses (TBC).
- in infanto-juvenile age.
- in case of delays of consolidation or pseudoarthrosis.
- on the anterolateral intercostal zone, on the abdomen and on the segments with means of endotissual synthesis.
- on the testicles.
- on the eyes.
- on the cartilages of growth.
- in the thromboses and thrombophlebitis.
- on the abdomen, in all the cases of calculous or inflammation.


## TECHNICAL CHARACTERISTICS

Supply Voltage	230 Vac
Frequency	50 Hz
Max absorbed power	40 VA
Internal Fuse	T 2.5A (TR5)
Ultrasound Frequency	1.0 / 3.0 MHz
Magneto Max	from 5 to 100 Gauss (from 0.5 to 10 mT)
Duty Cycle	25-50-75-100%
Modulation Waveform	continuous or mod. Imp.
Modulation Frequency	2 Hz
Output Intensity	2 W/cm <sup>2</sup>
Electric Safety Class (EN60601-1)	Class II Applied Part BF
EMC Classification (EN60601-1-2)	Group 2 Class B
Protection class (EN60529)	IP21
MDR 2017/745/UE Class	Ila
Physical characteristics	
- Dimensions (HxWxD)	110 x 270 x 220 mm
- Weight	1.4 Kg
Environmental operating conditions	
- Temperature	between +10°C and +40°C
- Relative humidity	between 30% and 75% non-condensing
- Atmospheric pressure	between 70 kPa and 106 kPa
Environmental storage conditions	
- Temperature	between -10°C and +50°C
- Relative humidity	between 10% and 85% non-condensing
- Atmospheric pressure	between 50 kPa and 106 kPa

## HARDWARE REQUIREMENTS

Microcontroller	ARM Cortex M4
Frequenza Clock	200 MHz
Flash	2048 KB
Ram	512 KB
SDRAM DDR2	512 MB
Nand Flash	1 Gb
Peripherals	UART, I2C, SPI, Watch-dog timer, USB2.0
Visual	Display Touchscreen 7" 800x480 px

## MEANING OF GRAPHICS SYMBOLS


### Meaning of the symbols:

The meaning of the graphic symbols printed on the device is as follows:


- |  |   |
|--|---|
| <ol style="list-style-type: none"> <li>1. Floating Neutral Electrode: patient circuit is isolated from earth at both high and low frequencies.</li> <li>2. The equipment is BF class, protected against electric shock than that provided BY TYPE B APPLIED PARTS.</li> <li>3. Not Ionizing Radiation emitted.</li> <li>4. Operating instructions.</li> <li>5. CE Mark (2017/745/UE) + Number of Notify Body 0051 = IMQ Italy</li> <li>6. The product mustn't be throwing in the containers for urban wastes, but it must be swallowed by a separate picking.</li> <li>7. Manufacturer</li> <li>8. Serial Number</li> <li>9. Date of manufacturer</li> <li>10. UDI Unique Device Identification</li> </ol> | <ol style="list-style-type: none"> <li>11. MD Medical Device</li> <li>12. Distributor</li> <li>13. Not to be serviced by users</li> <li>14. REF Reference Code Number</li> <li>15. Temperature limits</li> <li>16. Humidity limits</li> <li>17. Atmospheric Pressure limits</li> <li>18. This way up</li> <li>19. FRAGILE – Handle With Care</li> <li>20. KEEP AWAY FROM SUNLIGHT</li> <li>21. KEEP DRY</li> <li>22. Stacking limit by number</li> <li>23. Weight</li> <li>24. Size</li> <li>25. Number of pieces</li> <li>26. Recycle</li> <li>27. Model/Trade Name</li> <li>28. IP – Protection against harmful ingress of water or particulate matter</li> </ol> |
|--|---|

# BOX LABEL

# DESCRIPTION OF THE UNIT

1	Unit power supply socket
2	Power switch
3	Alphanumerical display
4	Control keyboard
4a	UP key to increase
4b	DWN key to decrease
4c	START/STOP key to start and interrupt
4d	ENTER key to confirm
4e	ESC key to go back
5	Warning/indications lights
5a	Output 1 supplying warning light (yellow)
5b	Unit switched-on warning light (green)
5c	Output connectors

## DESCRIPTION OF THE CONTROLS AND PARTS

### *Unit power supply socket*

1. The power supply socket is present on the rear of the unit. Introduce the supplied power cable into this socket. The cable must be connected to a power supply that has electric features that are compatible with the unit (see plate data present on the lower part of the unit).

### *Power switch*

2. The power switch is present on the left side of the unit. By placing the switch in position I, the unit switches on. Unit switch-on is indicated by the illumination of the alphanumerical display and the power supply indication light present on the central part of the unit.

### *Alphanumerical display*

3. There is an alphanumerical display in the centre of the unit which displays all indication or variable parameters in a determined procedure.

### *Control keyboard*

4. The control keyboard is made of:

**4a UP key** to increase the value, to select an option present on the previous line of the display or to pass to a previous page (indication on the display of the upward arrow “↑”).

**4b DWN key** to decrease the value, to select an option present on the following line of the display or to pass to a following page (indication on the display of the downward arrow “↓”).

**4c START/STOP key** to start/interrupt the supply of current.

**4d ENTER key** to confirm a data on the display.

**4e ESC key** to go back to the previous selection page.

### *Warning/indication lights*

5. Three luminous warning lights are present underneath the keys. They have the following meaning:

1	Output 1 supplying warning light (yellow)
2	Output 2 supplying warning light (yellow)
	Unit on indication light (green)

### *Output connectors*

The connection points of the supplied output cables are present in this area.

## PUTTING INTO SERVICE

- Unpack the appliance and the accessories, check that these have not undergone damage during transport. The carrier should be informed of any damage immediately. If damage has occurred, do not use the appliance and/or accessories and have them checked by qualified staff. If the appliance is returned, it is necessary to use the original wrapping or packaging that guarantees safety during transport.
- Carefully study the supplied documentation and operational instructions.
- Check that the power supply voltage available is in compliance with that requested by the appliance (see plate data positioned on the lower part of the unit). Position the unit on a hard surface in a way not to obstruct the cooling slots on the lower part (e.g., do not place on covers).
- Position the power switch, present on the left side, on O (off). Insert the power supply cable supplied into the socket present on the rear part.
- Insert the transducer hand piece connector to the front. For correct tightness, turn the lock nut.


**IMPORTANT:** a label is present on the transducer hand piece cable, in the vicinity of the connector, which has a 10-digit number. Check that this number corresponds to the serial number on the appliance (present on the lower part of the unit among the plate data), as the transducer hand piece and the unit are unmistakably associated.

If several units are available, do not invert the transducer hand pieces.

- Switch the appliance on, by positioning the power switch on I (on). Switch-on is indicated by the illumination of the green indication light as well as switch-on of the alphanumeric display. Your **Unisonic** is now ready for use.
- Use the indications given in the “Applications” chapter.
- Start the treatment by following the indications stated in the “Treatment” chapter.

## ACCESSORIES INSTALLATION

### ***POWER CABLE***

### ***HANDPIECE***

For the assembly of the accessories refer to the paragraph of the commissioning.

At the end of the treatment remove the accessories, turn off the device and remove the connection from the power outlet.

## APPLICATIONS

Any therapeutic protocols should be established by the attending practitioner. A table indicating several treatments is shown successively. Refer to this paragraph for application points. They must however be considered benchmark.

Generally, daily applications are initially recommended, passing to applications on alternate days, progressively increasing the pause days between the applications until the desired result is attained. If optimal results are obtained after only a few sittings it is however recommended to end the cycle of applications to obtain stability of the result.

## APPLICATIONS TABLE

Column P indicates the programs pre-set in the unit. For the others, use the free program mode (see Treatment chapter).

<i>P</i>	<i>Treatment</i>	<i>FREQUENCY(MHz)</i>	<i>MODULATION (%)</i>	<i>DURATION (min)</i>	<i>POWER (%)</i>	<i>Notes</i>
X	ANKLE	1	50	10	100	M.H.
	ANTERIOR TIBIALIS TENDONITIS	1	25	10	50	M.H.
X	ARTHROSIS	1	100	10	50	M.H.
	BICEPS FEMORIS MUSCLE	1	25	10	40	M.H.
	BICIPITAL TENDONITIS	1	100	10	40	M.H.
X	BURSITIS	1	50	10	100	M.H.
X	CAPSULITIS	1	50	10	50	M.H.
X	CARPAL TUNNEL	1	100	10	50	M.H.
X	COXO-ARTHROSIS	1	100	15	50	M.H.
X	DELTOID	1	100	12	40	M.H.
X	EPICONDYLITIS	1	100	10	30	M.H.
X	EPITROCLEYTIS	1	100	10	30	M.H.
	FLEXOR MUSCLES	1	25	10	30	M.H.
X	HANDENING	3	50	20	100	M.H.

<i>P</i>	<i>Treatment</i>	<i>FREQUENCY(MHz)</i>	<i>MODULATION (%)</i>	<i>DURATION (min)</i>	<i>POWER (%)</i>	<i>Notes</i>
X	LIGAMENT	1	75	10	50	M.H.
X	LUMBAGO	1	75	10	50	M.H.
X	METATARSUS PAIN	1	50	10	100	RFI.
X	MYALGIAS	1	75	10	30	M.H.
X	MYOTENDON	1	50	10	40	M.H.
	OCCIPITAL NEURALGIA	1	25	10	40	M.H.
X	PLANTAR FASCIITIS	1	50	10	100	RFI.
X	POLYARTHROSIS	1	100	10	80	M.H.
X	QUADRICEPS	1	75	10	80	M.H.
	SACRALE PAIN	1	100	10	40	M.H.
X	STYLOID	1	100	10	80	M.H.
	SUBDELTOID BURSITIS	1	100	10	30	M.H.
	SUBSCAPULAR PAIN	1	25	10	60	M.H.
	SURAL TRICEPS	1	50	10	50	M.H.
	TEMPOROMANDIBULAR PAIN	1	25	5	10	M.H.
X	TENDINITIS	1	75	10	40	M.H.
X	TONING	3	50	30	100	M.H.
X	TRICEPS	1	75	15	80	M.H.
	TRICEPS BRACCHI MUSCLE	1	100	10	40	M.H.



*RFI: Recommended for immersion (see Ultrasounds chapter – Methods of application)*

*M.H.: Mobile Head (see Ultrasounds chapter – Methods of application)*

## POINTS OF APPLICATION

Some application points are shown below.

Femoraltibial arthrosis	Subdeltoid bursitis	Pain in the sacral region	Subscapular pain
			
Coxoarthrosis pain	Temporomandibular pain	Epicondylitis	Plantar fasciitis
			

<p><b>Lumbago</b></p>	<p><b>Metatarsus pain</b></p>	<p><b>Myalgia of the trapezius muscle</b></p>	<p><b>Occipital neuralgia</b></p>
			
<p><b>Pathology of the flexor muscles</b></p>	<p><b>Pathology of the biceps femoris muscle</b></p>	<p><b>Pathology of the lateral ankle compartment</b></p>	<p><b>Pathology of the sural triceps muscle</b></p>
			
<p><b>Pathology of the popliteal</b></p>	<p><b>Muscular pathologies of the rotator cuff and the triceps brachii muscle</b></p>	<p><b>Bicipital tendonitis</b></p>	<p><b>Tendinitis of the anterior rectus muscle</b></p>
			
<p><b>Anterior tibialis tendonitis</b></p>			
			

## TREATMENTS

```

** MORETTI SpA **
UNISONIC
Re1 . 1.0.0
  
```

On switch-on using the power switch on the side of the unit, the initial screen will appear. The software version is displayed in the lower part. Press ENTER.

```

PROGRAMS 1/26
→FREE PROGRAM
ANKLE
ARTHROSIS ↓
  
```

The screen is displayed for the selection of the program to be used. The unit allows to work with a free program, where the parameters to be used are introduced from time to time or with pre-set programs. To select the program, press UP (to ascend) and DWN (to descend). Once the desired program has been selected, press ENTER to pass to the successive screen.

*NOTE:* If the treatment head is used in contact with the area to be treated, cover the area well with a layer of conductive gel specific for ultrasounds.

## FREE PROGRAM

```

WORK FREQUENCY:
→ 1 MHz
3 MHz
  
```

On selection of the “Free Program”, the screen is displayed for selection of the value of work frequency. With the scope of helping, the unit memorises all the latest settings used. To modify the value indicated, press the UP key (to increase the value) and DWN key (to decrease the value). Once the desired value has been selected, press the ENTER key. Press ESC to return to the previous screen.

```

MODULATION:
→ 50 %
  
```

After having selected the work frequency, the screen is visualised for selection of modulation. To modify the value indicated, press the UP and DWN keys.

Once the desired value has been selected, press the ENTER key. Press ESC to return to the previous screen.

TREATMENT TIME →10 min
---------------------------

After having selected the modulation, the screen is visualised for selection of the duration of the treatment. To modify the value indicated, press the UP and DWN keys. Press ENTER to confirm the value. Press ESC to go back to the previous screen.

PRESS START P. 0.0 W/cm <sup>2</sup> 0 W MODULATION      50 % FREQ 1MHZ      10:00 min	↔	<table> <tr> <td style="text-align: center;">P. 0.0 W/cm<sup>2</sup>      0 W MODULATION      50 % FREQ 1MHZ      10:00 min</td> </tr> </table>	P. 0.0 W/cm <sup>2</sup> 0 W MODULATION      50 % FREQ 1MHZ      10:00 min
P. 0.0 W/cm <sup>2</sup> 0 W MODULATION      50 % FREQ 1MHZ      10:00 min			

Once the duration of the treatment has been selected the unit is ready to work. "PRESS START" flashes on the screen with the indication of the previously selected settings. To start the treatment, press the Start/Stop button. Press ESC to go back to the previous treatment duration selection screen.

FREE PROGRAM P. 1.0 W/cm <sup>2</sup> 5 W ← MODULATION      50 % FREQ 1MHZ      9:50 min
---

In the Start status, the FREE PROGRAM name is displayed in the upper part. The arrow indicates the channel, which by pressing the UP and DWN keys raises and lowers the output level.

**IMPORTANT:** to prevent damage to the ultrasound transducer, activate emission after having dissipated the hand piece.

**NOTE:** Slowly raise the output level in a way to detect a slight sense of warmth without being uncomfortable. It is preferable to move the treatment head with a circular motion.


To pass from one alterable parameter to another, press ENTER.

The display shows the effective intensity value (W/cm<sup>3</sup>) and the power value (W). By modifying the power value, the effective intensity is automatically stated.

**IMPORTANT:** To prevent the danger of the formation of stationary waves it is necessary to move the head slowly and uniformly with circular movements. In the case of treatment with the static treatment head, set a low modulation level (25) and decrease the power (20-50%).


The treatment duration time only decreases with the output power equal to or exceeding 1 Watt. The yellow warning lights below are only switched-on when ultrasounds are emitted. The flashing speed depends on the modulation value (100% fixed light).

The treatment can be interrupted by pressing the Start/Stop key.

CONTINUE?  
→ YES  
NO

“Continue” appears on the screen. Use the UP and DWN keys to select the desired option. If the treatment is to be continued (YES), go back to the previous screen with channel levels at zero and the indication of the residual duration of the treatment. To re-start the treatment, press START again.

CONTINUE?  
YES  
→ NO

If the answer is NO or by pressing ESC, the treatment is finished.

TREATMENT  
ENDED

At the natural conclusion of the treatment time or after the negation in an interruption of the treatment, the indication of the end of the treatment appears on the screen. Press the ESC key to go back to the program selection screen.

## Pre-set program

Using the UP and DWN keys it is possible to select one of the programs present in the unit.

PRESS START	
P. 1.0 W/cm <sup>2</sup>	5 W
MODULATION	100 %
FREQ 1MHZ	10:00 min
↔	
P. 1.0 W/cm <sup>2</sup>	5 W
MODULATION	100 %
FREQ 1MHZ	10:00 min

Once a program has been selected, the unit is ready to work. "PRESS START" flashes on the screen with the indication of the parameters of the selected program.

To start the treatment press Start/Stop. Press the ESC key to return to the previous program selection screen.

ARTHROSIS	
P. 1.0 W/cm <sup>2</sup>	5 W
MODULATION	100 %
FREQ 1MHZ	9:58 min

In the Start status, the name of the program is displayed in the upper part. With pre-set programs the parameters cannot be varied.

The display shows the effective value of intensity (W/cm<sup>2</sup>) and power (W). On modifying the power value, the effective intensity is automatically stated.

The yellow warning lights below are switched-on when ultrasounds are being emitted. The intermittence of switch on depends on the modulation value (100% switched-on).

**IMPORTANT:** to prevent damage to the ultrasound transducer, activate emission after having dissipated the hand piece.

**NOTE:** In treatments not in immersion, move the treatment head with a circular motion.

The treatment can be interrupted by pressing Start/Stop.

CONTINUE?  
→YES  
NO

“Continue” appears on the screen. Use the UP and DWN keys to select the desired option. If the treatment is to be continued (YES) go back to the previous screen with channel levels at zero and the indication of the residual duration of the treatment. To re-start the treatment, press START again.

CONTINUE?  
YES  
→ NO

If the answer is NO or by pressing ESC, the treatment is finished.

TREATMENT  
ENDED

At the natural conclusion of the treatment time or after the negation in an interruption of treatment, the indication to end treatment appears on the screen. By pressing ESC, go back to the program selection screen.

## PREVENTIVE MAINTENANCE

**DISCONNECT ALL ELECTRICAL CONNECTIONS BEFORE CARRYING OUT ANY  
MAINTENANCE OPERATIONS.**

Preventive maintenance consists mainly in:

- The control before use of the status of the appliance and the accessories with the aim of examining the integrity of all cables, the insulations, the casing etc. designed to prevent access to live parts.
- Periodical cleaning of the appliance to be carried out on a regular basis with the scope of keeping it in the best working conditions and with a good appearance.
- The periodical replacement of the adhesive electrodes.

### *Cleaning the Appliance*

Clean the external container and control panels with a damp cloth: do not use abrasive products or solvents.

Do not place the appliance in liquids. Do not use the appliance if liquids have penetrated it and contact professional qualified staff.

### *Cleaning and Disinfection of the transducer handpiece*

Clean the transducer hand piece using a damp cloth: do not use abrasive products or solvents. Disinfect the parts in contact with the patient with an alcohol solution after each treatment.

Delicately remove any residues of gel present in the recesses using cotton wool soaked in alcohol.

**NOTE:** The appliance should be checked periodically (at least once a year) by qualified staff to verify the following electrical safety values:

- An earth leakage measurement.

## **CORRECTIVE MAINTENANCE**

Whenever malfunctioning is detected, it is suggested to check that there is no error in the layout of the controls.

If the appliance or transducer hand piece undergo external mechanical stress, for example after a serious fall or if the appliance or transducer hand piece have been dripped on by liquids, or if the appliance or transducer hand piece have been strongly overheated (e.g. direct sunlight, fire) or if the functionality of the

appliance or transducer hand piece seem to be altered or if parts of the casing of the appliance or cables are split, moved or missing, or if any connector or transducer hand piece shows signs of deterioration, the appliance and the relative accessories should be controlled by professionally qualified staff.

## TROUBLESHOOTING

Below are some simple solutions to problems that may occur on the Unisonic:

Even if the power switch is in position I (on) the appliance does not switch on.	Check that the power supply cable is introduced correctly. Ensure the effective functioning of your power supply socket.
It appears that the appliance does not function correctly when switched on.	Switch the appliance off, wait for 30-45 seconds and try and switch it on again.
Very painful treatment and/or evident redness	Decrease the output level and increase modulation (decrease value set on modulation).

If, despite the interventions, problems are detected, take the appliance and accessories, together with a description of the problem detected to an authorised after-sales centre.

**Do not open the appliance for any reason. Apart from being potentially dangerous, it will make all forms of warranty become void.**

## WARRANTY

- The appliance is guaranteed for 24 months, any electric accumulators are guaranteed for 3 months, accessories are excluded from the warranty.
- In the case of assistance under warranty, a copy of the receipt with purchase date must be sent with the appliance. Whenever the user cannot

show the documents proving the date of purchase, the repair will be performed under warranty if on the date of the intervention by our staff no more than 24 months have passed since production of the appliance.

- During the warranty period all parts that according to the unquestionable judgement of the manufacturer result as being faulty at a material or working level will be repaired or replaced free of charge. The intervention is carried out at the manufacturer's after-sales technical service and the appliance must be delivered by the user.
- The warranty does not cover mobile or removable plastic parts, bulbs, glass parts, external pipes, and any accessories.
- The warranty is not applied, and assistance will be carried out at a charge in the following cases:
  - For imperfect or incomplete installation.
  - For damage due to transport.
  - For damage due to accidents (lightening, earthquakes, fires, etc.).
  - If the serial number has been removed, deleted, or altered.
  - If the owner of the appliance is not the first purchaser.
  - For negligent use, carelessness, or bad maintenance.
  - For insufficient capacity or abnormality of electric plants.
  - For repairs, modifications or tampering carried out by staff not authorised by the manufacturer and for all cases that do not depend on the manufacturer.
- These conditions cannot be modified by other verbal or written agreements.
- The warranty excludes any compensation for the period of inefficiency of the appliance.

***Information about elimination of this product (Applicable in the European Union and other European countries with separate collection systems)***

On the end of the life, the present product mustn't be eliminated as urban refusal, but it must be eliminated in a separated collection.

If the product is eliminated in unsuitable way, it is possible that some parts of the product (for example some accumulators) could be negative for the environment and for the human health.

The symbol on the side (barred dustbin on wheel) denotes that the products mustn't throw into urban refuses container but it must be eliminated with separate collection.

In case of abusive elimination of this product, could be foreseen sanctions.

# EMC DECLARATIONS

## EMISSIONS

<b>Guide and manufacturer's declaration – electromagnetic emissions</b>		
The device is intended to operate in the electromagnetic environment specified below. The Customer or user of the device should ensure that it is used in this environment.		
<b>Emission Test</b>	<b>Conformity</b>	<b>Electromagnetic environment</b>
RF CISPR 11 Emissions	Group 1	The device uses RF energy only for its internal functioning. As a result, its RF emissions are very low and probably does not cause any interference in nearby electronic devices.
RF CISPR 11 Emissions	Class B	The device is suitable for use in all environments, including domestic and those connected directly to a low-voltage public mains power supply that powers buildings used for domestic purposes.
Emission of harmonics on the electricity network IEC 61000-3-2	Class A	
Voltage fluctuations and flicker generation IEC 61000-3-3	Conform	

## IMMUNITY

Electromagnetic immunity tests are intended to simulate disturbances typical of an electromagnetic environment. The device has been tested for immunity to disorders and their levels typical of a domestic, hospital and commercial environment.

**Guidance and manufacturer's declaration - electromagnetic immunity**

The device is intended to operate in the electromagnetic environment specified below. The Customer or the user of the device should ensure that it is used in this environment.

<b>Immunity Test</b>	<b>Conformity</b>	<b>Electromagnetic Environment and Precautions</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Floors should be made by wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Transient/sequence of fast electrical pulses IEC 61000-4-4	±2 kV for power supply lines	The mains power supply quality and the magnetic fields at mains frequency should be the
Impulse IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	

<p>Voltage dips, short interruptions, and voltage variations on power input lines IEC 61000-4-11</p>	<p>&lt;5 % of the nominal voltage (<math>U_T</math>) (Voltage dip &gt;95 %) for half a cycle 70 % <math>U_T</math> (Voltage dip 30 %) for 25 cycles &lt;5 % <math>U_T</math> (Voltage dip &gt;95 %) for 5 sec</p>	<p>typical one used in domestic, commercial or hospital environment.</p> <p>If the use of the device requires it to operate continuously during power outages, it is recommended to power it from an uninterruptible power supply (UPS) or batteries.</p>
<p>Magnetic fields at mains frequency (50/60 Hz) IEC 61000-4-8</p>	<p>30 A/m</p>	<p>The magnetic fields at mains frequency should have characteristic levels of a typical location in a commercial or hospital environment.</p>

<p>Conducted RF Fields IEC 61000-4-6</p>	<p>3 Vrms from 150 kHz to 80 MHz</p>			<p>Electromagnetic interference may occur near the devices bearing the symbol below.</p> 
<p>Radiated RF Fields IEC61000-4-3</p>	<p><b>Field (V/m)</b></p>	<p><b>Frequency</b></p>	<p><b>Modulation</b></p>	
	<p>3</p>	<p>80MHz÷2700MHz</p>	<p>1kHz AM 80%</p>	
	<p>27</p>	<p>380MHz÷390MHz</p>	<p>18Hz PM 50%</p>	
	<p>28</p>	<p>430MHz÷470MHz</p>	<p>18Hz PM 50%</p>	
	<p>9</p>	<p>704MHz÷787MHz</p>	<p>217Hz PM 50%</p>	
	<p>28</p>	<p>800MHz÷960MHz</p>	<p>18Hz PM 50%</p>	
	<p>28</p>	<p>1700MHz÷1990MHz</p>	<p>217Hz PM 50%</p>	
	<p>28</p>	<p>2400MHz÷2570MHz</p>	<p>217Hz PM 50%</p>	
	<p>9</p>	<p>5100MHz÷5800MHz</p>	<p>217Hz PM 50%</p>	

**Official Distributor**

Via Bruxelles, 3 - Località Meleto  
I-52022 Cavriglia (AR) ITALY  
Tel +39 055 9621111  
Fax +39 055 9621200  
info@morettispa.com