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
USER MANUAL

PLEASE READ THE USER MANUAL CAREFULLY AND FOLLOW THE INSTRUCTIONS BEFORE USE.

- 1. DEVICE USE AND APPLICANT SCOPE**
Intended use and application scope: This device is used to measure the forehead temperature of the human body. It is intended for use on people of all ages except pre-term babies or very small (small for gestational age) babies.
- Contraindications:** None.
- 2. DEVICE STRUCTURE AND COMPOSITION**
Device main structure and composition: The thermometer consists of a housing, sensor, a display, and a circuit board.

• LCD area.

3.SYMBOLS	Equipment with type BF applied parts
(A)	Warnings and precautions
(R)	Paper recycling
(M)	Manufacture
(E)	Authorized representative in the European community
(S)	Follow instructions for use
(X)	Symbol for the marking of electrical and electronics devices according to Directive 2012/19/EU
(Y)	Symbol for the marking according to Directive 2006/66/EC
(Z)	Safety and environmental protection use period for 10 years
(IP22)	Prevent the invasion of solid objects larger than 12mm. When the lamp is tilted 15 degrees from the vertical Angle, the dripping water will not cause harmful effects.
(L)	Temperature range
(U)	Humidity range
(P)	Atmospheric pressure range
(CE)	This device fulfills the provisions of EC directive 93/42/EEC. (Medical Device Directive).
(I)	See instruction for use
(LOT)	Lot Number
(REF)	Product catalogue reference code or catalog number

WARNINGS:

- The measurement result is only for reference, which is not a substitute for a physician's diagnosis. It is very dangerous to self-judge and treat only based on the measurement result. Please follow the doctor's instructions.
- Please put the battery out of the reach of children, otherwise it is dangerous.
- When the product is not used for a long time (more than 3 months), remove the battery from device to prevent the battery leakage.
- If there is a temperature difference compared environment of the device storage with that of the measuring environment, place the device in the measuring environment for more than 30 minutes, otherwise there may be error for the measurement result.
- It is forbidden to immerse the infrared thermometer in any liquid, and it is forbidden to use it for a long time under too high or low temperature condition. No collisions, drops and mixing with sharp objects.
- Do not put the battery close to the fire or into the fire to avoid the battery explosion. Do not use the battery when it leaks or molds; when discarding batteries or this product, it should be following local regulations to avoid contamination.
- This product contains sensitive electronic components, should avoid strong electromagnetic interference.
- Do not modify this equipment without the authorization of the manufacturer, which may lead to measurement error or machine malfunction.
- If the situations cannot be solved or unexpected problem happens, please consult the local distributor.

CAUTION:

- Do not try to measure when the device is wet, which may cause measurement results inaccurate.
- Before measurement, please make sure that there is no sweat, cosmetics or oil stains on the forehead of subject, please make sure that the subject does not take a bath, exercise or have a meal within 30 minutes, and the body is measured at a steady state.
- During measurement, do not let the subject directly face the air outlet of the fan, air conditioner or heater, which will change the temperature of the forehead. Please conduct measurement in a stable environment as far as possible.
- When the product gets wet due to contact with steam, use after it become dry or gently wipe with a soft dry cloth or cotton ball, otherwise it will cause measurement errors.
- For patients measuring their own temperature, it is recommended to measure close to the forehead.
- Please read this instruction carefully before use to confirm that the battery is installed.

This product contains batteries and recyclable electronic waste. To protect the environment, do not dispose of it in the household waste, but take it to appropriate local collection points.

4. MEASUREMENT AND TRANSPORT STORAGE ENVIRONMENT
Measurement environment: Temperature: 10°C (50°F) to 40°C (104°F)
 Relative humidity: 15% to 90%RH (non-condensing)
 Atmospheric pressure: 70kPa to 106kPa
Transportation and storage environment: Temperature: -20°C (-4°F) to +55°C (131°F)
 Relative humidity: 15% to 90%RH (non-condensing)
 Atmospheric pressure: 70kPa to 106kPa

5. INSTALLATION AND USAGE
Installing the battery.
 The device is supplied with 2 Section 7 alkaline batteries. Push the battery cover downward and load the battery into the battery compartment. At this time, the device will start self-inspection. Pay attention to the positive and negative poles instruction in the battery compartment and cover the battery cover. Refer to the picture at right.

Setting measurement unit.
 In the shutdown state, long press the starting measuring button for 10 seconds to enter the

temperature unit conversion state, and then short press the starting measuring button to select "C" and "F" temperature unit, then long press the starting measuring button to final confirmation.

Measurement of human body.
 1. Point the infrared thermometer probe at the center of eyebrow, within 0-5cm of the forehead.
 2. Press the Start measuring button . After about 1 second, the infrared thermometer will have vibration and sound prompts to display the measurement results.

- Notes:**
- When the measurement distance is more than about 5 cm, the word "0-5" will be displayed on the screen, please close to the subject, as the picture shown at right;
 - When the measurement distance is less than about 5 cm, the device will automatically measure after the measurement, it will vibrate and have "beep" sound, and measurement result will be displayed.
 - When temperature is 37.6°C (99.7°F) or higher, the device will make quick continuous three-time "beep-beep-beep" alarm sound.
 - When the measured result is beyond device measurement range 32°C ~ 43°C (89.6°F ~ 109.4°F), thermometer will make quick continuous three-time "hi-di-di" alarm sound, please make sure the measurement method and the external environment is normal at this time.
- Sound on/off setting.**
 When the prompt on/off button slides to , the sound is on.
 When the prompt on/off button slides to , the sound is off.
- Power off.**
 1. Long press Starting measurement button for 3-5 seconds to turn off the device;
 2. In the absence of any operation, the device will automatically shut down after about 30 seconds.

6. PRODUCT MAINTANCE AND UPKEEP METHOD
 As the product is reusable device, please pay attention to cleaning and disinfection after use. If the product is dirty, please keep the sensor and probe cavity clean, otherwise the measurement accuracy will be affected. How to clean sensor and probe cavity: gently wipe the inner cavity or sensor mirror with a clean soft cloth or cotton swab. Do not wipe it with other material, otherwise it may cause lens surface of sensor abrasion or machine malfunction. Disinfection of the product and its head: Gently wipe the surface of the product and the head used for measurement with a cotton swab dipped in a little medical alcohol (70%), it can be used until the alcohol has completely evaporated.

Precaution for storage.
 This product should be stored in a dust-free, dry place; please avoid direct sunlight; do not store in places with high temperature, humidity, dust and corrosive gas. This product is a high-precision device, please do not drop the device. Avoid drastic collision and jobs and other adverse possibilities for transport. If the probe or the product itself has been damaged, do not continue to use it. Please do not use this product for any purposes other than its intended use. When it is used for children, please observe the general safety precautions. We suggest to calibrate the monitor (at least once a year) according to local laws and regulations.

7. COMMON FAILURES AND TROUBLESHOOTING METHODS

Phenomenon of breakdown	Possible cause	Troubleshooting methods
The screen display: "Lo"	The measurement temperature is lower than 32°C (89.6°F), which is beyond the measurement range	Please re-measure following the product instruction
The screen display: "Hi"	The measurement temperature is higher than 43°C (109.4°F), which is beyond the measurement range	Please re-measure following the product instruction
The screen display: "E.F."	The environment temperature is too high or too low	Please measure under the 10°C (50°F) to 40°C (104°F) environment temperature
The screen display: "!"	The power of cell is shortage	Change the cell
No display or abnormal display	The positive and negative poles of the cell are reversed or the device is abnormal	Re-install the cell, or contact the dealer

9. PRODUCT TECHNICAL PARAMETERS

- Power source: DC 3V (2 section 7 alkaline batteries)
- Measurement range: 32.0°C ~ 43.0°F (89.6°C ~ 109.4°F)
- Display resolution: 0.1°C (0.1°F)
- Measurement accuracy: ±0.2°C (±0.4°F) in the range of 35.0°C ~ 42.0°C (95.0°F ~ 107.6°F) ±0.3°C (±0.5°F) in the range of 32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) and 42.1°C ~ 43.0°C (107.8°F ~ 109.4°F)
- Electric shock protection: the device that cannot be used in the presence of inflammable anesthetic gases mixed with air or oxygen or nitrous oxide
- Product size: 149x38x42(mm)
- Product weight: about 68g (not including cell)
- Measuring units: °C / °F
- Measuring site: forehead
- Time interval of each measurement: ≤1s.
- measurement time: ≤1s.
- Degrees of protection against ingress liquid: IP22
- Safety classification: the device that cannot be used in the presence of inflammable anesthetic gases mixed with air or oxygen or nitrous oxide
- Service life: 5 years
- The clinical thermometer is a direct mode
- Battery replacement cycle: use the new battery for no less than 3000 measurements

10. TEMPERATURE SENSE
 The normal body temperature of people is a range, different people's normal body temperature varies, and individual body temperature will change at different times.
 The normal body temperature of most people is provided as follows, just for reference. (The measurement of Chicco thermometer refers to axillary.)
 Axillary temperature: 36.0°C ~ 37.4°C / 96.8°F ~ 99.3°F
 Oral temperature: 36.3°C ~ 37.2°C / 97.34°F ~ 98.96°F
 Rectal temperature: 36.9°C ~ 37.9°C / 98.42°F ~ 100.22°F

11. REPLACING THE BATTERIES
 1. Battery replacement: When "!" symbol is displayed, Replace with 2 new AAA batteries, slide open the battery cover and remove old batteries. Replace the batteries being sure to align properly as indicated inside the battery compartment.
 2. Remove the battery from the product if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.
 3. To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

12. ELECTROMAGNETIC COMPATIBILITY INFORMATION
 a) Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
 b) This equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
 c) When the instrument is in use, never put it near other instruments or stack it on other instrument if you have to put it near other instrument or stacked other instruments, please inspect and verify if the instrument could run normally.

There is the potential risk of radio frequency interference between the device and other devices, if there is, please find out the problems and take the following measures:
 (1) Turn off the device, and turn on again.
 (2) Change the direction of the device.
 (3) Keep the product away from the interferences devices.

Table 1 - Compliance information for Emission test

Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B

Table 2 - Compliance information for Immunity test

Immunity test	Compliance level
Electrostatic discharge(ESD) IEC 61000-4-2	Contact: ±8 kV Air: ±2kV, ±4kV, ±8kV, ±15kV
Power frequency(50, 60Hz) magnetic field IEC 61000-4-8	30A/m 50Hz and60 Hz
Radiated RF IEC 61000-4-3	10V/m 80MHz ~ 2.7GHz 80% AM at 1 kHz

Table 3 - Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{b)}	Modulation ^{c)}
385	380 - 390	TETRA 400	Pulse modulation ^{d)} 18Hz
450	430 - 470	GMSR 460, FRS 460	FM (±) ±5kHz deviation 1kHz sine
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{e)} 217Hz
745	704 - 787	LTE Band 13, 17	Pulse modulation ^{f)} 217Hz
810		GSM 800/900, TETRA 800, CDMA 1900, DECT, LTE Band 1, 3, 4, 2S, UMTS	Pulse modulation ^{g)} 18Hz
870	800 - 960	GSM 800/900, TETRA 800, CDMA 1900, DECT, LTE Band 1, 3, 4, 2S, UMTS	Pulse modulation ^{h)} 18Hz
930		GSM 800/900, TETRA 800, CDMA 1900, DECT, LTE Band 1, 3, 4, 2S, UMTS	Pulse modulation ⁱ⁾ 217Hz
1.845	1.700 - 1.990	GSM 1800, CDMA 1900, DECT, LTE Band 1, 3, 4, 2S, UMTS	Pulse modulation ^{j)} 217Hz
1.970		Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{k)} 217Hz
2.450	2.400 - 2.570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{l)} 217Hz
5.240		Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{m)} 217Hz
5.500	5.100 - 5.800	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ⁿ⁾ 217Hz
5.785		Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{o)} 217Hz

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NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.
 b) The carrier shall be modulated using a 50% duty cycle square wave signal.
 c) As an alternative to FM modulation, 50% pulse modulation at 18Hz may be used because while it does not represent actual modulation, it would be worst case.

13. WARRANTY
Warranty: 12 months
 The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use. The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events.
 For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

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