MOPEDIA®

MA RP881 02 A_ITA_10 2020

IMBRACATURE

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ITALIANO

INDICE

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ENGLISH

CE | Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices**1. CODE****RP881** Standard sling with headrest**2. INTRODUCTION**

Thank you for purchasing a MOPEDIA sling by Moretti S.p.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

3. INTENDED USE

The MOPEDIA by MORETTI slings are designed to transport patient with mobility disease.

- Do not use the product for a purpose not indicated in this manual
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product
- The manufacturer reserves the right to change the information contained in this document without previous notice

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of STANDARD SLINGS complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

5. GENERAL WARNING

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual.
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the sling
- Disinfect the item after every usage
- SERVICE LIFE- the sling use limit is defined by the wear of the parts
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

- Carefully read the patient lift or stand up lifter manual and follow the instruction suggested on the same
- In case of use of harnesses with a different maximum load from the one provided for the lift or stand up lifter, always use the lowest maximum load.

6. SYMBOLS

Product code

Unique Device Identification

CE mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Maximum washing temperature 30°C; very mild process

Rinse very gently, do not wring

Do not bleach

Do not iron


Do not tumble dry

Do not dry clean

7. GENERAL DESCRIPTION

Sling for all kind of patient lift. Made in polyester (100%), especially suitable for water contact, strong and breathable, sewn with polyester thread. The border is made by a strong ribbon (100% nylon), which make three eyelets in every extremity to guarantee a perfect balance. The belt continuous and swen around the entire perimeter of the sling to ensure greater toughness. The seams are made using high strength polyester yarn.

1. Upper eyelet
2. Loops
3. Intermediate eyelet
4. Loops
5. Movement handle
6. Label
7. Movement handle
8. Lower eyelet
9. Loops

REF. RP881

8. BEFORE USE

Always verify the sling is in perfect condition, carefully check the seams and the loops.

9. USING WARNING

- Before use, make sure the device is not damaged: the tissue, the border, the sewing and suspension loops must be intact
- Periodically check for any damage or usury, every 6 month at least, in particular if you often clean the sling
- Ensure the loops are correctly hooked to the cradle
- Store and use away from flame or source of heat
- Do not exceed the maximum load

- Always consult the user manual of the patient lift or stand on which install the harness and strictly follow the instructions given

- Evaluate the risk, in order to guarantee the exact body support unit for the patient, according to the correct dimensions, typology and shape for him
- In case of use of sling with a different maximum load from the one provided for the lift or stand up lifter, always use the lowest maximum load

10. HOW TO USE

The MOPEDIA slings by Moretti can fit all kind of patient lift, fixed and movable.


Fig.1


Fig.2

Fig.3

Fig.4

10.1 From a sitting position

1. Ensure that the straps face outward. Pull the sling down until the user's back leaving the opening of the sling at the base of the spine.
2. Check that the sling is aligned properly at shoulder height. (Fig.1)
3. Lift the patient's leg and slip the strap under it and pull it up between their legs, ensuring that the sling is not creased under the thigh. (Fig.2)
4. Repeat the procedure with the other leg.
5. Cross the straps over the legs and attach them to the front hooks on the hoist spreader bar.
6. Raise the hoist slightly and attach the shorter, shoulder straps to the spreader bar. Move the hoist towards the patient while doing this. (Fig.3)

ATTENTION! The long straps are for the legs, and the short ones for the shoulders

10.2 From a lying position

1. Roll the patient aside and position the sling on the bed (Fig.4) (the sling has to be smooth without any wrinkle and it is aligned properly at shoulder height)
2. Make sure that the opening at the bottom of the sling is positioned at the base of the spinal column.
3. Roll the patient on his/her back again (supine position)
4. Cross the bands under the patient's legs and adjust the sling as shown for the sitting position (Fig.2), maintaining the shoulders's bands as short as possible.5. You are now ready to hang the straps to the cradle (Fig.3).

NOTE: Moretti recommends that two assistants (one on each side of the bed) be used when positioning the patient in a sling.

ATTENTION! Ensure that the patient's head is supported at all times

11. MAINTENANCE

The devices MOPEDIA by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 6 months.

12. CLEANING AND DISINFECTION

Maximum washing temperature 30°C; very mild process

Rinse very gently, do not wring

Do not bleach

Do not iron

Do not tumble dry

Do not dry clean

13. CONDITIONS OF DISPOSAL

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the sling in the appropriate disposal areas for recycling.


14. TECHNICAL FEATURES

REF. RP881

MATERIAL 100% Polyester and Webbing 100% nylon (Flame retardant fabric EN 1021)

MAX WEIGHT 180 Kg

Size	A	B	C	D
ONE SIZE	143 cm	80 cm	93 cm	35 cm

15. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to.

The correct intended use is specified in this manual.

Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use.

Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

16. EXEMPT CLAUSES

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

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