ITALIANO

MA LAD670-LAD672 03 A_ITA_10-2020

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ENGLISH

MA LAD670-LAD672 03_A_ENG_10-2020

DOMUS 3 ALTERNATING PRESSURE REDISTRIBUTION SYSTEM

INSTRUCTION MANUAL

INDEX

1. CODES	PAG.3
2. INTRODUCTION	PAG.3
3. INTENDED USE	PAG.3
4. DECLARATION OF CONFORMITY	PAG.4
5. STANDARDS	PAG.4
6. GENERAL WARNINGS	PAG.4
6.1 Electromagnetic compatibility warnings.....	pag. 4
7. SYMBOLS	PAG.5
8. KIT CONTENTS	PAG.6
9. GENERAL DESCRIPTION	PAG.6
9.1 Pump and mattress system	pag. 6
9.2 Pump unit.....	pag. 7
9.3 Opertions.....	pag. 7
9.3.1 Pressure Adjustment Knob	pag. 8
9.3.2 Normal Pressure Indicator	pag. 8
9.3.3 Low Pressure Indicator	pag. 8
9.3.4 Alternate / Static Switch	pag. 9
9.3.5 Alarm Mute.....	pag. 9
9.4 Emergency CPR Operations.....	pag. 9
10. INSTALLATION	PAG.9
10.1 Pump & Mattress Installation	pag. 10
11. BEFORE USE	PAG.11
12. USING WARNINGS	PAG.11
13. CLEANING	PAG.12
14. STORAGE	PAG.13
15. MAINTENANCE	PAG.13
15.1 General maintenance	pag. 13
15.2 Air filter replacement.....	pag. 13
16. CONDITIONS OF DISPOSAL	PAG.14
17. SPARE PARTS/ACCESSORIES	PAG.14
18. TROUBLESHOOTING	PAG.14
19. GUIDANCE AND MANUFACTURER'S DECLARATION	PAG.14
19.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions.....	pag. 14
19.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity.....	pag. 15
19.3 Recommended separation distances between portable and mobile RF communications equipment and this device	pag. 16
20. TECHICAL FEATURES	PAG.17
20.1 Pump.....	pag. 17
20.2 Mattress.....	pag. 18
21. WARRANTY	PAG.18
22. REPAIR	PAG.18
22.1 Warranty repairs	pag. 18
22.2 Repair not covered by warranty.....	pag. 18
22.3 Non-defective device.....	pag. 18
23. SPARE PARTS	PAG.18
24. EXEMPT CLAUSES	PAG.19

CE I Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices**1. CODES****2 OUTLET SYSTEM**

LAD670-2 Kit mattress and pump with pressure adjustment Domus 3

LAD672-2 Kit mattress and pump with pressure adjustment Domus 3

LAD675-2 Kit mattress and pump with pressure adjustment Domus 3

3 OUTLET SYSTEM

LAD670-3 Kit mattress and pump with pressure adjustment Domus 3

LAD675-3 Kit mattress and pump with pressure adjustment Domus 3

2. INTRODUCTION

Thank you for purchasing a LEVITAS mattress and pump Domus Kit by Moretti. The pump and mattress kits by Moretti are designed and manufactured to meet all your needs for a practical, correct and safe use. This manual contains small suggestions for proper use of the device you have chosen and valuable advice for your security. You may read the complete manual before using the kit pump and mattress, in case of doubts please contact your dealer, who will help and advise you properly.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

3. INTENDED USE

The kit pump and mattress have been designed for the bedsore treatment and prevention, which may occur during long hospital and/or house confinement.

**CAUTION!**

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bedsore, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of DOMUS ANTIDECUBITUS KIT complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

5. STANDARDS

The system has been tested and successfully approved to the following standards:

- EN 60601-1
- EN 60601-1-2
- EN 55011 Classe B
- EN 61000-3-2 Classe A
- EN 61000-3-3

6. GENERAL WARNINGS

- For proper use of the device, refer to the following manual carefully
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

6.1 Electromagnetic compatibility warnings

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

7. SYMBOLS

	Product code
	Unique Device Identification
	CE mark
	Manufacturer
	Batch Lot
	Read the instruction manual
	Medical Device
	Conditions of disposal
	Attention
	Production dates
	Waste Electrical & Electronic Equipment (WEEE) UE/19/2012
	Temperature limitation
	BF Type
	Class II equipment
IP21 	Protected against solid foreign objects of 12,5 mm and greater. Protection against vertically falling water drops.
	Professional dry cleaning; normal process
	Do not iron
	Tumble dry

Do not tumble dry

Do not bleach

Do not dry clean

Maximum washing temperature 60°C; normal process

Maximum washing temperature 95°C; normal process

8. KIT CONTENTS

The product can only be operated by personnels who are qualified to perform general nursing procedures and has received adequate training in knowledge of prevention and treatment of pressure ulcer.

This product has been designed for:

- Help reduce the incidence of pressure ulcers while optimizing the patient comfort
- Long-term home care for patients with pressure ulcers
- Pain management according to doctor's prescriptions

⚠ NOTE: This device is not suitable for use in the presence of flammable anesthetic mixture with air or with pure oxygen or nitrous oxide.

2 OUTLET SYSTEM

LAD670-2 Kit mattress and pump with pressure adjustment Domus 3

LAD672-2 Kit mattress and pump with pressure adjustment Domus 3

LAD675-2 Kit mattress and pump with pressure adjustment Domus 3


3 OUTLET SYSTEM

LAD670-3 Kit mattress and pump with pressure adjustment Domus 3

LAD675-3 Kit mattress and pump with pressure adjustment Domus 3


9. GENERAL DESCRIPTION

9.1 Pump and mattress system


1. CPR
2. Pump unit
3. Mattress
4. Quick Connector

9.2 Pump unit

Front


- 1. Quick connectors
- 2. Front panel


Rear

- 1. Power switch
- 2. Power cord
- 3. Mounting brackets
- 4. Air filter


9.3 Operations

2 OUTLET WORKING SYSTEM

3 OUTLET WORKING SYSTEM

9.3.1 Pressure Adjustment Knob

The pressure adjustment knob controls the air pressure in the mattress. Turning the knob clockwise will increase the pressure; counter-clockwise decreases the pressure. Higher pressures will support heavier patients. The pressure should be adjusted according to individual comfort preferences.

When the mattress is finished pressurizing to the desired level of firmness, the Normal Pressure LED indicator will illuminate.

⚠ NOTE: You can also check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. The air cells will alternately inflate and deflate. You should feel slight contact with the buttocks when the air cells beneath the buttocks deflate.

⚠ NOTE: Every time the mattress is initialized (inflated), it is recommended that you set the pressure knob to "Max" to hasten inflation. You can then adjust the air mattress to the desired firmness later.

9.3.2 Normal Pressure Indicator

When the green LED illuminates, the pressure inside of air mattress has reached the desired pressure setting (preset on the Pressure Adjustment Knob).

9.3.3 Low Pressure Indicator

When the mattress is initializing, the Low Pressure LED will illuminate until the appropriate pressure is reached (according to the Pressure Adjustment Knob); this is normal. Otherwise, the Low Pressure LED is a warning, indicating that the pressure in mattress is unusually low. Check to ensure all connections are secured and correctly installed as per the installation instructions.

⚠ NOTE: If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

9.3.4 Alternate / Static Switch


The ALTERNATE/STATIC switch selects between Alternate Pressure mode and Static Pressure mode. With Alternate Pressure mode, alternating air cells are partially deflated and inflated, avoiding prolonged pressure on any single point beneath the patient; this is to prevent pressure ulcers. With Static Pressure mode, all of the air cells are equally inflated.

9.3.5 Alarm Mute

When low-pressure situation occurred, both the LED light and buzzer will alarm to warn the patient. By pressing the button, it will mute the buzzer so the care giver may check for possible air leaks.

△ NOTE: To re-activate the audible alarm function, please press "Alarm Mute" 

9.4 Emergency CPR Operations

CPR must be performed on a firm surface. Therefore, if an emergency CPR situation occurs with the patient on the mattress, the mattress must be quickly deflated. This is done by quickly pulling on the CPR tag located at the head of the mattress on the patient's right-side. The quick connector from the pump unit can also be disconnected to fasten deflation.

10. INSTALLATION

Unpack the box and check the package contents for completeness.


Package content List:

- Mattress unit x 1 (may not be included if only pump unit is purchased)
- Pump Unit x 1
- User Manual x 1

Inspect equipment for damage which may have occurred during shipment. If there is damage, please contact your dealer immediately

10.1 Pump & Mattress Installation

ENGLISH


1. Place the mattress or pad on top of a bed frame. Orient the mattress so that the foot-markers on the mattress are situated towards the end of the bed (where pump will be located).


WARNING:

The overlay mattress must be applied on the underlying mattress.


2. Place the pump hangers over the bed railing at the foot of the bed. Support the pump from the bottom and turn the hangers outward (or inward) to secure the pump against the railing. If no bed railing is available, place the pump on a flat, sturdy surface or on the floor, beneath the bed.

3. Connect the Quick Connector on the air hose (from air mattress) to the pump unit Quick Connector Slot. When you hear a "click", the connection will be secured.

NOTE: Check and ensure the air hoses are not kinked or tucked under the mattress.

4. Plug the power cord into a electrical outlet with the appropriate voltage

NOTE: The pump unit must be suitable for the local power voltage.

5. Turn the main power switch (on pump) to the ON position. After you turn on power to the pump, you will have to wait at least 20 minutes for the mattress to become fully inflated.

The low pressure indicator (yellow LED) will illuminate while the mattress is pressurizing.

When the appropriate pressure is reached, the Low Pressure indicator (yellow LED) will go off and the Normal Pressure indicator (green LED) will illuminate. The bed is now ready for use.

 NOTE: You can unplug the unit to disconnect the device. Do not position the equipment so that it is difficult to operate the disconnecting device.

CAUTION: The pump supplied with the mattress can only be used for mattresses recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)

 NOTE: (For models WITHOUT the low-air-loss function) During power outages, you can cover the Quick Connector with the transport cap to maintain air pressure inside the air cells.

 NOTE: After installation, make sure any excess cord or tubing is removed from possible foot traffic to avoid accidental tripping. All EQUIPMENT should be positioned to always allow unhindered patient access by physicians and caregivers.

11. BEFORE USE

Check the correct assembly, paying particular attention to air tubing connection between pump & mattress

12. USING WARNINGS

DANGER - To reduce risk of electrocutione

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

WARNINGS - To reduce risk of burns, electrocution, fine or injury to people

1. Never leave this product unsupervised when connected to the mains
2. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
3. This system is not intended to be used with patients who have spinal cord injury.
4. Close supervision is necessary when this product is used on or near children.
Electrical burns or shocking accident may occur if a child swallow a small part detached from the device
5. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer
6. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
7. Keep the cord away from heated surfaces
8. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
9. Never drop or insert any object into any opening or hose
10. Do not modify this equipment without authorization of the manufacturer
11. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
12. Do not leave excessively long tubes in the upper area of the bed, to avoid the possibility of strangulation

**CAUTION!**

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3,3m) between devices or turn off the mobile phone.

13. CLEANING

It is important to follow the cleaning procedures before using the equipment on human bodies.

Wipe the pump unit with a damp cloth and a mild detergent. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit. Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

 CAUTION: Do not immerse or soak pump unit.

 CAUTION: Do not use phenolic based product for cleaning.

 CAUTION: After cleaning, dry the mattress without direct exposure sunlight.

**WARNING!**

Unplug the pump before starting the cleaning procedure

Cover fabric: Nylon/PU
Cover fabric: Silver + Nylon/PU

Do not iron

Tumble dry

Do not bleach

Do not dry clean

Maximum washing temperature 60°C; normal process

Cover fabric: Stretch

Do not iron

-  Tumble dry
-  Do not bleach
-  Professional dry cleaning; normal process
-  Maximum washing temperature 95°C; normal process

14. STORAGE

1. Lay the mattress on a flat surface and upside down
2. Roll-up the mattress from the head end towards the foot end
3. Foot-end strap can then be stretched around the rolled pad/ mattress to prevent unrolling

 **CAUTION:** Do not fold, crease or stack the mattresses.

15. MAINTENANCE


The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by the manufacturer, who recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function. Otherwise, service and inspection of the devices generally should be done every 2 years.

15.1 General maintenance

1. Check main power cord and do not plug it if there is an abrasion or excessive wear.
2. Check mattress cover for signs of wear or damage.
3. Disconnect the air tube from the mattress. Then check the airflow coming from the two air outlets on the pump. They should be alternately delivering air when the pump is set at "alternate" mode.
4. Check the air hoses if there is kink or breaks. For replacement, please local distributor.

15.2 Air filter replacement

1. Open the air filter cap located at the back of pump.
2. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
3. Check and replace air filter regularly if environment is dirty.

16. CONDITIONS OF DISPOSAL

(Directive 2012/19/UE) 

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin

17. SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to Moretti main catalogue.

18. TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem.

PROBLEM	SOLUTION
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains
Alarm is ON	<ul style="list-style-type: none"> • Check if the power is suddenly shut down. • Check if the CPR is sealed. • Check if the connection between air tube connector to pump unit is tightly secured. • Check if all tubing connections along mattress are secured
Patient is bottoming out	<ul style="list-style-type: none"> • Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort.
The mattress doesn't inflate	<ul style="list-style-type: none"> • Check if the connection between mattress and pump unit is tightly secured • Check if there is any leakage on air cells
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened • Check if the mattress is fixed to the bed frame by straps
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> • This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time

19. GUIDANCE AND MANUFACTURER'S DECLARATION

19.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<p>⚠ Warning:</p> <p>1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.</p> <p>2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.</p> <p>3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p>		

19.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	+8kV contact +15kV air		+8kV contact +15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	+2kV for power supply line +1kV for input/output line		+2kV for power supply line +1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	+ 1 kV line(s) to line(s) + 2 kV line(s) to earth	+ 1 kV line(s) to line(s)	+ 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

ENGLISH

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $=\sqrt{P}$ 150kHz to 80MHz $=0.6\sqrt{P}$ 80MHz to 800MHz $=1.2\sqrt{P}$ 800 MHz to 2.7G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10V/m	
NOTE 1: UT is the a.c. mains voltage prior to the application of the test level NOTE2: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people				
a)Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.				

19.3 Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz =√P	80 MHz to 800 MHz =0.6√P	800 MHz to 2,7 GHz d =1.2√P
0,01	0,1	0,06	0,12
0,1	0,31	0,19	0,38
1	1	0,6	1,2
10	3,1	1,9	3,8
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
100	10	6	12

20. TECHICAL FEATURES

20.1 Pump

FEATURES	LAD670-2 / LAD670-3 / LAD672-2 / LAD675-2 / LAD675-3
Power Supply	230V/50Hz
Fuse T1AL	250V
Power cord	4,8 mt
Pressure output	20~30~40mmHg
Air release	>8LPM
Cycle time fixed	9.6 min
Dimensions	29x18,5x12,6 cm
Weight	2,2 Kg
Environment Temperature	Operation: +10°C / +40°C Storage: -15°C/+50°C Shipping: -15°C/+70°C
Environment Humidity	Operation: 10%/90%non-condensing Storage: 10%/90% non-condensing Shipping: 10%/90% non-condensing
Environment Atmospheric Pressure	Operation: 700/1013,25 hPa Storage: 50-106 kPa Shipping: 50-106 kPa
Classification	Class II, Type BF, IP21 Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

20.2 Mattress

SPECIFICATION	LAD670-2/ LAD670-3	LAD672-2	LAD675-2/ LAD675-3
Model	12.7 cm Mattress	5" Mattress +3" Foam pocket	20.3 cm Mattress
Dimension	200x90x12.7 cm	200x90x12.7 cm	200x90x20.3 cm
Material element	Nylon/TPU		
Material cover	Nylon/PU		
Weight	5.2 kg	5.6 Kg	6.9 Kg
Max user's weight	180 Kg	200 Kg	200 Kg
Pressure range	SOFT to FIRM setting		

21. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for a period of (two) 2 years from the purchasing date, subject to any limitations and exclusions specified below. The warranty does not apply in case of damages caused by improper use, abuse or modification of the product, and the warranty will not be valid if the instructions for use are not strictly followed. The correct intended use of the device is specified in this manual. Moretti is not responsible for damages resulting from personal injury or otherwise caused by or in relation to wrong installation or improper use. Moretti does not guarantee the products against damage or defects in the following conditions: natural disaster, not authorized maintenance or repair, damage due to power supply problems (where provided) use of spare parts not supplied by Moretti, improper use, not authorized alteration, shipping damage (other than original shipment by Moretti), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

22. REPAIR**22.1 Warranty repairs**

In the event that a product Moretti to be defective in material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, in its sole discretion, can repair or replace the item under warranty at a Moretti dealer or at Moretti location . The labor cost for the repair may be charged to Moretti if the repair is covered by warranty. Repair or replacement does not renew or extend the warranty

22.2 Repair not covered by warranty

A product not covered by warranty can be returned for repairing only after receiving prior authorization from Moretti customer service. The labor and shipping costs for a repair not covered by warranty will be fully borne by the customer or dealer. Repairs of products not covered by warranty are warranted for six (6) months from the date you receive the product repaired.

22.3 Non-defective device

The customer will be informed if, after having examined and tested returned product, Moretti claims that the product is not defective. In this case the good will be returned to the customer, all shipping costs will be paid by customer.

23. SPARE PARTS

Moretti original spare parts are guaranteed for 6 months from the day you receive the replacement.

24. EXEMPT CLAUSES

Except as expressly specified in this warranty and within the limits of the law, Moretti does not offer any representation, warranty or condition, express or implied, including any representations, warranties or conditions of merchantability, fitness for a particular purpose, non-infringement and interference. Moretti does not warrant that the use of the product Moretti will be uninterrupted or error free. The duration of any implied warranties that may be imposed by law are limited to the warranty period, within the limits of the law. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitation of incidental or consequential damages in connection with products for consumers. In such states or countries, certain disclaimers and limitations of this warranty may not apply to you. This warranty is subject to change without notice.

ENGLISH**WARRANTY CERTIFICATE****Product** _____**Purchased on (date)** _____**Retailer** _____**Address** _____ **Town/city** _____**Sold to** _____**Address** _____ **Town/city** _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: **info@morettispa.com****MADE IN P.R.C.**

MORETTI S.P.A.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com