**LEM**

ITALIANO

MA LTM560 01F\_ITA\_11-2020

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Venduto a \_\_\_\_\_

Via \_\_\_\_\_ Località \_\_\_\_\_

Xiamen Senyang Co Ltd Factory: 4/F, Xingbei Industry, No.95-99,  
West 2 Road, Jiu Tian Hu, Xing Lin, Xiamen, 361022, China

EC

REP

Dromos Europe S.L.

San Francisco 44-1, 48003 Bilbao, Vizcaya,

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**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

[www.morettispa.com](http://www.morettispa.com) email: [info@morettispa.com](mailto:info@morettispa.com)**MADE IN P.R.C.**

**LEM**

ENGLISH

MA LTM560 01F\_ENG\_11-2020

**COMPRESSION THERAPY WITH MODULATING  
COMPRESSOR AND LCD DISPLAY**

---

**INSTRUCTION MANUAL**

## INDEX

1. CODES.....	PAG.3
2. INTRODUCTION.....	PAG.3
3. INTENDED USE.....	PAG.3
4. EC DECLARATION OF CONFORMITY .....	PAG.3
5. GENERAL CONSIDERATIONS.....	PAG.4
5.1 Lymphatic system.....	pag. 5
5.2 Field of application .....	pag. 5
6. GENERAL WARNING.....	PAG.6
7. SYMBOLS USED.....	PAG.7
8. COMPRESSION THERAPY SYSTEM PARTS.....	PAG.8
9. GENERAL DESCRIPTION .....	PAG.9
10. CONTENTS.....	PAG.10
11. BEFORE USE .....	PAG.11
12. HOW TO USE.....	PAG.11
13. MAINTENANCE .....	PAG.12
14. CLEANING AND DISINFECTION.....	PAG.12
15. CONDITION OF DISPOSAL .....	PAG.13
15.1 General conditions of disposal .....	pag. 13
15.2 Information for the correct disposal of the product in accordance with European Directive 2012/19/UE: .....	pag. 13
16. TROUBLESHOOTING.....	PAG.13
17. SPARE PARTS AND ACCESSORIES .....	PAG.13
18. ELECTROMAGNETIC COMPATIBILITY DECLARATION .....	PAG.14
19. TECHNICAL FEATURES .....	PAG.17
20. WARRANTY .....	PAG.18
21. REPAIRING .....	PAG.18
21.1 Warranty repair .....	pag. 18
21.2 Repair not covered by warranty.....	pag. 18
21.3 Non defective device .....	pag. 18
22. SPARE PARTS.....	PAG.18
23. EXEMPT CLAUSES .....	PAG.19

0197 IIa Class Medical Device

According to directive 93/42/EEC and further modifications

## 1. CODES

LTM560	Modulable therapy pump with LCD display
LTM563-S	Sleeve S/M - Ambidextrous
LTM563-L	Sleeve L/XL - Ambidextrous
LTM566-S	Boot S/M - Ambidextrous
LTM566-L	Boot L/XL - Ambidextrous
LTM566-XXL	Boot XXL - Ambidextrous
LTM566-3XL	Boot 3XL - Ambidextrous
LTM568	Extension
LTM569-S	Abdominal band S
LTM569-L	Abdominal band L

## 2. INTRODUCTION

Thank you for purchasing a compression therapy Advance LEM by Moretti S.p A. designed and built to meet all your needs for a correct, safe and practical use. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, please contact the dealer, who will help and advise you properly.

## 3. INTENDED USE

Advance compression therapy is a device intended for patients who need a lymphatic drainage in the limbs, i.e. it increases the lymph drainage from congested areas facilitating the body's natural skill to eliminate toxins.

### WARNING!

**The pressure therapy is a medical device, in the presence of any health problem use the device only after consulting your doctor.**

- Don't use the device for a purpose not indicated in the present manual.
- Moretti S.p.A. declines all responsibilities for any consequences resulting from improper use of the device or from an unauthorized alteration of the device.
- The manufacturer reserves the right to make changes to the device and to this manual without notice in order to improve its features.

## 4. EC DECLARATION OF CONFORMITY

Moretti S.p.A. Declares under his sole responsibility that the product marketed by the same Moretti S.p.A. Is in accordance to the relevant provisions of Directive 93/42/EEC on medical devices and subsequent amendments and are produced by xiamen senyang co ltd according to the standards set out in the below standards:

### Standards and directives

The pressure therapy Advance has been tested and approved according to the following guidelines and standards:

**Safety standard EN 60601-1** Medical electrical equipment part 1: General requirements for safety

**EMC standard EN 60601-1-2** - Medical electrical equipment part 1: General requirements for safety- Collateral standard: Electromagnetic compatibility - Requirements and tests

**EN 60601-1-11** Medical electrical equipment part 1: General requirements for basic safety and essential performance

**Electromagnetic compatibility warnings**

This device has been tested and proved in conformity with the limits for medical devices of EN 60601-1-2: 2007. These limits are designed to provide reasonable protection against harmful interference in a typical medical application. This equipment generates, uses, and can radiate radio frequency energy and if not installed and used in accordance with the instructions, it may cause harmful interference to other devices in the vicinity. However there are no total guarantees that the electromagnetic interference may occur in special circumstances. If this device does cause harmful interference to other devices, which can be determined by tuning off, the user is encouraged to follow the following steps:

- Reorient or relocate the receiving device
- Increase the distance between the devices
- Connect the device to a plug or circuit different from other devices
- Consult the manufacturer or an authorized service center

- Do not use mobile phones and other devices in the vicinity of the device. They can generate strong electrical or electromagnetic fields as they can interfere with the operation of the device.
- Please do not use in an area with strong electromagnetic radiation, and during the use of this to keep at least a distance of 40cm from the user's cardiac pacemaker.

**5. GENERAL CONSIDERATIONS**

The pressotherapy is a lymphatic massage which operates on the venous and lymphatic circulation. Its purpose is to normalize the venous circulation, avoiding lymphatic stasis and increasing the peripheral microcirculation.

The pressure sequence improves the lymph flow from the interstitial spaces to the lymph nets, helping to promote the body's natural toxin clearing functions. This increases the venous flow (blood flow to the heart).

Although its beneficial effects, the pressotherapy is not indicated for everyone. In certain states of disease there may be an increased risk, for this reason it is recommended to consult a doctor before undergoing the therapy.


Among the possible uses please note:

- Normalize the lymphatic system.
- Body shaping and profiling
- To give support to the doctor in the treatment and prevention of vascular diseases
- Increase Endurance and Oxygenation for Sport Performance
- To relieve discomforts such as: swollen and heavy feet, ankles and legs, water retention, cramps, tired legs, muscle fatigue.
- Re-absorption of edema.

**CAUTION!**

While representing a seemingly simple therapeutic aid and often proposed and practiced by non-medical professionals (beauty salons), please note that it is essential that this treatment is prescribed and implemented on the basis of precise instructions and after a careful clinical examination (visit) to identify the frequency, the type, the duration of treatment and the possible contraindications. It is therefore recommended to always check with your doctor's physical and health condition before undergoing sessions of pressure therapy.

### 5.1 Lymphatic system

The lymphatic system is a complex network composed of lymphatic vessels (the most important of which is the thoracic duct), and lymphoid tissue. The lymphoid tissue is found in the lymph nodes, the so-called MALT (mucose-associated lymphoid tissue) and in the structures dedicated to the production and putting into circulation of lymphocytes, such as bone marrow and thymus.

In lymph vessels flows the sap, which is drained by the peripheral tissues and poured into the venae cavae.

The lymphatic system is an important component of the immune system and has many functions in relation with each other:

- Remove excess water from the tissues;
- Absorption of fatty acids and subsequent transport of fat into the circulatory system;
- Production and transport of immune cells such as lymphocytes, monocytes;
- Transport of Antigen-Presenting Cell from the tissues to the lymph nodes, the site of stimulation of the immune response;
- Production of antibodies;

### 5.2 Field of application

The pressure therapy is a treatment based on the application of external pressure on the body according to a specified gradient and with a sequence programmed using an electrical device and accessories worn by the patient. The inflatable garments (boots/sleeves/waist) are made of air chambers. The pump fills the air chambers in order to pressurize the tissues in the limbs and force fluids out of the pressurized area. The pressure sequence starts from the periphery of the lymphatic and blood vessels and progressively activates the circulation.

It is very important to consult a doctor before undergoing the pressotherapy, because its use could damage your health. Please note that more benefits do not always correspond to a higher pressure. The maximum pressure used should not exceed the minimal arterial pressure of the patient in order not to stop the oxygenated blood flow to the tissues.

It is important to unblock the lymph nodes manually before each treatment (starting draining from the superclavicular, then to the inguinal and finally to the popliteal).

**The use of pressotherapy is indicated in many diseases such as:**

- Lymphedema either congenital or postoperative
- Edema by water retention
- Edema for venous insufficiency
- Primary and secondary lymphedema
- Post-mastectomy lymphedema
- Lipedema
- Cellulite (adiposis edematosa)
- Venous or lymphatic edema in pregnancy
- Post-thrombotic syndrome or post-phlebitis
- Prevention of venous thrombosis before and after surgery

Here are some examples of the above indications where the pressure therapy is recommended (indicative list):

- Lymphatic drainage: legs - arms - swollen feet or legs
- Lymphedemas
- Edema because of immobility
- Hematoma after cosmetic surgery treatment
- Vascularisation
- Drainage
- Cellulite (adiposis edematosa)
- Vascularisation: legs - buttocks - abdomen - arm
- Water retention
- Relax and preventing cramps
- Fast recovery
- Slow deep recovery
- Heavy legs after exercise
- Muscle fatigue

These are the cases in which the pressure therapy is contraindicated:

- Acute vein inflammation (phlebitis)
- Acute lymphangitis
- Unstabilized deep vein thrombosis
- Dermatitis and acute dermatocellulitis
- Erysipelas
- Arterial
- Cardiac and/or respiratory failure
- Not controlled high blood pressure
- Chronic diseases
- Malignancies
- Women who are pregnant (for waist garment )

#### CAUTION !

- Consult your doctor before use.
- Before use read the manual carefully in all its parts, keep this manual near the pressure therapy, so you can consult it at any time
- Do not use the product on people suffering from inflammation or wounds at risk of bleeding both internal and external, do not use the product on people suffering from purulent inflammation.
- Failing to follow instructions can cause serious injury or damage to the person, up to the danger of death.

#### 6. GENERAL WARNING

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual;
- For proper use of the device always consult a doctor. Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use;
- Keep the packed product away from any heat source; the packaging is made of cardboard;
- Avoid direct exposure of the device to direct sunlight or any other heat source;
- SERVICE LIFE of the device is defined by the wear of the parts;
- Always pay attention to the presence of children;
- Product is designed for its intended use only; Do not abuse it in any way;
- Product is not intended for infants or individuals who cannot express their intentions;
- DO NOT disassemble or attempt to repair;
- The device should only be used by qualified medical personnel;

- This manual has to be totally read and understood before using the device on a patient;
- The device is not suitable for use in presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide;
- Unplug the unit before cleaning. Use a damp cloth only. Do not pour any liquid in the compressor;
- Always place the pump unit on a flat, stable surface before use. If placed on an inclined surface, the unit may not work properly;
- Do not cover the unit during operation; use the power indicated on the back label of the compressor only;
- Do not use a damaged power cord;
- Unplug the unit and contact qualified personnel if any of the below conditions will occur:
  - damaged power cord or plug;
  - liquid has entered in the unit;
  - the external case has been damaged;
  - the unit is not working properly, it operates continuously, it works in a noisy way or with low performance;
- Always turn off the unit before deflating the garments;

**CAUTION!**

- If during use you feel ill or have pain, do not use the equipment. Please immediately stop using.
- When wearing the garments does not move and/or walk.
- Do not use the garments directly on the skin. Wear comfortable dress under the garments.

**7. SYMBOLS USED**

Product code

Unique Device Identification

CE Mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Attention

Production dates

Waste Electrical & Electronic Equipment (WEEE) UE/19/2012

BF Type

Class II equipment

**IPX2**

Protection against solid objects greater than 12 mm- Water drops Protection inclined up to 15°

Temperature limitation

ON

OFF

## 8. COMPRESSION THERAPY SYSTEM PARTS

### Component part

**LTM560**

### Description

4-OUTLET PUMP

**LTM563-S**

SLEEVE SIZE - S/M (70x25 cm) - AMBIDEXTROUS

**LTM563-L**

SLEEVE SIZE L/XL (90x25 cm) - AMBIDEXTROUS

**LTM566-S**

BOOT SIZE - S/M (90x30 cm) - AMBIDEXTROUS

**LTM566-L**

BOOT SIZE - L/XL (110x30 cm) - AMBIDEXTROUS

**LTM566-XXL**

BOOT SIZE - XXL (90x40 cm) - AMBIDEXTROUS

**LTM566-3XL**

BOOT SIZE - 3XL (110x36 cm) - AMBIDEXTROUS

**LTM568**

EXTENSION ONLY FOR BOOT L/XL (100x10 cm)

**LTM569-S**

ABDOMINAL BAND S (120x 40 cm)

**LTM569-L**

ABDOMINAL BAND L (145x 40 cm)

## 9. GENERAL DESCRIPTION

The garments sections are partially overlapping to eliminate any opening between the cells and reduce the possibility of stagnant liquids. The sections inflate and deflate sequentially with the possibility of excluding some chambers. The compression begins from the section the most distant from the heart and moves towards the closest section to help the drainage of the blood. The system pressure and the timer for the treatment can be adjusted using the controls on the display of the unit. Using the buttons on the display you can select **two compression mode A or B** depending on the treatment required.

### Mode A

The compression starts from the end part of the limb and goes to the heart. Once a section inflates, the previous deflates. And so on, repeating the cycle. This mode simulates the massage of a therapist that makes a pressure starting from your feet to the groin with his circle closed hands. The simulation of manual lymphatic drainage is useful for fluid retention and swelling.

Fig.2

### Mode B

The compression starts from the end part of the limb and goes to the heart. In this case the air chambers will deflate all together at the end of the inflation sequence. And so on, repeating the cycle. This kind of massage is excellent to drain stagnant liquids from limbs. Thanks to the homogeneous and continuous pressure stagnation zones are avoided.

Fig.3

#### CAUTION!

- Check that the data indicated on the unit label are always corresponding to those of your mains supply.
- Do not place the device where it could be difficult to disconnect it.


#### DANGER - To reduce the risk of electrocution:

- Disconnect the device from the mains supply immediately after use.
- Do not place or store the product in a place where it could come into contact with water or liquid.
- Do not attempt to remove the product if it fell into the water

## 10. CONTENTS

The contents may differ depending on the components you have chosen

- A. 2 treatment modes: "A" and "B"
- B. 4 compression sections/outputs selectable in accordance with the patient's treatment: 1-2-3-4, each figure represents a section. If that section is switched ON, it means that part will be inflated.
- C. 3 selectable treatment durations: patients can choose 10, 20, 30 minutes depending on their conditions
- D. "ON /OFF" button
- E. "START" / "STOP" button: Press the button to start or stop treatment.
- F. Range / pressure value: there are 10 types of pressure, use "+/-" button for the pressure regulation. If you press "+", the pressure will increase, if you press "-", the pressure will be reduced.

Boot


Fig.5

Abdominal band


Fig.6

Sleeve

Fig.7

Extension


Fig.8

## 11. BEFORE USE

- Check the correct assembly and operation of the device with particular attention to the connecting tubes between the pump and the garments.
- Always check the wear and operation status of the parts to ensure a safe use of the device. Replace them if necessary.
- Make sure the zipper closures of garments are ok.

## 12. HOW TO USE

1. Insert the plug into a standard mains socket 220 V.
2. Connect the air tubing to the unit air output and to the garment. Make sure that all connections are well connected.

**NOTE:** During the treatment it is better to wear a lightweight dress to reduce perspiration. It should be without buttons or zippers that can go in contact with the garments.

3. a. Boot/sleeve: connect the air tubing to the garments according to their length: the short tube must be connected to the extremity of the boot or sleeve next to the pump and the longest tube to the most distant. **So connect the ends of the tubes with the same color connectors (pink, yellow, purple, light blue).**  
b. Abdominal band: connect the tubes to the chambers. Pay attention not to cross the tubes. Connect the connector to the tubes and to the pump.
4. Do not wear: watch, rings and bracelets.
5. Place the garment on the limb to be treated. The patient should feel comfortable. Firmly close the zip or the Velcro (abdominal band). The patient could feel more comfortable not tightening the garment too much, always selecting the pump pressure level indicated by the doctor.
6. Turn the unit on with the ON / OFF button. The display will light up.
7. Select the pressure according to the value recommended by the doctor.
8. Set the time of treatment (10 min - 20 min - 30 min).
9. Select the desired treatment: "A" or "B".
10. If needed, deselect some sections.
11. Press START / STOP button to start the treatment.

### ATTENTION!

Pay attention to the inflation of the chambers: it has to begin from the most distant chamber from the heart to the closest.

12. Press START / STOP button to stop the treatment.

### Boot - AMBIDEXTROUS

- Wear the boot, close the zipper.
- Connect all the chambers to the tubes (the shortest tube on the foot and the rest up to the groin).  
So connect the ends of the tubes with the same color connectors (pink, yellow, purple, light blue).
- Connect the connector to the compressor.
- The inflation must start from the foot and move on to the groin.

### Sleeve - AMBIDEXTROUS

- Wear the sleeve, close the zipper.
- Connect all the chambers to the tubes (the shortest tube on the foot and the rest up to the groin).  
So connect the ends of the tubes with the same color connectors (pink, yellow, purple, light blue).
- Connect the connector to the compressor.
- The inflation must start from the hand and move on towards the shoulder.

**Abdominal band/gluteus**

- Wear the abdominal band, fasten the Velcro belt.
- Connect all the chambers to the tubes. So connect the ends of the tubes with the same color connectors (pink, yellow, purple, light blue).
- Connect the connector to the compressor.
- The inflation must start from the distal chamber and move on up to the proximal.

**NOTE:** Always consult a doctor to decide the suitable treatment for the patient.

**13. MAINTENANCE**

The Advance compression therapy LEM by Moretti is carefully checked and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by approved laboratory every 2 years. In case of repair, use only original spare parts and accessories.

**14. CLEANING AND DISINFECTION****Pump unit**

For the pump unit use only a damp cloth. Pay attention to the display.

**NOTE: DO NOT immerse the device in water!**

**NOTE:** For a deeper cleaning you can use a small quantity of neutral cleanser. Do not let any liquid seep into the device and dry it carefully after cleaning.

- Never use acids, alkalis or solvents such as acetone or thinner
- Avoid dropping the device and slamming it against other objects
- Avoid strong temperature excursions. Keep away from heat sources and direct sunlight.
- Do not disassemble the product

**Garments**

For a better preservation it is recommended to fully deflate the garments after every use. Do not fold the garments too much. To clean the aid and the structure use a solution of alcohol and neutral soap (in equal parts) mixed with water.

**WARNING!**

Do not use abrasive substances, acids, alcohol, chlorine-based cleaners, disinfectants and acetone when cleaning, because these substances cause the abrasion of the plastic surfaces and rusting of metal surfaces

The manufacture company cannot be held responsible for the damage that is caused by using materials that could damage the product surface or corrosive chemicals during cleaning.

If you need to disinfect the product use a common disinfectant neutral detergent.

**To remember:**

1. The garments life depends on their maintenance, as well as on the type and intensity of the treatment.
2. The garments should be replaced when they do not inflate properly
3. Garments must be used on intact skin
4. You must stop treatment and consult your doctor in case of skin irritation.
5. Garments are designed for personal use

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### 15.2 Information for the correct disposal of the product in accordance with European Directive 2012/19/UE:

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin the product bears the barred bin symbol.

## 16. TROUBLESHOOTING

If the following answers to your problems are not sufficient, please contact your dealer for a qualified technical assistance.

### PROBLEM

The pump cannot turn ON

### SOLUTION

- Check whether the plug is connected to the power outlet
- Check that there are no blown fuses

The garments cannot inflate

- Check whether the tube connecting compressor and accessories is well fixed
- Wait a few minutes after power to evaluate the cyclic operation of the compressor

## 17. SPARE PARTS AND ACCESSORIES

For spare parts and accessories only refer to Moretti general catalogue.

<b>LTM563-S</b>	<b>LTM563-L</b>	<b>LTM566-S</b>	<b>LTM566-L</b>	<b>LTM566-XXL</b>	<b>LTM566-3XL</b>	<b>LTM568</b>
Sleeve S/M	Sleeve L/XL	Boot S/M	Boot L/XL	Boot XXL	Boot 3XL	Extension
Ambidextrous	Ambidextrous	Ambidextrous	Ambidextrous	Ambidextrous	Ambidextrous	

<b>LTM569-S</b>	<b>LTM569-L</b>
Abdominal band S	Abdominal band L

### 18. ELECTROMAGNETIC COMPATIBILITY DECLARATION

The Advance devices is intended for use in the electromagnetic environment specified below. The customer or user of the Advance devices must ensure that it is used in such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT- GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC61000-3-3	Complies	

#### Guidance and manufacturer's declaration - Electromagnetic immunity -

The Advance devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	⊕ 6kV contact ⊕ 8kV air	⊕ 6kV contact ⊕ 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	⊕ 2kV for power supply line ⊕ 1kV for input/out line	⊕ 2kV for power supply line ⊕ 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	⊕ 1 kV line(s) to line(s) ⊕ 2 kV line(s) to earth	⊕ 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60% dip in UT)for 5 cycles 70 % UT (30% dip in UT)for 25 cycles <5 % UT (>95 % dip in UT)for 5 sec	<5 % UT (>95% dip in UT)for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
<b>NOTE UT is the a.c. mains voltage prior to the application of the test level</b>			

**Guidance and Manufacturer's Declaration - Electromagnetic Immunity -**

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bandsa	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz $d = 1.2\sqrt{P}$ 150 kHz to 80 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device

NOTE 4 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	Da 150 KHz to a 80MHz d=1.2√P	Da 80MHz to a 800MHz d=1.2√P	800MHz fino a 2.5 GHz d=2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

**19. TECHNICAL FEATURES****Pump unit**

Power source	LTM560 AC 220-240V~50 Hz 0,05A
Power	30 W
Fuse	T2A 250V
Pressure output	20-250 mmHg
Pump unit dimension	26x7x3 cm
Cycle time	10-20-30 minutes
Outlets	4
Weight	2,6 Kg
Ambient temperature	Operating environment +10°C / +40°C Storage environment -10°C/+50°C Transport environment -10°C/+60°C

**Ambient humidity**

Operating environment 20%-90% no condensa  
Storage environment 10%-95% no condensa  
Transport environment 10%-90% no condensa

**Ambient pressure**

Operating environment 86-106KPa

**Classification**

Class II, BF type, IP2X (keep dry)  
Applied part: garments the device is not suitable for use in presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide

	<b>Arm</b>	<b>Arm</b>
	LTM563-S	LTM563- L
Size	S/M	L/XL
Material	TPU	TPU
Number of chambers	4	4
Dimension	70x25 cm	90x25 cm

	<b>Leg</b>	<b>Leg</b>	<b>Leg</b>	<b>Leg</b>
	LTM566-S	LTM566-L	LTM566-XXL	LTM566-3XL
Size	S/M	L/XL	XXL	3XL
Material	TPU	TPU	TPU	TPU
Number of chambers	4	4	4	4
Dimensions	90x30 cm	110x30 cm	90x40 cm	110x36 cm

	<b>Extension</b>	<b>Waist</b>	<b>Waist</b>
	LTM568	LTM569-S	LTM569-L
Material	TPU	TPU	TPU
Number of chambers	-	4	4
Dimensions	100x10 cm	120x40 cm	145x40 cm

## 20. WARRANTY

Every Moretti product is guaranteed for 2 (two) YEARS from the purchase date, except the restrictions specified below. The warranty shall not apply in case of damage caused by improper use, abuse, alterations of the product or any other use not in compliance with the user manual. The proper use of the product is indicated in the user manual. Moretti S.p.A. shall not be liable for any damages or personal injuries due to incorrect assembly/improper use of the device, not compliant with the user manual. Moretti does not guarantee its products against damages or faults in the following cases: natural disasters, non authorized maintenance or reparations, damages caused by electrical supply problems (if provided), use of parts or components not provided by Moretti, failure in accomplishing the guidelines and instructions for use, non authorized alterations, damages due to the transport (different from Moretti original delivery) or failure in maintaining the device as indicated in the manual.

The components subject to wear are not covered by this warranty if the damage is due to the normal use of the product.

### Rechargeable batteries warranty (If provided)

The original batteries and the spare batteries are covered by a 90 days warranty during a period of 6 months, with regards to manufacturing defect or to what concern the law regulation. If charge batteries are unused for more than 3 consecutive months, the warranty is void. If exhausted batteries are unused for more than 3 consecutive days, the warranty is void.

## 21. REPAIRING

### 21.1 Warranty repair

If a Moretti product shows material or manufacturing faults during the warranty period, Moretti S.p.A., together with the client, will consider whether the product fault is covered by warranty. Moretti S.p.A., at its unquestionable discretion, may replace or repair the item under warranty in a Moretti reseller point or at its own qualified offices. Labour costs for repairing the product may be at Moretti's expense if it is determined the reparation is under warranty. A reparation or replacement does not renew nor extend the warranty.

### 21.2 Repair not covered by warranty

A product, not covered by warranty, may be returned for reparation only after receiving prior authorization from Moretti client service. Labour and delivery costs for an out-of-warranty reparation will be fully at the client's or reseller's expense. Reparations of out-of-warranty products are guaranteed for 6 (six) months, from the day the client receives the repaired product.

### 21.3 Non defective device

The client is informed if, after examining and testing a returning product, Moretti decides the product is not defective. The product will be sent back to the client and all the return delivery costs will be at his expense.

## 22. SPARE PARTS

The original Moretti spare parts are guaranteed for 6 (six) months from the day you receive them.

**23. EXEMPT CLAUSES**

Except what clearly specified in this warranty and in accordance with the law, Moretti S.p.A. does not offer its clients any other declaration, warranty or condition, expressed or implicit, including any declaration, warranty or condition of merchantability, of non violation or non interference, of suitability for a specific aim. Moretti S.p.A. does not guarantee that the use of Moretti products will be continuous and faultless. The duration of possible implicit warranties that may be applied in accordance with the law is limited to the warranty period.

Some states or countries do not allow limitations to the duration of an implicit warranty, or the exoneration, or the limitations of accidental or indirect damages related to products for consumers. In such states or countries some exonerations or limitation of this warranty may not be applied to the user. This warranty is subject to modifications without any prior notice.

**WARRANTY CERTIFICATE**
**Product** \_\_\_\_\_

**Purchased on (date)** \_\_\_\_\_

**Retailer** \_\_\_\_\_

**Address** \_\_\_\_\_ **Town/city** \_\_\_\_\_

**Sold to** \_\_\_\_\_

**Address** \_\_\_\_\_ **Town/city** \_\_\_\_\_

Xiamen Senyang Co Ltd Factory: 4/F, Xingbei Industry, No.95-99,  
West 2 Road, Jiu Tian Hu, Xing Lin, Xiamen, 361022, China

EC

REP

Dromos Europe S.L.

San Francisco 44-1, 48003 Bilbao, Vizcaya,

product distributed by:

**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

**www.morettispa.com** email: **info@morettispa.com****MADE IN P.R.C.**

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