MOPEDIA®

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Vendido a _____

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Via Bruxelles, 3 - Melegnano 20022 Cavriaglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com - email: info@morettispa.com**MADE IN CHINA.**

MORETTI S.p.A.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com

Hospital beds three-jointed electric beds Gardenia 2.0
INSTRUCTION MANUAL

INDICE

1. CODES.....	3
2. INTRODUCTION.....	3
3. INTENDED USE	3
3.1 Intended use	3
3.2 Application environment	3
3.3 Warnings.....	3
3.4 Unauthorized use	4
4. DECLARATION OF CONFORMITY UE.....	4
4.1 Reference standards and directives.....	4
5. GENERAL WARNINGS	4
5.1 Warnings and instructions for correct use of the device	5
6. SYMBOLS	5
7. GENERAL DESCRIPTION.....	6
7.1 View of the parts.....	6
8. ASSEMBLY	6
8.1 Bed assembly.....	6
8.2 Motor assembly.....	7
8.3 Cables connection.....	8
9. BEFORE EVERY USE	8
10. WARNINGS ON USE.....	8
11. HOW TO USE.....	9
12. MAINTENANCE	9
13. CLEANING AND DISINFECTION	9
13.1 Cleaning.....	9
13.2 Disinfection.....	10
14. CONDITIONS OF DISPOSAL	10
14.1 General conditions of disposal	10
14.2 Avvertenze per il corretto smaltimento del prodotto ai sensi della direttiva europea 2012/19/UE:	10
15. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY	10
16. REPLACEMENT PARTS AND ACCESSORIES	13
17. TROUBLESHOOTING	14
18. TECHNICAL SPECIFICATIONS.....	15
18.1 Technical specifications for all models	15
18.2 Dimensions	15
19. WARRANTY	16
20. REPAIRS	16
21. NON-LIABILITY CLAUSE	16

CE Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices

1. CODES

- MP392** Three-jointed electric dual "Trendelenburg" hospital bed, variable height
Gardenia 2.0
- MP393** Three-jointed electric dual "Trendelenburg" hospital bed, variable height
Gardenia 2.0
- MP394** Three-jointed electric "Trendelenburg" hospital bed with wooden net, variable
height Gardenia 2.0

EN

2. INTRODUCTION

Thank you for choosing a device from the MOPEDIA line manufactured by MORETTI S.p.A. Carefully designed and manufactured to the highest quality standards, this device guarantees optimal comfort, advanced safety, and long-term reliability, making it ideal for effectively and practically meeting the needs of users and operators and for use in the home and healthcare settings.

This manual provides:

- Detailed instructions for the correct use of the device;
- Useful information for proper maintenance;
- Essential recommendations to ensure the safety of the user and staff. Important Please read the entire manual carefully before using the device. If you have any questions or need further clarification, please contact your trusted dealer or the authorized MORETTI S.p.A. technical support service.

3. INTENDED USE

3.1 Intended use

MOPEDIA hospital beds are intended for the hospitalization of people admitted for health reasons to institutions, healthcare facilities, or in a home environment.

3.2 Application environment

Gardenia beds can be used in medical facilities such as nursing homes, rehabilitation centers, and geriatric centers and/or in the home, as indicated in points "Application environment 3 and 4" of standard EN60601-2-52.

3.3 Warnings

- The Gardenia bed is designed for adults with a height of 146 cm and a body weight of 40 kg to max. 145 kg.
 - The Gardenia bed is not designed for use in hospitals.
 - The Gardenia bed is not explosion-proof and must not be used in hazardous areas.
 - The Gardenia bed may only be used in dry indoor rooms.
- They are only suitable for transporting patients within the patient's room and with the support surface adjusted to the lowest horizontal position.

- The Gardenia bed cannot be connected to an equipotential bonding system.

This must be taken into account when combining the Gardenia bed with other electro-medical devices or other mains-powered products.

- If cables from other devices are present, precautions must be taken to prevent them from being crushed between parts of the bed. Please refer to the safety information and instructions in the manuals of the electrical appliances you wish to combine with the Gardenia bed.

WARNING

It is forbidden to use this device for purposes other than those defined in this manual.

- Moretti S.p.A. declines all responsibility for damage resulting from improper use of the device or use other than that indicated in this manual.
- The manufacturer reserves the right to make changes to the device and this manual without prior notice in order to improve its characteristics.

3.4 Unauthorized use

Any use that deviates from the intended use and may therefore pose a hazard.

For example:

- The load on the medical bed exceeds the safe working load.
- Attempting to move the medical bed when the brakes are applied.
- Using the medical bed on a non-horizontal surface (max. inclination 5°).

4. DECLARATION OF CONFORMIT UE

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of MEDICAL BEDS - MOPEDIA complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017. For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

4.1 Reference standards and directives

To ensure safety standards for users of the product, Moretti S.P.A tests and repeats the legislation: UNI CEI 60601-2-52.

5. GENERAL WARNINGS

WARNING

Take extreme care in the presence of moving parts that could trap limbs and cause injury.

5.1 Warnings and instructions for correct use of the device

To ensure safe and effective use of the device, it is recommended that you strictly adhere to the following instructions:

Consulting the user manual

- Before each use, carefully read this manual, paying particular attention to the sections on routine maintenance and safety warnings.

Medical advice

- The device must be used after consulting a doctor or qualified therapist to ensure that it is fully compatible with the user's therapeutic needs.

Storage of packaging

- Store the original packaging away from heat sources or open flames, as it is made of flammable paper material.

Device service life

The service life of the device is determined by the wear and tear of components that cannot be repaired or replaced and, in any case, shall not exceed 10 years from the date of manufacture indicated on the label.

Child safety

- Pay close attention to children in the vicinity of the device in order to avoid accidents.

Load limits

- Never exceed the safe working load specified by the manufacturer. Exceeding this limit may compromise the structural integrity of the device and the safety of the user.
- Reporting adverse events

The user or patient is required to promptly report any serious incident related to the use of the device to both the manufacturer and the competent authority of the Member State in which they reside.

6. SYMBOLS

 REF	Product code		Attention, should read the instructions
 UDI	Unique Device Identification		Production dates
	CE mark		Waste Electrical & Electronic Equipment (WEEE)
	Manufacturer		B Type
 LOT	Batch Lot		Degree of protection against dusts and liquids
	Read the instruction manual		Insulation class
 MD	Medical Device		! Maximum permissible load
	Conditions for Waste Disposal		Maximun Load

7. GENERAL DESCRIPTION

MP392, MP393 and MP394 Trendelemburg electric dual 3-articulation, height-adjustable hospital bed.

EN

7.1 View of the parts


Fig. 1

List of parts

- A N°1 mesh frame head raiser part
- B N°1 mesh frame leg raiser part
- C N°2 frame connection rectangular tubes
- D N°4 knobs to lock the frame
- E N°1 headboard
- F N°1 footboard
- G N°4 wooden slats for the sides
- H N°4 plastic slides for the sides

8. ASSEMBLY

8.1 Bed assembly


Fig.2

1. Insert the rectangular tubes into the headrest frame and push them up to the stop (Fig.2)

2. Join the headrest frame with the leg frame by pushing the frame all the way down and tighten the 4 locking knobs (Fig. 3)(Fig. 3a)


Fig. 3


Fig. 3a

3. Disengage the black pin and lock it in the open position by turning it 90°, then join the frame to the headboard and footboard (identical elements), taking care to insert the metal pins correctly into the appropriate slots and that the castors are all locked during this process. Reinsert the black pin by releasing it (Fig. 4)

4. Insert the circular tubes of the slides into the holes provided on the sides of the tail lifts, making sure that the shaped and pointed parts of the tail lifts and the slide are facing upwards. Mount the slides on both sides (Fig. 5).

Fig.4

Fig.5

5. Insert the two pairs of side rails into the bottom of the guides provided on the headboard and on the wooden footboard. Carry out this operation first one side and then the other by pressing the locking button for insertion of the slide into the guide, then release the button to lock them (Fig. 6).

6. Position the safety pins at the base of the rails, after inserting the side rails (Fig. 7).

Fig. 6

Fig. 7

8.2 Motor assembly

WARNING!

Provide a soft surface underneath the motor to cushion a possible fall, and do not remove it until the end of the process described below

1. Place the actuator below the grid, ensuring that each motor moves the correct section (Fig.8, Fig.9, Fig.10)


Fig.8

Fig.9

Fig.10

2. Insert the bar tabs inside the motor housing until they both snap in. (Fig.11)
3. Close the mechanism compartment with the covers, sliding them to the side, taking care not to detach the motor from the sections (Fig.12)


Fig.11


Fig.12

8.3 Cables connection

The double motor acts as both an actuator and a control unit, so it will be necessary to connect the remote control and the cables of the headboard and footboard actuators directly to it in the appropriate connections.

9. BEFORE EVERY USE

- Always check for wear and tear of moving mechanical parts, in order to ensure that the wheelchair can be used in total safety for both people and property
- Check that the device is correctly assembled, taking particular care over the locking screws of the headboard, the footboard and the mechanical joints in general

10. WARNINGS ON USE

- Ensure that all the wheels are locked with their brake before making any adjustment to the bed
- Always disconnect the power cable from the mains electricity before performing any cleaning or maintenance operations or simply moving the bed
- Do not use additional mobile multiple plug sockets or extension cables
- Do not connect any electrical devices not intended for normal functioning to the bed

WARNING!!! The bed includes electrically controlled, moving parts. Operating the moving parts of the bed inappropriately could result in crushing of parts of the patient's limbs. Extreme care must be taken in using the electrical movements, making sure that the patient is positioned with all limbs stably resting on the top of the bed

11. HOW TO USE

The MP393 (**Gardenia 2.0**) hospital bed offers the possibility of adjusting the inclination of the head raiser and the leg raiser, using the relative buttons on the remote control.

The device is equipped with a magnetic key that unlocks the functions of the remote control by bringing it close to the base near the symbol shown below.


WARNING

Always lock the bed wheels with the brakes on each wheel before Any adjustment

WARNING

Connect the power cable to a 220-230 VAC mains plug

- A Head raiser adjustment button (UP)
- B Head raiser adjustment button (DOWN)
- C Leg raiser adjustment button (UP)
- D Leg raiser adjustment button (DOWN)
- E Combined head raiser/leg raiser adjustment button (UP)
- F Combined head raiser/leg raiser adjustment button (DOWN)
- G Mesh lifting adjustment button (UP)
- H Mesh lifting adjustment button (DOWN)
- I "Reverse Trendelenburg" adjustment button (head UP, legs DOWN)
- L "Trendelenburg" adjustment button (head DOWN, legs UP)
- X Magnetic zone for unlocking functions

N.B. The device in question might not be installed on this remote control model. The manufacturer reserves the right to decide independently, based on the type of electronics with which the bed is equipped

12. MAINTENANCE

MOPEDIA devices manufactured by Moretti S.p.A. are CE marked and carefully checked before being commercialised. For patient and doctor safety, have the product checked at least once a year by the manufacturer or by an authorised laboratory to verify its suitability for use. Periodic inspection means examination of the following parts:

- The bearing structure of the bed and the lifting mechanism with relative attachments, controls, safety devices and patient support devices
- The sides of the bed (check attachment and movement)
- Check wheel attachment and brake functioning
- Check that the electric cables and plugs are undamaged
- Check all welded points

Use only original replacement parts and accessories for repairs.

13. CLEANING AND DISINFECTION

13.1 Cleaning

Use only a damp cloth and mild soap. Then dry thoroughly before use. Use water at a maximum temperature of 30 °C. Do not use water jet and/or steam cleaning machines.

13.2 Disinfection

Use a common disinfectant cleaner if you need to disinfect the wheelchair.

N.B. Never use acid, alkali or solvents (acetone or diluent) to clean it.

14. CONDITIONS OF DISPOSAL

14.1 General conditions of disposal

Never dispose of the product as normal domestic waste. Dispose of the product at a sorted waste collection centre for recycling.

14.2 Instructions for correct disposal in accordance with Europe an Directive 2002/96/EC:

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment.

15. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY

The device is designed for use in the electromagnetic conditions specified in the following tables. The user of the device must make sure that these conditions obtain when he is using it.

Table 1

Manufacturer's guide and declaration - electromagnetic emissions		
The device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.		
Emission tests	Conformity	Electromagnetic environment - guide
RF emissions CISPR 11	Group 1	The device uses RF energy only for internal functioning and for charging the battery. Its RF emissions are therefore very low and in all likelihood do not interfere with nearby electronic equipment.
Emissions RF CISPR 11	Class B	The device is intended to be used in all buildings, including the home and buildings connected directly to the low voltage public power grid for residential applications.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage flutter/flicker IEC 61000-3-3	Conforming	

EN

Table 2

Manufacturer's guide and declaration - electromagnetic emissions			
The device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Information for using the device and charging the battery. The floor must be wood, concrete or ceramic. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst immunity IEC 61000-4-4	±2KV for power lines ±1 kV for I/O lines	2KV for power lines ±1 kV for I/O lines	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery
Over voltages IEC 61000-4-5	1 kV differential mode ±2 KV common mode	1 kV differential mode ±2 KV common mode	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery
Voltage drop outs, brief interruptions and variations of voltage on the power supply lines. IEC 61000-4-11	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	The quality of the mains voltage must be that of a typical business environment or hospital when charging the battery.
Magnetic field Mains frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the mains frequency should have the typical characteristics of a commercial or hospital environment. (Applies both to using the device and charging the battery).
Note: UT is the AC mains voltage prior to application of the test level.			

Table 3

Manufacturer's guide and declaration - electromagnetic immunity			
The device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions. Make sure it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
conducted RF IEC 61000-4-6	3 V eff. From 150 KHz to 80 MHz	3 V eff.	Portable and mobile RF communications devices should not be used near any part of the device, including the cables, unless the recommended separation distances are observed, calculated from the equation applicable to the transmitter frequency. Recommended distances $d = 1.2\sqrt{P}$ From 80 MHz to 800 Mhz $d = 2.3\sqrt{P}$ from 800 MHz to 2.5 GHz $d = 7.7$ m during charging operations where "P" indicates the maximum rated output of the transmitter in Watt (W) according to the transmitter manufacturer and "d" is the recommended separation distance in metres (m). The intensity of the field generated by fixed RF transmitters, as determined in an on-site electromagnetic inspection, might be below the level of conformity in each frequency band. Interference may be experienced in the proximity of equipment marked with this symbol: 
radiated RF IEC 61000-4-3	3 V/m From 80 MHz to 2.5 GHz	3 V/m	
NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a The intensity of the field for fixed transmitters, such as base stations for radio-telephones (cellulars and cordless) and terrestrial mobile radios, amateur radio equipment, AM and FM radio transmitters and TV transmitters may not be known precisely in theory. To establish an electromagnetic environment caused by fixed RF transmitters, one must run an on-site electromagnetic inspection. If the field intensity measured in the place where the DEVICE is used exceeds the applicable RF level of conformity indicated above, normal functioning of the DEVICE must be monitored. If abnormalities are observed, supplementary measures may be necessary, such as re-orientation or repositioning of the DEVICE. ^b The field intensity over the frequency range from 150 kHz to 80 MHz should be less than 3 V/m.			

Table 4

Recommended separation distances between portable and mobile radio communications equipment and the DEVICE						
The device is intended to operate in an electromagnetic environment in which RF disturbances are under control. The client or user of the DEVICE can contribute to preventing EM interference by ensuring a minimum distance between the mobile and portable RF communications devices (transmitters) and the DEVICE battery charger during charging and the entire device during its operation as recommended below, in relation to the maximum output power of the radio communication equipment						
Tests	Separation distance as a function of the frequency of the transmitter m					
Maximum rated output power of transmitter W	During battery charging	During use OF THE DEVICE	During battery charging	During use OF THE DEVICE	During battery charging	During use OF THE DEVICE
	from 150 kHz to 80 MHz $d = 1.2\sqrt{P}$	from 150 kHz to 80 MHz $d = 1.2\sqrt{P}$	from 80 kHz to 800 MHz $d = 1.2\sqrt{P}$	from 80 kHz to 800 MHz $d = 1.2\sqrt{P}$	from 800 kHz to 2.5 GHz $d = 1.2\sqrt{P}$	from 800 kHz to 2.5 GHz $d = 1.2\sqrt{P}$
0.01	0.12	0.12	0.12	0.12	0.23	0.23
0.1	0.38	0.38	0.38	0.38	0.73	0.73
1	1.2	1.2	1.2	1.2	2.3	2.3
10	3.8	3.8	3.8	3.8	7.3	7.3
100	12	12	12	12	23	23
For transmitters whose maximum rated output power is not listed above, the recommended separation distance d, in metres (m), can be calculated using the equation applicable to the transmitter frequency, where P is the maximum rated output power of the transmitter, in Watt (W), as given by the transmitter's manufacturer.						
NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.						
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.						

16. REPLACEMENT PARTS AND ACCESSORIES

Use only replacement parts and accessories listed in the Moretti general catalogue. Available accessories:

MPA301 Patient lifter

MPA303 Drip rod

MPA305 Hospital bed table

MPA307 Extension for sides

17. TROUBLESHOOTING

PROBLEMS	PROBABLE CAUSES	SOLUTION
No control is functioning	1. The electric plug is not connected to the mains 2. The remote control is not connected to the control unit 3. Probable fault on the remote control or control unit	1. Connect the electric plug 2. Check and, if necessary, connect the remote control to the control unit 3. Contact the assistance centre for a precise diagnosis of the problem
The electric head raiser is not working	1. The head raiser actuator is not connected to the control unit 2. The actuator or the remote control or the control unit may be faulty	1. Check and, if necessary, connect the actuator to the control unit 2. Contact the assistance centre for a precise diagnosis of the problem
The electric leg raiser is not working	1. The leg raiser actuator is not connected to the control unit 2. The actuator or the remote control or the control unit may be faulty	1. Check and, if necessary, connect the actuator to the control unit 2. Contact the assistance centre for a precise diagnosis of the problem
Electric lifting of the bed is not working	1. The lifting actuator is not connected to the control unit 2. The actuator or the remote control or the control unit may be faulty	1. Check and, if necessary, connect the actuator to the control unit 2. Contact the assistance centre for a more precise diagnosis of the problem
Trendelenburg/reverse Trendelenburg lifting of the bed is not working	1. One or both of the actuators (headboard/footboard) are not connected to the control unit 2. The actuator(s) or the remote control or the control unit may be faulty	1. Check and, if necessary, connect the actuators to the control unit 2. Contact the assistance centre for a more precise diagnosis of the problem

18. TECHNICAL SPECIFICATIONS

18.1 Technical specifications for all models

Units of force needed to control the commands on the remote control	5N
Power supply	220 VAC 50-60 Hz
Voltage output from the control unit	24 VDC
Protection class	IP 54
Noise level	< 45 dB
Insulation class	Class 2
Materials used	Powder-coated steel tube structure, welded steel mesh. Wooden headboard, footboard and aluminium sides

18.2 Dimensions

Bed dimensions: 216x102x127h cm

Mesh dimensions: 187x80x40 cm÷81h

Maximum height at the headboard/footboard: 127 cm

Minimum height at the headboard/footboard: 25 cm

Maximum angle of the head raiser: 70°

Maximum angle of the leg raiser: 30°

Maximum "Trendelenburg" angle (head down, legs up): 12°

Maximum "reverse Trendelenburg" angle (head up, legs down): 12°

Maximum weight	135 kg
Total weight	94 kg only for MP392 85 kg
Headboard weight	22 kg only for MP392 18 kg
Footboard weight	22 kg only for MP392 18 kg
Frame weight headboard	18 kg
Frame weight footboard	21 kg
Electrical kit weight	9 kg
Peso máximo de trabajo	170 kg

19. WARRANTY

Moretti products are guaranteed for 2 (two) years from the date of sale against material and manufacturing defects, subject the following limitations. The warranty is voided by improper use, abuse, modifications to the product and failure to follow the instructions.

The intended use of the product is given in the user manual

Moretti is not liable for damage, injury or any other consequences resulting from installation or use which are not scrupulously conforming with the instructions given in the installation, assembly and user manual.

Moretti does not guarantee its products against damage or defects in the following circumstances:

natural disasters, unauthorised repair or maintenance, improper electric power supply (as applicable), use of parts or components not supplied by Moretti, failure to follow the guidelines and instructions for use, tampering, shipping damage (other than the original shipping by Moretti), or failure to perform maintenance as indicated in the manual.

Components subject to wear and tear are not covered by this warranty if the damage is caused by normal use of the product.

20. REPAIRS

Repairs under warranty

If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti reseller or its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

Repairing a product not covered by the warranty

A product not covered by warranty may be returned for repair only if authorised in advance by Moretti customer service. The costs of labour and shipping incurred by repairs not covered by the warranty are borne by the client or reseller in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for 6 (six) months from the day of reception of the repaired product.

Non-defective products

The client will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the client at his expense.

Replacement parts

Moretti original replacement parts are guaranteed for 6 (six) months from the day of delivery.

21. NON-LIABILITY CLAUSE

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti makes no declaration, guarantee or condition, express or implicit, including any future declaration, guarantee or condition of sale, suitability for a given purpose, non violation and non interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit guarantee under the law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit guarantee or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice.

**WARRANTY CERTIFICATE**

Product _____

Purchased on (date) _____

Retailer _____

Address _____ **Town/city** _____

Sold to _____

Address _____ **Town/city** _____

**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com - email: info@morettispa.com

MADE IN CHINA

MORETTI S.p.A.

Via Bruxelles, 3 - Meleto
52022 Cavriglia (Arezzo)

Tel. +39 055 96 21 11
Fax. +39 055 96 21 200

www.morettispa.com
info@morettispa.com