size: 110x145mm

LOGIKODIGIT

MA DM490-491-492 01 D_ITA_03-2022

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 CE⁰¹⁹⁷ IIa Class Medical Device

According to directive 93/42/EEC and further modifications

1. CODES

DM490	Digital sphygmomanometer desk type, LCD display 3"
DM491	Digital sphygmomanometer desk type, LCD display 4"
DM492S	Digital sphygmomanometer desk type, LCD display 4,8"

2. INTRODUCTION

Thank you for purchasing LOGIKO sphygmomanometer by Moretti S.P.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

WARNING!

- Do not use the product for a purpose not indicated in this manual
- Moretti S.P.A. declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.
- This device is intended for non-invasive measuring of an adult individual's systolic, diastolic blood pressure and heart rate using the oscillometric method.
- The device is not intended to use on infants and children.
- The device is designed for home or clinical use.
- All values can be read out in one LCD Display.
- Measurement position is on adult wrist only.

3. INTENDED USE

The digital sphygmomanometer wrist type is measurement device for systolic and diastolic blood pressure.

4. REGULATIONS AND DIRECTIVE OF REFERENCE

The LOGIKO sphygmomanometer by Moretti S.p.A are in accordance to the following regulations: This Blood Pressure Monitor complies with the European regulations. This blood pressure monitor also complies with mainly following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2 Medical electrical equipment part 1-2: General requirements for safety- Collateral standard: Electromagnetic compatibility- Requirements and tests

Performance standards:

EN ISO 81060 -1 Non-invasive sphygmomanometers - General requirements

EN 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

EN ISO 81060 - 2 Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers

5. GENERAL WARNINGS

- DO NOT use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment, otherwise injury or damage may occur.
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the sphygmomanometer.
- Keep the packed sphygmomanometer away from heat sources.
- SERVICE LIFE- the sphygmomanometer use limit is defined by the wear of the parts.
- DO NOT allow children to play on or operate with the sphygmomanometer.
- Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.
- Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.
- Product is designed for its intended use only.
- DO NOT misuse in any way.
- Product is not intended for infants or individuals who cannot express their intentions.
- DO NOT disassemble or attempt to repair.
- DO NOT use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.

6. SYMBOLS USED

 serial number
  warning
  EC REP
  EC REP

 type BF
  instruction
  LOT
  lot
  date of production

 disposal
  CE Mark
  manufacturer

7. GENERAL DESCRIPTION


What is Blood Pressure?

Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure. Blood pressure is measured in mmHg. It is submitted to many fluctuations and it depends on many elements. It is usually lower in summer and higher in winter. Blood pressure can change according to atmospheric pressure and it is affected by physical efforts, temperament, stress, meals, medicines, alcohol, smoking, age...etc... It is recommended to note down your daily measurements and to consult with your physician who will tell you your normal blood pressure range. For the above reasons do not worry if you encounter some measurements at high level, consult with your physician who will tell you the point at which you will be considered at risk.

NOTE: Blood pressure can change with the age, it is necessary to consult with a physician to know your normal blood pressure range. Never change drug dosage without consulting your physician.

Why measure your blood pressure?

Among the various health problems afflicting modern people, those associated with blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressure as a necessity in identifying those at risk.

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7.1 Questions and answers

Q: What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

A: Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

NOTE: Abnormal test results may be caused by:

1. Improper cuff placement Make sure cuff is snug-not too tight or too loose.
2. Improper body position Make sure to keep your body in an upright position.
3. Feeling anxious or nervous Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different readings?

A: Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.

Q: Should I apply the cuff to the left or right wrist? What is the difference?


A: Either wrist can be used when testing, however, when comparing results, the same wrist should be used. Testing on your left wrist may provide more accurate results as it is located closer to your heart.

Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.

8. CONTENTS


8.1 Contents DM490


Monitor Unit

Owner's Manual


Arm Cuff

Soft bag


8.2 Contents DM491-DM492S


Monitor Unit


Owner's Manual

Arm Cuff


Soft pouch
(DM491)


Soft bag
(DM492)

9. PARTS VIEW


9.1 Parts view DM490

9.2 Parts view DM491-DM492S


9.3 Display


1. Groups of memories (1/2) only mod. DM491-492
2. Low battery indicator
3. WHO blood pressure classification indicator
4. Last 3 measurements average
5. Irregular heartbeat indicator
6. Time
7. Date
8. Systolic blood pressure
9. Diastolic blood pressure
10. Pulse rate/min


NOTE: The arrow indicator will show your current measurement value according to WHO classification. See the scale here below:

Blood Pressure Classification Indication

The following scheme has been developed by The World Health Organization (WHO), and it identifies the low- and high-risk blood pressure areas. This standard, however, is a general guideline as individuals' blood pressure varies among different people and different age groups, it is always subjective.

10. CAUTIONS BEFORE USE

- Avoid any type of physical exercise for at least 30 min before the measurement
- Sit in a calm environment for at least 5 minutes before the measurement
- Avoid speaking during the measurement
- Do not move during the measurement
- Do not stand during the measurement. Sit in a relaxed position, rest your elbow and be sure the cuff is at heart level
- During the measurement, avoid electromagnetic interference with other electronic devices
- Wait 3 min before a new measurement
- If you need to measure your blood pressure every day, do it at the same time of the day
- Results comparisons should be made only when the monitor is used on the same arm, in the same position and at the same time of the day
- This blood pressure monitor is not recommended for people with strong arrhythmia.
- Do not use the device when damage and/or malfunction has occurred
- Remove batteries when the unit is not in use for an extended period of time


10.1 Battery precautions (⊗)

- Do not mix new and old batteries simultaneously.
- Replace batteries when Low Battery Indicator appears on screen (⊗).
- Be sure battery polarity is correct.
- Do not mix battery types. Long-life alkaline batteries are recommended.
- Remove batteries from device when not in operation for more than 3 months.
- Dispose batteries properly; observe local laws and regulations.


11. HOW TO USE

11.1 Quick start



1. Install batteries.

2. Insert cuff air plug into the left side of monitor unit.

3. Rest for several minutes prior to testing. Sit down in a quiet place, preferably at a desk or table, with your arm resting on a firm surface and your feet flat on the floor


4. Apply cuff to your left arm and keep level with your heart. Bottom of cuff should be placed approximately 1-2cm (0.4-0.8") above elbow joint.

5. Press "START/STOP" button to start testing

11.2 Battery installation


- Slide battery cover off as indicated by arrow.
- Install 2 new AAA alkaline batteries according to polarity.
- Close battery cover.

NOTES: Replace batteries when Low Battery Indicator appears on screen. ☒
2. Batteries should be removed from device when not in operation for an extended period of time.

11.3 Using the AC/DC adapter (optional)

AC Adapter jack is on the back side of the monitor. Medical AC adapter (DC 6.0 V,600mA) can be used with the device (recommended, not provided). The adapter connect pin should be positive inside and negative outside with a 2.1mm coaxial joint. Do not use any other type of AC adapter as it may harm the unit.

12. START


Applying The Wrist Monitor

Do not apply over clothing. If wearing a long sleeved shirt, be sure to roll sleeve back to forearm. Apply monitor to wrist as illustrated. Tighten cuff firmly as not to wiggle. Do not stand while testing. Sit in a comfortable position and place wrist level with heart.

1. Power On


Press "START/STOP" button until a beep sounds. The LCD screen will appear for one second as unit performs a quick diagnosis. A long tone indicates device is ready for testing.

NOTE: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash "↓" until pressure is stabilized.

2. Pressurization

Initial pressure is first pumped to 190mmHg. If the current user's systolic blood pressure is over 190mmHg, the unit will automatically re-inflate to the proper shelf

NOTE: Pressurization will gradually subside and ultimately stop when cuff is not properly applied to the wrist. If this occurs, pressing the "START/STOP" button to turn the unit off.

3. Testing

After cuff inflation, air will slowly subside as indicated by the corresponding cuff pressure value. A flashing "♥" will appear simultaneously on screen signaling heart beat detection.

NOTE: Remain relaxed during testing. Avoid speaking or moving body parts


4. Result Display

Three short beeps sound when testing is complete.

The screen will display measurements for systolic and diastolic blood pressure. An indicator representing the current measurement will appear next to the corresponding WHO Classification.

When the measurement is finished, the following results will be shown on the display:

- Systolic pressure (SYS)
- Diastolic pressure (DIA)
- Beats per minute (PULSE)
- Classification of values according to standard WHO/OMS (World Health Organization)
- Indicator of irregular beats (Arrhythmia) "♥"

NOTES: This symbol will be shown only if the device detects two or more irregular beats during measurement.

Do not be alarmed if an abnormal reading occurs. each day over an extended period of time. Consult your physician if test results remain abnormal. A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time.

Irregular Heartbeat Indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol "♥" appears on screen along with measurement results. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult your physician if the Irregular Heartbeat Symbol.

Consult your physician if the Irregular Heartbeat Symbol "♥" frequently appears with your test results.

5. Storing Test Results

User may delete their current test result due to unfavorable testing conditions or for any other reason.

- To delete the last test result, press the "SET" button after result is displayed. If result is not deleted,
- It will automatically store by date.

NOTE: Wait 3 minutes or longer prior to re-testing

If the number of tests surpasses the allotted 120 memories, the most recent tests will appear first, thus eliminating the oldest readings.

6. Power Off


- The "START/STOP" button can be pressed to turn off the unit in any mode.
- The unit can turn off the power itself about 3 minutes no operation in any mode.

Safety Precaution: If pressure in cuff becomes too extreme while testing:

- press the "START /STOP" button to turn power off. The cuff pressure will rapidly dissipate once the unit is off.

7. Memory Check

- With power off, you may check past test results by using the "M" button.
- Upon activating test results, you can press the "M" button to scroll through all test results stored in memory. The LCD will display the last measurement memory as NO:1 reading.

8. Memory Deletion

- Memory for a selected group may be deleted while in Memory Check mode.
- Press and hold the "SET" button for approximately 3 seconds to delete all memory records.
- The monitor will beep indicating successful deletion and then transfer into testing mode.

9. Last 3 Tests Average AVG (just for DM491 - DM492S)

- With power off, press the "M" button to activate screen display.
- After the unit performs a self diagnosis, the screen will display the average test results from the last 3 readings of the last group used.
- The "AVE" symbol will appear along with the corresponding WHO Blood Pressure Indicator.


10. Low Battery Indicator

4 short warning beeps sound when battery life is depleting and unable to inflate cuff for testing. The "🔋" appears simultaneously for approximately 5 seconds prior to shutting down. No memory loss will occur throughout this process. Replace batteries at this time.

13. SELECT MEMORY GROUP (ONLY DM491 AND DM492S)

While in the System Setting mode, you may accumulate test results into 2 different groups. This allows multiple users to save individual test results (up to 60 memories per group.) Press "M" button to choose a group setting. Test results will automatically store in each selected group.

13.1 Time/Date Settings

DM490

- When power off, press "START/STOP" button for 3 seconds to set the Time/Date mode.
- "M" (stands for month) is flashing
- Press "MEM" button for set up
- Press "START/STOP" button to confirm it
- "D" (stands for day) is flashing
- Press "MEM" button for set up
- Set the day ("D"), hour (in AM/PM format) and minute in the same way.
- Press "START/STOP" button for 3 seconds to turn the unit off. All information will be saved

DM491

- When power off, press "SET" button to start settings
- The Memory Group icon flashes in the upper left corner (see Memory Group settings)
- Press "SET" button to set the "Time/Date" mode
- "M" (stands for month) is flashing
- Set the month by pressing "M" button
- Press SET again to confirm current month
- "D" (stands for day) is flashing
- Set the day ("D"), hour (in AM/PM format) and minute in the same way
- Press "START/STOP" button for 3 seconds to turn the unit off. All information will be saved

DM492S

- When power off, press "SET" button to start settings
- The Memory Group icon flashes (see Memory Group settings)
- Press "SET" button to set the "Time/Date" mode
- When the "Yr" icon flashes select M button to set current year
- Press "SET" button again to confirm it. (Every time the "SET" button is pressed, it will lock in your selection and continue in succession)
- Set the month, day, hour and minutes in the same way
- Press "SET" button again to set the time format setting mode. Set the time format by adjusting the "M" mode. EU means European time while, US means U.S. time

NOTE: While in any setting mode, press "START/STOP" button to turn the unit off. All information will be saved. Unit will automatically save all information and shut off if left it for 3 minutes.

14. MAINTENANCE

The LOGIKO devices by Moretti S.P.A are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

15. CLEANING AND DISINFECTION

- DO NOT immerse the device!
- Cleaning: use only a damp cloth. Pay attention to the display
- Disinfection: use a cleanser, avoiding the liquid seep into the device

NOTE: Never use acids, alkalis or solvents such as acetone

- Avoid dropping, slamming, or throwing the unit
- Avoid extreme temperatures. Do not use outdoors
- Do not use petrol, thinners or similar solvents
- Remove batteries when not in operation for an extended period of time
- Do not disassemble product
- It is recommended the performance should be checked every 2 years

16. CONDITIONS OF DISPOSAL

General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the sphygmomanometer in the appropriate disposal areas for recycling.

TREATMENT OF BATTERIES (Directive 2002/96/EC):

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

DISPOSAL OF WASTE BATTERIES - Directive 2006/66/EC):

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

17. TROUBLESHOOTING

PROBLEM 1

Blood pressure results are not reliable

POSSIBLE CAUSE

1. Cuff is too tight or not properly positioned on the arm
2. Body or device movements during measurement

SOLUTION

1. Firmly reposition cuff as illustrated in point 5, paragraph Pressure Measurement
2. Stand still and relaxed during measurement, avoiding speaking

PROBLEM 2

"ERR" displayed

POSSIBLE CAUSE

1. Cuff fails to inflate properly
2. Improper operation
3. Pressurization is over 300mmHg

SOLUTION

1. Check connections between cuff and device
2. Read user manual carefully and do the operation properly
3. Read user manual carefully and do the operation properly

18. ACCESSORIES

- DR205 - Charger
 - DR215 - Large size cuff
 - DR211 - Small size cuff
 - DR475 - USB cable + CD software (only for DM492P - DM492S)
- For spare parts and accessories refer to main catalogue.

19. ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturer's declaration- Electromagnetic emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT-GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR11	Class B	The device is suitable for use in all establishment. Including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and manufacturer's declaration - Electromagnetic immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to the application of the test level			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity -

The devices are intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bands	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey c, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device

NOTE 4 Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER M		
	150 KHz to a 80MHz d=1.2√P	80MHz to 800MHz d=1.2√P	800MHz to 2.5 GHz d=2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

20. SPECIFICATIONS

PRODUCT	Fully automatic blood pressure monitor
MODELS	Cod.DM490 Mod.BSP11 - Cod.DM491 Mod.BSP12 Cod.DM492S/DM492P Mod.BSP13
DISPLAY	DM490: 62,3 x 46 mm (3") DM491: 84,1 x 55,1 mm (4") DM492S: 102,1 x 68,9 mm backlit (4,8")
MEASUREMENT METHOD	Oscillometric
MEASUREMENT ACCURACY	Pressure: 0mmHg / 300 mmHg Pulse: from 30 to 180 Puls/min
MEASUREMENT ACCURACY	Pressure: +/- 3mmHg Pulse: +/- 5%
PRESSURIZATION	Automatic pressurization
MEMORY	DM490: 120 memories DM491: 120 memories in two groups DM492S: 120 memories in two groups
FUNCTION	Irregular beat detection (arrhythmia) WHO/OMS classification indicator Last 3 result average - AVG Automatic Power off
POWER SOURCE	4 AA batteries - 1,5V Power Supply AC-DC 6.0VDC 600mA (OPTIONAL)
BATTERY LIFE	Approximately 2 months at 3 Tests day

UNIT WEIGHT	DM490: Approx. 340 g (excluding batteries) DM491: Approx. 382 g (excluding batteries) DM492S: Approx. 430 g (excluding batteries)
UNIT DIMENSION	DM490: 140 x 98 x 48h mm DM491: 148 x 100 x 56h mm DM492S: 164 x 111 x60h mm
CUFF CIRCUMFERENCE	From 22 cm to 36 cm
OPERATING ENVIRONMENT	Temperature: from 10°C to 40°C Humidity: < 85% RH Pressure: Atmospheric
STORAGE ENVIRONMENT	Temperature: from -20°C to 55°C Humidity: < 95% RH Pressure: from 525mmHg to 795 mmHg
CLASSIFICATION	Internal Powered Equipment type BF 

21. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual. Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use. Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

21.1 Rechargeable batteries warranty (if provided)

The original batteries and the spare batteries are covered by a 90 days warranty during a period of 6 month, with regards to manufacturing defeat or to what concern the law regulation. If charge batteries are unused for more than 3 consecutive months, the warranty is void. If exhausted batteries are unused for more than 3 consecutive days, the warranty is void.

22. REPAIRING

22.1 Warranty repair

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

22.2 Repair not covered by warranty

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

22.3 Non-defective products

The customer will be notified if, after examining and tried a returned product, Moretti concludes that the product is not defective. The product will be returned to the customer and will be liable to pay the costs caused by the return.

23. SPARE-PARTS

Moretti's genuine parts are guaranteed for 6 months with effect from the day you receive the replacement.

24. EXEMPT CLAUSES

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user.

LOGIKODIGIT

WARRANTY CERTIFICATE

Product _____

Purchasing date _____

Authorized dealer _____

Street _____ **Place** _____

Sold to _____

Street _____ **Place** _____

Joytech Healthcare Co., Ltd
No.365, Wuzhou Road, Yuhang Economic Development Zone,
311100, Hangzhou, China

Shanghai International Holding Corp. GmbH (Hamburg)
Eiffestrasse 80, 20537 Hamburg, Germany

Distributed by: MORETTI S.P.A. Via Bruxelles, 3 - Melegnano 20122 Cavriaglia (Arezzo)

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