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N° DI LOTTO: _____

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GB

WARRANTY CERTIFICATE Valid for 36 months from the date of purchase

WARRANTY CONDITIONS

- The device is guaranteed for 36 months from the date of purchase against all defects arising from materials or manufacture, provided it has not been tampered with by the customer or unauthorised persons.
- The warranty covers the replacement or repair of construction components.
- For health reasons, the device will not be replaced under any circumstances, since it is a strictly personal medical device.
- This warranty does not cover parts subject to normal wear and tear and any damage arising from improper use, dropping, transport, failure to carry out proper maintenance or any circumstances beyond the manufacturer's control.
- 3A Health Care declines all responsibility for any direct or indirect damage arising from improper or reckless use of the product.
- In the event of a malfunction, the device, suitably cleaned and packed, must be sent to our dealer, enclosing this warranty certificate, properly filled in, and the tax receipt or purchase invoice; otherwise the warranty shall not be considered valid and the services rendered shall be charged for.
- Outward and return shipping costs shall be paid by the customer.
- 3A Health Care does not accept liability for further extensions of the warranty period offered by third parties.

WARNING: THE WARRANTY IS ONLY VALID IF COMPLETELY FILLED IN AND ACCOMPANIED BY THE TAX RECEIPT/INVOICE.

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USE AND MAINTENANCE MANUAL

Manufactured in accordance with the latest regulations, the **LTK185** is a quick effective nasal irrigation device for rinsing the nasal and rhino-pharyngeal cavities (Risk Class IIa - Medical Device Directive 93/42/EEC and subsequent amendments).

SAFETY GUIDELINES

- The device should be treated as mono-patient.
- If there are risks of infection or microbe contamination, the device should be used by one person only (always consult your physician).
- The device must only be used with aerosol therapy compressor devices compliant with the latest regulations, having the use of nasal devices of this type envisaged in their instruction manuals.
- This device is not suitable for anaesthesia or pulmonary ventilation.
- The device must only be used with the original accessories described in this instruction manual.
- The information relating to the connection to and use in combination with aerosol treatment devices is provided in the instructions for use of the devices concerned. Always comply with the instructions and safety guidelines provided in the instruction manuals of the aerosol devices which envisage use with nasal devices of this type.
- Always use the device and its accessories in accordance with your physician's instructions. Use only medicinal specialities prescribed by your physician and administer them in accordance with his or her instructions.
- Never use the device for any purpose except for its intended use as a nasal irrigation device. The manufacturer is not responsible for improper use.
- The device is not supplied in a sterile pack; always clean and disinfect it before and after each use.
- The device contains small components which may be removed and easily swallowed. The device must therefore only be used by juveniles and the differently able in the presence of an adult in full possession of their mental faculties. Never leave the device unattended in places easily accessible to juveniles and the differently able.
- Store in a clean, dry place, protected from light, heat and weather. Dispose of the device in accordance with the relevant legal requirements.

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INSTRUCTIONS FOR CLEANING THE DEVICE

At the end of each application, clean all the components to remove residues of drug and any dirt. Clean the spray device and its parts as described below, except for the air hose. If the device is to be used on another patient, or if it is dirty, change the hose. Clean the components (except for the air hose) in running hot tap water (max. 60° C) for about 5 minutes, adding a very small amount of detergent if necessary (following instructions concerning quantities and restrictions on use provided by the detergent manufacturer). Rinse with plenty of water, remove the excess water and leave to dry in the air in a clean place. The spray device and accessories (except for the air hose) can be disinfected with chemical disinfectants easily available from pharmacies, suitable for use on plastics (Milton®, Amuchina®, etc.). Follow the instructions and details of quantities and restrictions on use provided by the disinfectant manufacturer. All parts of the spray device may be sterilised in steam (except for the air hose and masks) at 121° C (20 min.) or 134° C (7 min.) (EN554/ISO11134). Sterile packaging must be compliant with EN868/ISO11607 and must be suitable for sterilisation with steam. After sterilisation, always allow the device and its parts to cool to room temperature before use. Never repeat the sterilisation cycle on the hot device without allowing it to cool.

INSTRUCTIONS FOR PREPARING THE DEVICE AND PUTTING IT INTO OPERATION

The device is **not sterile**, so before use, perform the cleaning and disinfection procedures described in the specific section of this manual. Connect the system as shown in figure 2. To fill with solution, squeeze the bottom of the spray tightly and turn the top anti-clockwise until it is disconnected. Remove the top and place the quantity of solution required in the bottom section (3); for quantities, refer to the graduated scale marked on the bottom section (3) with the symbol:

DEVICE WORKING LIFETIME

The **LTK185** spray device should be replaced after max. 100-120 applications or after about 20 sterilisation cycles.

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