LOGIKODIGIT

MA DM590 DM591 DM592S DM592P 01 C*_ITA_03-2023

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Via _____ Località _____

Joytech Healthcare Co., Ltd

No.365, Wuzhou Road, Yuhang Economic Development
Zone, Hangzhou City, 311100 Zhejiang P.R. China

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2. INTRODUCTION

Thank you for purchasing LOGIKO sphygmomanometer by Moretti S.P.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

WARNING!

- Do not use the product for a purpose not indicated in this manual
- Moretti S.P.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.
- This device is intended for non-invasive measuring of an adult individual's systolic, diastolic blood pressure and heart rate using the oscillometric method.
- The device is not intended to use on infants and children.
- The device is designed for home or clinical use.
- All values can be read out in one LCD Display.
- Measurement position is on adult wrist only.

3. INTENDED USE

The digital sphygmomanometer wrist type is measurement device for systolic and diastolic blood pressure.

4. REGULATIONS AND DIRECTIVE OF REFERENCE

Safety standard:

EN 60601-1

Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2

Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard: Electromagnetic Disturbances - Requirements And Tests

Performance standards:

IEC 80601-2-30

Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

EN 1060-3

Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

ISO 81060-2

Non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type.

5. GENERAL WARNINGS

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.

2. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
3. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
4. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
5. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.
6. Too frequent measurements can cause injury to the patient due to blood flow interference.
7. The cuff should not be applied over a wound as this can cause further injury.
8. DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
9. The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
10. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
11. A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
12. Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.
13. Product is not intended for infants or individuals who cannot express their intentions.
14. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.
15. Use only the approved arm cuff for this unit. Use of other arm cuffs may result in incorrect measurement results.
16. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Keep the device away from childrens, pets and pests.
17. Advising operator that Instruction manual/ Booklet must be consulted.
18. Do not use the device during patient transport outside healthcare facility for interference source existing as well.
19. Essential performance:

| | |
|---|-----------------------------------|
| Electrosurgery interference recovery | Refer 202.6.2.101 IEC 80601-2-3 |
| Limits of the error of the manometer | Refer 202.12.1.102 IEC 80601-2-30 |
| Reproducibility of the BLOOD PRESSURE DETERMINATION | Refer 201.12.1.107 IEC 80601-2-30 |

| | |
|---|---|
|  CAUTION! | |
| Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use. |  |
| Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first. | |
| Product is designed for its intended use only. Do not misuse in any way. |  |
| Product is not intended for infants or individuals who cannot express their intentions. | |
| Do not disassemble or attempt to repair. | |
| Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device. | |
| Only use a recommended AC adaptor double-insulated complying with EN 60601-1 and EN 60601-1-2(see page 6). An unauthorized adaptor may cause fire and electric shock. | |
|  ADVICE ON THE USB CABLE | |
| Please use the USB cable after thoroughly reading the User Manual. |  |
| Do not move the cable during data transfer. Replace the batteries after removing the cable. Do not touch the cable with wet hands. Do not leave the USB connected to the device when not in use. |  |

5.1 Battery Precautions

- Do not mix new and old batteries simultaneously.
- Replace batteries when Low Battery Indicator “” appears on screen.
- Be sure battery polarity is correct.
- Do not mix battery types. Long-life alkaline batteries are recommended.
- Remove batteries from device when not in operation for more than 3 months.
- Dispose batteries properly; observe local laws and regulations.

6. SYMBOLS

Caution

Mandatory

Prohibited

Type BF Equipment

 Instructions For Use
 MUST be Consulted

The product conforms to the requirements of the EC Directive MDD(93/42/EEC) on medical devices

Manufacturer

Authorised Representative in the European Community

Keep Dry

Keep off Sunlight

SN

Serial Number

LOT

Batch Lot

Discard the used product to the recycling collection point according to local regulations

Manufacturing Date

IP21

Ingress protection code against dusts and liquids

7. GENERAL DESCRIPTION

7.1 What is blood pressure?


Pressure measures the action of the force exerted by blood on the walls of blood vessels. Systolic or maximum pressure indicates the maximum force exerted against the vessel walls during cardiac contraction. DIASTOLIC or minimum pressure indicates the force exerted on the vessels during cardiac relaxation. It is measured in mmHg (millimetres of mercury).

It is necessary to know that an individual's blood pressure fluctuates greatly and depends on many factors. It is generally lower in summer and higher in winter. Pressure can vary with atmospheric pressure and is greatly influenced by physical loads, emotional excitability, stress, meals, medication, alcohol, smoking, age, etc. It is recommended to note down your daily measurements and to consult with your physician who will tell you your normal blood pressure range. For the above reasons do not worry if you encounter some measurements at high level, consult with your physician who will tell you the point at which you will be considered at risk

NOTE: Blood pressure can change with the age, it is necessary to consult with a physician to know your normal blood pressure range. Never change drug dosage without consulting your physician.


7.2 Why measure your blood pressure?

Among the various health problems afflicting modern people, those associated with blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressure as a necessity in identifying those at risk.

**8. CONTENTS**


8.1 Contents


1. Monitor Unit

2. Owner's
Manual

3. Arm Cuff


4. Storage
Case5. MOPP Medical AC
Adapter (DC5.0
V,1000mA)
(recommended,
not provided)**9. PARTS VIEW**9.1 Monitor unit
DM590

DM591



DM592S/DM592P


9.2 Arm cuff


9.3 Display
DM590

DM591


DM592S/DM592P


9.3.1 Blood pressure WHO classification

The LOGIKO by Moretti monitors are equipped with a classification indicator based on established guidelines from the World Health Organization. The chart below (colour coded on monitor unit) indicates test results.

Classification of values according to the standard

DM590


DM591

DM592S / DM592P

The following scheme has been developed by The World Health Organization (WHO), and it identifies the low- and high-risk blood pressure areas. This standard, however, is a general guideline as individuals' blood pressure varies among different people and different age groups, it is always subjective.

10. IMPORTANT TESING GUIDELINES

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position.
4. Avoid speaking or moving body parts while testing.

5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

11. HOW TO USE

11.1 Quick start

1. Install batteries (See figure A)
2. Insert cuff air plug into the left side of monitor unit. (See Figure B)

DM590


Figure A


Figure B

DM591


Figure A

Figure B

DM592S/DM592P


Figure A

Figure B

- Remove thick clothing from the arm area.
- Rest for several minutes prior to testing. Sit down in a quiet place, preferably at a deskor table, with your arm resting on a firm surface and your feet flat on the floor.(See Figure C)

Figure C

- Apply cuff to your left arm and keep level with your heart. Bottom of cuff should be placed approximately 1-2cm (0.4-0.8") above elbow joint. (See Figures D&E)


Figure D



Figure E

- Press " START/STOP " Button to start testing.


11.2 Battery installation

- Slide battery cover off as indicated by arrow.
- Install 3 new AAA alkaline batteries according to polarity.
- Close battery cover.


DM590

DM591


DM592S/DM592P

Medical AC adapter (DC 5.0 V, 1A, 5.0 W) can be used with the device (recommended, not provided). The adapter connect pin should be positive inside and negative outside with a 2.1mm coaxial joint. Do not use any other type of AC adapter as it may harm the unit.

DM590

USB Adapter
Jack

DM591

USB Adapter Jack

DM592S/DM592P


USB Adapter Jack

Note: Power supply is specified as part of MEDICAL ELECTRONIC EQUIPMENT.


**12.START****12.1 Cuff positioning**

1. Firmly insert air plug into opening located on left side of monitor unit.

DM590


DM591

DM592S/DM592P

2. With sticky nylon section facing outward, insert end of cuff underneath metal ring of cuff
3. Fasten cuff about 1-2cm (0.4-0.8") above the elbow joint. For best results apply cuff to bare arm and keep level with heart while testing.

Note: Do not insert air plug into opening located on right side of monitor unit. This opening is designed for an optional power supply only.


12.2 Power on

Press and hold "START/STOP" button to turn the unit on. The LCD screen will appear for one second as unit performs a quick diagnosis. A voice tone will indicate when unit is ready for testing.

DM590

DM591

DM592S/DM592P

Note: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash "▼" until pressure is stabilized.

12.3 Testing

After cuff inflation, air will slowly rise as indicated by the corresponding cuff pressure value. A flashing  will appear simultaneously on screen signaling heart beat detection.

DM590

DM591

DM592S/DM592P

Note: Keep relaxed during testing. Avoid speaking or moving body parts.

12.4 Result

The screen will display measurements for systolic and diastolic blood pressure with voice broadcast (DM592P only). An indicator representing the current measurement will appear next to the corresponding WHO Classification

DM590

DM591

DM592S/DM592P

Note: Refer to "Blood Pressure WHO Classification" for detailed information

12.5 Irregular heartbeat indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol "IHB" appears on screen along with measurement results. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult your physician if the Irregular Heartbeat Symbol "IHB" frequently appears with your test results.

12.6 Arm shake indicator

If there is arm movement during the measurement, the "" icon may flash. Indicates that the measurement results may be inaccurate, and the situation will be recorded at the end of the measurement as a reminder.

12.6.1 Cuff loose indicator

When starting the measurement, "C" will be displayed when the cuff is properly wound. When the cuff is too loose, "C" will be displayed. At this time, please wear the cuff correctly and start measuring again.

12.7 Power off


The "START/STOP" button can be pressed to turn off the unit in any mode. The unit can turn off the power itself about 3 minutes no operation in any mode.

Safety Precaution: If pressure in arm cuff becomes too extreme while testing, press the "START/STOP" button to turn power off. The cuff pressure will rapidly dissipate once the unit is off.

12.8 Last 3 tests average

With the power off, press the "MEM" button to activate the screen. After self-diagnosis, the screen will display the average of the last 3 readings of the last group used. The "AVG" symbol will appear together with the corresponding WHO/WHO classification indicator.

To check the average results of other groups, select the desired group before activating the "MEM" button. (see "Select memory group")

12.9 Memory check

You may check past test results by using the "MEM" button. The most recent test result and oldest test result in memory can be viewed by pressing and holding the "MEM" button. Upon activating test results, you can press the "MEM" button to scroll through all test results stored in memory. To check the average results from other groups, please select desired group first.

12.10 Memory deletion

Memory for a selected group may be deleted while in Memory Check mode. Press and hold "SET" ("Group1/Group2") button for approximately 3 seconds to delete all memory records from the selected group with voice broadcast "Memory Clear" and then transfer into testing mode. Press the "START/STOP" button to turn the unit off.

Note: Memory cannot be recovered once it has been deleted.


12.11 Low battery indicator

The unit will broadcast "Low Battery" when battery life is depleting and unable to inflate cuff for testing. The "⊗" appears simultaneously for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.

12.12 Static pressure measurement

In the power down state, press and hold the "START/STOP" button, and then install the batteries. Until the LCD screen is full, release the "START/STOP" button.

When the LCD screen displays the double zero, the bloodpressure meter is in static state. Software version is displayed: 10 is a software version in the figure.

Note: Only Service personnel can have access to this mode.

13. SYSTEM SETTINGS DM590

With power off, press "SET" button to activate System Settings. The Memory Group icon flashes.

13.1 Select memory group

While in the System Setting mode, you may accumulate test results into 2 different groups. This allows multiple users to save individual test results (up to 150 memories per group.) Press "MEM" button to choose a group setting. Test results will automatically store in each selected group.

13.2 Date/Time setting

Press the SET button to set the Time/Date mode press the 'MEM' button to change the current datum. Press the SET button to confirm the current datum. Continue setting date, hour and minute as described. Press the SET button to confirm the selection and continue in succession (month, day, hour, minute, 12/24 hour).

13.3 Time format

Press the 'SET' button again to continue to the time format selection, press the 'MEM' button to choose 12/24 hours.

13.4 Saving

In any mode, press the " START/STOP " button to switch off the unit. All information will be saved.

Note: If the unit is switched on and not used for 3 minutes, it will automatically save all information and switch off..

14. SYSTEM SETTINGS DM591

With the power off, press the "  " button to activate the system settings. The memory group icon flashes.

14.1 Select memory group

In system settings mode, it is possible to choose which group to store test results in from 2 different groups. This allows multiple users to save individual test results (up to 150 memories per group). Press the " MEM " button to choose the group. Test results will automatically be saved in the selected group.

14.2 Time/Date setting

Press the "  " button to set the Time/Date mode press the 'MEM' button to change the current datum. Press the "  " button to confirm the current datum. Continue setting date, hour and minute as described. Press the "  " button to confirm the selection and continue in succession (month, day, hour, minute, 12/24 hour).


14.3 Time format

Press the "  " button again to continue to the date format selection, press the 'MEM' button to choose either EU (European format) or US (American format).

14.4 Unit setting

Press the “” button to access the unit settings. Set the format by pressing the “MEM” button.

14.5 Saving

In any mode, press the “START/STOP” button to switch off the unit. All information will be saved.

Note: If the unit is switched on and not used for 3 minutes, it will automatically save all information and switch off.

15. SYSTEM SETTINGS DM592S-DM592P

With the power off, press the ‘Group 1’ or ‘Group 2’ button to activate the system settings. The memory group icon flashes.


15.1 Date/Time setting

Press the “Group1” or “Group2” button to set the Time/Date mode, press the “MEM” button to change the current data.

Press the “Group1” or “Group2” button again to confirm the current datum. Continue setting date, hour and minute as described. Press the ‘Group1’ or ‘Group2’ button to confirm the selection and continue in succession (month, day, hour, minute, 12/24 hour).


15.2 Date Format

Press the ‘Group1’ or ‘Group2’ button again to continue to the date format selection, press the ‘MEM’ button to choose either EU (European format) or US (American format).

15.3 Unit setting

Press the 'Group1' or 'Group2' button to access the unit settings. Set the format by pressing the "MEM" button.

15.4 Voice settings (DM592P only)

Press the "Group1" or "Group2" button to access the voice settings. Set the voice format to ON or OFF by pressing the "MEM" button.

15.5 Saving

In any mode, press the " START/STOP " button to switch off the unit. All information will be saved.

Note: If the unit is switched on and not used for 3 minutes, it will automatically save all information and switch off.

16. MAINTENANCE

The LOGIKO devices by Moretti S.P.A are checked carefully and supplied with EC mark, once launched on the market. For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

17. CLEANING AND DISINFECTION

- DO NOT immerse the device!
- Cleaning: use only a damp cloth. Pay attention to the display
- Disinfection: use a cleanser, avoiding the liquid seep into the device

NOTE: Never use acids, alkalis or solvents such as acetone

- Avoid dropping, slamming, or throwing the unit
- Avoid extreme temperatures. Do not use outdoors
- Do not use petrol, thinners or similar solvents
- Remove batteries when not in operation for an extended period of time
- Do not disassemble product
- It is recommended the performance should be checked every 2 years

18. CONDITIONS OF DISPOSAL

General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the sphygmomanometer in the appropriate disposal areas for recycling.

**18.1 TREATMENT OF BATTERIES (Directive 2002/96/EC):**

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

**18.2 DISPOSAL OF WASTE BATTERIES - Directive 2006/66/EC):**

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

19. TROUBLESHOOTING

| Abnormal phenomenon | Cause analysis | Processing method |
|---------------------------|--|---|
| Abnormal sphygmomanometer | The armband is tied too tight or too loose, Or the arm strap is tied incorrectly; | Roll the armband correctly |
| | Move the arm during measurement or Electronic sphygmomanometer | Stay quiet, keep your arm steady, and do not move the monitor |
| | Speaking, nervous or emotional during measurement | Instead of talking, take deep breaths to calm your mood and relax your body |
| | Incorrect measurement posture | Adjust posture, see "Blood pressure gauge Wearing" |
| | There is interference in charging process or improper operation in measuring process | See operation Instructions. |

The following table shows the error signs that may occur during measurement, possible causes and handling methods. Please measure again using the correct method

| Error display | The cause of the problem | The solution |
|---|---|---|
| Er1 | Can't detect high and low pressure | Please fasten the cuff before measuring |
| Er2 | Cuff too loose or loose | Please fasten the cuff before measuring |
| Er3 | Improper compression caused by arm or body movement | Hold the arm or body still and measure again |
| Er4 | The pressure exceeds 300mmHg | Please fasten the cuff before measuring |
| Er5 | The pressure exceeds 15mmHg for 3 minutes | Check whether the cuff is knotted or the vent valve is blocked. If the problem persists, contact the manufacturer |
|  | Battery dead | Replace the battery or connect the power adapter (if any). |

Note: If you cannot solve the abnormal situation by yourself, you can consult the manufacturer or the manufacturer's designated unit by phone. It is forbidden to disassemble and repair without permission. If necessary, professional maintenance personnel can ask the manufacturer for the list of components and circuit schematic diagram.

20.Q&A

Q: What is the difference between measuring blood pressure at home or at a professional health-care clinic?

A: Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement
Make sure cuff is snug-not too tight or too loose.
2. Improper body position
Make sure to keep your body in an upright position.
3. Feeling anxious or nervous
Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different readings?

A: Blood pressure varies throughout the course of a day. Many factors may affect an individual's blood pressure.

Q: Should I apply the cuff to the left or right arm? What is the difference?

A: Either arm can be used when testing, however, when comparing results, the same arm should be used. Testing on your left arm may provide more accurate results as it is located closer to your heart.

Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.

21. SPECIFICATIONS

| | | |
|---|--|--------------------------|
| Product Description | Digital sphygmomanometer desk type | |
| Model | DM590 - DM591 - DM592S - DM592P | |
| Display | Display LCD DM590: 66.4mm x 43.1mm (2.61" x 1.70") DM591: 83.1mm x 53.1mm (3.27" x 2.09") DM592S - DM592P: 94mm x 90.6mm (3.7" x 3.57") | |
| Measurement Method | Oscillometric | |
| Pressure measurement range | Systolic | 60mmHg~260 mmHg |
| | Diastolic | 30mmHg~200 mmHg |
| | Pressure | 0mmHg~299mmHg |
| | Accuracy | ± 3mmHg |
| | Pulse | 30 ~ 180 Beats/Minute |
| | Accuracy | ± 5% |
| Pressurization | Automatic | |
| Memory | 2x150 memories divided into two groups | |
| Functions | Detection of irregular heartbeat | |
| | WHO classification | |
| | Last 3 average | |
| | Low battery indicator | |
| | Automatic power off | |
| | Backlight (solo DM592S/DM592P) | |
| | Talking (solo DM592P) | |
| Power Source | 3 AAA alkaline batteries or external power adapter (Optional) | |
| Battery Life | Around 2 month with 3 measures per day | |
| Unit Weight | DM590: circa 155 g (without battery) DM591: circa 213 g (without battery) DM592S - DM592P: circa 319 g (without battery) | |
| Unit Dimensions (L x W x H) | DM590 - DM591: 149.8 x 80.4 x 45.7mm (5.90" x 3.17" x 1.80") DM592S - DM592P: 150mm X 108mm X 65mm (5.90" x 4.25" x 2.56") | |
| Cuff Circumference | Around:140mmX568mm (arm diameter: 22-36cm) | |
| Operating environment | Temperature | 10°C ~ 40°C (50°F~104°F) |
| | Humidity | 15%~93%RH |
| | Pressure | 800hPa~1060hPa |
| Storage environment | Temperature | -25°C~55°C (-13°F~131°F) |
| | Humidity | ≤93% RH |
| Degree of protection against dusts and liquids | IP21, indore use only | |
| Classification | Internal powered equipment, Type BF _{II} , cuff in the applied part | |

Specifications are subject to change without notice.

**22. ELECTROMAGNETIC COMPATIBILITY DECLARATION**

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

| Guidance and declaration of manufacturer-electromagnetic emissions | | |
|--|-------------------|---|
| The device is intended for use in the electromagnetic environment specified below. | | |
| The customer or the user of the device should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment -guidance |
| Radiated emission CISPR 11 | Group 1, class B. | The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| Conducted emission CISPR 11 | Group 1, class B. | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies | |

Table 2

| Guidance and declaration of manufacturer-electromagnetic immunity | | | |
|---|--|--|--|
| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | | |
| IMMUNITY test | IEC 60601 test level | Compliance level | Electromagnetic environment -guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | ± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrostatic transient/burst IEC 61000-4-4 | ± 2 kV , 100kHz, for AC power port | ± 2 kV , 100kHz, for AC power port | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 0.5kV, ± 1kV (differenziale) | ± 0.5kV, ± 1kV (differenziale) | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle | 0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle | Mains power quality should be that of a typical commercial or hospital environment. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m; 50Hz or 60Hz | 30 A/m; 50Hz or 60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

| | | | |
|--|--|--|---|
| <p>Radiated RF EM fields IEC 61000-4-3</p> | <p>3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz</p> | <p>3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80MHz to 800 MHz 800 MHz to 2.7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p> |
| <p>Conducted disturbances Induced by RF fields IEC 61000-4-6</p> | <p>3 V in 0.15 MHz- 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz</p> | <p>3 V in 0.15 MHz- 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p> |

Table 3

| Guidance and declaration of manufacturer-electromagnetic immunity | | | | | | |
|--|------------|--|--|-------------------|--------------|---------------------------|
| <p>Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-type Fully Automatic Digital Blood Pressure Monitor has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/ or systems as recommended below.</p> | | | | | | |
| Test frequency (MHz) | Band (MHz) | Service | Modulation | Maximum power (W) | Distance (m) | Immunity test level (V/m) |
| 385 | 380-390 | TETRA 400 | Pulse modulation 18Hz | 1.8 | 0.3 | 27 |
| 450 | 430-470 | GMRS 460 FRS 460 | FM ± 5 kHz deviation 1 kHz sine | 2 | 0.3 | 28 |
| 710 | 704-787 | LTE Band 13, 17 | Pulse modulation 217Hz | 0.2 | 0.3 | 9 |
| 745 | | | | | | |
| 780 | | | | | | |
| 810 | 800-960 | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 | Pulse modulation 18Hz | 2 | 0.3 | 28 |
| 870 | | | | | | |
| 930 | | | | | | |
| 1720 | 1700-1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation 217Hz | 2 | 0.3 | 28 |
| 1845 | | | | | | |
| 1970 | | | | | | |
| 2450 | 2400-2570 | Bluetooth, WLAN, 802.11 b/g/n,RFID 2450,LTE Band 7 | Pulse modulation 217Hz | 2 | 0.3 | 28 |
| 5240 | 5100-5800 | WLAN 802.11 a/n | Pulse modulation 217Hz | 0.2 | 0.3 | 9 |
| 5500 | | | | | | |
| 5785 | | | | | | |

Table 4

| Recommended separation distances between portable and mobile RF communications equipment and the device |
|--|
| <p>The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.</p> |

| Rated maximum output power of transmitter | Separation distance according to frequency of transmitter | |
|---|---|---|
| | 80 MHz to 800 MHz $d=[3.5/E1]\sqrt{P}$ | 80 MHz to 2.7 GHz $d=[7/E1]\sqrt{P}$ |
| 0.01 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.73 |
| 1 | 1.2 | 2.3 |
| 10 | 3.8 | 7.3 |
| 100 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

24. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual. Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use. Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

24.1 Repairing

- Warranty repair

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

- Repair not covered by warranty

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

- Non-defective products

The customer will be notified if, after examining and tried a returned product, Moretti concludes that the product is not defective. The product will be returned to the customer and will be liable to pay the costs caused by the return.

24.2 Spare-parts

Moretti's genuine parts are guaranteed for 6 months with effect from the day you receive the replacement.

**24.3 Exempt clauses**

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user.

**LOGIKODIGIT****WARRANTY CERTIFICATE**

Product _____

Purchasing date _____

Authorized dealer _____

Street _____ Place _____

Sold to _____

Street _____ Place _____

Joytech Healthcare Co., Ltd
No.365, Wuzhou Road, Yuhang Economic Development
Zone, Hangzhou City, 311100 Zhejiang P.R. China

Shanghai International Holding Corp. GmbH(Europe)
Eiffelstrasse 80, 20537 Hamburg, Germany

Distributed by: MORETTI S.P.A. Via Bruxelles, 3 - Melegnano 20122 Cavriaglia (Arezzo)

Tel. +39 055 96 21 11 - www.morettispa.com email: info@morettispa.com

MADE in P.R.C.

** Please consult our website for the latest available version of the user manual*

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MORETTI S.p.A.