ITALIANO

MA LAD462-LAD463-LAD464 Piuma up 02 A\_ITA\_09 2020

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MADE IN P.R.C.

MA LAD462-LAD463-LAD464 Piuma up 02 A\_ENG\_09 2020

## **PIUMA UP ALTERNATING PRESSURE REDISTRIBUTION SYSTEM**

---

### **INSTRUCTION MANUAL**

ENGLISH

## INDEX

<b>1. CODES</b> .....	<b>PAG.3</b>
<b>2. INTRODUCTION</b> .....	<b>PAG.3</b>
<b>3. INTENDED USE</b> .....	<b>PAG.3</b>
<b>4. DECLARATION OF CONFORMITY</b> .....	<b>PAG.3</b>
<b>5. STANDARDS</b> .....	<b>PAG.4</b>
<b>6. GENERAL WARNINGS</b> .....	<b>PAG.4</b>
6.1 Electromagnetic compatibility warnings.....	pag. 4
<b>7. SYMBOLS</b> .....	<b>PAG.5</b>
<b>8. KIT CONTENTS</b> .....	<b>PAG.6</b>
<b>9 GENERAL DESCRIPTION</b> .....	<b>PAG.6</b>
9.1 Pump unit .....	pag. 6
9.2 Front panel .....	pag. 6
9.2.1 Pressure Adjust Knob .....	pag. 6
9.2.2 Main Power Switch .....	pag. 6
<b>10. INSTALLATION</b> .....	<b>PAG.6</b>
<b>11. BEFORE USE</b> .....	<b>PAG.7</b>
<b>12. USING WARNINGS</b> .....	<b>PAG.7</b>
<b>13. OPERATION</b> .....	<b>PAG.8</b>
13.1 General operation .....	pag. 8
<b>14. CLEANING</b> .....	<b>PAG.9</b>
<b>15. STORAGE</b> .....	<b>PAG.9</b>
<b>16. MAINTENANCE</b> .....	<b>PAG.10</b>
<b>17. CONDITIONS OF DISPOSAL</b> .....	<b>PAG.10</b>
17.1 General conditions of disposal .....	pag. 10
17.2 Information for the correct disposal of the product in accordance with European Directive 2012/19/UE: .....	pag. 10
<b>18. SPARE PARTS/ACCESSORIES</b> .....	<b>PAG.10</b>
<b>19. TROUBLESHOOTING</b> .....	<b>PAG.10</b>
<b>20. GUIDANCE AND MANUFACTURER'S DECLARATION</b> .....	<b>PAG.11</b>
20.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions.....	pag. 11
20.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity .....	pag. 11
21.2 Mattress.....	pag. 14
<b>22. WARRANTY</b> .....	<b>PAG.14</b>
<b>23. REPAIR</b> .....	<b>PAG.14</b>
23.1 Warranty repairs .....	pag. 14
23.2 Repair not covered by warranty.....	pag. 14
23.3 Non-defective device .....	pag. 15
<b>24. SPARE PARTS</b> .....	<b>PAG.15</b>
<b>25. EXEMPT CLAUSES</b> .....	<b>PAG.15</b>

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL  
of 5 April 2017 concerning medical devices

## 1. CODES

<b>LAD462</b>	Kit mattress and pump with pressure adjustment Piuma UP (bubble pad)
<b>LAD463</b>	Kit mattress and pump with pressure adjustment Piuma UP (single cells pad)
<b>LAD464</b>	Kit mattress and pump with pressure adjustment Piuma UP (bubble pad increased height)

## 2. INTRODUCTION

Thank you for purchasing a LEVITAS mattress and pump Piuma UP Kit (single cell pad) by Moretti. The pump and mattress kits by Moretti are designed and manufactured to meet all your needs for a practical, correct and safe use. This manual contains small suggestions for proper use of the device you have chosen and valuable advice for your security. You may read the complete manual before using the kit pump and mattress, in case of doubts please contact your dealer, who will help and advise you properly.

### NOTE, CAUTION AND WARNING STATEMENTS:

NOTE: Indicate some tips.

CAUTION: Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING: Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

## 3. INTENDED USE

The kit pump and mattress have been designed for the bedsore treatment and prevention, which may occur during long hospital and/or house confinement.

### CAUTION!

- Do not use the product for a purpose not indicated in this manual
- Only qualified personnel or trained for the treatment and prevention of bedsore, can operate with the device
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.

## 4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of ANTIDECUBITUS PIUMA UP complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

**Note:** Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

## 5. STANDARDS

The system has been tested and successfully approved to the following standards: EN 60601-1:2005+A1:2012  
 EN 60601-1-2:2007  
 IEC 60601-1-11:2010  
 IEC 60601-1-6:2010+A1:2013  
 IEC 62366:2007+A1:2014  
 EN ISO 13485:2012+AC:2012  
 EN ISO 14971:2012  
 EN ISO 15223-1:2012  
 EN 1041:2008  
 EN 71-3:2013+A1:2014

## 6. GENERAL WARNINGS

- For proper use of the device, refer to the following manual carefully
- Advise with physician or therapist to determinate the correct adjustment and the correct use of the device
- Keep the packed kit away from heat sources
- SERVICE LIFE- the device use limit is defined by the wear of the parts
- DO NOT allow children to play on or operate the pump
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

### 6.1 Electromagnetic compatibility warnings

This equipment has been tested and it complies with medical devices limits belonging to EN 60601-1-2:2007. These limits have been designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guaranteed that the interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference following one or more of the suggestions below:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other devices are connected
- Consult the manufacturer or field service technician for help

**7. SYMBOLS**

	Product code
	Unique Device Identification
	CE mark
	Manufacturer
	Batch Lot
	Read the instruction manual
	Medical Device
	Conditions of disposal
	Attention, should read the instructions
	Production dates
	Temperature Limitation
	Waste Electrical & Electronic Equipment (WEEE)
<b>IP21</b> 	Protected against solid foreign objects of 12.5 mm and greater. Protection against vertically falling water drops

 BF Type

 Class II Equipment

| ON

○ OFF

## 8. KIT CONTENTS


The Piuma UP series is composed by 2 kind of kits which can satisfy all the requirement about treatment and prevention of complications caused by long hospital and/or house confinement, in particular bedsores.

- LAD462** Kit mattress and pump with pressure adjustment Piuma UP (bubble pad)  
**LAD463** Kit mattress and pump with pressure adjustment Piuma UP (single cells pad)  
**LAD464** Kit mattress and pump with pressure adjustment Piuma UP (bubble pad increased height)

ENGLISH

## 9 GENERAL DESCRIPTION

### 9.1 Pump unit

### 9.2 Front panel

#### 9.2.1 Pressure Adjust Knob

Pressure adjust knob controls the air pressure output. Please consult the physician for a suitable setting

#### 9.2.2 Main Power Switch

To turn the pump unit on/off

## 10. INSTALLATION

Unpack the box and check the package contents for completeness. If there are any damages, please immediately contact your retailer.

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard

**ATTENTION!**

We recommend to place the mattress with interchangeable elements over a mattress with a height of at least 5 cm

2. Hang the pump onto the bed rail (footboard side), or place it on a flat surface
3. Connect air hose connectors from air mattress to the pump unit

**ATTENTION!**

Ensure the air hoses are not kinked or tucked under mattress

4. Plug the power cord into the electrical outlet

- ⚠ NOTE: 1. Make sure the pump unit is suitable for the local power voltage  
2. The plug is also served to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device

**ATTENTION!**

The pump should only be used with mattress recommended by the manufacturer. Do not use it for any other purpose

5. Turn the switch ON

Several installation tips are listed below:

After installation, the extra length of the power cord, if any, should be neatly arranged to avoid any tripping accidents. The EQUIPMENT should be firmly placed at position where users/doctors can access easily.

## 11. BEFORE USE

Check the correct assembly, paying particular attention to air hoses connector.

## 12. USING WARNINGS

### DANGER - To reduce risk of electrocution

1. Always unplug this product immediately after using
2. Do not use while bathing
3. Do not place or store this product where it can fall or be pulled into a tub or sink
4. Do not place in or drop into water or other liquid
5. Do not reach for a product that has fallen into water. Unplug immediately

### WARNING - To reduce risk of burns, electrocution, fine or injury to people

1. Evaluate patients for entrapment risk according to facility protocol and monitor patient appropriately
2. This system is not intended to be used with patients who have spinal cord injury.
3. Close supervision is necessary when this product is used on or near children.  
Electrical burns or choking accident may occur if a child swallow a small part detached from the device
4. Use this product only for its intended use as described in this manual. Do not use attachment not recommended by the manufacturer




5. Never operate this product if it has a damaged cord or plug, if it is not working properly, or if it has been dropped, damaged, or dropped into water. Return the product to a service centre for examination and repair
6. Keep the cord away from heated surfaces
7. Never block the air openings of this product or place it on a soft surface, such as a bed or couch, where their openings may be blocked. Keep the air openings free of lint, hairs, and other similar particles
8. Never drop or insert any object into any opening or hose
9. Do not modify this equipment without authorization of the manufacturer
10. Mattress cover has passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately
11. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance to 3,3m between devices or turn off the mobile phone

### 13. OPERATION


⚠ NOTE: Always read the operating instruction before use

#### 13.1 General operation

1. Open the package: it is recommended to use a cutter with short blade or a tool not to damage the mattress.
2. The overlay mattress is rolled up: take it out of the box paying attention to the compressor. Unroll it on the bed mattress.
3. Fold and put the flaps of the overlay mattress between the mattress and the metal net, to secure it and prevent unintentional movements.
4. Take out and unpack the compressor: then hang it to the footboard of the bed by the hooks on the back of the compressor, or place it on a flat surface.


5. Make sure that the switch "ON/OFF" of the compressor is OFF.
6. Connect the tubes placed at the bottom of the overlay mattress to the compressor, checking they are connected correctly.

8. Turn the knob clockwise, putting the narrow's knob in proximity of "+" symbol. In this way the compressor will inflate the mattress more quickly.
9. Connect the compressor plug to the 220V outlet.
10. Switch on the compressor with the "ON/OFF" switch and wait about 25 minutes, so that the compressor inflate the mattress completely.

⚠ NOTE: every time the mattress is set up for use, it is recommended that the pressure first to be set to the max. The user/career can then adjust air mattress weight levels to the desired softness after set up has been completed

#### 14. CLEANING

It is important to follow the cleaning procedures before using the equipment on human bodies; otherwise, patients and/or doctors may have the possibility of getting infection.

⚠ CAUTION: Do not immerse or soak pump unit.

Wipe the pump unit with a damp cloth and a mild detergent. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.

⚠ CAUTION: Do not use phenolic based product for cleaning.

⚠ CAUTION: After cleaning, dry the mattress without direct exposure sunlight.

#### **WARNING!**

**Unplug the pump before starting the cleaning procedure**

#### 15. STORAGE

- Fold the mattress accurately to prevent creases
- Roll from the head end towards the foot end with the valve open

For single cells mattress:

1. Lay the bubble pad or mattress on a flat surface and upside down
2. Roll-up the mattress from the head end towards the foot end
3. Foot-end strap can then be stretched around the rolled pad/ mattress to prevent unrolling

⚠ CAUTION: Do not fold, crease or stack the mattresses.

## 16. MAINTENANCE

The LEVITAS devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market.

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function. Otherwise, service and inspection of the devices generally should not be required.

1. Check main power cord and do not plug it if there is an abrasion or excessive wear
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly
3. Check the airflow from the air hose. The airflow should alternate between each connector every half-cycle time
4. Check the air hoses if there is kink or breaks. For replacement, please contact Apex Medical Corp. or your suppliers

## 17. CONDITIONS OF DISPOSAL

### 17.1 General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the thermometer in the appropriate disposal areas for recycling.

### 17.2 Information for the correct disposal of the product in accordance with European Directive 2012/19/UE:

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin the product bears the barred bin symbol.

## 18. SPARE PARTS/ACCESSORIES

For spare parts and accessories refer to main catalogue

## 19. TROUBLESHOOTING

If your questions can't be answered with above information, please contact your local agent directly. They might require a technician to take care of the problem

PROBLEM	SOLUTION
<b>Power is not ON</b>	<ul style="list-style-type: none"> <li>• Check if the plug is connected to mains</li> <li>• Turn "ON" the compressor switch</li> </ul>
<b>The mattress doesn't inflate</b>	<ul style="list-style-type: none"> <li>• Check if the connection between mattress and pump unit is tightly secured</li> <li>• Check if there is any leakage on air cells</li> </ul>
<b>Mattress form is loose</b>	<ul style="list-style-type: none"> <li>• Check if all the snap buttons or straps of mattress are all securely fastened</li> <li>• Check if the mattress is fixed to the bed frame by straps</li> </ul>
<b>No air produced from some air outlets of the air tube connector</b>	<ul style="list-style-type: none"> <li>• This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time</li> </ul>

 **NOTE:** If the pressure level is consistently low, check for any leakage (tubes or air hoses). If necessary, replace any damaged tubes or hoses or contact your local qualified dealer for repair.

## 20. GUIDANCE AND MANUFACTURER'S DECLARATION

### 20.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.


EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

### 20.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	± 2kV for power supply line ± 1kV for input/out line	± 2kV for power supply line ± 1kV for input/out line	Mains power quality should be that of atypical commercial or hospital environment
Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT)for 0,5 cycle  40 % UT (60 % dip in UT)for 5 cycles  70 % UT (30 % dip in UT)for 25 cycles  <5 % UT (>95 % dip in UT)for 5 sec	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of atypical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment
 <b>NOTE:</b> UT is the A.C. mains voltage prior to the application of the test level			

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bandsa	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 150kHz to 80MHz 150kHz to 80MHz 80 MHz to 2.5G MHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). b Field strengths from fixed RF transmitters, as determined by an electromagnetic site surveyc, should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

 NOTE: At 80 MHz and 800 MHz, the higher frequency range applies

 NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

 NOTE: The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land / mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m  
 Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2,5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

⚠ NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

⚠ NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

## 21. TECHNICAL FEATURES

### 21.1 Pump

MODEL	LAD462	LAD463 - LAD464
<b>Power supply</b>	~ 220V 50/60Hz	~ 220V 50/60Hz
<b>Fuse T1AL</b>	250V	250V
<b>Power cord</b>	2,20 m	2,20 m
<b>Pressure output</b>	40~100mmHg	40~100mmHg
<b>Air release</b>	6/7LPM	6/7LPM
<b>Cycle time</b>	6 min	12 min
<b>Noise</b>	32-34 dB	32-34 dB
<b>Dimensions</b>	27x13x10,5 cm	27x13x10,5 cm
<b>Weight</b>	~ 1,1 Kg	~ 1,1 Kg
<b>Environment temperature</b>	Operation: +10°C/+40°C Storage: -18°C/+43°C Shipping: -18°C/+43°C	Operation: +10°C/+40°C Storage: -18°C/+43°C Shipping: -18°C/+43°C
<b>Environment humidity</b>	Operation: 10%/75% non condensing Storage: 10%/95% non condensing Shipping: 10%/95% non condensing	Operation: 10%/75% non condensing Storage: 10%/95% non condensing Shipping: 10%/95% non condensing
<b>Atmospheric pressure</b>	70/101,325 KPa	70/101,325 KPa
<b>Classification</b>	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)	Class II, Type BF, IP21 Applied part: air mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)

## 21.2 Mattress

MATTRESS	MATTRESS BUBBLE PAD - LAD462 MATTRESS BUBBLE PAD INCREASED HEIGHT - LAD464	ELEMENT MATTRESS LAD463
Model	Bubble pad	Overlay mattress
Components	130 bubbles (LAD462) 260 bubbles (LAD464)	18 cells
Material	Medical-grade PVC	Medical-grade PVC Cell: Nylon/PVC
Dimension	~ 200x90x7,6 cm (LAD462) ~ 200x90x15,2 cm (LAD464)	~ 190x80x11 cm
Cell high	-	~ 10,2 cm
Flaps	50 cm	50 cm
Air tube	1,5 m	1,5 m
Weight	2,3 Kg (LAD462) 4,6 Kg (LAD464)	3,8 Kg
Max user's weight	135 Kg (LAD462) 145 Kg (LAD464)	145 Kg
Thickness	0,8 mm	-

## 22. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for a period of (two) 2 years from the purchasing date, subject to any limitations and exclusions specified below. The warranty does not apply in case of damages caused by improper use, abuse or modification of the product, and the warranty will not be valid if the instructions for use are not strictly followed. The correct intended use of the device is specified in this manual.

Moretti is not responsible for damages resulting from personal injury or otherwise caused by or in relation to wrong installation or improper use. Moretti does not guarantee the products against damage or defects in the following conditions: natural disaster, not authorized maintenance or repair, damage due to power supply problems (where provided) use of spare parts not supplied by Moretti, improper use, not authorized alteration, shipping damage (other than original shipment by Moretti), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

## 23. REPAIR

## 23.1 Warranty repairs

In the event that a product Moretti to be defective in material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, in its sole discretion, can repair or replace the item under warranty at a Moretti dealer or at Moretti location. The labor cost for the repair may be charged to Moretti if the repair is covered by warranty. Repair or replacement does not renew or extend the warranty.

## 23.2 Repair not covered by warranty

A product not covered by warranty can be returned for repairing only after receiving prior authorization from Moretti customer service. The labor and shipping costs for a repair not covered by warranty will be fully borne by the customer or dealer. Repairs of products not covered by warranty are warranted for six (6) months from the date you receive the product repaired.

### 23.3 Non-defective device

The customer will be informed if, after having examined and tested returned product, Moretti claims that the product is not defective. In this case the good will be returned to the customer, all shipping costs will be paid by customer.

### 24. SPARE PARTS

For spare parts only refer to MORETTI master catalogue.

### 25. EXEMPT CLAUSES

Except as expressly specified in this warranty and within the limits of the law, Moretti does not offer any representation, warranty or condition, express or implied, including any representations, warranties or conditions of merchantability, fitness for a particular purpose, non-infringement and interference. Moretti does not warrant that the use of the product Moretti will be uninterrupted or error free. The duration of any implied warranties that may be imposed by law are limited to the warranty period, within the limits of the law. Some states do not allow limitations on how long an implied warranty lasts or the exclusion or limitation of incidental or consequential damages in connection with products for consumers. In such states or countries, certain disclaimers and limitations of this warranty may not apply to you. This warranty is subject to change without notice

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