ITALIANO

MA CM900x_910x_02 C_ITA_04 2021

CARROZZINA ELETTRICA "ESCAPE"

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INDICE

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**MORETTI S.P.A.**

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com**MADE IN TAIWAN**

ENGLISH

MA CM900x_910x_02 C_ENG_04 2021

"ESCAPE" ELECTRIC WHEELCHAIR

INSTRUCTION MANUAL

INDEX

1. CODES	PAG.3
2. INTRODUCTION	PAG.3
3. INTENDED USE	PAG.3
4. EU DECLARATION OF CONFORMITY	pag.3
4.1 Applicable regulations and directives.....	pag. 4
5. GENERAL PRECAUTIONS	PAG.4
6. GENERAL INFORMATION ON ELECTROMAGNETIC DISTURBANCE	PAG.4
7. SYMBOLS USED	PAG.5
8. SAFETY	PAG.6
9. WARNINGS ON USE	PAG.6
9.1 General warnings	pag. 6
9.2 Changes	pag. 7
9.3 Before every use.....	pag. 7
9.4 Weight limits.....	pag. 7
9.5 Inflating the tyres	pag. 7
9.6 Temperature.....	pag. 8
10. COMPONENTS	PAG.8
11. ASSEMBLY INSTRUCTIONS	PAG.9
11.1 Installing the controller and the connection cables	pag. 9
11.2 Installing the footplates.....	pag. 10
11.3 Installing the batteries and connecting the cables	pag. 10
11.4 Installing the rear push bar (only for CM900-xx)	pag. 11
11.5 Adjusting comfort	pag. 11
12. OPERATION	PAG.13
12.1 Functioning of the lights (CM900L/CM910L)	pag. 17
12.2 Reclining backrest	pag. 18
13. BATTERY & CHARGING	PAG.18
13.1 Charging	pag. 19
13.2 Charging instructions	pag. 19
13.3 Attention.....	pag. 19
14. REPLACING THE WHEELS	PAG.20
15. MAINTENANCE AND REPAIR OF THE ELECTRIC WHEELCHAIR	PAG.21
15.1 Daily checks.....	pag. 22
15.2 Weekly checks.....	pag. 22
15.3 Six-monthly checks	pag. 22
16. CARE AND MAINTENANCE	PAG.24
16.1 Tyre pressure.....	pag. 24
16.2 Cleaning and disinfection.....	pag. 24
16.3 Battery terminal connections	pag. 24
16.4 Cables.....	pag. 24
16.5 Axle supports and motor/differential.....	pag. 24
16.6 Control panels, battery charger and rear electronics	pag. 24
16.7 Storing the electric wheelchair	pag. 25
17. TROUBLESHOOTING	PAG.25
18. DISPOSAL CONDITIONS	PAG.26
18.1 General conditions of disposal	pag. 26
18.2 Instructions for correct disposal in accordance with European directive 2012/19/EU:.....	pag. 26
18.3 Disposal of exhausted batteries - (Directive 2006/66/EC):	pag. 26
19. REPLACEMENT PARTS AND ACCESSORIES	PAG.26
20. TECHNICAL SPECIFICATIONS	PAG.27
21. WARRANTY	PAG.29
22. REPLACEMENT PARTS	PAG.29
23. NON-LIABILITY CLAUSE	PAG.29

CE Class I medical device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 on medical devices

1. CODES

- CM900-41 "Escape DX" electric wheelchair with 41 cm seat
- CM900-46 "Escape DX" electric wheelchair with 46 cm seat
- CM900L-41 "Escape DX" electric wheelchair with lights and 41 cm seat
- CM900L-46 "Escape DX" electric wheelchair with lights and 46 cm seat
- CM910-46 "Escape LX" electric wheelchair with reclining backrest and 46 cm seat
- CM910L-46 "Escape LX" electric wheelchair with reclining backrest, lights and 46 cm seat

2. INTRODUCTION

Thank you for choosing an **ARDEA MOBILITY** electric wheelchair by Moretti. Moretti electric wheelchairs have been designed and manufactured to satisfy all your needs for practical, correct and safe use. This manual contains useful suggestions for using your device properly and safely. Please read this manual in full before you use the electric wheelchair. Should you have any queries, please contact your retailer for suitable advice and assistance.

Note: Check to ensure that no parts of the wheelchair have been damaged during shipment. Do not use the product if it is damaged and contact the retailer for further instructions.

3. INTENDED USE

The electric wheelchair is intended for people with mobility impairments.

- Do not use the product for any purpose other than that specified in this manual
- The electric wheelchair must always be used by people of sound body and mind and not under the influence of drugs or alcohol.
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than that specified in this manual
- The manufacturer reserves the right to modify the device and its relevant manual without prior notice, in order to improve the product's features

4. EU DECLARATION OF CONFORMITY

MORETTI SpA declares, under its own exclusive responsibility, that the products manufactured and sold by MORETTI SpA in the ELECTRIC WHEELCHAIRS product family conform with the applicable provisions of Regulation (EU) 2017/745 of 5 April 2017 on MEDICAL DEVICES.

To this end, MORETTI SpA guarantees and declares as follows, under its own exclusive responsibility:

1. The devices in question satisfy the general safety and performance requirements set out in Annex I to Regulation (EU) 2017/745, as prescribed by Annex IV to said regulation.
2. The devices in question ARE NOT MEASUREMENT TOOLS.
3. The devices in question ARE NOT INTENDED FOR CLINICAL INVESTIGATIONS.
4. The devices in question are sold in NON STERILE packing.
5. The devices in question should be considered as belonging to class I, in accordance with the rules set out in Annex VIII to said regulation.
6. MORETTI SpA maintains and provides to the competent authorities, for at least ten years from the date of manufacture of the last production lot, the technical documentation proving conformity with Regulation (EU) 2017/745.

Note: The complete product codes, the manufacturer's single registration number (SRN), the basic UDI-DI code and any references to standards used are indicated in the EU Declaration of Conformity that MORETTI SPA issues and provides through its own channels.

4.1 Applicable regulations and directives

- EN 12184:2014
- EN 60601-1:2007;
- EN 60601-1-2:2006.

5. GENERAL PRECAUTIONS

- Please consult this manual carefully for correct use of the device
- Always consult your doctor or therapist for correct use of the device.
- Keep the packaged product away from sources of heat, as the packaging is cardboard
- Avoid exposure of the device to direct sunlight or other sources of heat
- The lifetime of the device depends on wear and tear of non-replaceable and/or non-replaceable parts
- Always take care when children are present
- The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

6. GENERAL INFORMATION ON ELECTROMAGNETIC DISTURBANCE

The rapid growth of electronics, especially in communications, has saturated the environment with electromagnetic radiation (EM) emitted by television, radio and communication signals. These EM waves are invisible and are stronger the closer you are to their source. All electrical conductors act as antennas for EM signals and, to varying extents, all electric wheelchairs and scooters may be affected by electromagnetic interference (EMI). EMI may cause unintentional movements and/or irregular control of the vehicle.

The electromagnetic energy emitted by sources such as:

- radio stations;
- TV stations;
- amateur radio stations;
- shop alarm systems;
- mobile telephones and cordless telephones can interfere with electric wheelchairs and electric scooters;

The interference can cause the electromagnetic brake to release, unintentional movement of the device and movement in the wrong direction. It can also permanently damage the device's electronic control unit. In our everyday environment, there are a certain number of intense sources of EMI.

The intensity of the EM energy can be measured in Volts per metre (V/m). Your electric wheelchair is equipped with a protection against EMI which has been tested and certified as compliant with current international regulations, up to a certain intensity. This is called the "level of immunity" The greater the level of immunity, the better the device is protected. At this time, technology can provide a level of immunity of at least 20 V/m; this is sufficient protection against the most common sources of EMI. However, we consider that the following precautions can further minimise the risk.

- A) Do not use portable radio transmitters, CB radios and avoid switching on mobile telephones or other personal communications equipment when the electric wheelchair is switched on.
- B) Do not approach powerful radio and TV transmitters too closely.
- C) If an unintended movement or unintended release of the brake occurs, turn the key to the OFF position as soon as you can do so safely.
- D) Installing additional electrical equipment or accessories and modifying the motor of your electric wheelchair may make it more susceptible to EMI and hence make its use hazardous. Note that doing so also voids the warranty.
- E) In any case, report any unintentional movement of the vehicle to your retailer or to MORETTI SPA, note it down in this manual and note down the fact that there is a source of electromagnetic radiation in your vicinity, if any.

7. SYMBOLS USED

-  Product code
-  Serial number
-  Unique device identification
-  CE marking
-  Manufacturer
-  Batch Lot
-  Read the instruction manual
-  Medical Device
-  Disposal conditions
-  Warning
-  Production date
-  Disposal of the product according to directive EC/19/2012
- IP54** Protection against dust and water spray
-  Type B applied part

8. SAFETY

- Always use the safety belt and keep your feet resting on the footplates at all times.
- Do not use the electric wheelchair under the influence of alcohol.
- Do not use electronic radio transmitters such as walkie-talkies or mobile telephones.
- Make sure there are no obstacles nearby when you are driving the electric wheelchair.
- Do not make a sharp turn or brake suddenly when you are driving the electric wheelchair.
- Do not drive the electric wheelchair in traffic.
- Do not attempt to mount curbs higher than those indicated in the technical specifications.
- Do not put your hands or legs outside the electric wheelchair while driving.
- Do not drive the electric wheelchair on slippery roads or when it is snowing.
- Do not allow unsupervised children to play near the electric wheelchair when the batteries are charging.

Do not use the electric wheelchair without first reading and understanding this manual completely.

9. WARNINGS ON USE**9.1 General warnings**

1. Do not use the electric wheelchair on public highways. Please be aware that it may be difficult for drivers to see you when you are sitting on the electric wheelchair. Always use pedestrian routes. Comply with the rules of pedestrian traffic. Wait for the way to be clear, then proceed with extreme caution.
2. Before getting into or out of your electric wheelchair, always check that it is switched off and that the battery charger is disconnected, to avoid injury to yourself or to others.
3. Always check that the drive wheels are in drive mode before driving. Do not switch off the electric wheelchair when it is still moving, as it would stop very suddenly.
4. Do not use this product or any other optional accessory available without first reading and understanding this manual completely. If you do not understand the warnings, the precautions or the instructions, contact the retailer or a healthcare operator before attempting to use the device, to avoid injury or damage.
5. In certain cases, such as certain medical conditions, the user of the electric wheelchair will need practice in its use in the presence of a qualified assistant, who may be a relative or a professional qualified to assist the user, in use of the electric wheelchair during various everyday activities.
6. Avoid lifting or moving the electric wheelchair by gripping one of its removable parts, including the armrests, seat or fairing. This could cause personal injury and/or damage the wheelchair.
7. Avoid using the electric wheelchair beyond its limits, as described in this manual.
8. Do not sit on the electric wheelchair when it is inside a moving vehicle.
9. Keep your hands away from the wheels (tyres) while driving the electric wheelchair. Please remember that loose clothing could get caught in the wheels.
10. Consult your doctor if you are taking drugs or if you have specific physical limitations. Certain drugs and certain limitations could affect the ability to use the electric wheelchair safely.
11. Always check that the wheel locking/unlocking lever is locked or unlocked.
12. Do not remove the roll-over protection wheels if they are present on the electric wheelchair.
13. Contact with tools could cause an electric shock; do not connect an electrical extension cord to the AC / DC converter or to the battery charger.
14. Although the electric wheelchair is designed for slopes with a gradient up to 10° (17.5%), try to avoid slopes with a gradient over 6° (10%) and take extreme care when turning on a slope.

15. Do not drive on a slope over the limit of the electric wheelchair.
16. Avoid reversing off steps, pavements or other obstacles. This could cause the electric wheelchair to fall or tip over.
17. Always reduce speed and maintain a stable centre of gravity on sharp bends. Do not use the maximum power on sharp bends.
18. Using the electric wheelchair in the rain, snow or foggy conditions and on icy or slippery surfaces could damage the electrical wiring.
19. Never sit in the electric wheelchair when it is being used with any type of lift or lifting device. The electric wheelchair is not designed for this use. The manufacturer is not liable for any damage or injury caused by improper use.

9.2 Changes

The electric wheelchairs in the ARDEA MOBILITY range by Moretti have been designed and manufactured to satisfy all your needs for practical, correct and safe use. Do not alter, add, remove or disable any part or function of the electric wheelchair under any circumstances. This could cause damage to the electric wheelchair and/or personal injury.

1. The control parameter must only be altered by authorised technicians, for safety reasons.
2. It is recommended to carry out a safety check before every use of the electric wheelchair, in order to ensure that it is functioning safely.

9.3 Before every use

1. Check the tyres are correctly inflated, if present.
2. Check all the electrical connections and make sure they are properly connected and not corroded.
3. Check all the cables and make sure they are connected properly.
4. Check the brakes.

9.4 Weight limits

1. Check the information on maximum capacity in the table of the technical specifications in this manual or on the label on the device. The power of the electric wheelchair is assessed according to the capacity.
2. Comply with the weight limits indicated for your electric wheelchair. Exceeding the maximum capacity invalidates the warranty. The manufacturer cannot be held liable for injury or damage to property caused by failure to comply with the weight limit.
3. Do not carry passengers on the electric wheelchair, as this could alter the centre of gravity and cause tipping over or falls.

9.5 Inflating the tyres


1. If the electric wheelchair has tyres, check the air pressure at least once a week.
2. Correct inflation pressure prolongs the duration of the tyres and guarantees better performance while driving.
3. Do not underinflate or overinflate the tyres. It is extremely important for the tyre pressure to remain between 30-25 psi (2-2.4bar) at all times.
4. Inflating the tyres with an unregulated air source can cause them to be overinflated and burst.

9.6 Temperature

1. Certain parts of the electric wheelchair are affected by changes in temperature. The controller only functions at a temperature ranging between -8 C° and 50 C° .
2. The batteries could freeze at extremely low temperatures and prevent your electric wheelchair from operating. At extremely high temperatures, the electric wheelchair might only function at lower speeds, due to a safety feature of the controller that prevents damage to the motor and other electrical components.
3. Exposing your electric wheelchair to water could cause damage to the electronic or mechanical components. Water can cause an electronic malfunction and/or premature corrosion of the electric components and the frame.

10. COMPONENTS

Your electric wheelchair is shipped partially disassembled, for greater protection during transport. After unpacking the wheelchair, check that you have received all the main components, as listed below. (fig.1)



1. Main structure with motor
2. Battery charger
3. Controller
4. Controller attachment bar
5. Left footplate
6. Right footplate
7. Battery attachment plate
8. Cable
9. Batteries

11. ASSEMBLY INSTRUCTIONS


11.1 Installing the controller and the connection cables

Slide the fixed bar of the controller (# 4) into the tube of the frame sufficiently (fig.2), then tighten the knob (fig.3).

(fig.2) Fixed bar of the controller (fig.3) Knob


Mount the controller (# 3) (fig.4) onto the bracket of the right or left armrest, as needed, with two screws (fig.4A). Then connect the attachment of the plug to the jack of the control unit. (fig.5)

(fig.4)

(fig.4A)

(fig.5) Control unit connector

Proceed in the same manner with the lighting system connector (if present) (fig.6).

(fig.6) Lighting connector

ENGLISH

11.2 Installing the footplates

Insert the footplates (left # 5 and right # 6) into the structure, by rotating from the outside inwards (fig.7) (fig.8).


(fig.7)

(fig.8)

11.3 Installing the batteries and connecting the cables

Slide the grooves of the battery box into the support bar. (fig.9) Connect the black plug of the cable (ref.8) to the control unit and the red plug to the red connector of the first battery (fig.10). Check that the connections are secure.

(fig.9) Support bar

(fig.10)

Repeat the same operation with the second battery, connecting the two black connectors. (fig.11) Attach the attachment bar (#7) to the outside of the battery container and tighten the knob. (fig.12)


(fig.11)


(fig.12) Attachment bar

11.4 Installing the rear push bar (only for CM900-xx)

Raise the rear bar positioned behind the backrest (fig.13). Slide the end of the support tube into the attachment and insert the pin with relative iron. (fig.14)

(fig.13) Rear support bar

(fig.14) Pin with iron


11.5 Adjusting comfort

Before starting to use your **ARDEA MOBILITY** electric wheelchair by Moretti, carry out all the adjustments indicated below to obtain a comfortable driving and seating position suited to the user.


1. Adjusting the armrest height
2. Tipping and extracting the armrests
3. Adjusting the controller position
4. Adjusting the height of the roll-over protection wheels

1. Adjusting the armrest height

The height of the armrests on your electric wheelchair is adjustable. To adjust the height, pull the knob outwards (fig.15), slide the armrest until it reaches the desired height (fig.16), then lock the armrest into the preferred position.


(fig.15) Knob


(fig.16)

2. Tip-up and extractable armrests


Your folding electric wheelchair has tip-up armrests. To tip up the armrests, press the pin and lift the armrest (fig.16-17), then press the rear pin to extract the armrest. (fig.19)

(fig.17) Pin

(fig.18)

(fig.19)

3. Adjusting the controller position

The controller of your folding electric wheelchair can be positioned on both the right and the left. Insert the support bar into the specific guide positioned underneath the armrest on the chosen side. Check that the distance of the drive knob is the best one for use.

- Loosen the knob
- Slide the controller bar into the desired position.
- Tighten the knob

4. Adjusting the roll-over protection height

- Raise the plastic disc
- Slide the tube of the roll-over protection wheels up to the desired height
- Block the tube in the corresponding position (fig.20)

(fig.20)

Note: There are five height adjustment positions, with one every 20 mm (fig.36)

12. OPERATION

The electric wheelchair is simple to use. However, it is recommended to read the following instructions carefully, to familiarise yourself with your new vehicle.

WARNING

Before switching on the power supply, always be aware of your surroundings when selecting the desired speed. It is recommended to select the slowest speed for use indoors. It is recommended to select a speed that allows you to control the wheelchair safely for use outdoors.

The following steps are necessary to ensure that the wheelchair operates safely.

A. CONTROLLER A JOYSTICK

The controller has a multi-purpose visual display. It can provide four types of information:

Controller (fig.21)

- | | |
|------------------------------|-------------------------------|
| 1. Battery charge status LED | 7. Reduce speed |
| 2. ON / OFF status | 8. Increase speed |
| 3. Horn | 9. Left direction indicator |
| 4. Emergency indicators | 10. Right direction indicator |
| 5. Switch on lights | 11. Joystick |
| 6. Speed set LED indicator | |

The joystick controls the direction and speed of the wheelchair. Operate the wheelchair by pressing gently on the joystick in the direction you intend to take. Returning the joystick to the neutral position (central) reduces the speed and stops the wheelchair by automatically applying the electromechanical brakes.

1. ON / OFF status

When the wheelchair is switched on, the controller LED is lit up. When the LED is not lit up, the controller is switched off.

Notes:

- Press the ON / OFF button, then wait two seconds before using the joystick. This is a safety feature that prevents a sudden start-up.
- Gentle use of the joystick results in smoother movements in speed and direction, whereas brusque use of the joystick results in sudden movements in speed and direction.
- In an emergency, release the joystick. The wheelchair will stop immediately.

2. Battery status indicator

The battery status indicator is formed of ten segments (three red, four yellow and three green). It allows you to monitor the battery charge. The battery indicator indicates the approximate quantity of the remaining charge.

- The red, yellow and green LEDs lit up together indicate that the battery is completely charged.
- The red and yellow LEDs lit up together indicate that the battery needs charging as soon as possible.
- The red LEDs lit up on their own indicate that the battery needs charging as soon as possible, as it is almost flat. This could cause a loss of power in the wheelchair.

Note: When the batteries start to discharge, the first red LED will start flashing slowly, to remind you that the batteries need to be charged immediately.

B. USING THE ELECTRIC WHEELCHAIR

Before you get into your folding electric wheelchair, make sure you are on a flat surface with no obstacles. Although your wheelchair is designed for use on slopes, it is safer to practise using it on a flat surface.

Several checks to make before every use are indicated below:

- Check functioning of the horn
- Check that the tyres are in good condition
- Check the battery charge status of your wheelchair
- The motor unlocking levers must be in position D (DRIVE).
- It is recommended to select a low speed for the first few times you use it.

Follow these steps for use of the electric wheelchair:

1. Sit on your wheelchair and press the ON/OFF button. All the LED indicator lights of the battery charger must be lit up. The control and self-diagnosis lights must not be flashing.
2. Rest your forearm firmly on the armrest and place your hand on the controller. After switching on the wheelchair, wait two seconds before activating the joystick. This is a safety feature to prevent sudden start-ups.
3. The joystick central lever controls all movements of the wheelchair. If you move the joystick forwards, the wheelchair moves forwards. If you move the joystick backwards, the wheelchair moves backwards. If you move the joystick to the left, the wheelchair moves to the left. If you move the joystick to the right, the wheelchair moves to the right. When you release the joystick in the vertical or central position, the automatic brake is activated and the wheelchair stops automatically. The joystick also operates as an accelerator. The further you push the lever, the faster your wheelchair will move. The speed depends on the position of the speed control button.
4. Move the joystick central lever gently, applying stable and controlled pressure.
5. In an emergency, release the joystick central lever and the wheelchair will stop automatically.
6. Never drive at a speed beyond your ability to control the wheelchair.
7. Practise driving in a place where there are no obstacles. Start at the lowest speed and move forwards and backwards. Do this several times. When you feel secure, increase the speed gradually by pressing the speed button.
8. The number of LEDs lit up on the speed control indicator identifies the level of speed set.
9. If the battery charge indicator only has two LEDs lit up, the batteries must be charged as soon as possible.

10. If the wheelchair stops and does not function, identify the circuit breakers in the battery containers underneath the seat to reset the circuit. Use the ON/OFF button to switch off the wheelchair, press the reset button and try the electric wheelchair again.
11. When you have finished travelling or before putting the wheelchair away, always switch off the controller by pressing the ON/OFF button.

C. RECOMMENDATIONS FOR SAFE USE

- Before changing direction from forwards to backwards or backwards to forwards, bring the control lever into the central position and wait for the wheelchair to come to a complete stop.
- Only use the wheelchair where it would also be safe to walk.
- If you need to go up or down a slope, check that the gradient does not exceed 6°/10% and that the surface has a good grip and is not slippery.

Use a slow speed to reverse, go up a slope, down a slope, turn a corner or on uneven surfaces.

NEVER use the leg rests to lift the wheelchair. The leg rests are removable and therefore must not be used for lifting.

D. OTHER INFORMATION FOR USE

DRIVING UP A SLOPE

International standards set the maximum safe gradient on which an electric wheelchair can be driven without any difficulty. This value is set as 6°/10%. The Moretti folding electric wheelchair can travel on slopes with a gradient up to 10°/17.5%. If you need to use it on the maximum gradient, you may need to travel at a higher speed. Press the speed button to increase the speed. Be extremely careful when travelling on steep slopes.

DRIVING DOWN A SLOPE


Proceed slowly down steep descents. Press the speed button to reduce the speed. This allows better control of driving. However, your wheelchair will not accelerate on its own down a slope, due to the automatic braking system, which starts functioning if you attempt to drive too fast. In general, you should drive in a linear manner on slopes, without sudden changes in direction, and never stop the wheelchair on a slope that exceeds the safety limit.

GOING OVER OBSTACLES

Go over obstacles slowly and in a forwards direction. Do not attempt to go over obstacles from the side. This will allow you to go over obstacles up to 5 cm. Do not go over obstacles higher than 5 cm. Your wheelchair has a roll-over protection system. These devices may interfere when you are ascending onto or descending from a sidewalk. Take maximum care when going over these obstacles.

E. FREE WHEEL - NEUTRAL

The electromechanical brakes are activated when the wheelchair is not in use or is switched off. The wheelchair has a "free wheel" (neutral) function that is activated by turning the lever from the wheel in drive position to the free wheel position. (fig.22)

(fig.22) Free wheel "neutral" lever

WARNING!

- Never use the free wheel position on a slope.
- Never use the free wheel position when driving the wheelchair.
- Always remember to place the lever in "wheel in drive" position before switching on the wheelchair

F. ELECTROMECHANICAL BRAKES

The electric wheelchair has electromechanical brakes. This is an automatic magnetic disc, also known as a safety brake. The electromechanical brakes are automatic and are engaged when the electric wheelchair is switched on but stationary, even when the wheelchair is on a slope. The electromechanical brakes are also engaged when the electric wheelchair is switched off, but the motor levers are in the drive position (vertical).


Note: Please consult the section on maintenance and repair to ensure that the brakes are in good condition.

G. THERMAL PROTECTION

The controller on your electric wheelchair has a thermal safety system. An incorporated circuit controls the temperature of the controller and the motors and the controller reduces the voltage of the motor and the speed of the wheelchair to allow the electrical components to cool down. Although your wheelchair will return to its normal speed when the temperature returns to safe levels, it is recommended to cut off the power supply and wait five minutes before restarting it, to allow the components to cool down if the speed has dropped suddenly.

H. MAIN CIRCUIT BREAKER:

The reset button is close to the seat. (fig.23)

(fig.23) RESET button

The main circuit breaker controls the electric current absorbed by the battery. This is a safety function of the wheelchair that is included for your safety. When the motor batteries are heavily stressed (such as by excess charging), the main circuit breaker trips to prevent damage to the motors and to the electronics. If the automatic circuit breaker trips, wait around one minute and then press the reset button. Then switch on the control unit and continue with normal functioning. If the main circuit breaker continues to trip, contact your authorised retailer.

12.1 Functioning of the lights (CM900L/CM910L)

The lights on the electric wheelchair are activated and controlled by the light command initially positioned on the left side. Press the ON/OFF button to switch on the light command and activate any type of light. Use the light-blue button to switch on the front and rear lights and press the same button again to switch them off. Press one of the two green arrow buttons to switch on the direction indicators. Press the same button again to switch off the direction indicator. Press the triangle button to switch on the emergency lights (hazard lights). Press the same button again to switch them off.

WARNING!

Use of the lights affects the battery and reduces the duration of the batteries, although to a limited extent.

12.2 Reclining backrest

- Position the plastic shell of the backrest as shown in (photo 24)

(photo 24)

- Insert the aluminium spacer inside the guide on the back of the backrest and tighten the screws as shown in (photo 25)

(photo 25)

- Repeat this step for both sides of the backrest
- When all the screws (6) have been screwed into the relative holes, tighten completely
- Cover the screw heads with the plastic caps provided, as shown in (photo 26)

(photo 26)

- Finish by securing the padding of the backrest in place with the Velcro strap

13. BATTERY & CHARGING

It is recommended to use sealed and maintenance-free deep cycle batteries for your wheelchair. Both lead-acid batteries (SLA) and gel batteries are deep cycle batteries and are similar in terms of performance. Deep cycle batteries are designed to supply current and, once depleted, can be charged again relatively quickly.

Lead batteries must be charged as often as possible.

Specifications of the batteries we recommend:

- Type: **sealed deep cycle or gel**
- Size: **34AH(U1)**
- Voltage: **12V each**
- Amp-hours: **34**

Depending on the conditions of use, the terrain and the manner of driving, the batteries supply around 32 kilometres of travel. However, even if the wheelchair is not in use, you should still charge the batteries periodically.

Notes:

- Do not use car batteries. They are not designed for managing a lengthy and deep discharge and are also not safe for use on electric wheelchairs.
- The lifetime of a battery often depends on its care.

13.1 Charging

The battery charger functions at the normal voltage of a socket (alternating current) and converts it into VDC (voltage direct current). The batteries use direct current to operate the wheelchair. When the batteries are completely charged, the amperage of the battery charger is almost zero. This allows the battery charger to maintain the charge but not overcharge the battery.

Note: The batteries cannot be charged after they have discharged almost to zero voltage

13.2 Charging instructions


To charge the batteries, follow the steps below to charge the batteries:

- Position the wheelchair close to a standard wall socket
- Switch on the controller
- Connect the power cable of the battery charger (round head) to the controller
- Connect the power cable of the battery charger to a standard wall socket.
A yellow LED and a red LED will switch on (if the battery charge is low)
- When the battery is completely charged, the yellow LED turns green
- Disconnect the battery charger cable from the wall socket when the batteries are completely charged. (fig.27)

WARNING!


Do not use a different voltage input to the one indicated. Check the current input voltage (110V or 220V) and adjust manually (fig.28)

(fig.27)

Controller socket

Wall socket

(fig.28)

13.3 Attention

1. Read all the instructions and warnings before using the battery charger.
2. Use the battery charger in a well-ventilated area.
3. To avoid injury, only charge lead-acid or gel batteries.
4. Disconnect the power after charging.

NB: only use the battery charger provided by the supplier of the electric wheelchair. Use of any other type of battery charger can be dangerous and requires the manufacturer's approval.

- Always charge the batteries in a well-ventilated place.
- The battery charger is only for internal use. Protect against the damp.
- For better performance, it is recommended to replace both batteries at the same time if they are flat
- If the wheelchair is not used for a lengthy period of time, it is advisable to charge the batteries at least once a month to avoid their deterioration.
- Can a different battery charger be used? Please note that the battery chargers are selected specifically for special applications and combined with the type and dimensions of the batteries required. In order to charge your wheelchair safely and efficiently, it is recommended to use the battery charger provided. Any type of individual charging of the batteries is prohibited.

Depending on the battery type and state, charging usually takes around 4-10 hours. Charging is complete when the status light on the side panel of the battery charger turns green. Charging the battery for longer than necessary does not damage the battery. It is recommended to charge the batteries for eight to ten hours after everyday use. Do not charge the batteries for more than twenty-four hours.

14. REPLACING THE WHEELS

WARNING! The wheels must be replaced by an authorised retailer or a qualified technician in a workshop

15. MAINTENANCE AND REPAIR OF THE ELECTRIC WHEELCHAIR

The ARDEA MOBILITY range of devices by Moretti are CE marked and carefully checked before being placed on the market. For the safety of the patient, your device should be checked at least once every year by the manufacturer or an authorised laboratory to verify its suitability for use. Use only original replacement parts and accessories for repairs.

MAINTENANCE	Daily	Weekly	Monthly	Semi-annual
ELECTRICAL				
Battery indicator - Check the battery indicator to see whether the battery needs charging	•			
Controller / Display - Check for damage and exposed cables			•	
Check that the sockets and cables are not loose			•	
Check that the batteries are completely charged before daily use	•			
Check that all bolts and tight	•			
Check that the lighting system is working properly (if present)	•			
WHEELS AND TYRES				
Check the tyre pressure	•	•		
The front and rear wheels must rotate smoothly without any obstruction		•		
The front and rear wheels must rotate without wobbling			•	
Inspect the tread of the tyre. If it is less than 1 mm (1/32"), replace it with a new tyre (available from your local retailer)			•	
OTHER				
The motor brushes. It is recommended to have them checked by an authorised retailer every six months, if the electric wheelchair does not function properly. If the inspection shows that the brushes are excessively worn, they must be replaced to avoid damage to the motor				•

ENGLISH

15.1 Daily checks

1. With the controller switched off, check the joystick. Check to ensure that it is not bent or damaged and that it returns to the centre on being released. Visually inspect the rubber sleeve at the base of the joystick to check for any damage. Do not touch it or try to repair it
2. Visually inspect the controller cables. Check to ensure that there are no damaged or exposed wires
3. Check the battery status indicator on the controller to see whether the batteries need to be charged

15.2 Weekly checks

1. Disconnect the controller and inspect the connections to check for any signs of corrosion
2. Ensure that all parts of the controller are securely attached to the electric wheelchair. Do not over-tighten the screws
3. Check the tyre pressure The wheelchair has solid tyres. If your electric wheelchair is equipped with optional tyred wheels, ensure that the tyre pressure remains between 30-35 psi
4. Check the brakes. This test must be carried out on a flat surface with at least one metre of space around your electric wheelchair

To check the brakes (the wheelchair should move slightly during performance of this test):

- Switch on the controller and lower the speed and the regulation knob
- After one second, check that the battery level indicator remains on
- Push the joystick forwards slowly until you hear the parking brakes go click

Immediately release the joystick. You should hear each parking brake operate within a few seconds of releasing the joystick.

- Repeat the check on the brake for the backwards, left and right positions of the joystick

15.3 Six-monthly checks

1. Check the motor brushes. We recommend having the wheelchair checked every six months by an authorised retailer, or sooner if it is not functioning properly. If the inspection identifies excessive wear of the brushes, they must be replaced, to avoid damage to the motor. Check the state of the battery terminals every six months. Check that they are not corroded and that the connections are secure. Periodically apply a thin layer of Vaseline on the surface of the terminals to protect them against corrosion.


Check the state of the battery terminals every six months. Check that they are not corroded and that the connections are secure. Periodically apply a thin layer of Vaseline on the surface of the terminals to protect them against corrosion.

- Failure to perform maintenance on the brushes could invalidate the warranty on your wheelchair.

To check and replace the motor brushes:

- Unscrew the motor brush cap;
- Remove the brushes;
- Check wear of the brushes;
- Replace the brushes if necessary

New brush

Worn brush

Brush cap

ENGLISH

CHECKS & REMINDERS:

- Keep the controller clean and protect it from rain or water. Never expose the electric wheelchair to direct contact with water.
- Keep the wheel clean of residue, hair, sand and carpet fibres.
- All the upholstery can be washed with warm water and neutral soap. Check the seat and backrest for rips and/or tears every now and then. Replace them if necessary. Do not store your electric wheelchair in damp environments, as this can cause mould to form and quickly damage the upholstery.
- Inspect and lubricate all moving parts. Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains. Always check that all nuts and bolts are tight.

PERIODIC CHECKS

1. Keep the controller clean and protect it from rain or water. Never expose the electric wheelchair to direct contact with water.
2. Keep the wheels free of fluff, hair, sand and textile fibres.
3. Visually inspect the tread of the tyres. If it is less than 1/32", replace the tyres
4. All the upholstery can be washed with warm water and neutral soap. Check the seat and backrest for rips and/or tears every now and then, and replace them if necessary. Do not store your electric wheelchair in damp environments, as this can cause mould to form and quickly damage the upholstery.
5. Inspect and lubricate all moving parts. Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains. Always check that all nuts and bolts are tight.

16. CARE AND MAINTENANCE

Your electric wheelchair requires a minimum amount of care and maintenance. The following parts require checking and/or care and maintenance

16.1 Tyre pressure

- If the wheelchair has tyres, maintain the air pressure as indicated on the tyres.
- It is important to maintain the specified tyre pressure at all times. Do not over-inflate the tyres. Low pressure can result in loss of control and excessively high pressure can cause the tyres to burst. Failure to maintain the specified tyre pressure at all times can damage the tyre and/or wheel. Inspect the tyres regularly for wear.

16.2 Cleaning and disinfection

- Use a damp cloth and non-abrasive neutral detergent to clean the plastic and metal parts of the electric wheelchair. Avoid using products that could scratch the surface of your electric wheelchair.
- DO NOT immerse the device in water for any reason!
- If necessary, clean the product with an approved disinfectant. Make sure the disinfectant is safe for the product before you use it.
- Read all the safety instructions for correct use of the disinfectant and/or detergent before using it on your product. Failure to do so may result in irritation of the skin or premature damage to the upholstery and/or finishes of the electric wheelchair.
- Never use acids, alkalis or solvents (acetone or diluent) to clean it.

16.3 Battery terminal connections

- Make sure the connectors are tightly connected and not corroded.
- The batteries must be in their compartments
- The battery terminals must be turned inwards towards the electric wheelchair

16.4 Cables

- Regularly check the electrical connections.
- Regularly check the electrical insulation, including the battery charger cable, for signs of wear or damage.
- Have any damaged connector, connection or insulation replaced or repaired by an authorised retailer before using the electric wheelchair.

16.5 Axle supports and motor/differential

These components are pre-lubricated and sealed and do not require further lubrication.

REPLACING A WHEEL

If the electric wheelchair has tyred wheels and one of the tyres is flat, replace the inner tube.

If the electric wheelchair has a solid rubber tyre, replace the entire wheel.

Contact an authorised retailer for details on spare wheels for your electric wheelchair.

MOTOR BRUSHES

The motor brushes are housed inside the differential/motor unit. They must be inspected periodically. Have them replaced by an authorised retailer if they are worn.

16.6 Control panels, battery charger and rear electronics

- Protect these components from damp.
- If exposed to damp, leave them to dry completely before you use the electric wheelchair again.

16.7 Storing the electric wheelchair

If you do not plan to use the electric wheelchair for a lengthy period of time, it is advisable to:

- Completely charge the batteries before you store it.
- Disconnect the batteries from the electric wheelchair.
- Keep the electric wheelchair in a warm and dry place.
- Do not store the electric wheelchair in places where it could be exposed to extreme temperatures.
- Operating conditions (-25° ~ +50°).
Storage conditions (-40° ~ +65°).

Batteries that are completely discharged, charged rarely, stored at extreme temperatures or stored without being completely charged may suffer permanent damage, resulting in their unreliability and limited durability. It is advisable to charge the batteries of your electric wheelchair periodically in the case of lengthy storage, to guarantee proper performance.

17. TROUBLESHOOTING

P&G S-Drive 45 Amp Controller: The electric wheelchair has a controller that continually controls the operating conditions of your electric wheelchair. If it detects a problem, the ON / OFF button flashes. You must count the number of flashes and consult the list to check which type of error has been detected

Number of flashes	Description
1 ■ □ □ □ □ □ □ □ □ □	The battery needs charging or there is a battery connection fault Check the battery connections If the connections are correct, try charging the battery
2 ■ ■ □ □ □ □ □ □ □ □	The left motor is not properly connected. Check the connections
3 ■ ■ ■ □ □ □ □ □ □ □	The left motor has a short circuit to the battery connection. Contact the assistance centre
4 ■ ■ ■ ■ □ □ □ □ □ □	The right motor is not properly connected. Check the connections
5 ■ ■ ■ ■ ■ □ □ □ □ □	The right motor has a short circuit to the battery connection. Contact the assistance centre.
6 ■ ■ ■ ■ ■ ■ □ □ □ □	The wheelchair has been prevented from driving by an external signal. Contact the assistance centre (for numerous reasons and causes)
7 ■ ■ ■ ■ ■ ■ ■ □ □ □	A joystick error is indicated. Ensure that the joystick is in the central position before switching on the control system
8 ■ ■ ■ ■ ■ ■ ■ ■ □ □	A possible fault of the control system is indicated. Ensure that all the connections are secure
9 ■ ■ ■ ■ ■ ■ ■ ■ ■ □	The parking brakes are not properly connected. Check the brake and motor connections
10 ■ ■ ■ ■ ■ ■ ■ ■ ■ ■	An excessively high voltage has been applied to the control system This is usually caused by a poor battery connection. Check the battery connections

Note: If technical problems occur, it is recommended to check the device with a local retailer before trying to solve the problems yourself.

The following symptoms may indicate that your electric wheelchair has a serious problem. Contact your local retailer if any one of the following cases occurs:

1. Motor noise
2. Frayed harnesses
3. Bent or broken connectors
4. Uneven wear on any one of the tyres
5. Sudden movements
6. The electric wheelchair pulls more to one side
7. Bent or broken wheel units
8. The electric wheelchair does not switch on
9. The electric wheelchair switches on, but does not move

18. DISPOSAL CONDITIONS

18.1 General conditions of disposal

Never dispose of the product as normal domestic waste. Dispose of the product at a sorted waste collection centre for recycling.

18.2 Instructions for correct disposal in accordance with European directive 2012/19/EU:

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment.

18.3 Disposal of exhausted batteries - (Directive 2006/66/EC):

This symbol on the product indicates that the batteries must not be considered as normal domestic waste. Ensuring that the batteries are discarded correctly contributes to preventing potentially negative consequences for health and the environment that would otherwise be caused by their inappropriate disposal. Recycling materials helps to conserve natural resources. Take depleted batteries to the collection points indicated for recycling. You can contact your local council, the local waste disposal service or the shop where you purchased the device for more detailed information on disposal of depleted batteries or the product.

19. REPLACEMENT PARTS AND ACCESSORIES

Use only replacement parts and accessories listed in the Moretti general catalogue

20. TECHNICAL SPECIFICATIONS

CODE- MODEL	CM900 - HP5 CM900L - HP5L
MAXIMUM CAPACITY	135kg
SEAT	41 cm (CM900-41) 46 cm (CM900-46)
REAR WHEELS	320 mm x 60 mm
FRONT WHEELS	200 mm x 50 mm
ROLL-OVER PROTECTION WHEELS	35 mm x 17 mm x 2pcs
MAXIMUM SPEED	6Km
BATTERY SPECIFICATIONS	12V 34Ah x 2pz
BATTERY DURATION	32km
BATTERY CHARGER TYPE	5Amp, external, 120/240 Volt, 50/60 Hz
CONTROLLER	VSI 50Amp
MOTOR	3900rpm 200W x 2pcs
WEIGHT WITH BATTERY	62kg
WEIGHT WITHOUT BATTERY	38kg
STEERING RADIUS	865 mm
SUSPENSION	NO
LENGTH	1100 mm
DEPTH	685 mm
HEIGHT	900 mm
SEAT WIDTH	460 mm
SEAT HEIGHT	508 mm
SEAT DEPTH	406 mm
BACKREST HEIGHT	460 mm
WHEELBASE	520 mm
LEG REST	400 mm~500 mm
HEIGHT FROM THE GROUND	80 mm
MAXIMUM SLOPE SAFETY	10°
MAXIMUM SLOPE	12°
SURMOUNTING OF AN OBSTACLE	50 mm

CODE- MODEL	CM910 - HP8 CM910L - HP8L
MAXIMUM CAPACITY	135kg
SEAT	46 cm
REAR WHEELS	320 mm x 60 mm
FRONT WHEELS	200 mm x 50 mm
ROLL-OVER PROTECTION WHEELS	35 mm x17 mm x 2pcs
MAXIMUM SPEED	6Km
BATTERY SPECIFICATIONS	12V 34Ah x 2pz
BATTERY DURATION	32km
BATTERY CHARGER TYPE	5Amp, external, 120/240 Volt, 50/60 Hz
CONTROLLER	VSI 50Amp
MOTOR	3900rpm 200W x 2pcs
WEIGHT WITH BATTERY	66kg
WEIGHT WITHOUT BATTERY	42kg
STEERING RADIUS	865 mm
SUSPENSION	NO
LENGTH	1085 mm
DEPTH	640 mm
HEIGHT	950 mm
SEAT WIDTH	460 mm
SEAT HEIGHT	508 mm
SEAT DEPTH	406 mm
BACKREST HEIGHT	406 mm
WHEELBASE	520 mm
LEG REST	400 mm~500 mm
HEIGHT FROM THE GROUND	80 mm
MAXIMUM SLOPE SAFETY	10°
MAXIMUM SLOPE	12°
SURMOUNTING OF AN OBSTACLE	50 mm

21. WARRANTY

Moretti products are guaranteed for 2 (two) years from the date of sale against material and manufacturing defects, subject the following limitations. The warranty is voided by improper use, abuse, modifications to the product and failure to follow the instructions. The intended use of the product is given in the user manual. Moretti is not liable for damage, injury or any other consequences resulting from installation or use which are not scrupulously conforming with the instructions given in the installation, assembly and user manual. Moretti does not guarantee its products against damage or defects in the following circumstances: natural disasters, unauthorised repair or maintenance, improper electric power supply (as applicable), use of parts or components not supplied by Moretti, failure to follow the guidelines and instructions for use, tampering, shipping damage (other than the original shipping by Moretti), or failure to perform maintenance as indicated in the manual. Components subject to wear and tear are not covered by this warranty if the damage is caused by normal use of the product.

21.1 Warranty on rechargeable batteries (if provided)

The original and replacement batteries are covered by a ninety-day warranty on performance and a six-month warranty on manufacturing defects or as required by law. If completely charged batteries are left unused for more than three consecutive months, the warranty is automatically invalidated. If completely discharged batteries are left unused for more than three consecutive days, the warranty is automatically invalidated.

21.2 Repairs

Repairs under warranty If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti reseller or its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

Repairing a product not covered by the warranty Product not covered by warranty may be returned for repair only if authorised in advance by Moretti customer service. The costs of labour and shipping incurred by repairs not covered by the warranty are borne by the client or reseller in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for 6 (six) months from the day of reception of the repaired product.

Non-defective products The client will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the client at his expense.

22. REPLACEMENT PARTS

Moretti original replacement parts are guaranteed for 6 (six) months from the day of delivery.

23. NON-LIABILITY CLAUSE

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti makes no declaration, guarantee or condition, express or implicit, including any future declaration, guarantee or condition of sale, suitability for a given purpose, non violation and non interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit guarantee under the law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit guarantee or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice.

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