M**ARDEA**bility

ITALIANO

MA TB030X-040X-075X-100X 03 A_ITA_10 2020

PROPULSORE ANTERIORE

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INDICE

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INSTRUCTION MANUAL

INDEX

1. CODES	PAG.5
2. INTRODUCTION.....	PAG.5
3. INTENDED USE.....	PAG.5
4. DECLARATION OF CONFORMITY	PAG.5
4.1 Applicable regulations and directives.....	pag. 6
5. GENERAL WARNINGS	PAG.6
5.1 General warnings.....	pag. 6
5.2 General information on electromagnetic disturbance	pag. 8
6. SYMBOLS	PAG.9
7. DESCRIPTION OF THE DEVICE	PAG.10
7.1 Description of parts.....	pag. 11
8. BEFORE EVERY USE	PAG.11
9. INSTRUCTIONS FOR USE	PAG.11
9.1 Fitting the quick connector.....	pag. 12
9.2 Fitting the battery	pag. 12
9.3 Coupling.....	pag. 12
9.4 Uncoupling	pag. 13
9.5 Using and recharging the battery	pag. 14
9.6 Battery recharging procedure.....	pag. 15
9.6.1 Guidelines for the use of lithium batteries.....	pag. 16
9.6.2 Initialising a lithium battery for the first time.....	pag. 16
9.6.3 Charging and discharging procedures	pag. 16
9.6.4 How to keep a lithium battery in case of prolonged disuse.....	pag. 16
9.6.5 Recommendations	pag. 16
9.6.6 How to extend the lithium battery's life.....	pag. 16
9.6.7 Transporting the battery by plane.....	pag. 17
9.7 Using the throttle.....	pag. 17
9.8 Using the brake	pag. 17
9.9 Using cruise control	pag. 17
9.10 Reverse button	pag. 18
9.11 Adjusting the handlebar	pag. 18
10. MAINTENANCE	PAG.18
10.1 Ordinary maintenance.....	pag. 18
10.2 Annual maintenance.....	pag. 20
10.3 Wheel maintenance.....	pag. 20
11. CONTROLS AND REMINDERS	PAG.21
11.1 Steering tube.....	pag. 22
11.2 Tyre pressure	pag. 22
11.3 Battery terminal connections	pag. 22
11.4 Cabling, controller, charger, electronics.....	pag. 22
11.5 Storing the device.....	pag. 22
11.6 Transport	pag. 22

12. SPARE PARTS AND ACCESSORIES..... PAG.23

- 12.1 TBA102 Digital display - TIBODA pag. 24
 - 12.1.1 Installation pag. 24
 - 12.1.2 Display interface pag. 25
 - 12.1.3 Reading the information pag. 25
 - 12.1.4 Battery charge indicator pag. 26
 - 12.1.5 Power output indicator pag. 26
 - 12.1.6 Turning the LCD display's backlighting on and off pag. 27
 - 12.1.7 Units pag. 27
 - 12.1.8 General settings pag. 28
 - 12.1.9 Saving the settings pag. 28
- 12.2 TBA120 Led light with horn 36V
 - TIBODA - TBA121 Led light with horn 48V - TIBODA pag. 29
 - 12.2.1 Installation pag. 29
- 12.3 TBA130 Handlebar throttle 36V
 - TIBODA; TBA131 Handlebar throttle 48V pag. 30
 - 12.3.1 Installation pag. 30

13. CLEANING AND DISINFECTION PAG.30

14. TROUBLESHOOTING..... PAG.31

- 14.1 Error codes pag. 31

15. CONDITIONS OF DISPOSAL PAG.32

- 15.1 General terms of disposal pag. 32
- 15.2 Instructions for correct disposal in accordance with European directive 2012/19/UE pag. 32
- 15.3 Treatment of exhausted batteries (Directive 2006/66/EC) pag. 32

16. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY PAG.32

17. TECHNICAL SPECIFICATIONS PAG.35

- 17.1 Dimensions pag. 35
- 17.2 Technical specifications pag. 36

18. WARRANTY PAG.38

- 18.1 Rechargeable battery warranty pag. 38

19. REPAIRS..... PAG.39

- 19.1 Repairs under warranty pag. 39
- 19.2 Repairing a product not covered by the warranty pag. 39
- 19.3 Non-defective products pag. 39

20. SPARE PARTS PAG.39

21. NON-LIABILITY CLAUSE PAG.39

ENGLISH

TIBODA

CE | Class Medical Device

REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices

1. CODES

TB030X	TIBODA front drive for wheelchairs, power 300W
TB040X	TIBODA front drive for wheelchairs, power 400W
TB075X	TIBODA front drive for wheelchairs, power 750W
TB100X	TIBODA front drive for wheelchairs, power 1000W

2. INTRODUCTION

Thank you for choosing a **TIBODA by ARDEA MOBILITY di MORETTI S.p.A.** front drive. It is designed and specified for unbeatable comfort, safety and reliability. This user manual contains useful information for using your product properly and safely. Please read this manual in full before you use the product. Should you have any queries, please contact your dealer for suitable advice and assistance. The manual is part of the product and must be kept with care and included with it if it is transferred to another user. The manual is intended for operators, owners, users and maintenance staff. The manual provides information on: use, technical specifications, transport, keeping, maintenance, disposal and safety. Any modifications to the manufacturer's instructions which affect the safety of the patient or operator will be communicated to the owner of the product by Moretti S.p.A. as it deems suitable. Any non-substantial modifications and additions to the product will not be notified by the manufacturer.

3. INTENDED USE

The **TIBODA by ARDEA MOBILITY di MORETTI S.p.A.** front drive is an electric drive for wheelchairs, fitted with a quick mounting, which will revolutionise your urban mobility by eliminating the need for you to push the wheelchair yourself.

WARNING!

- Do not use the product for any purpose other than that specified in this manual.
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than that specified in this manual.
- The manufacturer reserves the right to modify the device and its relevant manual without prior notice, in order to improve the product's features.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of FRONT WHEELCHAIR DRIVE complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

4.1 Applicable regulations and directives

The TIBODA front wheelchair drive is conforming and has been tested and approved pursuant to the following directives and standards:

- EN 14971:2012 - Medical devices - Application of risk management to medical devices
- EN ISO 15223-1:2012 - Medical devices - Symbols for use in medical device labels, and in obligatory labelling and information - Part 1: General requisites
- EN 10993-1:2010 - Biological assessment of medical devices - Part 1: Assessment and testing as part of the risk management process
- EN 12182:2012 - Products for the assistance of disabled persons - General requisites and test methods
- ¹EN 12183:2014 - Hand powered wheelchairs - Requisites and test methods
- ²EN 12184:2014 - Electric wheelchairs, scooters and charging systems - Requisites and test methods
- EN 60601-1:2007 - Electric medical equipment, part 1: General safety and performance prescriptions
- EN 60601-1-2:2003 - Electric medical equipment, part 1: General safety prescriptions - Collateral standard: Electromagnetic compatibility - Prescriptions and tests
- EN 60601-1-11 - Electric medical equipment, part 1: General safety and performance prescriptions for electric medical equipment
- The "TIBODA" is also conforming with Directive EU2012/19/EC relating to the scrapping of waste electrical and electronic equipment

¹EN 12183:2014 and ²EN 12184:2014: the following standards have been used as guidelines; where possible, the tests specified in these standards have been applied

5. GENERAL WARNINGS

5.1 General warnings

Take great care when using the device the first few times. Every manoeuvre should be considered potentially hazardous. It is essential that you familiarise yourself with every manoeuvre and gain confidence in using the "TIBODA". Take special care when steering and accelerating. Choose a clear, safe and level area to try out the "TIBODA" the first time, and have someone with you who has read this manual and is able to help you if necessary.

Using the "TIBODA" in a way that poses no danger to yourself or other people takes training and experience, and familiarity with the rules governing foot and road traffic in your country.

WARNING!

- Do not use the product for any purpose other than that specified in this manual
- Moretti s.p.a. declines all liability for damage or sanctions resulting from the use of the front drive in violation of traffic regulations in the country/region in which Tiboda is used
- The device may only be used by persons in possession of their full intellectual faculties
- Keep the device away from children under the age of 12
- Make sure the package is undamaged before you open it
- Do not exceed a speed of 6 kph on irregular ground
- Do not exceed a speed of 10 kph on cycle paths
- Do not exceed a speed of 12 kph on private roads
- Do not exceed a speed of 3 kph on bends. Failure to observe these limits is dangerous and may cause the wheelchair to tip over
- Do not use TIBODA on lengthy uphill. Doing so may overheat the battery
- Do not use the battery charging cable as a tourniquet or rope, use it only as intended by the manufacturer

- Always slow down when approaching steps, corners and narrow passages. Do not drive over steps of more than 5 cm, always take them orthogonally to the direction of travel.
- Moderate your speed on wet and icy ground, gravel or irregular terrain
- Keep your hands on the handlebar when using the product
- TIBODA may only be used for transporting a single person at a time
- The weight of the passenger should not exceed 120 kg for reasons of safety
- In case of battery leaks, avoid contact with the fluid.
Contact your reseller or service centre
- If you replace any components yourself, Moretti is not liable for any resulting damage or injury.
- The device has the following parts: anticorrosive 6082 aluminium + S355 steel + polypropylene + stainless steel fasteners and lithium battery. Do not use the device if you are allergic to any of its materials
- Dispose of replaced or broken parts as required by local regulations
- If TIBODA is subjected to direct sunlight or low temperatures over a long period of time, the parts may overheat. The operating and storage temperatures are: 0°C - +40° C
- Service and clean the product as indicated under "MAINTENANCE"; accumulated fluff, dust and parasites can compromise the operation of the device
- Do not use in the rain, the product is rated IP21
- Do not use together with other equipment which generates EM interference
- Do not operate the power button and throttle lever at the same time.
This will make the product move suddenly and in an uncontrolled manner
- Check the condition of the cabling and connectors and that the latter are properly connected
- Check that the rear wheels of your wheelchair are securely installed
- Use ramps and avoid steps wherever possible. Take care when crossing roads
- We recommend wearing a helmet
- Always turn the lights on when using the product at night
- Do not attempt to drive up or down slopes in excess of the maximum slope indicated in "technical specifications"
- When driving up steep slopes, move your body forwards to increase the adherence of the drive wheel, have an assistant aid you if possible
- Adjust your speed and driving on slippery ground (for instance, wet ground or loose soil and gravel); never use TIBODA on mud or ice
- Always take great care when crossing roads
- When travelling a long distance, keep a spare battery with you
- Keep your hands away from moving parts
- Never use TIBODA with pets or children on your lap
- Use only the included battery charger for charging the batteries
- Use only the accessories provided by the manufacturer
- Do not tamper with the device, doing so compromises its performance and safety
- Always remove the battery when transporting TIBODA
- TIBODA devices are lower than a normal bicycle and are thus less visible to a motorist
- After use, store the device on a flat surface and away from potential collisions. Failure to observe this instruction can cause the product to tip over
- If the device is transported in a car, place it in the luggage compartment in a stable position and drive at moderate speed to prevent collisions and excess vibration
- Do not place objects weighing more than 5 kg on the product
- Do not place objects on the battery mount when the battery has been removed. The mount is designed only for the battery
- Moretti S.p.A. declines any and all liability for damages resulting from improper use of the device or any use other than that specified in this manual.
- The manufacturer reserves the right to modify the device and its relevant manual without prior notice, in order to improve the product's features.

- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.
- It is forbidden to modify the parts of the TIBODA including the setting parameters of the TIBODA engine control unit. MORETTI S.p.A. declines any responsibility for damage caused by a non-original setting of the control unit.
- TIBODA is not suitable to be mounted on pushchairs with transit wheels.

5.2 General information on electromagnetic disturbance

The rapid growth of electronics, especially in communications, has saturated the environment with electromagnetic radiation (EM) emitted by TV, radios and communications equipment. These EM waves are invisible and are stronger the closer you are to their source. All electrical conductors act as antennas for EM signals and, to varying extents, wheelchair drives may also be affected by electromagnetic interference (EMI). EMI may cause unintentional movements and irregular control of the drive.

The electromagnetic energy emitted by sources such as:

- radio stations
- TV stations
- amateur radio
- shop alarm systems
- cellphones and cordless phone can interfere with electric wheelchairs and electric scooters

The interference can cause the regenerative brake to release, cause unintentional movement of the device, or movement in the wrong direction. It can also permanently damage the device's controller. In our everyday environment, there are a certain number of intense sources of EMI. The intensity of the EM energy can be measured in Volts per metre (V/m). Your drive is equipped with protection against EMI which has been tested and certified to be conforming with international regulations, up to a certain intensity of the signal. This is called the "level of immunity". The greater the level of immunity, the better the device is protected. At this time, technology can provide a level of immunity of at least 20 V/m; this is sufficient protection against the most common sources of EMI. However, we consider that the following precautions can further minimise the risk.

- A) Do not use portable radio transmitters, CB radios, cellphones or other personal communications equipment when the drive is powered on.
- B) Do not approach powerful radio and TV transmitters too closely.
- C) If you feel that your scooter has moved or the brake has released in a way that you do not intend, press the OFF button as soon as you can do so safely.
- D) Installing additional electrical equipment and modifying your scooter's motor may make it more susceptible to EMI and hence make using the equipment hazardous. Note that doing so also voids the warranty.
- E) In any case, report any unintentional movement of the vehicle to your reseller or Moretti S.p.A., note it down in this manual and note down the fact that there is a source of electromagnetic radiation in your vicinity, if any.

6. SYMBOLS

All the warnings in this user manual are highlighted by symbols. The text employs symbols and words which indicate the gravity of the hazard and the need to take particular care.

	Product code
	Unique Device Identification
	CE mark
	Manufacturer
	Batch Lot
	Read the instruction manual
	Medical Device
	Serial number
	Condition of disposal
	Attention
	Danger of electric current
	Production dates
	Waste Electrical & Electronic Equipment (WEEE)
	Notes
	B Type
	It is mandatory to consult the instruction manual in case of maintenance

7. DESCRIPTION OF THE DEVICE

TIBODA is a universal wheelchair drive with very fast coupling and rechargeable battery.

It is very easy to couple TIBODA to a hand-powered wheelchair so that the front caster wheels are raised off the ground and the device's steering wheel is used to manoeuvre the wheelchair.

Standard features of all models:

- Brushless motor with internal sine wave controller for smooth, silent running and increased comfort
- 160mm mechanical disk brake for effective braking
- Regenerative brake for recharging the battery and increased braking power (intervenes automatically when the speed exceeds 8 kph)
- Reverse gear for manoeuvring even in the most restricted spaces
- Cruise control maintains constant speed over long distances
- Parking stand
- Return-to-centre steering
- USB port for recharging battery devices such as smartphones, tablets, cameras, etc.

1
Cruise Control
(Optional horn)

2
Reverse gear


3
Accelerator and battery
charge indicator

4
Ultra quick coupling

5
Parking stands with
dia. 40mm wheels

6 7
16" wheel
Regenerative brake
(EBS)

8
Mechanical disk brake
160mm

9
Return-to-centre
steering

TIBODA also has a vast range of accessories, so that you can adapt it to your particular requirements. Please go to www.tiboda.it for further information.

7.1 Description of parts

Fig.1

ENGLISH

8. BEFORE EVERY USE

WARNING!

- Always make sure the device's components are present and properly located and assembled. If they are not, DO NOT USE THE DEVICE; contact your service centre
- Make sure the tyre is properly inflated and in good condition. Refer to "MAINTENANCE" for the inflation pressure
- Make sure the battery has sufficient charge for the route you want to cover
- Make sure that the front drive is coupled properly to the wheelchair; refer to "COUPLING"
- Make sure the brake is working properly; test the brake before you set off
- When you first use the Tiboda, set off slowly so as to gain familiarity with it. Every new manoeuvre is potentially hazardous. Pay special attention to your speed around corners
- Service and clean the product as indicated under "CLEANING"; accumulated dust and dirt can compromise the operation of the device
- Check that the rear wheels of your wheelchair are securely installed
- We recommend wearing a helmet and other protective equipment
- Always check that the throttle lever and reverse button are working properly


9. INSTRUCTIONS FOR USE

WARNING!


- To proceed with the quick coupling/user instructions, the device must be set up by the reseller or an authorised Moretti S.p.A. workshop using the "TIBODA installation manual"
- Only couple or uncouple the device on level, clear ground

9.1 Fitting the quick connector

- A Couple the quick coupling unit (2 - Fig.2) to the drive unit (1 - Fig.1)
- B Turn the handle on the quick coupling clockwise until it is secure (there should be no play between the parts)

9.2 Fitting the battery

- A Make sure the battery power button (A - Fig.3) is set to OFF
- B Fit the battery into its mount on the drive unit (1 - Fig.1)

9.3 Coupling

Fig.4

- 1 Before you couple the device to the wheelchair, make sure that the battery's power button is set to OFF (Ref. A - Fig. 3)
- 2 Align your wheelchair with the drive (Fig.4)
- 3 Grasp the handlebar and pull the Tiboda towards yourself, until its mounting is resting on the top pin of the bracket mounted to the wheelchair

Fig.5



Fig.6


- 4 Press the handlebar downwards and pull the TIBODA towards your (Fig.5), until the first pin (Ref. 1 - Fig. 6) fits into the quick mounting
- 5 Push the handlebar forwards (to rotate it as shown in Fig.6) while moving your body back, until the quick mounting engages with the second pin (Ref.2 - Fig.6). The front wheels of the wheelchair are now raised off the ground
- 6 Set the battery's power button to ON
- 7 The Tiboda is now ready to go

9.4 Uncoupling

Fig.7

Fig.8


- 1 Before you uncouple the device from the wheelchair, make sure that the battery's power button is set to OFF (Ref. A - Fig. 3)
- 2 Release and lower the quick mounting levers by pulling them towards you (Fig.7)
- 3 Grasp the handlebar and push it forwards while moving your body back so that the second pin disengages (Rif.2 - Fig.6); when the stand is resting on the ground, push the handlebar downwards, to allow the first pin to disengage, then push the TIBODA forwards and away from you (Fig.8)
- 4 The Tiboda is now uncoupled from the wheelchair

9.5 Using and recharging the battery

The TIBODA front drive is equipped with compact 36/48V batteries (depending on the model of drive) with Samsung lithium cells. The battery has a lock and a USB port for recharging portable devices such as smartphones, tablets, etc.

It also has an internal electronic controller (BMS = Battery Management System) which shuts off power delivery when the voltage drops below or rises above a certain range. The controller also protects the battery against current overloads and short circuits, thanks to its internal fuse.

9.6 Battery recharging procedure

- 1 Remove the battery from its mount and take it to a cool, dry place
- 2 Make sure its power button is set to OFF
- 3 Insert the charger connector into the socket under the cover
- 4 Hook the charger up to the mains (220-240V); the charger's led turns RED while it is charging the battery;
- 5 A full battery charge cycle takes 4-5 hours (depending on the model of battery and charger)
- 6 Once the battery is fully charged, the led turns GREEN
- 7 You can now disconnect the charger from the battery and close the cover on the battery socket

This battery is rated for over 400 charging cycles and has a range of around 50 km. The performance of the battery is affected by a variety of factors: performance requirements, weight of passenger, running of the wheelchair, route, weather. The data indicate its average performance on mixed routes with 50% power assistance

WARNING!

- Use only official Tiboda batteries and chargers
- You may see a spark when connecting the battery to the charger. Make sure no flammable matter or fluids are in the vicinity
- Charge the battery under adult supervision and in a well ventilated area
- Charge battery far away from flammable matter
- Do not leave or charge the battery in the vicinity of a heat source
- Do not charge the battery at temperatures below 0°C
- Keep the battery charged at all times
- Do not leave the charger powered at 230 V for too long after it has finished charging the battery, as indicated by the green light. (BMS FUNCTION)
- Do not short circuit the battery (by touching the two terminals to a metal object)
- Do not immerse the battery in water
- Do not attempt to open the battery pack
- Do not allow children to play with the battery pack
- Do not dispose of the battery as domestic waste; have it disposed of by a sorted waste centre
- Do not use the battery to power other equipment
- Do not leave the battery in the sun for a long time, and never leave it under the rain
- Avoid contact with the battery fluid if it leaks out of the battery
- Never run two full charge cycles in a row
- Always combine partial and full charge cycles as described above
- Never discharge the battery completely, i.e. below 20% (when the battery indicates that it is discharged). Fully discharging the battery can damage it beyond repair

9.6.1 Guidelines for the use of lithium batteries

Lithium ion batteries are generally composed of one or more cells depending on their voltage rating (3.7 V per cell) and their charge capacity (in mAh - milliamp hours) generally determines their size. The memory effect problem characteristic of old style nickel (NiCd or NiMH) batteries no longer affects lithium batteries, or at least only to a very much smaller extent. A lithium battery can be recharged even if it is only partly discharged, without affecting its capacity.

However, it is important to allow it to discharge completely every 20/30 partial recharge cycles (so that the indicator is red, around 25%) and then fully recharge it. NEVER completely discharge the battery to 0%.

9.6.2 Initialising a lithium battery for the first time

When you acquire a tablet or cellphone with a rechargeable battery, you never know how long the battery has been in it; it is therefore essential to recharge it completely before using the device. You do not need to recharge it for a long time; 5/6 hours are sufficient (even if the devices indicates that the battery is charged after only 1 or 2 hours). This indication is invalid when the battery has not yet been initialised. The next step is to run 2 to 4 discharge and recharge cycles (as noted above, only discharge the battery to around 20% charge, not completely). After this, the battery can be used normally, even with partial recharge cycles.

9.6.3 Charging and discharging procedures

As already mentioned, partially discharging and recharging a lithium battery does not create the memory effect; indeed it is very convenient to recharge the battery frequently even when it is not fully discharged. After the first initialisation (see above), always disconnect the charger when the charge is reading 100%. After 20/30 partial charging cycles, run a complete charging cycle (recalibration) by allowing the battery to discharge until the battery discharged indication displays, then fully recharge it. If you do not recalibrate the battery in this way, the charge sensor will become increasing inaccurate and the battery will lose its efficiency.

9.6.4 How to keep a lithium battery in case of prolonged disuse

If you do not intend to use the product for a long time, certain precautions are in order. First of all, set the power switch to OFF, remove the battery from the product and store it in a cool, dry place. The battery must be kept at 40/50% charge (never completely discharged) and fully recharged before it is used again. Finally, note that no matter how carefully you treat it, a battery will only keep its full capacity for 2 or 3 years or 300/600 charging cycles!

9.6.5 Recommendations

- Keep the battery in a cool place at all times. It is sensitive to heat; for instance, do not keep it in the car under the summer sun!
- Remove the battery from the device when it is charging
- Clean the contacts from time to time, both on the battery and the product, with some cotton wool and alcohol
- If you have two batteries, always discharge one fully and keep the other in a cool place

9.6.6 How to extend the lithium battery's life

We give below the information and recommendations from Battery University for extending your battery's life and duration:

1. Do not discharge the battery completely and frequently; this stresses the battery
2. Discharge the battery partially and charge it frequently. Recharging a partially charged lithium battery does not cause the memory effect, in contrast with nickel or nickel cadmium batteries
3. If a battery lasts for a short time only, this more likely due to heat than discharge/charge cycles.
4. Lithium batteries must be recalibrated with a full charge cycle every 30 charge/discharge cycles
5. If you do not recalibrate the battery, its charge sensor will become increasingly inaccurate and the battery may no longer work to full effect
6. Keep your lithium battery in a cool place, do not leave it in a car exposed to the sun
7. If you do not use the battery for a long time, keep in a cold place (0°C is ideal, at 25°C the battery loses 16% more of its charge than at 0°C) and with 40/50% of its charge (so that its internal

circuits have enough power to keep working)

8. Do not buy spare lithium batteries for future use; always check the date of manufacture, and do not buy old batteries even if they are discounted
9. If you have a backup battery, use one to full capacity and keep the backup in a cool place with 40% charge

9.6.7 Transporting the battery by plane

Disabled mobility aids are usually accepted by air carriers, with limitations depending on local regulations, which are always changing. We recommend that you ALWAYS contact your carrier and indicate the specifications of your TIBODA in terms of battery power, to prevent it being refused at check-in. The battery power can be calculated by multiply its Ah rating by its Voltage rating. $[Wh]=[Ah]*[V]$ Check the ET Label. Contact Moretti Spa for further details.

9.7 Using the throttle

Fig.10

The throttle is located on the side of the RH handlebar handle. Press the lever shown in Fig.10 to accelerate.

9.8 Using the brake


Fig.11

Pull the lever on the LH handle (Fig.11). If the speed is higher than 8 kph when you use the mechanical brake, the regenerative brake will cut in to recharge the battery. Below that speed, the brake lever only actuates the mechanical disk brake. It is in any case impossible to accelerate until the brake lever has been completely released.

9.9 Using cruise control

Fig.12

Press the throttle lever until you reach your desired speed, then press the red button shown in Fig.12 at the same time to enable cruise control. When you release the two controls, the drive will hold the speed constant.

Using the brake or throttle disables cruise control.

9.10 Reverse button

The reverse button is on the left handle next to the brake lever and cruise control button (Fig.13). Set the button to “I” to enable reverse drive, and return it to “O” to resume normal forwards drive.

9.11 Adjusting the handlebar

To adjust the height and angle of the handlebar, undo the indicated bolts with a 5 mm Allen key.

10. MAINTENANCE

ARDEA devices manufactured by Moretti S.p.A. Are CE marked and carefully checked before being commercialised. For the safety of the patient, your device should be checked at least once every year by the manufacturer or an authorised laboratory to verify its suitability for use. Use only original spare parts and accessories for repairs and have the system checked by an electrician afterwards.

10.1 Ordinary maintenance

- Clean the controller and buttons
- Clean the frame to prevent any build up of dirt and fluid
- Check that the electrical connections are hooked up properly
- Keep the wheels clean of dirt and residue of any kind
- Lightly lubricate the wheel axles and springs every 3 months if necessary
- Inspect the tyres. The tread should not wear more than 1.5 mm
- Charge the batteries regularly
- Check that the throttle works properly, and that it returns to centre when released. In case of doubt, contact an authorised service centre
- Check the cable covers for damage and wear. In case of doubt, contact an authorised service centre

MAINTENANCE	DAILY	WEEKLY	MONTHLY	ANNUAL
MECHANICAL				
LUBRICATE THE MOUNTING SYSTEM			✓	
LUBRICATE THE STEERING RETURN MECHANISM			✓	
CHECK THE BRAKE	✓			
CHECK THE TIGHTNESS OF THE FASTENERS	✓			
CLEANING				
HAVE THE TIBODA CHECKED BY AN AUTHORISED SERVICE CENTRE				✓
ELECTRICAL				
BATTERY INDICATOR - CHECK THE BATTERY INDICATOR TO SEE WHETHER THE BATTERY NEEDS RECHARGING	✓			
CONTROLLER / DISPLAY - CHECK FOR DAMAGE AND EXPOSED CABLES			✓	
CHECK THAT THE SOCKETS AND CABLES ARE NOT LOOSE			✓	
CHECK THAT THE BATTERIES ARE COMPLETELY CHARGED BEFORE DAILY USE	✓			
CHECK THAT ALL BOLTS AND TIGHT	✓			
CHECK THAT THE LIGHTING SYSTEM IS WORKING PROPERLY (IF PRESENT)	✓			
WHEELS AND TYRES				
CHECK THE TYRE PRESSURE	✓	✓		
THE WHEEL MUST RUN SMOOTHLY WITHOUT OBSTRUCTION		✓		
THE WHEEL MUST RUN WITHOUT WOBBLING			✓	
INSPECT THE TREAD OF THE TYRE. IF IT IS LESS THAN 1MM (1/32"), REPLACE IT WITH A NEW TYRE (AVAILABLE FROM YOUR LOCAL RESELLER)			✓	

10.2 Annual maintenance

We recommend that you have your product serviced completely by an authorised service centre at least once a year. This ensures that it works properly without future problems.

The annual service should include the following checks:

- Wheels and tyres
- Battery, cable and connections
- The controller's programming in relation to the user's requirements
- Frame, bracket and fasteners (nuts and bolts)

10.3 Wheel maintenance

Replacing the TIBODA tyre is a seemingly simple operation that does not require any special skills, however, given the importance of the part in question it is always advisable to contact your trusted dealer.

WARNING!

Always use original parts TIBODA, MORETTI S.p.A. is not responsible for incorrect assembly or assembly with unauthorized components including bolts..

A **torque wrench** is required for correct mounting of the wheel so that the correct tightening of the bolt can be adjusted. Avoid tightening without proper measurement.

750 & 1000 W MOTORIZATION - Exact sequence of spacers and screws

Right view

The cable coming out of the engine must be placed in the slot created specifically on the fork.

Left view

300 & 400 W MOTORIZATION - Exact sequence of spacers and screws

Left view

Right view

The cable exiting the motor must be positioned at the front and make a bend to the lower fixing point.

BOLT TIGHTENING

Tightening of the wheel seal bolt must be done with a torque wrench set to 70 N*m. Avoid making non-measurable clamps this could compromise the tightness of the structure.

11. CONTROLS AND REMINDERS

WARNING!

Never expose the device to direct contact with water; it is rated IP21

- Keep the controller clean and protect it from water and rain
- Keep the wheel clean of residue, hair, sand and carpet fibres
- Do not store the product in a damp area; this can cause mould to form and quickly damage its parts
- Inspect and lubricate all moving parts.
Lubricate with Vaseline or light oil. Do not use too much lubricant, the small drops can cause stains
- Always check that all nuts and bolts are tight
- Do not drop the product or bang it against property or persons
- Avoid extreme temperature differences
- Do not expose the product to heat or direct sunlight for long periods
- Do not disassemble the product

11.1 Steering tube

- These components are lubricated and sealed and do not require further lubrication

11.2 Tyre pressure

- Maintain the tyre pressure marked on the tyre itself - **(recommended: 2.5 bar)**
- It is important to maintain the specified tyre pressure at all times. Do not over-inflate the tyres. Low pressure can result in loss of control, over-inflation can cause the tyre to explode
- Failure to maintain the specified tyre pressure can damage the tyre and wheel. Inspect the tyres regularly for wear
- REPLACING THE WHEEL: If the tyre is punctured, replace its inner tube and contact your authorised reseller

11.3 Battery terminal connections

- Make sure the connectors are tightly connected and not corroded
- The batteries must be in their compartments

11.4 Cabling, controller, charger, electronics

- Regularly check the electrical connections
- Regularly check the cable insulation, including the charge cable, for wear and damage
- Have any connectors, connections or damaged cables replaced or repaired by an authorised service centre before you use the product
- Protect these components from damp
- if exposed to damp, dry them off completely before you use the product

11.5 Storing the device

If you do not intend to use the product for a long time, certain precautions are in order:

- Completely recharge the battery before you store them
- Disconnect the batteries from the device
- Keep the device in a warm, dry place
- Do not keep the TIBODA in a place exposed to extreme temperatures
- Operating conditions (0°* ~ +40°) - *With display (2° ~ +40°) - disconnect the display if necessary
- Storage conditions (-10° ~ +60°)

11.6 Transport

1. Remove the battery from its bracket
2. Remove the stainless steel mounting
3. The TIBODA can be charged

WARNING!

Turn the device off before carrying it onto public and private transport. Always turn "TIBODA" off during transport. During transport, make sure that the cables and handlebar controls are in no danger of being crushed. This can damage the device and render it inoperable

12. SPARE PARTS AND ACCESSORIES

Use only spare parts and accessories listed in the Moretti general catalogue:

TBA001	CLAMP PROGEO TECNA ADVANCE
TBA002	CLAMP QUICKIE HELIUM
TBA003	CLAMP D20.5
TBA004	CLAMP D22.5
TBA005	CLAMP RH D25.4 EXELLE
TBA006	CLAMP LH D25.4 EXELLE
TBA007	CLAMP D26
TBA008	CLAMP D28
TBA009	CLAMP D29
TBA010	CLAMP D30
TBA020	PAIR OF SIDE PLATES STANDARD
TBA021	PAIR OF SIDE PLATES L.36
TBA022	PAIR OF SIDE PLATES EXELLE
TBA023	PAIR OF SIDE PLATES EXELLE L.36
TBA030	TETRA ACCELERATOR
TBA031	TETRA BRAKE
TBA100	CUSTOM PAINT - TIBODA
TBA101B	PAIR OF ROTARY STANDS, WHITE - TIBODA
TBA101N	PAIR OF ROTARY STANDS, BLACK - TIBODA
TBA102	DIGITAL DISPLAY - TIBODA
TBA110D	REAR VIEW MIRROR, RH - TIBODA
TBA110S	REAR VIEW MIRROR, LH - TIBODA
TBA120	LED LIGHT WITH HORN 36V - TIBODA
TBA121	LED LIGHT WITH HORN 48V - TIBODA
TBA125	FRONT LED LIGHT, USB - TIBODA
TBA127	REAR LED LIGHT, BATTERY - TIBODA
TBA129	GLOBE LED LIGHT KIT, BATTERY - TIBODA
TBA130	HANDLEBAR ACCELERATOR 36V - TIBODA
TBA131	HANDLEBAR ACCELERATOR 48V - TIBODA
TBA150	KIT FOR INSTALLATION. PLATES + COMPUTER - TIBODA
TBA151	TILT METER FOR INSTALLATION - TIBODA
TBA204*	AIRPLANE BATTERY 36V 7.8Ah PER TIBODA 300W
TBA205*	AIRPLANE BATTERY 48V 5.2Ah PER TIBODA 400W
TBA206*	AIRPLANE BATTERY 36V 5.8Ah PER TIBODA 750W
TBA207*	AIRPLANE BATTERY 48V 5.8Ah PER TIBODA 1000W


* Read paragraph 8.6.7 in this manual

12.1 TBA102 Digital display - TIBODA

TIBODA can be equipped with an LCD display which provides the following information:

- speed
- total distance (km)
- battery charge
- motor power (current)
- error codes

All display functions can be controlled with the keypad next to the handlebar. The accessory can be retrofit to the drive, to do so proceed as indicated below.

Display control buttons

- 1 UP
- 2 MODE
- 3 DOWN

Fig.15

The accessory can be retrofit to the drive, to do so proceed as indicated below.

12.1.1 Installation

1. Install the display mount and control buttons as shown. (Fig.16)
2. Open the cable compartment by undoing the 3 screws and route the display cable through the hole at the top of the compartment. (Fig.17)
3. Connect the 5 pin socket in the TIBODA cable compartment (Fig.16). Close the compartment. (Fig.18)

Fig.17 - A

Fig.17 - B

Fig.17 - C

Fig.18

12.1.2 Display interface


Fig.19

12.1.3 Reading the information




Fig.20

Press "MODE" briefly to display, in sequence:


- SPEED (current speed in kph)
- AVG SPEED (average trip speed in kph)
- SPEED MAX (maximum trip speed in kph)
- TRIP (Trip distance)
- TRIP TIME (Trip time)
- ODO (Total distance)


Current speed in kph


Average trip speed in kph


Maximum trip speed in kph

Trip distance

Trip time

Total distance (km)

*PAS, BMS, ECO, CITY, TOUR, POWER, BOOST are not shown on the display

12.1.4 Battery charge indicator

When the battery is fully charged, the display shows 5 bars.

When the battery is undercharged, the display shows a single flashing bar.

Battery fully charged


Battery almost discharged

Battery discharged
The indicator's last bar flashes

12.1.5 Power output indicator

While you are driving, the LCD display shows the motor output power at the bottom right. The value is shown in Watt.

12.1.6 Turning the LCD display's backlighting on and off

Hold down the "UP + MODE" buttons to turn the backlighting on

Fig.21

Hold down the "UP + MODE" buttons again for 3 seconds to turn the backlighting off



Fig.22


*PAS, BMS, ECO, CITY, TOUR, POWER, BOOST are not shown on the display

12.1.7 Units


The speed and total distance can be expressed in km or miles


Total distance in km

Total distance in miles

Speed in kph

Speed in mph

*PAS, BMS, ECO, CITY, TOUR, POWER, BOOST are not shown on the display

12.1.8 General settings

WARNING!

Read the following instructions carefully before you make any changes

- Turn the device on with the battery button. The display will turn on
- When the display is on, hold the "UP + DOWN" buttons down for 3 seconds to enter the general settings menu

Press the "MODE" button for one second to scroll through the settings: **SE1, SE2, SE3.**

SE1 - wheel diameter setting

Use the "UP" and "DOWN" buttons to modify the wheel diameter and make sure the speed is correct. The default setting is 17.5 inches. (This includes the tyre)

SE2 - wheel diameter factor adjustment

Set the wheel diameter factor for more precise speed indication with the "UP" and "DOWN" buttons. The default setting is 1.00, with a range of 1 to 3

SE3 - display backlighting setting

The brightness of the display can be set in a range of 1-3

Use the "UP" and "DOWN" buttons to set the brightness

- 0: none
- 1: high
- 2: moderate
- 3: low

The default setting is 0.

12.1.9 Saving the settings

After you have made the adjustments, hold down the "MODE" button to save them and quit the menu.

12.2 TBA120 Led light with horn 36V - TIBODA - TBA121 Led light with horn 48V - TIBODA
 The TIBODA is equipped with a led light with horn.

Fig.23a

Fig.23b

- Press the button under the battery indicator on the handlebar (Fig.23a) to turn on the led light
- Press the green button next to the cruise control button (Fig.23b) to sound the horn

ENGLISH

The accessory can be retrofit to the drive, to do so proceed as indicated below.

12.2.1 Installation

Fig.24

1. Mount the led light to the bracket on the RH rearview mirror as shown in the photograph. (Fig.24)
2. Open the cable compartment by undoing the 3 screws and route the light cable through the hole at the top of the compartment. (Fig.25)

Fig.25



Fig.26 (illustrative image)

3. Connect the socket in the TIBODA cable compartment (Fig.26). Close the compartment.

12.3 TBA130 Handlebar throttle 36V - TIBODA; TBA131 Handlebar throttle 48V

The TIBODA also has a handlebar throttle.

This accessory can be retrofit to the drive, to do so proceed as indicated below.

12.3.1 Installation


Fig.27

1. Fit the handle throttle on the handlebar as shown in figure (Fig.27)
2. If the lever throttle is already present, disconnect and remove it.
 - a) Now remove the handle using a compressed air gun to blow into the hole on the side of the handle itself. Remove the other handle too if necessary
 - b) Undo the grubscrow on the throttle mount and remove the throttle
3. Fit the throttle handle to the handlebar and secure it with the grubscrow

Fig.28


Fig.29

4. Open the cable compartment by undoing the 3 screws and route the throttle cable through the hole at the top of the compartment (Fig. 28)
5. Connect the 5 pin socket (Fig.29) in the TIBODA cable compartment.
Close the compartment

13. CLEANING AND DISINFECTION

Use a damp cloth and non-abrasive neutral detergent to clean the device's plastic and metal parts. Do not use products which can scratch the product.

WARNING!

- DO NOT immerse the device in water for any reason!
- Sea water and sand can damage parts of the "Tiboda" if they remain in contact with it for a long time. Always clean the "Tiboda" after using it at the seaside
- Never use acid, alkali or solvents (acetone or diluent) to clean it


NOTES: For more thorough disinfection, use a small amount of neutral detergent and dry it off very carefully, making sure that no fluid enters the product

- If necessary, clean the product with an approved disinfectant. Make sure the disinfectant is safe for the product before you use it
- Read the disinfectant's and detergent's safety instructions before using them. Failure to do so may result in irritation of the skin and damage to the device's controls and finish

14. TROUBLESHOOTING

FAULT	SOLUTION
TIBODA DOES NOT TURN ON	CHECK THE BATTERY CHARGE - RECHARGE THE BATTERY CHECK THAT THE BATTERY IS PROPERLY INSTALLED AND TURN THE KEY CLOCKWISE. IF THE BATTERY IS CHARGED AND PROPERLY INSTALLED, CONTACT YOUR SERVICE CENTRE
TOO MUCH VIBRATION IN USE	CHECK THAT THE CLAMPS, PLATES AND STRUT ARE PROPERLY FASTENED TO THE WHEELCHAIR
THE BATTERY SWITCH AND STATUS INDICATOR LIGHT UP BUT TIBODA DOES NOT MOVE	CHECK THAT THE PLUG/CABLE CONNECTING THE MOTOR TO THE CONTROLLER, ON THE FORK, IS PROPERLY CONNECTED; IF IT IS NOT, CONTACT YOUR SERVICE CENTRE
THE DISPLAY IS LOCKED AND AN ERROR IS DISPLAYING	SHUT THE PRODUCT OFF AND ON AGAIN. IF THE ERROR PERSISTS, CONTACT YOUR SERVICE CENTRE

14.1 Error codes

When the TIBODA's electrical system has problems, the display automatically reports the error code at the bottom right.

The error codes are listed in the following table

CODE	DEFINITION
0	NO PROBLEM
1	BRAKE LEVER FAULT
2	ACCELERATOR FAULT
3	CONTROLLER FAULT
4	UNDER VOLTAGE
5	MOSFET FAULT
6	OVER TEMPERATURE
7	OVER CURRENT
8	OVER VOLTAGE
9	CONTROLLER COMMUNICATIONS FAULT

WARNING!

The display stops working if a fault occurs. Contact a service centre to resolve the problem

15. CONDITIONS OF DISPOSAL

15.1 General terms of disposal

Never dispose of the product as normal domestic waste. Dispose of the product at a sorted waste collection centre for recycling

15.2 Instructions for correct disposal in accordance with European directive 2012/19/UE

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment

15.3 Treatment of exhausted batteries (Directive 2006/66/EC)

At the end of its working life, the product must not be disposed of together with normal urban waste. It must be delivered to municipal separated waste collection facilities, or to appropriate dealers that provide this service. Separated waste disposal helps to reduce possible negative effects on the environment and health deriving from improper disposal and allows for recycling the materials comprising the product, which translates into significant energy and resource savings. The product bears the barred bin symbol to underline the obligation of disposing of electro-medical equipment

16. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY

The TIBODA device is designed for use in the electromagnetic conditions specified in the following tables. The user of the TIBODA device must make sure that these conditions obtain when he is using it.

Table 1

Manufacturer's guide and declaration - electromagnetic emissions		
The TIBODA device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.		
Emission tests	Conformity	Electromagnetic environment - guide
RF emissions CISPR 11	Group 1	The TIBODA device uses RF energy only for charging the battery. Its RF emissions are therefore very low and in all likelihood do not interfere with nearby electronic equipment.
RF emissions CISPR 11	Class B	The TIBODA device is intended to be used in all buildings, including the home and buildings connected directly to the low voltage public power grid for residential applications.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuation/ flicker emissions IEC 61000-3-3	Conforming	

Table 2

Manufacturer's guide and declaration - electromagnetic emissions			
The TIBODA device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Information for using the device and charging the battery. The floor must be wood, concrete or ceramic. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2KV for power lines ±1 kV for I/O lines	2KV for power lines ±1 kV for I/O lines	The quality of the mains supply must be that of a typical commercial hospital when charging the battery
Surge IEC 61000-4-5	1 kV differential mode ±2 KV common mode	1 kV differential mode ±2 KV common mode	The quality of the mains supply must be that of a typical commercial or hospital when charging the battery.
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
Voltage drop outs, brief interruptions and variations of voltage on the power supply lines. IEC 61000-4-11	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	<5% UT (>95% UT drop out) for 0.5 cycles 40% UT (60% UT drop out) for 5 cycles 70% UT (30% UT drop out) for 25 cycles <5% UT (>95% UT drop out) for 5 s.	The quality of the mains supply must be that of a typical commercial hospital when charging the battery.
Mains frequency magnetic field (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. (Applies both to using the device and charging the battery).
Note: UT is the AC mains voltage prior to application of the test level.			

ENGLISH

Table 3

Manufacturer's guide and declaration - electromagnetic immunity			
The TIBODA device is intended to be used in the following electromagnetic conditions. The client or user must ensure that it is used in these conditions. Make sure it is used in these conditions.			
Immunity tests	Test level IEC 60601	Conformity level	Electromagnetic environment - guide
conducted RF IEC 61000-4-6	3 V eff. From 150 KHz to 80 MHz	3 V eff.	<p>Portable and mobile RF communications equipment should be used no closer to any part of the [equipment or system], including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d = 1,2\sqrt{P}$</p> <p>$d = 1,2\sqrt{P}$ from 80 MHz to 800 Mhz</p> <p>$d = 2,3\sqrt{P}$ from 800 MHz to 2,5 GHz</p> <p>$d = 9,3$ m during recharge operation</p> <p>where "P" is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and "d" is the recommended separation distance in metres (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
radiated RF IEC 61000-4-3	3 V/m From 80 MHz to 2.5 GHz	3 V/m	
<p>NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a The intensity of the field for fixed transmitters, such as base stations for radio-telephones (cellulars and cordless) and terrestrial mobile radios, amateur radio equipment, AM and FM radio transmitters and TV transmitters may not be known precisely in theory. To establish an electromagnetic environment caused by fixed RF transmitters, one must run an on-site electromagnetic inspection. If the field intensity measured in the place in which the Tiboda is used exceeds the RF level of conformity indicated above, the function of the Tiboda must be monitored. If abnormalities are observed, supplementary measures may be necessary, such as re-orientation or repositioning.</p> <p>^b The field intensity over the frequency range from 150 kHz to 80 MHz should be less than 3 V/m.</p>			

Table 4


Recommended separation distances between portable and mobile radio communications equipment and TIBODA						
<p>The TIBODA device is intended to operate in an electromagnetic environment in which RF disturbances are under control. The client or user of TIBODA can contribute to preventing EM interference by ensuring a minimum distance between the mobile and portable RF communications devices (transmitters) and the TIBODA battery charger during charging and the internal TIBODA device during its operation as indicated below, in relation to the maximum output power of the transmitters</p>						
Tests	Separation distance as a function of the frequency of the transmitter m					
Maximum rated output power of transmitter W	During battery charging	During the use of TIBODA	During battery charging	During the use of TIBODA	During battery charging	During the use of TIBODA
	from 150 kHz to 80 MHz d = 1,2√P	from 150 kHz to 80 MHz d = 1,2√P	from 80 kHz to 800 MHz d = 1,2√P	from 80 kHz to 800 MHz d = 1,2√P	from 800 kHz to 2.5 GHz d = 1,2√P	from 800 kHz to 2.5 GHz d = 1,2√P
0.01	0.12	0.12	0.12	0.12	0.23	0.23
0.1	0.38	0.38	0.38	0.38	0.73	0.73
1	1.2	1.2	1.2	1.2	2.3	2.3
10	3.8	3.8	3.8	3.8	7.3	7.3
100	12	12	12	12	23	23
<p>For transmitters whose maximum rated output power is not listed above, the recommended separation distance d, in metres (m), can be calculated using the equation applicable to the transmitter frequency, where P is the maximum rated output power of the transmitter, in Watt (W), as given by the transmitter's manufacturer.</p> <p>NOTE 1: At 80 MHz and 800 MHz, apply the separation distance for the higher frequency band.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>						

ENGLISH

17. TECHNICAL SPECIFICATIONS

17.1 Dimensions

TB030X - TB040X - TB075X - TB100X TIBODA front drive for wheelchairs

Weight of frame, TB030X and TB040X:	13.2 Kg
Weight of frame, TB075X and TB100X:	15.2 Kg
Weight of battery, 36V 10.4Ah:	2.5 Kg
Weight of battery, 36V 13Ah:	2.9 Kg
Weight of battery, 48V 10.4Ah:	3.1 Kg
Weight of battery, 48V 11.6Ah:	3.2 Kg
Weight of mounting equipment:	2.5 Kg
Total weight, TB030X:	18.2 Kg
Total weight, TB040X:	18.8 Kg
Total weight, TB075X:	20.2 Kg
Total weight, TB100X:	20.8 Kg

*The weights and dimensions may vary by $\pm 2\%$

17.2 Technical specifications

TB030X: TIBODA front drive for wheelchairs

- ¹Motor 300W
- Reinforced aluminium rim 16"x1.95 with tyre and inner tube
- Mechanical disk brake 160mm
- ²Motor brake (EBS) with regeneration, actuated by brake lever
- Samsung lithium ion cell battery: 36V - 10.4Ah
- USB socket on battery
- Powder coated aluminium frame (RAL9010 and RAL9005)
- Stand-alone frame with dia. 40mm wheels (stand)
- ³Stainless steel mountings with quick coupling, for use by tetraplegic users
- Automatic steering return
- Reverse gear
- Cruise control
- Battery charge indicator
- ⁴Autonomy 50km (mixed route with 15% uphill - tested with 70kg patient - Tb040x)
- Speed control with quad-style throttle lever, homologated (15 kph) pursuant to EN12184:2014 - Electric wheelchairs, scooters and their charging systems
- Battery charger (standard supply); Input 220-240V 50hz; Output 42V - 2A DC
- ⁵Charging time 4/5 hours (depending on wear of battery)

TB040X: TIBODA front drive for wheelchairs

- ¹Motor 400W
- Reinforced aluminium rim 16"x1.95 with tyre and inner tube
- Mechanical disk brake 160mm
- ²Motor brake (EBS) with regeneration, actuated by brake lever
- Samsung lithium ion cell battery: 48V - 10.4Ah
- USB socket on battery
- Powder coated aluminium frame (RAL9010 and RAL9005)
- Stand-alone frame with dia. 40mm wheels (stand)
- ³Stainless steel mountings with quick coupling, for use by tetraplegic users
- Automatic steering return
- Reverse gear
- Cruise control
- Battery charge indicator
- ⁴Autonomy 50km (mixed route with 15% uphill - tested with 70kg patient - Tb040x)
- Speed control with quad-style throttle lever, homologated (15 kph) pursuant to EN12184:2014 - Electric wheelchairs, scooters and their charging systems
- Battery charger (standard supply); Input 220-240V 50hz; Output 54.6V - 2A DC
- ⁵Charging time 4/5 hours (depending on wear of battery)

TB075X: TIBODA front drive for wheelchairs

- ¹Motor 750W
- Reinforced aluminium rim 16"x1.95 with tyre and inner tube
- Mechanical disk brake 160mm
- ²Motor brake (EBS) with regeneration, actuated by brake lever
- Samsung lithium ion cell battery: 36V - 13Ah
- USB socket on battery
- Powder coated aluminium frame (RAL9010 and RAL9005)
- Stand-alone frame with dia. 40mm wheels (stand)
- ³Stainless steel mountings with quick coupling, for use by tetraplegic users
- Automatic steering return
- Reverse gear
- Cruise control
- Battery charge indicator
- ⁴Autonomy 50/55km (mixed route with 15% uphill - tested with 70kg patient - Tb040x)
- Speed control with quad-style throttle lever, homologated (15 kph) pursuant to EN12184:2014 - Electric wheelchairs, scooters and their charging systems
- Battery charger (standard supply); Input 220-240V 50hz; Output 42V - 2A DC
- ⁵Charging time 4/5 hours (depending on wear of battery)

TB100X: TIBODA front drive for wheelchairs

- ¹Motor 1000W
- Reinforced aluminium rim 16"x1.95 with tyre and inner tube
- Mechanical disk brake 160mm
- ²Motor brake (EBS) with regeneration, actuated by brake lever
- Samsung lithium ion cell battery: 48V - 16,6Ah
- USB socket on battery
- Powder coated aluminium frame (RAL9010 and RAL9005)
- Stand-alone frame with dia. 40mm wheels (stand)
- ³Stainless steel mountings with quick coupling, for use by tetraplegic users
- Automatic steering return
- Reverse gear
- Cruise control
- Battery charge indicator
- ⁴Autonomy 50/55km (mixed route with 15% uphill - tested with 70kg patient - Tb040x)
- Speed control with quad-style throttle lever, homologated (15 kph) pursuant to EN12184:2014 - Electric wheelchairs, scooters and their charging systems
- Battery charger (standard supply); Input 220-240V 50hz; Output 54.6V - 2A DC
- ⁵Charging time 4/5 hours (depending on wear of battery)

¹Motor

- Effective power output to the wheel
- Motor with integrated sine wave controller. Eliminates vibration for improved comfort and control. Can handle slopes of around 15%.

²Motor brake (EBS: Electronic Brake System)

- Regenerative brake provides additional braking power (activated automatically when the brake lever is applied at over 8 kph). Switches the motor to generator mode, thus recharging the battery.
- The EBS increases the safety of the vehicle and its circulation on the road thanks to reduction of braking, braking stability and control braking system.

³Stainless steel mountings with quick coupling

- Just a few moments and motions are enough to couple the TIBODA to your wheelchair, even for persons with reduced manual ability (TETRAPLEGICS)

4The autonomy of the TIBODA varies in relation to:

- Weight of the patient
- Slopes along the route
- Battery wear
- Climate, operating temperature, atmospheric agents

5The battery recharge time can vary in relation to:

- Battery wear
- Wear of the battery charger
- Output of the battery charger
- Climate, operating temperature, atmospheric agents

Operating conditions:	Temperature: min 0°C - max +40 °C Humidity: max 90 %, non condensing Atmospheric pressure: 800 -1060 hPa
Transport conditions:	Temperature: min -10°C - max +60 °C Humidity: max 90 %, non condensing Atmospheric pressure: 500 -1060 hPa
Storage conditions:	Temperature: min 0°C - max +60 °C Humidity: max 90 %, non condensing Atmospheric pressure: 700 -1060 hPa

18. WARRANTY

Moretti products are guaranteed for 2 (two) years from the date of sale against material and manufacturing defects, subject the following limitations. The warranty is voided by improper use, abuse, modifications to the product and failure to follow the instructions. The intended use of the product is given in the user manual.

Any modification to the product, including reprogramming of the TIBODA controller to change its factory settings, voids the warranty. Moretti is not liable for damage, injury or any other consequences resulting from installation or use which are not scrupulously conforming with the instructions given in the installation, assembly and user manual.

Moretti does not guarantee its products against damage or defects in the following circumstances: natural disaster, unauthorised repair or maintenance, improper electric power supply (as applicable), use of parts or components not supplied by Moretti, failure to follow the guidelines and instructions, tampering, shipping damage (other than the original shipping by Moretti), or failure to run maintenance as indicated in the manual.

The following are not covered by the warranty:

- Parts subject to normal wear and tear, such as the tyre
- Mechanical or structural failures caused by improper assembly
- Damage due to improper use resulting in accidents of any kind
- Loss of the product
- Damage due to improper or inadequate maintenance
- Damage caused during shipping and return by third parties
- Damage to products whose serial number has been removed or tampered with
- Damage due to use in competition of any kind
- Damage caused by use on irregular ground and similar activities
- Damage caused by transporting more than one person at a time
- Damage caused by use not envisaged in this manual

18.1 Rechargeable battery warranty

The original and spare batteries are guaranteed for 90 days over a period of 6 months as regards their production and the law. If the battery charger is not used for more than 3 consecutive months, the warranty is void. If the exhausted batteries are not used for more than 3 consecutive days, the warranty is void.

19. REPAIRS**19.1 Repairs under warranty**

If a Moretti product has material or manufacturing defects during the warranty period, Moretti will agree with the client whether the defect is covered by the warranty. Moretti, at its sole discretion, may replace or repair the article at a specified Moretti reseller or its own premises. The costs of labour incurred in repairing the product will be borne by Moretti if it determines that the repair is covered by the warranty. Repair and replacement do not renew the warranty period.

19.2 Repairing a product not covered by the warranty

A product not covered by warranty may be returned for repair only if authorised in advance by Moretti customer service. The costs of labour and shipping incurred by repairs not covered by the warranty are borne by the client or reseller in their entirety. Repairs on products not covered by the warranty are themselves guaranteed for 6 (six) months from the day of reception of the repaired product.

19.3 Non-defective products

The client will be notified if Moretti concludes that the product is not defective after having received and examined it. The product will be returned to the client at his expense.

20. SPARE PARTS

Moretti original spare parts are guaranteed for 6 (six) months from the day of delivery.

21. NON-LIABILITY CLAUSE

Unless otherwise expressly specified in this warranty and within the limits of the law, Moretti makes no declaration, guarantee or condition, express or implicit, including any future declaration, guarantee or condition of sale, suitability for a given purpose, non violation and non interference. Moretti does not guarantee that the use of its product will be uninterrupted and problem-free. The duration of any implicit guarantee under the law is limited to the warranty period, within the limits of the law. Certain states and countries do not permit limitations on the duration of an implicit guarantee or the exclusion of limitation of accidental or indirect damages in relation to consumer products. In said states and countries, certain exclusions and limitations of this warranty may not apply to the user. This warranty is subject to modification without notice.

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