

Aspimed 1.9

- [IT] Aspiratore chirurgico professionale
- [DE] Professionelles chirurgischen Absauggerät
- [EN] Professional surgical suction unit
- [FR] Aspirateur chirurgical professionnel
- [ES] Aspirador cirúrgico profesional

**VERSIONE CON VASO DA 1000 ml / AUSFÜHRUNG MIT 1000 ml-BEHÄLTER / 1000 ml VESSEL
VERSION / VERSION AVEC BOCAL DE 1000 ml / VERSIÓN CON FRASCO DE 1000 ml**


Fig. - Abb. 1a

I

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GB

1. 1000 ml Vessel
2. Vessel plug
3. Protection device
4. ON-OFF button
5. Vacuum regulator
6. Vacuum gauge
7. Multi-voltage power supply unit and battery charger socket
8. 24 cm silicon tube
9. Antibacterial filter
10. Working light
11. Light ON-OFF button
12. Accessory compartment
13. Air inlet connector INLET

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GB

1. 2000 ml Vessel
2. Vessel plug
3. Protection device
4. ON-OFF button
5. Vacuum regulator
6. Vacuum gauge
7. Multi-voltage power supply unit and battery charger socket
8. 20 cm silicon tube
- 8a. 10 cm silicon tube
9. Antibacterial filter
10. Working light
11. Light ON-OFF button
12. Accessory compartment
13. Air inlet connector INLET

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
Fig. - Abb. 2

GB

14. Bag with shoulder strap
15. Disposable manual regulator
16. Multi-voltage power supply unit
17. Power supply cord
18. Disposable sterile cannula
19. Cord with cigarette lighter plug

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- GB**
- Flashing red Led:** the battery is no longer capable of supplying maximum charge duration; replace it.
 - Red Led permanently on:** the battery is flat, use the appliance for a maximum of 3 minutes and re-charge the battery
 - Yellow Led permanently on:** the battery is in reserve (about 10 minutes' operating time left); recharge it if possible.
 - Green Led permanently on:** appliance operating, battery charged.
Green Led flashing: battery under charge.

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The LTA420 aspirator, specific for secretion removal, is a medical and surgical device suitable for use in the home, in clinics, and also in field facilities and on vehicles. It has a vacuum regulator, a vacuum gauge, a 1000 ml sterilisable vessel, or (depending on the version) a 2000 ml sterilisable vessel; these vessels are equipped with a protective device to prevent liquid from entering the suction pump and interrupting the aspiration flow. It is lifetime lubricated, easy to handle, simple to use, reliable, strong and silent. This model is equipped with a 12V rechargeable battery and can also operate directly through the cigarette lighter in car, boat or other vehicles with suitable socket and voltage.

The LTA420 aspirator, in all its versions, has the following 3A accessories: a 1000 ml vessel or 2000 ml vessel (depending on the version) with protection device, cigarette lighter power supply cord, power supply cord with multi-voltage power supply unit, 24 cm long sterilisable silicone connection tube (only for 1000 ml vessel version), 10 cm and 20 cm long sterilisable silicone connection tubes (only for 2000 ml vessel version), 130 cm sterilisable silicone connection tube, **sterile disposable** cannula, **sterile disposable** manual flow regulator, **disposable** bag (only for 1000 ml vessel version) and **disposable** antibacterial filter.
N.B.: Only use original 3A accessories.

IMPORTANT WARNINGS

This is a medical device and must be used by qualified staff. It must be operated as indicated in this user instruction manual. It is important for the operator to read and understand the information on use and maintenance of the unit. If you have any questions, contact your stockist. MICROBIAL CONTAMINATION: in the case of illnesses with a risk of infection or microbial contamination, the accessories should be thoroughly cleaned and sterilised after each use.

The manufacturer has made every effort to ensure that all its products are of the highest quality and are reliable. Nevertheless, as for all electrical appliances, fundamental safety standards must be observed:

- Children and persons who are not self-sufficient may only use the unit under the strict supervision of a responsible adult who has read this manual.
- The device must always be used by specifically trained staff who have read this manual.
- Never use adapters for supply voltages different to the voltage shown on the data plate on the back of the unit. Keep the cord away from hot surfaces.
- The unit should not be used in the presence of inflammable anaesthetic mixtures with air, oxygen or nitrous oxide.
-  Never handle the supply cord plug with wet hands or use the unit when taking a bath or a shower. Never leave the unit close to water, do not submerge it in any liquids. If by chance it falls into water, pull the plug out of the socket immediately before recovering it. Do not use the unit if the plug or supply cord are worn or wet (send it immediately to your stockist).
- Unit casing not protected against liquid penetration.
- Only authorised personnel may perform maintenance and/or repair work. Unauthorised repairs annul the guarantee.
- Ensure that the connections and container closing are performed carefully to prevent suction losses.
- Do not tip the container over while connected to the unit when operating, as liquid may be sucked inside the appliance causing pump damage. If this happens, switch the aspirator off immediately and empty and clean the container. Send the appliance to your stockist.
- When the protection device intervenes aspiration is interrupted; empty the container and perform the cleaning operations.
- The cannula and the manual control of the aspirated flow are **sterile, disposable** products: They must be replaced after every application.
- Check the use-by-date on the original packaging of the cannula and check the integrity of the sterile packaging. If expired and/or deteriorated, replace it.
- The **disposable** antibacterial filter must be replaced after every application.
- Never use the battery charger with other appliances or for uses other than that established by this manual. Never use LTA420 with other power supply units.
- The power supply cord with cigar lighter plug has a safety fuse for inspection in the event of a fault.
- In view of their length, the power cord and connection hose could constitute a strangulation hazard.

- Use of the device in ambient conditions other than those specified in the manual may seriously impair its safety and technical characteristics.
- In the event of aspiration without the container and/or antibacterial filter, or if it is suspected that substances have entered the aspiration circuit, contact your stockist at once.
- Always use the device in the vertical position on an unobstructed, stable, flat surface.
- The lifetime of the device itself is 5 years, and that of the container and silicone tubes 1 year or 30 sterilisation cycles. However, the lifetime may vary depending on the operating environment.

USE INSTRUCTIONS

Each time before use, ensure that all the accessories are perfectly clean according to the instructions indicated in the “CLEANING AND DISINFECTION OPERATIONS”.

1. OPERATION WITH THE CIGARETTE LIGHTER POWER SUPPLY CORD

- 1.1 Before using the device, check the charge level of the vehicle battery.
- 1.2 For use in a car, camper, boat, etc....., connect the cigarette lighter power supply cord (19) to the device socket (7).
- 1.3 Connect the device as indicated in figure 1a of page 1 for the 1000 ml vessel version or figure 1b of page 2 if the 2000 ml vessel is used.
- 1.4 Start the device by placing the switch to the “I” (ON) position (4) (green pilot light on).
- 1.5 The vacuum regulator (5) can be used to set the level of vacuum required (bar/KPa). Turn the knob clockwise, in the “+” direction, to obtain a higher vacuum or anticlockwise, in the “-” direction, for a lower vacuum; these values are read on the vacuum gauge (6). **Important: the vacuum values on the control decal are purely for guidance; always refer to the vacuum gauge reading.**
- 1.6 After the application, switch off the appliance, detach the cigarette lighter power supply cord from the power socket (19) and perform the cleaning operations as illustrated in the “CLEANING AND DISINFECTION” paragraph.
- 1.7 In the event of a blown fuse in the cigarette lighter plug, replace with a 6.3A – 250V delayed fuse – dimensions Ø 6.3 x 30 mm, by unscrewing the end of the cigarette lighter plug.

2. OPERATION WITH A BATTERY AND/OR WITH THE MULTI-VOLTAGE POWER SUPPLY UNIT

- 2.1 The device is supplied with the battery partially charged. We therefore recommend charging it before use.
- 2.2 To charge the battery, with the device switched off connect the multi-voltage power supply unit (16) to the device's socket (7) and to the mains power supply by means of the cord (17). Charging time: about 4 hours. Charge duration: about 40/45 minutes at maximum suction power.
- 2.3 Operation with internal battery only:
Switch the device on by pressing the “ON-OFF button (4) (green light comes on). If the green light goes out during use, and the yellow light comes on, there is about 10/15 minutes of battery charge left (reserve level). Therefore: terminate the application if possible. If the application cannot be terminated, it can be continued until the red (battery flat) light comes on, but if this occurs do not use the appliance for more than 3 minutes to avoid damaging the battery. To continue the application with the battery flat (red light on), connect the multi-voltage power supply unit (16) (as described in point 2.2).
- 2.4 At the end of the application, switch the device off by pressing the button (4). Connect the power supply unit (16) to charge the battery (as described in point 2.2). When the appliance is not in use, we recommend leaving the power supply unit (16) connected so that the battery charge level is always optimum.
- 2.5 For aspiration procedures see points 1.3; 1.4; 1.5.

N.B.: when the battery power drops below a set level the aspirator shuts down to avoid damage to the battery. If absolutely necessary, the user can restart the aspirator for 1 minute by pressing the ON/OFF button again.

3. SECRETION COLLECTION VESSEL - 1000 ml

The 1000 ml collection bottle supplied with the aspirator can be used in two ways: as a collection vessel which can be sterilised as shown in Figure 3a or as a collection vessel with disposable bag (9) as illustrated in Figure 3b.

- 3.1 Sterilisable secretion vessel (1): the secretion vessel set consists of an overflow valve, a vessel (1) in clear material (polycarbonate) and a blue plug (2). Fit the antibacterial filter (6) straight into the plug (2); it will only fit into the hole marked VACUUM/VUOTO, which is outlined in yellow. The antibacterial filter also protects the aspiration circuit from any contaminating agents sucked in during use. **Do not use the aspirator without**

the antibacterial filter, because from a bacteriological point of view, it becomes dangerous for the patient. Keep the device vertical to allow the overflow to function correctly.

All the components of the vessel can be sterilised using a conventional system in an autoclave at a temperature of 121° C, or by boiling for 10 minutes. We recommend replacing the complete vessel at every 30 sterilisation cycles. Do not overturn the vessel during use, in order to prevent the intervention of the non-return valve (3-4-5); should this occur, switch the aspirator off and detach the tube connected to the antibacterial filter. Never use the aspirator without the secretion collection vessel and/or without the antibacterial filter.

3.1.1 Connection: connect one end of the short, sterilisable silicon tube (7) to the antibacterial filter connector (6) and insert this latter into the "VACUUM" hole of the blue top (2); connect the other end to the "INLET" connector (13) of the aspirator. Connect one end of the long sterilisable silicon tube (8) to the "PATIENT/PAZIENTE" connector of the blue top (2); to the other end, connect the disposable, sterile manual regulator (11) and connect the disposable, sterile cannula (10) to this latter.

- 3.2 Secretion collection vessel with disposable bag (9). The aspirator can be used with the re-usable 1000 ml. transparent secretion collection vessel (1) and with the disposable bag (9) supplied. In this case, the antibacterial filter is integrated in the disposable bag and therefore the antibacterial filter (6) and blue top with valve (2) are not used. The filter integrated in the bag also performs the function of preventing the return of aspirated liquids towards the aspirator. **In this case, to restore device operation, it is necessary to replace the disposable bag.** For the cleaning and disinfecting operations of the tubes (7-8) and vessel (1), sterilise the single parts in an autoclave at a maximum temperature of 121° C, or by boiling for 10 minutes. The bag is disposable and must always be replaced after every use. **The bag must be completely inserted in the vessel in order to prevent any vacuum losses.**

N.B.: maximum disposable bag usage vacuum: – 0.75 bar (75 kPa).

3.2.1 Connection: connect one end of the short sterilisable silicon tube (7) to the yellow connector (VACUUM) of the top (photograph 3b) and the other end to the "INLET" connector (13) of the aspirator. Connect one end of the long sterilisable silicon tube to the red connector (PATIENT) and connect the disposable sterile manual regulator (11) and the disposable sterile cannula (10) to the other.

N.B.: only use the disposable bags supplied by 3A - Code 3A1687.

4. SECRETION VESSEL - 2000 ml (Fig. 4)

The secretion vessel set consists of an overflow valve, a vessel (1) in clear material (polycarbonate) and a blue plug (2).

Connection: connect one end of the 10 cm tube (8a) to "VACUUM" connector of the blue top (2) and the other end connect it to the antibacterial filter connector (7). Then connect one end of the 20 cm tube (8) to the free connector of the antibacterial filter (7) and the other end to the "INLET" connector (13) of the aspirator. The antibacterial filter also protects the aspiration circuit from any contaminating agents sucked in during use. Do not use the aspirator without the antibacterial filter, because from a bacteriological point of view, it becomes dangerous for the patient. Keep the device vertical to allow the overflow to function correctly. Connect one end of the 130 cm tube (9) to the "PATIENT" connector of the blue cap (2), while at the other end to connect the disposable, sterile manual regulator (11) and connect the disposable, sterile cannula (10) to this latter.

All the components of the vessel can be sterilised using a conventional system in an autoclave at a temperature of 121°C, or by boiling for 10 minutes. We recommend replacing the complete vessel at every 30 sterilisation cycles. Do not overturn the vessel during use, in order to prevent the intervention of the non-return valve (3-4-5); should this occur, switch the aspirator off and detach the tube connected to the antibacterial filter. Never use the aspirator without the secretion collection vessel and/or without the antibacterial filter.

DESCRIPTION OF USE OF WORKING LIGHT

The LTA420 is fitted with a LED working light (10) which allows use of the device in dark conditions and/or at night.

Press the button (11) to switch on the light, as shown in Figure 6 of page 5; press the button (11) again to switch it off. The light switches off automatically after about 60 minutes.

The light is of LED type, so the amount of current it consumes is minimal and it can be used even during treatment; however, it should be switched off whenever it is not required and/or at the end of the treatment.

CLEANING AND DISINFECTION OPERATIONS

N.B.: If using chemical disinfectants, follow the manufacturer's instructions exactly.

- The cannula and the aspirated flow manual command are **sterile, disposable** products and must be replaced after every application.
- The antibacterial disposable filter must be replaced after every application.
- Never leave the appliance in water or submerged; clean the external casing of the appliance using only a damp cloth with detergent (non abrasive).

PROBLEMS, CAUSES AND SOLUTIONS

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Excessive noise.	Damaged pump or blockages in the internal aspiration circuit.	Send to the assistance circuit.
The unit switches on but does not aspirate.	- Damaged pump. - Vacuum regulator fully open. Connection tubes disconnected and/or badly connected, broken connection tubes. Container not in a vertical position, full, or defective overflow valve. Possible blockage of the hydraulic circuit inside the unit.	- Send to the assistance circuit. - Check the position of the vacuum regulator. Check the connections and the integrity of the tubes. Position the container in a vertical position, check the overflow valve (blocked) and/or replace the silicon tubes.
The vacuum rate cannot be regulated.	Damage to the internal hydraulic system or blockage of the connection tubes to the aspiration unit.	Send to the assistance circuit.
When the appliance is switched on, the protection fuse always trips.	Pump probably damaged or in shortcircuit.	Send to the assistance circuit.
The vacuum gauge does not work.	Liquid penetrating the pneumatic circuit.	Send to the assistance circuit.

Note: if you experience faults or malfunctioning problems different to those listed above, always and exclusively contact authorised assistance centres.

TECHNICAL CHARACTERISTICS

Professional portable medical device for secretion removal, suitable for fixed and transportable use. Complete with flame retardant ABS V0 casing and piston-cylinder electrical compressor lubricated for life.

Risk class under Directive 93/42/EEC:	Ila
Suction class:	High vacuum / High flow
Voltage:	12VDC
Absorption:	3A
Adjustable vacuum level:	from 0 to -0.85 bar (-85KPa)
Vacuum gauge precision class:	2.5% under UNI EN 837
Air flow:	30 lt/min operating without restrictions +/- 10 %
Temporary use:	max. 45 minutes
Dimensions:	38 (length) x 13.5 (width) x 26 (height) cm (1000ml Vessel) 38 (length) x 13.5 (width) x 35 (height) cm (2000ml Vessel)
Weight:	4.5 kg approx (1000ml Vessel) 4.7 kg approx (2000ml Vessel)
Noise level:	55dBA
Multi-voltage switching battery charger:	PRI: 100 ÷ 240V~ 50 – 60Hz SEC: 14VDC 4,28A
Internal battery:	12VDC 4Ah Lead, hermetic
Battery charge duration:	45 minutes at maximum suction power
Fuse for cigarette lighter power supply cord:	F6.3A - 250V delayed fuse - dimensions Ø 6.3 x 30 mm
Operating conditions:	
Temperature:	min. 0° C; max 40° C
Air humidity:	min. 10 %; max 95 %
Storage conditions:	
Temperature:	min. -10° C; max 50° C 
Air humidity:	min. 10 %; max 95 % 
Operating-storage conditions:	min. 690 hPa; max 1060 hPa 

SYMBOLS USED

-  Type BF appliance
-  It is compulsory to carefully read the instructions before using this device
- IP21** The device's casing is protected against solids 12.5 mm or more in diameter, against vertically falling drops of water and against finger access to hazardous parts.
-  Switch On
-  Switch off (or battery on charge)
-  Direct current
-  Alternate current
-  Do not use the unit when taking a bath or a shower
-  0051 Conforms to Directive 93/42/EEC for medical devices
-  The device contains a hermetic lead battery. It must be disposed of in accordance with current regulations on the disposal of toxic-harmful waste
-  STERILE EO Sterilisation by ethyl oxide
-  Disposable

Electromagnetic Compatibility Compliance levels according to EN 60601-1-2:2015 standard

- ESD immunity: 15 kV air, 8 kV contact (EN 61000-4-2)
- Burst immunity: 2 kV/100 kHz (EN 61000-4-4)
- Surge immunity (EN 61000-4-5): 1 kV common mode /2 kV differential mode
- Magnetic field (EN 61000-4-8): 30 A/m
- Immunity to rf currents in the 150 kHz-80 MHz range (EN 61000-4-6) 3 V modulation 80% 1 kHz
- RF emissions, CISPR 11: Class B
- Harmonics emissions, EN 61000-3-2: Class A

Rf field immunity (EN 61000-4-3):		
Field (V/m)	Frequency	Modulation
3	80MHz 2700MHz	1kHz AM 80%
27	380MHz 390MHz	18Hz PM 50%
28	430MHz 470MHz	18Hz PM 50%
9	704MHz 787MHz	217Hz PM 50%
28	800MHz 960MHz	18Hz PM 50%
28	1700MHz 1990MHz	217Hz PM 50%
28	2400MHz 2570MHz	217Hz PM 50%
9	5100MHz 5800MHz	217Hz PM 50%

Warnings:

Although compliant with the EN 60601-1-2 standard, the LTA420 medical device may interfere with other devices in the vicinity. The device must not be used in proximity to or stacked on top of other equipment. Install the device well away from other equipment that emits high frequencies (short waves, microwaves, electric scalpels, cell phones).

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are under control. The customer or user can help prevent electromagnetic interference by maintaining a minimum distances between mobile and portable RF communication equipment (transmitters) and the medical device as recommended below, according to the maximum output power of the radio communication equipment.

Rated maximum output power of transmitter (W)	Separation distance (m) in relation to transmitter frequency		
	from 150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	from 80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	from 800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with rated maximum output power not listed above, the recommended separation distance d in metres (m) may be determined using the equation adopted for the transmitter frequency, where P is the maximum rated output power of the transmitter in Watts (W) stated by the transmitter manufacturer.

Notes:

- (1) At 80 MHz and 800 MHz the highest frequency range applies.
- (2) These guidelines might not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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DISPOSAL PROCEDURE (Dir. 2012/19/Ue-Weee) The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. **Caution!** The wrong disposal of electric and electronic equipment may involve sanctions.

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3A HEALTH CARE S.r.l.
MORETTI S.p.A.