
KYARA®

MA LTK170 02 A_02-2017

[IT] Aerosolterapia
[EN] Aerosol therapy
[FR] Aérosolthérapie
[DE] Aerosoltherapie
[ES] Aerosolterapia

SIMBOLOGIA - SYMBOLS - SYMBOLIQUE - VERWENDETE SYMBOLE - SIMBOLOGÍA

-  Apparecchio di classe II/Class II/Appareil de classe II/Schutzklasse II/Aparato de clase II.
-  Apparecchio di tipo BF/Type BF device/Appareil de type BF/Geräteyp BF/Aparato de tipo BF.
-  Fusibile di protezione/Safety fuse/Fusible de protection/Sicherung/Fusible de protección.
-  Interruttore spento/Switch "OFF"/Interrupteur éteint/Ein-Ausschalter ausgeschaltet/Interruptor apagado.
-  Interruttore acceso/Switch "ON"/Interrupteur allumé/Ein-Ausschalter eingeschaltet/Interruptor encendido.
-  Corrente alternata/Alternating current/Courant alternatif/Wechselstrom/Corriente alterna.
-  Non utilizzare l'apparecchio mentre si fa il bagno o la doccia/Never use the unit while having a bath or shower/Né pas utiliser l'appareil pendant le bain ou la douche/Das Gerät darf nicht während des Badens oder in der Dusche/No utilizar el aparato mientras se hace el baño o la ducha.
-  È obbligatorio leggere attentamente le istruzioni prima di utilizzare questo dispositivo/It is compulsory to carefully read the instructions before using this device/Il est obligatoire de lire attentivement les instructions avant d'utiliser ce dispositif/Vor der Verwendung des Gerätes ist die Bedienungsanleitung aufmerksam zu lesen/Es obligatorio leer atentamente las instrucciones antes de usar este dispositivo.
-  Tenere all'asciutto/Keep dry/Garder au sec/Trocken halten/Mantener seco.

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CE 0434

- | | | | | | |
|---|-----------------------------|----|---------------------|----|--------------------------------|
| 1 | Piston compressor | 7 | Filter lid | 13 | Nosepiece |
| 2 | Power cord | 8 | Air hose | 14 | Assembling nebuliser kit |
| 3 | ON/OFF Switch | 9 | Fasterjet Nebulizer | 15 | Filling the nebulizer |
| 4 | Vial pressure display gauge | 10 | Mouthpiece | 16 | "breath enhanced" nebulization |
| 5 | Vial pressure regulator | 11 | Adult face mask | 17 | Replacing air filter |
| 6 | Air connector | 12 | Child face mask | | |

Dear Customer, this nebuliser is a high quality device for inhalation therapy for asthma, chronic bronchitis, and other respiratory illnesses. This medical device generates aerosols with variable characteristics to adapt to the patient's respiratory pathology and generating a medicinal mist of particles that are small enough to reach even the deepest regions of your lungs and be of maximal benefit. The operation and usage of the device is very straightforward. All common liquid medication can be used for inhalation therapy. Please read through these instructions carefully so that you understand all functions and safety information. We want you to be happy with your LTK170 product. For any further questions you may have, please ask your chemist or usual retailer.

IMPORTANT SAFETY INSTRUCTIONS

This device may only be used for the purposes described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application. This device is not suitable for anaesthesia and lung ventilation. This device should only be used with original accessories as shown in these instructions. Do not use the device if you think it is damaged or notice anything unusual. Never open the device. This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section. Protect it from: - water and moisture - extreme temperatures - impact and dropping - contamination and dust - direct sunlight - heat and cold. Only use the medication prescribed for you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the therapy. All accessories should only be used by one person. Never bend the nebuliser over 60°.

- The unit should not be used in the presence of anesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Device not protected against splashes

Ensure that children do not use the device unsupervised some parts are small enough to be swallowed.

 Use of this device is not intended as a substitute for consultation with your physician.

 In view of their length, the power cord and connection hose could constitute a strangulation hazard.

PREPARATION AND USAGE OF THIS DEVICE

Prior to using the device for the first time, we recommend cleaning it as described in the section «Cleaning and Disinfecting».

1. Assemble the nebuliser kit 14. Ensure that all parts are complete.
2. Fill the nebuliser with the inhalation solution as per your doctor's instructions. Ensure that you do not exceed the maximum level.
3. Connect the nebuliser 9 with the air hose 8 to the compressor 6 and plug the power cord 2 into the socket.
4. Switch the ON/OFF switch 3 to position «I» to turn on the device and place the mouthpiece 10 in the mouth or fit one of the face masks over mouth and nose.
 - The mouthpiece 10 gives you a better drug delivery to the lungs.
 - Choose between adult 11 or child face mask 12 or the nosepiece 13 and make sure that it encloses the mouth and nose area completely.
5. Breathe in and out calmly during the therapy. Sit in a relaxed position with the upper body upright. **Do not lie down while inhaling.** Stop inhalation if you feel unwell.
6. After completing the inhalation period recommended by your doctor, switch the ON/OFF switch 3 to position «O» to turn off the device and unplug it from the socket.
7. Empty the remaining medication from the nebuliser and clean the device as described in the section «Cleaning and Disinfecting».

FASTERJET NEBULIZER

This compressed air system for aerosoltherapy is equipped with the innovative Fasterjet nebuliser with a breath-enhanced valve system 16. This enables the correct quantity of medication to be featured to the respiratory capacity of each patient. The valve controls the flow of medication during inhalation (16/A), thus reducing medication loss during exhalation (16/B). Fasterjet increases the quantity of inhaled medication.

Fasterjet is capable of nebulising efficiently all medications for inhalation therapy. In order to use the nebulizer synchronized to the breathing, it is necessary to mount on the valve-holder the membrane and use the mouthpiece with expiratory valve (14/B). A synchronized functioning of the nebuliser is possible only by using the mouthpiece with expiratory valve.

CLEANING AND DISINFECTING

Thoroughly clean all components to remove medication residuals and possible impurities after each treatment. The compressor 1 and the air hose 8 should be cleaned with a clean, moist cloth.

-  Always wash your hands well before cleaning and disinfecting the accessories.
-  Do not expose the compressor to water or heat.
-  Replace the air hose for each treatment with a new patient or in case of impurities.

 Unplug the device before cleaning.

Cleaning with water

Wash all nebuliser components (except air hose) under warm tap water (max 60 °C) for about 5 minutes adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer. Rinse thoroughly making sure that all deposits are washed away and leave to dry.

Disinfecting

All nebuliser components (except air hose) can be disinfected with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer. Disinfectants are usually available at pharmacies.

Sterilizing with steam

All nebuliser components (except air hose and masks) can be heat steam sterilized up to 121 °C (20 min.) or 134 °C (7 min.). EN554/ISO11134. The sterilization packaging must conform to EN868/ISO11607 and be suitable for steam sterilization. After sterilization always let all components cool down to ambient temperature before further use. Do not repeat sterilization cycle when components are still warm.

MAINTENANCE, CARE, AND SERVICE

Order all spare parts from your dealer or pharmacist.

It is advisable to replace nebulizer after some 100 to 120 treatments on single patient or after about 20 sterilization cycles. Check the filter continually for cleanliness and replace it if dirty, or after a maximum of 70 applications. Spare filters are provided with the device. To replace the filter 17, lift the filter lid 7 on the bottom of the device, applying leverage with a flathead screwdriver; remove the filter you wish to replace and put the new one. Then place the filter lid in its housing, making sure it adheres to the surface.

MALFUNCTIONS AND ACTIONS TO TAKE

The device cannot be switched on

Ensure the power cord 2 is correctly plugged into the socket. Ensure the ON/OFF switch 3 is in the position «I». The device has automatically switched off due to overheating. Wait until it has cooled down and try again.

The device is on but it does not nebulise

Make sure that the nebulisation nozzle (pisper) 9a is inserted in the vial. Ensure the air hose 8 is correctly connected at both ends. Ensure the air hose is not squashed, bent, dirty or blocked. If necessary, replace with a new one. Ensure the nebuliser 9 is correctly assembled 14. Make sure that the air intake filter 17 is not occluded or dirty. In this case, replace it. Make sure that the drug is inside the vial 9.

The device does not work

Thermoprotection intervened for the following reasons: the device was working near heat sources or in environments with temperature higher than 40°C. Let the device cool off for at least 60 minutes, then turn it on again.

GUARANTEE

This device is covered by a **3 year guarantee** from the date of purchase. The guarantee is valid only on presentation of the guarantee card completed by the dealer (see back) conforming date of purchase or the receipt. The guarantee covers only the compressor. The replaceable components like nebuliser, masks, mouthpiece, air hose, and filters are not included. Opening or altering the device invalidates the guarantee. The guarantee does not cover damage caused by improper handling, accidents or noncompliance with the operating instructions.

TECHNICAL DATA

Piston electrocompressor with thermal cut

Nominal voltage: see data label

Safety fuse: see data label

Dimensions: 196x357x185H mm

Weight: 2.5 Kg approx.

Professional aspiration filter: 50/70 µm

Continuous use compressor

MDD 93/42/EEC - Class of risk IIa

Maximum pressure: 250 kPa

Air flow: 15 lt/min

Noise level: 55 dBA (EN13544-1)

Neb rate: 0.60 ml/min (tap closed)

Nebulizer capacity: min. 2 ml - max. 16 ml

Operating pressure: from 0 to 110 kPa approx. (suggested range: from 70 to 110 kPa)

Operating flow: from 0 to 6.8 l/min

MWAD: 2.7 µm GSD: 2.98 Aerosol Output: 314 µl Output Rate: 170 µl/min

Operating Temperature/Humidity: +10° C to +40° C/10% to 95% RH

Storage and Transport Temperature/Humidity/Air Pressure:  -25° C to +70° C /  10% to 95% RH /  690 - 1060hPa

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WARRANTY CERTIFICATE

Validity 36 months from date of purchase

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IMPORTANT INFORMATION REGARDING ELECTRO MAGNETIC COMPATIBILITY (EMC)

Electronic devices (PC's, mobile phones, televisions, etc) that generate electromagnetic interference can affect the correct functioning of medical devices. This situation can lead to incorrect functioning creating unsafe conditions.

In order to show compliance with the Directives on Electromagnetic Compatibility (EMC), all the tests introduced by the standard EN60601-1-2:2007 have been implemented. This standard defines the levels of immunity to electromagnetic interference as well as maximum levels of electromagnetic emissions for medical devices: **3A HEALTH CARE** products comply with the levels specified in the standard.

However, when the device is in use, it is best practise not to use mobile phones or other electronic devices which generate strong electromagnetic fields near the medical device. This may create dangerous conditions that put the successful outcome of treatment at risk. On the basis of experience and tests conducted, it is recommended that a minimum distance of 7 m be kept: if the distance is less, the correct functioning of the device should be verified.

Documentation regarding the standard EN60601-1-2:2007 is available from **3A HEALTH CARE** (see the address in the instruction manual) and on the website www.3-a.it.

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**3A HEALTH CARE S.r.l.**
Moretti S.p.A. Cavriglia (AR)