MOPEDIA®

MA M100 - M110 - M120 - M130 - M140 - M140 03 B_ITA_03 2021

LETTI DEGENZA

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MOD.

Tulipano 

MADE IN ITALY

ITALIANO

INDICE

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REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices

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2. INTRODUCTION

Thank you for choosing a MOPEDIA hospital bed by Moretti S.p.A. MOPEDIA hospital beds manufactured by Moretti S.p.A. are designed to satisfy all your requirements for a practical, correct and safe use. This manual provides some useful suggestions to correctly use the device you have selected and offers valuable advices for your safety. We recommend carefully reading the entire manual before using the hospital bed you have purchased. Should you have any queries, please contact your retailer for proper advice or assistance.

Note: Check all parts of the product have not been damaged during shipment. In case of damage do not use the product and contact your dealer for further instructions

CAUTION!

- Do not use the following device for purposes other than those described in this manual.
- Moretti S.p.A. disclaims any liability for damage resulting from improper use of the device or use other than that described in this manual.
- We reserve the right to make changes to the device and this manual without notice in order to improve the characteristics.

3. INTENDED USE

MOPEDIA hospital beds are designed for people that have been hospitalized for health reasons in institutes, hospitals and/or health assistance facilities or within a family context.

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of BEDS - TULIPANO complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017.

For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

4.1 Standards and guidelines reference

All safety tests required by UNI CEI 60601-2 -52 were executed with a foam mattress.
Density 30 kg / m³, size L: (3x) 63 tot.195 P: 85 H: 12.5 cm.

5. GENERAL WARNINGS

ATTENTION!

- Always pay close attention to the presence of moving parts that could cause limb entrapment and injury.

- For proper use of the device refer to this manual.
- For proper use of the device consult your doctor or therapist.
- Keep the packaged product away from any source of heat because the packaging is made of cardboard.
- The useful life of the accessory is determined by wear of parts cannot be repaired and / or replaced.
- Be careful when children are nearby.
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

6. SYMBOLS

Product code

Unique Device Identification

CE mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Attention, should read the instructions

Production dates

Waste Electrical & Electronic Equipment (WEEE)

IP54

Protection against dust and water splashes

B Type

Insulation class


Maximum permissible load

Maximun Load

7. GENERAL DESCRIPTION

7.1 View and List of Components

(Fig.1) (Picture with illustrative purpose only)


1. Base frame - Single-Crank (MI100), - Double-Crank (MI110 - MI130)
2. Head frame
3. Foot frame
4. Screws and nuts Kit
5. Crank-Kit MI100 (1 crank pre-installed), MI110 - MI130 (2 cranks pre-installed)
6. Electric kit actuators, control unit, remote control and power supply MI120 - MI140 (pre-installed)

8. ASSEMBLING/DISASSEMBLING


8.1 Assembling

Note: Illustrations for the installation and adjustments refers to the hospital bed 2 cranks. Also assembling and adjusting operations are valid for the single crank bed hospital.

- Remove the product from its packaging and check that contents comply with the illustrated component description.
- Hold the head frame and the foot frame while directly inserting the base frame bayonet hooks into the relative slots (Fig.2). (M1130 - M1140 head frame and foot frame are equal).

(Fig.2) (Picture with illustrative purpose only)


(Fig.3) (Picture with illustrative purpose only)

- Secure the head frame and foot frame using the 4 screws and nuts (Fig 3)
- Your bed is fully assembled and ready to use

8.2 Electrical assembling kit (only MI120 - MI140)

Once get like the Fig. 4

- A Remote Controller Timotion Pc.1
- B Power supply Timotion Pc.1
- C Control unit Timotion Pc.1
- D/E Actuator Timotion Pc.1 Pz.2

(Fig.4) (Picture with illustrative purpose only)

(Fig.5) (Picture with illustrative purpose only)

- Plug the end connector with the electric actuator, paying attention to the tooth position, then secure the cables with the plastic cap. (Figure n.5)
- Repeat this operation for the other components

8.3 Disassembling

- Follow the procedure described above in reverse.

9. BEFORE TO USE

- Always check the status of the mechanical parts in order to guarantee they are used in total safety for persons and environment.
- Check the proper assembly of the device with attention to the correct insertion and locking screws and nuts of the headrest and footrest.
- Check the correct raising/lowering of foot rest frame and head frame.

10. PRECAUTIONS FOR USE

- Do not exceed the maximum load capacity for any reason.
- When lowering the head or foot section, make sure that nobody inserts hands, legs or other body parts under the bed base.
- Make sure that the floor surface supporting the bed is uniformly level.
- Before any operation of cleaning, maintenance, or simply move the bed, disconnect the plug of the mains voltage.
- Do not use multiple plugs or extension cables.
- Do not connect to the bed, some electrical devices not due by the normal operation.

11. CHOICE OF MATTRESS


Use mattresses conforming to requirements of reference norms, with size is not higher or lower than those of the bed. All mattresses Levitas range by MORETTI are suitable for TULIPANO beds.

Note: The usage of incompatible mattresses can be dangerous. Refer to the instructions to use.


12. HOW TO USE

12.1 Raising and lowering the head section (MI100, MI110 e MI130)

- Pull the crank outwards
- Turn the handle clockwise to raise the head section (Fig. 6), and anti-clockwise to lower it (Fig. 7)

(Fig.6) (Picture with illustrative purpose only)


(Fig.7) (Picture with illustrative purpose only)

- Push the crank inwards after adjusting the bed to avoid accidental knocks.


Note: As the head of the bed nears the end of its range of motion (both when raising or lowering), the crank will become harder to turn. Do not turn the crank further as this may damage mechanical components!

12.2 Raising and lowering the foot section (MI110 e MI130)

- Pull the crank outwards
- Turn the handle clockwise to raise the foot section (Fig. 8), and anti-clockwise to lower it (Fig. 9)

(Fig.8) (Picture with illustrative purpose only)


(Fig.9) (Picture with illustrative purpose only)

- Push the crank inwards after adjusting the bed to avoid accidental knocks

Note: As the foot of the bed nears the end of its range of motion (both when raising or lowering), the crank will become harder to turn. Do not turn the crank further as this may damage mechanical components!

12.3 Use of remote control (MI120- MI140)

Using the remote control (Fig. 10) act on the buttons to move the headrest and legrest

1. Button 1 to raise the headrest;
2. Button 2 to lower the headrest;
3. Button 3 to raise the footrest;
4. Button 4 to lower the footrest;
5. Button 5 to lift simultaneously headrest / foot rest;
6. Button 6 to lower simultaneously headrest / foot rest;

(Fig.10) (Picture with illustrative purpose only)

Nota: In point 7 shows the security key on the remote control. If you want to prevent the operation of the bed for safety reasons, we have to remove the key, otherwise it back in to enable the operation just put it back.

13. MAINTENANCE

Accessories MOPEDIA by Moretti at the moment placing on the market are carefully checked and carry the CE mark. For the safety of the patient and the doctor are recommended to check by the manufacturer or by an approved laboratory, at least every two years, the suitability for use of your product. In case of repair should be used only original parts and accessories.

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14. CLEANING AND SANITISATION

14.1 Cleaning

Clean the components exclusively with the aid of a damp cloth and neutral detergents.

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ATTENTION!

- Sanitising agents are corrosive. Avoid contact with the eyes and skin. Wear protective clothing, gloves, protective glasses and safety shoes. In case of accidental contact, rinse the affected part thoroughly using warm water and consult a doctor immediately.

15. CONDITIONS OF DISPOSAL

When disposing of the device never use the normal return systems of municipal solid waste. It is recommended instead to disposal through the joint municipal ecological islands for the expected operations of recycling the materials used.

15.1 warnings to correct disposal of the product in accordance with UE Directive 2012/19/UE:

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

16. SPARE PARTS AND ACCESSORIES

MI150 - Lift up patient pole

MP160 - Lift cowith floor stand

MP170 - Universal bedside rail

MP175 - Telescopic bedside rail

MI180 - Folding bedside rail

MIRO50 - 4 wheels with brakes kit

For spare parts and accessories, refer to the MORETTI main catalogue.



Note: Use Only original accessories by Moretti. Otherwise we decline all responsibility for accidents caused by improper use.

17. TECHNICAL FEATURES

17.1 - Dimension



MI100

ENGLISH

- A - Maximum height of the head frame 104 cm
- B - Maximum height of the foot frame 96 cm
- C - Height of base frame from the ground 50 cm
- D - Base frame Width 90 cm
- E - Max bed width 93 cm
- F - Max base frame length 200 cm
- G - Max bed length 206 cm
- H - Over hang handle extracted 10 cm
- I - Head frame length 65 cm
- L - Base fixed frame Length 123 cm
- M - Maximum angle headrest 72 °

MI110-MI120



- A - Maximum height of the head frame 104 cm
- B - Maximum height of the foot frame 96 cm
- C - Height of base frame from the ground 50 cm
- D - Base frame Width 90 cm
- E - Max bed width 93 cm
- F - Max base frame length 200 cm
- G - Max bed length 206 cm
- H - Over hang handle extracted 10 cm
(only for MI110)
- I - Head frame length 65 cm
- M - Maximum angle headrest 72°
- N - Maximum angle leg rest 54°
- O - Femoral support length 38 cm
- P - Foot frame length 66 cm

ENGLISH

MI130 - MI140

ENGLISH

- A - Maximum height of the head frame 104 cm
- B - Maximum height of the foot frame 104 cm
- C - Height of base frame from the ground 50 cm
- D - Base frame width 120 cm
- E - Max bed width 123 cm
- F - Max base frame length 200 cm
- G - Max bed length 206 cm
- H - Over hang handle extracted 10 cm
(only for MI130)
- I - Head frame length 65 cm
- M - Maximum angle headrest 72 °
- N - Maximum angle leg rest 54 °
- O - Femoral support length 38 cm
- P - Foot frame length 66 cm

17.2 Technical specification

MI100

Base frame:	powder coated rectangular 50x20 mm tubular steel structure
Bed Rete:	electro-welded steel mesh (4 mm thick wire)
Head frame/Foot frame:	powder coated round tubular steel structure with 30 mm diameter
Panels:	made of laminate with natural shade
Backrest inclination:	from 0° to 72° max
Head frame weight:	roughly 5.5 Kg
Foot frame weight:	roughly 4.5 Kg
Base frame weight:	roughly 33 Kg
Tot weight:	roughly 43 Kg
Max. load capacity:	135 Kg

MI110-MI120

Base frame:	powder coated rectangular 50x20 mm tubular steel structure.
Bed Rete:	electro-welded steel mesh (4 mm thick wire).
Head frame/Foot frame:	powder coated round tubular steel structure with 30 mm diameter
Panels:	made of laminate with natural shade.
Backrest inclination:	from 0° to 72°max
Footrest inclination:	from 0° to 54° max
Head frame weight:	roughly 5.5 Kg
Foot frame weight:	roughly 4.5 Kg
Base frame weight:	roughly 35 Kg
Tot weight:	roughly 45 Kg
Max. load capacity:	135 Kg

MI130-MI140

Base frame:	powder coated rectangular 50x20 mm tubular steel structure.
Bed Rete:	electro-welded steel mesh (4 mm thick wire).
Head frame/Foot frame:	powder coated round tubular steel structure with 30 mm diameter.
Panels:	made of laminate with natural shade.
Backrest inclination:	from 0° to 72°max
Footrest inclination:	from 0° to 54° max.
Head frame weight:	roughly 6.5 Kg
Foot frame weight:	roughly 5.5 Kg
Base frame weight:	roughly 37 Kg
Tot weight:	roughly 49 Kg
Max. load capacity:	135 Kg (MI130) - 180 Kg (MI140)

17.3 Electrical kit (only for MI120 e MI140)

Weight:	roughly 5 Kg
Force required to operate the controls on the remote control:	5N
Power supply:	220 vac 50-60 hz
Output voltage from the control unit:	24 Vdc
Class of protection:	IP 54
Noise level:	<45 dB
Insulation class:	Class 2

18. TROUBLE SHOOTING

PROBLEM

The command does not works

LIKELY CAUSES

1. The electrical plug is not connected to the power supply voltage
2. The remote controller is not connected to the control unit
3. Likely defect to the remote controller or the control unit

SOLUTION

1. Connect the plug
2. Check and connect the remote control to the control unit
3. Contact the service center for a more accurate diagnosis of the problem

PROBLEM

The electric head rest and / or leg rest does not work

LIKELY CAUSES

1. The actuator backrest / leg section is not connected to the control unit
2. The actuator or remote controller or control unit may be defective

SOLUTION

1. Check and connect the actuator to the control unit
2. Contact your service representative for a more precise diagnosis of the problem

19. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual. Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use. Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

19.1 Repairing

-Warranty repair-

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be cover from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

-Repair not cover by warranty-

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

-Non defective device-

The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be to be paid by the customer.

20. SPARE PARTS

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

21. EXEMPT CLAUSES

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

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