SKEMA®FORM

MA PM100X - PM120X - PM140X 01 A _ ITA_03-2022

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
Via _____ Località _____

Via Bruxelles, 3 - Meleto 52022 Cavriglia (Arezzo) Tel. +39 055 96 21 11

www.morettispa.com email: info@morettispa.com**MADE IN ITALY**

SKEMA®FORM

MA PM100X - PM120X - PM140X 01 A _ ENG_03-2022

**MULTI-FUNCTION
PROFESSIONAL CHAIR**

INSTRUCTIONS MANUAL

Serie  SKEMA FORM
RUGY
MADE IN ITALY

ENGLISH

INDEX

1. CODES.....	PAG.3
2. INTRODUCTION.....	PAG.3
3. INTENDED USE.....	PAG.3
4. DECLARATION OF CONFORMITY	PAG.3
4.1 Applicable regulations and directives.....	pag. 4
5. GENERAL WARNINGS	PAG.4
6. SYMBOLS	PAG.4
7. GENERAL DESCRIPTION.....	PAG.5
7.1 View and description of main parts model PM100	pag. 5
7.2 View and description of main parts model PM120	pag. 6
7.3 View and description of main parts model PM140	pag. 7
8. ASSEMBLING.....	PAG.8
8.1 Assembling of the medical chair	pag. 8
8.2 Assembly of the PMA200 IV pole accessory.....	pag. 8
8.3 Assembly of the accessory for the PMA210 electric pedal control kit	pag. 8
9. WARNINGS FOR USE	PAG.9
10. MODE OF USE	PAG.9
10.1 Chair use and handling - for all models	pag. 9
10.2 Adjustment of thigh supports (for PM12X models).....	pag. 11
10.3 Adjusting the armrests (for PM14X models).....	pag. 12
10.4 Use of the PMA200 2-hooks IV pole accessory.....	pag. 12
10.5 Accessory use of the PMA210 electric pedal control kit.....	pag. 13
11. MAINTENANCE.....	PAG.13
12. CLEANING AND DISINFECTION	PAG.13
12.1 Cleaning.....	pag. 13
12.2 Disinfezione	pag. 13
13. CONDITIONS OF DISPOSAL	PAG.13
13.1 General conditions of disposal.....	pag. 13
13.2 Correct treatment of electrical parts (Directive 2012/19/UE):	pag. 14
14. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY	PAG.14
14.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions.....	pag. 14
14.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity	pag. 15
14.3 Recommended separation distances between portable and mobile RF communications equipment and this device	pag. 17
15. SPARE PARTS AND ACCESSORIES.....	PAG.17
16. TECHNICAL FEATURES.....	PAG.18
16.1 Dimension and weight PM100	pag. 18
16.2 Dimension and weight PM120	pag. 19
16.3 Dimension and weight PM140	pag. 20
16.4 Technical specifications.....	pag. 21
17. TROUBLE SHOOTING	PAG.21
18. WARRANTY	PAG.22
19. REPAIRING	PAG.22
19.1 Warranty repair	pag. 22
19.2 Repair not covered by warranty	pag. 22
20. SPARE PARTS.....	PAG.22
21. EXEMPT CLAUSES	PAG.22

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1. CODES

PM100X Standard Ruggy 3-section examination chair with Trendelenburg

PM120X Gynecological Ruggy 3-section examination chair with Trendelenburg

PM140X Blood sampling and infusion Ruggy 3-section chair with Trendelenburg

X: Indicates the color of the upholstery

2. INTRODUCTION

Thank you for choosing a Ruggy professional multifunctional chair belonging to the SKEMA line by Moretti. The examination chairs have been designed and built to satisfy all your needs for practical and safe use. This manual contains some suggestions for proper use of the device you have chosen and valuable advice for your safety.

It is advisable to carefully read the entirety of this manual before using the chair for medical examination. In case of doubts please contact your dealer, who will help and advise you properly.

3. INTENDED USE

Ruggy by Moretti professional multifunctional medical chairs are intended to facilitate the positioning of the patient in the most suitable situation during medical examinations or medical-therapeutic treatments in general, during visits or gynecological examinations, during operations of blood sampling and infusion.

WARNING!

- Do not use the product for a purpose not indicated in this manual
- Moretti S.p.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product
- The manufacturer reserves the right to change the information contained in this document without previous notice

4. DECLARATION OF CONFORMITY

Moretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of RUGY Multi-function professional chairs complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017. For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

Note: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any

references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

4.1 Applicable regulations and directives

In order to satisfy safety standards for users, Moretti S.p.A. complies with the following standards:
 - UNI CEI EN 60601-2-52:2016 Electro-medical equipment. Part 2-52. Special requirements for the basic safety and essential performance of medical beds
 - UNI EN 12182:2012 Products intended for the assistance of people with disabilities

5. GENERAL WARNINGS

For a correct use of the product, please refer to the present manual

- Keep the packed examination bed away from heat sources
- **SERVICE LIFE**- the examination bed service life is determined by wear of parts not repairable and/or replaceable and in any case not exceed 10 years
- **ALWAYS** pay close attention to the presence of moving parts that could cause entrapment limbs and injuries
- **DO NOT** allow children to play on or operate the bed
- The user and/or the patient will have to report any serious accident that have occurred related the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

6. SYMBOLS

Product code

Unique Device Identification

CE mark

Manufacturer

Batch Lot

Read the instruction manual

Medical Device

Conditions of disposal

Production date

Waste Electrical & Electronic Equipment (WEEE) according to CE/19/2012

Protection against dust and water splashes

B Type applied part

Maximum permissible load

Double insulation class

7. GENERAL DESCRIPTION

7.1 View and description of main parts model PM100


Fig.1

1. Base with ABS casing and pivoting castors with locking brake
2. Leg section
3. Seat Section
- 3a. Bushing for tightening padded armrest with locking handwheel
4. Backrest section with padded adjustable headrest
5. Lifting/lowering beam
6. Roll holder with knobs provided (for disposable paper rolls of maximum width 60cm)
7. Padded armrest
8. Remote control

7.2 View and description of main parts model PM120


Fig.2

1. Base with ABS casing and pivoting castors with locking brake
2. Leg section
3. Seat Section
- 3a. Bushing for tightening padded armrest with locking handwheel
4. Backrest section with padded adjustable headrest
5. Lifting/lowering beam
6. Roll holder with knobs provided (for disposable paper rolls of maximum width 60cm)
7. Padded armrest
8. Slide and housing for thigh support
9. Thigh support with joint for 360° rotation
10. Stainless steel liquid collection tray
11. Remote control

7.3 View and description of main parts model PM140


Fig.3

1. Base with ABS casing and pivoting castors with locking brake
2. Leg section
3. Seat Section
- 3a. Bushing for tightening padded armrest with locking handwheel
4. Backrest section with padded adjustable headrest
5. Lifting/lowering beam
6. Roll holder with knobs provided (for disposable paper rolls of maximum width 60cm)
7. Armrest for blood sampling with joint for 360° rotation
8. Remote control

8. ASSEMBLING

8.1 Assembling of the medical chair

- Open the packing box of the bed and check that the content corresponds to the model you requested and that the parts have not been damaged during transport. If not, contact your dealer immediately.
- Install the roll holder (ref.6 in fig. 1, 2, 3) on the back of the backrest, in the 4 holes provided (2 on each side) using the screws supplied. (Note: for some models the roll holder is already pre-installed).
- For PM100 and PM120 models: install the padded armrests (ref.7 in fig. 1 and 2) with the provided handwheel on the appropriate support bushing (ref.3a in fig.1 and 2).
- For PM100 and PM120 models: install the padded armrests (ref.7 in fig. 1 and 2) with the provided handwheel on the appropriate support bushing (ref.3a in fig.1 and 2).
- For the PM140 model: insert the armrests (ref.7 in fig.3) in the appropriate support bush (ref.3a), adjusting the height using already-installed handwheel. (ref.3a) and adjust the height using the provided handwheel.

8.2 Assembly of the PMA200 IV pole accessory

The IV pole is supplied complete with bracket and fixing details for application to each model of RUGY medical chair. We recommend the application in correspondence of one of the two rear housings.

- Remove the screw in one of the 2 rear attachments of the casing to secure it to the frame
- Place the plastic spacer and the steel bracket in correspondence with the attachment of the casing, aligning all the holes and fixing everything using the screw (with washer) supplied in the accessory set tightening the screw well.
- Fix the lower part of the IV pole to the steel bracket using the screw (with washer) supplied in the accessory set.
- Insert the support with 2 plastic hooks in the upper part of the IV pole

Fig.4


Fig.4a

8.3 Assembly of the accessory for the PMA210 electric pedal control kit

- Install the fixing bracket of the 'female' sockets of the two Y-gears supplied in the kit by removing one of the screws that fixes the base cover to the frame
- Install the bracket using the screw supplied with the kit, tightening it firmly to the base.
- Install the two 'female' sockets of the 2 Y cables supplied in the kit on the top of the special bracket going to use the special screw provided in the kit. Overrun using the screw supplied, the two female sockets and the appropriate hole in the bracket and complete tightening with the nut provided.
- Disconnect the remote control electrical connector
- Connect the 'female' cable that goes into the control unit to the 'male' connector of one of the two Y cables installed; then connect the 'male' connector of the unused Y cable to one of the two 'female' sockets of the newly connected Y-cable. Then connect the connector of the

remote control to one of the inputs of one of the two 'female' sockets just fixed. Then connect the connectors of the 2 pedal controls to the other two free female sockets. See fig.5 for a suggestion for the installation.

Fig.5

9. WARNINGS FOR USE

- Before making any adjustments to the chair, make sure it is in position stable on the legs
- Before carrying out any cleaning, maintenance or simply moving the chair disconnect the mains voltage plug
- Do not use additional mobile multiple sockets or extension cables
- Do not connect any electrical device to the chair that is not intended for normal use of the chair
- Do not move the chair with patients on board
- Before seating the patient, check that the brakes are applied and the castors locked

ATTENTION!

Do not sit with your weight resting entirely on the back section or section leg as the maximum load declared for the bed is to be understood as "uniformly distributed on its support surface"; otherwise the frame structure could be seriously damaged. Moretti Spa declines any responsibility on damage resulting from improper use and / or use different from that indicated in this manual of use.

10. MODE OF USE

10.1 Chair use and handling - for all models

When used with patients, the chair must always be locked by activating the brakes of the 4 wheels provided by means of the appropriate lever. When it is necessary to move it or When it needs to be moved or handled, the brakes on each wheel can be released (see fig.6)


Fig.6

Before each use connect the electrical power cord to the 220 vac power outlet. Use the supplied remote control to adjust the height of the bed and the inclination of the back and leg sections as well as the trendelenburg angle. (see figs. 7 e 8).

List of buttons and remote control functions:

Fig.7

Button and Function	Description of functionality
A - Elevation of the armchair	Adjust the height of the chair by raising it
B - Lowering of the chair	Adjust the height of the chair by lowering it
C - Backrest section lifting	Adjusts the position of the backrest section lifting it
D - Back section Lowering	Adjusts the position of the back section lowering it
E - Leg section lifting	Adjusts the position of the upper section lifting it
F - Leg section lowering	Adjusts the position of the upper section lowering it
G - CPR button	Brings the chair to a minimum height from ground, with the backrest completely lowered and the leg totally raised, 0 ° trendelenburg tilt (the chair forms a horizontal surface)
H - 'Chair' function	Brings the chair to the minimum height from the ground, with the backrest completely raised (position 85 °), the leg totally lowered, 0 ° trendelenburg tilt
I - Restore of horizontal position of sitting floor	Restore the chair top to the horizontal position starting from a condition of inclination, bringing it back to the horizontal position (Trendelenburg 0°)
L - Seat surface inclination (Trendelenburg function)	Bring the chair to be inclined, towards the Trendelenburg position by adjusting the inclination rate
M - Safety nut	Safety key, without which the operation of the remote control and therefore of the chair is not allowed. During the operation it must always be inserted.


Fig.8
(picture for illustration purposes)

10.2 Adjustment of thigh supports (for PM12X models)


Fig.9
(photo for illustration purposes)

Movement A: adjustment in rotation and orientation of the PU thigh support, going to loosen the knob below it and tighten it in correspondence with the orientation desired

Movement B: raising or lowering the thigh support, loosening the knob of the stem adjustment bush and tightening it at the desired height

Movement C: adjustment of the depth of the thigh support, by pulling or pushing the thigh support itself horizontally moving through the slide

10.3 Adjusting the armrests (for PM14X models)


Fig.10
(picture for illustrative purposes)

Movement A: adjustment in rotation and orientation of the armrest support in PU, going to loosen the knob below it and tighten it in correspondence with the orientation desired

Movement B: raising or lowering the armrest support, loosening the knob of the stem adjustment bush and tightening it at the height desired

10.4 Use of the PMA200 2-hooks IV pole accessory

Adjust the height by unlocking the rod by turning the appropriate handle counterclockwise; once the height has been adjusted, tighten the handle clockwise to lock the rod.

Attention: max. 2.5 kg for each hook.


Fig.11

10.5 Accessory use of the PMA210 electric pedal control kit

Fig.12

Fig.13

Fig. 12 - Electric pedal control to activate the chair head lift

Fig. 13 - Electric pedal control to activate the elevation and lowering of the chair

The two models of electric foot controls differ only in the function they activate

1. To lift, press the arrow at the top (of the headboard or of the entire chair) with the foot
2. To lower, press the down arrow (on the head or the whole armchair) with the foot

11. MAINTENANCE

The SKEMA devices by MORETTI are checked carefully and supplied with EC mark, once launched on the market. For the safety of the patient and the doctor we recommend to check the suitability of your product at least once a year.

For periodic inspection we refer to the check of the following parts:

- The supporting structure of the armchair, of the leg and backrest movement mechanisms and of the lifting mechanism with the relative attachments, controls, safety devices
- Check installation , integrity and operation of the castors and their brake
- Check the integrity of cables and plugs
- Check all welding points

If a repair is needed, please use only approved parts and accessories. If damage is found during inspection or maintenance, immediately put it out service the product until repaired or replaced

12. CLEANING AND DISINFECTION

12.1 Cleaning

Use a damp cloth and mild soap. Then dry thoroughly before use Use water at a temperature not exceeding 30° C. Do not to use washing machines jet of water and steam.

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If you need to disinfect the device use a common disinfectant cleaner.

NOTE: Never use acids, alkalis or solvents such as acetone

13. CONDITIONS OF DISPOSAL

13.1 General conditions of disposal

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the beds in the appropriate disposal areas for recycling.

13.2 Correct treatment of electrical parts (Directive 2012/19/UE):

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

14. DECLARATION OF ELECTROMAGNETIC COMPATIBILITY

14.1 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<p>⚠ Warning:</p> <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

14.2 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz - 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz - 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $= \sqrt{P} \begin{matrix} 150\text{kHz to } 80\text{MHz} \\ = 0.6\sqrt{P} \end{matrix}$ $= 1.2\sqrt{P} \begin{matrix} 800\text{MHz to } 2.7\text{GHz} \end{matrix}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency ranged. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz 385-6000 MHz, 9-28V/m, 80% AM(1kHz) pulse mode and other modulation	10V/m	

NOTE 1: UT is the a.c. mains voltage prior to the application of the test level
 NOTE2: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

14.3 Recommended separation distances between portable and mobile RF communications equipment and this device

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz = \sqrt{P}	80 MHz to 800 MHz = $0.6\sqrt{P}$	800 MHz to 2,7 GHz d = $1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters specified for a maximum output not listed above, the recommended separation distance "d" in meters (m) can be calculated using the equation applicable to the frequency of the transmitter, where "P" is the maximum rated power of transmitter output in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz if the high frequency range is applied

NOTE 2: Guideline question for not applying in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

15. SPARE PARTS AND ACCESSORIES


For spare parts and accessories refer to Moretti main catalogue

- PMA200** Aluminum IV pole kit with 2 hooks (max load 2,5 Kg per hook)
PMA210 Electric pedal control kit for chair elevation and backrest movement

16. TECHNICAL FEATURES


16.1 Dimension and weight PM100

ENGLISH

PM100 external dimensions:	mm 1800x880
Minimum height:	mm 520
Maximum height:	mm 920
Maximum head lift angle:	85°
Maximum Trendelenburg angle:	15°
Maximum leg rest angle:	107°

 Pack size PM100	 L: 1160 mm	 W: 790 mm	 H: 1370 mm
 Pack weight PM100	92 kg		
 Chair weight PM100	80 kg		
 Max weight	200 kg		


16.2 Dimension and weight PM120

ENGLISH

PM120 external dimensions:: mm 1800x880
 Minimum height: mm 520
 Maximum height mm 920
 Maximum head lift angle: 85°
 Maximum Trendelenburg angle: 15°
 Maximum leg rest angle: 107°

 Pack size PM120	 L: 1160 mm	 W: 790 mm	 H: 1370 mm
 Pack weight PM120	107 kg		
 Chair weight PM120	95 kg		
 Max weight	200 kg		

16.3 Dimension and weight PM140

ENGLISH

External dimensions PM140: mm 1800x785
 Minimum height: mm 520
 Maximum height: mm 920
 Maximum head lift angle: 85°
 Maximum Trendelenburg angle: 15°
 Maximum leg rest angle: 107°

 Pack size PM140	 L: 1160 mm	 W: 790 mm	 H: 1370 mm
 Pack weight PM140	97 kg		
 Chair weight PM140	85 kg		
 Max weight	200 kg		

16.4 Technical specifications

Required force to operate the controls on the remote control	5 N
Required force to push the elevation pedal (MI39X)	30N
Input	100-240V AC 50/60Hz 145VA
Output	24 VDC
Protection Class	IP 54
Noise level	< 45 DB
Insulation Class	2 Class
Used materials for bed mainframe	Varnished steel frame, foam bed cover in PVC

17. TROUBLE SHOOTING

SYMPTOMS	POSSIBLE CAUSES	SOLUTION
No command working	<ol style="list-style-type: none"> 1. The power plug is not connected to the mains voltage socket 2. The remote control is not connected to the control unit 3. Probable defect in the remote control or control unit 	<ol style="list-style-type: none"> 1. Connect the power plug 2. Check and if necessary connect the remote control to the control unit 3. Contact the service center for a more precise diagnosis of the problem
Backrest, legrest and trendelenburg electric movements do not work	<ol style="list-style-type: none"> 1. The backrest, legrest or trendelenburg actuators are not connected to the control unit 2. The actuator or remote control or control unit may be faulty 	<ol style="list-style-type: none"> 1. Check and if necessary connect the actuators to the control unit 2. Contact the service center for a more precise diagnosis of the problem
The elevation of the chair does not work	<ol style="list-style-type: none"> 1. The elevation actuator is not connected to the control unit 2. The actuator or remote control or control unit may be faulty 	<ol style="list-style-type: none"> 1. Check and if necessary connect the actuator to the control unit 2. Contact the service center for a more precise diagnosis of the problem
The electric foot controls (pedals) does not work	<ol style="list-style-type: none"> 1. The foot control does not connected to the Y cable 2. Y cable is not connected to the control unit 3. The electronic control unit and / or the actuators may be defective 	<ol style="list-style-type: none"> 1. Check and if necessary connect the foot control to the Y cable 2. Connect the Y-cable to the controller 3. Contact the service center for a more precise diagnosis of the problem

18. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual.

Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use.

Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual.

The warranty doesn't cover components subject to wear and tear during the correct use of the device.

19. REPAIRING

19.1 Warranty repair

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty

19.2 Repair not covered by warranty

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received

19.3 Non-defective devices

The customer will be informed if, after the device return and examination, Moretti declares that the device is not faulty. In this case the good will be sent back to customer, the shipping charge will be paid by the customer

20. SPARE PARTS

The original Moretti spare parts are guaranteed for 6 months from the receiving date.

21. EXEMPT CLAUSES

Moretti does not offer any other declarations, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage.


The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

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