LEM

laserplus

MORETTI

MA271c_IT

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CE
0051

SOMMARIO

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TREATMENT	36
FREE PROGRAM	37
MAINTENANCE.....	42
PREVENTIVE MAINTENANCE	42
EQUIPMENT CLEANING	42
LASER HANDPIECES CLEANING AND DISINFECTION	42
CORRECTIVE MAINTENANCE	43
TROUBLESHOOTING GUIDE	43
WARRANTY	44
EMC DECLARATIONS	46
EMISSIONS	46
IMMUNITY	46

IMPORTANT

No part of this document may be photocopied, reproduced or translated into another language without the written consent of the manufacturer.

If you need Technical Assistance, please contact your dealer.

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www.led.it

WARNINGS AND PRECAUTIONS

THESE WARNINGS HAVE BEEN PREPARED FOR YOUR OWN SAFETY AND THAT OF OTHERS, SO PLEASE READ THEM CAREFULLY BEFORE INSTALLING AND USING THE EQUIPMENT.

- After removing the packaging make sure of the integrity of the equipment in case of doubt do not use and contact professionally qualified personnel.
- The electrical safety of this equipment is ensured only when it is correctly connected to a system that complies with current electrical safety regulations. It is necessary to verify this fundamental safety requirement and, in case of doubt, request a thorough control of the system by qualified personnel.
- Before connecting the equipment, make sure that the license plate data (on the rear panel) corresponds to those of the electrical distribution network.
- In case of incompatibility between the socket and the power cord of the equipment, replace it with another suitable type. In general, the use of adapters, power strips and/or extension cords is not recommended. If their use is essential, it is necessary to use only simple or multiple adapters and extensions that comply with current safety standards.
- The use of any electrical equipment involves the observance of some basic rules. Especially:
 - Do not touch the equipment with wet or damp hands or feet
 - Do not use the equipment barefoot
- Do not leave the equipment exposed to atmospheric agents (rain, sun, etc.)
- Do not leave the equipment unnecessarily inserted. Turn off the main switch of the system when it is not in use.
- This equipment shall be intended only for the use for which it was expressly designed. Any other use is to be considered improper and

therefore dangerous. The manufacturer cannot be held responsible for any damage resulting from improper, erroneous and unreasonable use.

- It is dangerous to modify or attempt to modify the characteristics of this equipment.
- Before carrying out any cleaning or maintenance operation, disconnect the equipment from the power supply network, either by unplugging the plug, or by turning off the main switch of the system.
- In case of failure and/or malfunction of the equipment, turn it off. For any repair, contact only an authorized technical assistance center and request the use of original spare parts. Failure to comply with the above may compromise the safety of the equipment and the user.
- It is advisable to carry out the treatments after medical consultation.
- The equipment has been designed to meet current requirements for electromagnetic compatibility. If there is a suspicion that the operation of the equipment is interfering or interfering with the normal operation of other electrical and electronic equipment, it is recommended to power the equipment with a different electrical outlet and / or try to arrange the equipment differently until the interference ceases.
- Avoid using mobile phones near the equipment to avoid possible interference with it.
- Many of the components used within this equipment are sensitive to electrostatic discharge. In the need to handle any electronic component, eliminate any static electricity that you could be carrying, touching a ground plane to discharge it to the ground. If possible, wear a bracelet attached to the ground. Failure to observe these precautions may result in permanent damage to the electronic components of the equipment.
- To avoid damage to the eyes (retinal detachment) the use of goggles is mandatory.
- The beam (which is invisible) must never be directed at the eye.

- Before each use, check the status of the equipment and cables to identify uncovered conductive parts or other electrical hazard conditions.
- The equipment is not suitable for use in explosive environments or environments saturated with oxygen and/or gas.
- Do not leave the appliance unattended with the presence of children and / or incapacitated people.
- Use the contact treatment technique (tip of the handpiece in contact with the part/area to be treated).
- Users of this product should be instructed about potential eye risks.
- Do not look at the target area without eye protection (glasses) during treatment.
- Any serious accident occurring in relation to the device must be reported to LED SpA (via Selciatella n.40, 04011 Aprilia (LT) ITALY) and to the competent authority:
Ministero della salute – Direzione generale dei dispositivi medici e del servizio farmaceutico
Viale Giorgio Ribotta, 5 – Roma
E-mail: segr.dgfdm@sanita.it
Tel.: +39 06 5994 3199 / +39 06 5994 3207
- The equipment, periodically (at least once a year) must be subjected to control by qualified personnel for the verification of the following electrical safety values:
 - a measurement of leakage currents
 - a measure of the resistance between the ground clamp and each accessible conductive part for which grounding is planned.
- When used in a residential environment, this equipment may not offer adequate protection for radio frequency communication services. You may need to take mitigation measures, such as transferring or reorienting equipment.

GENERAL DESCRIPTION

Laserplus is a control unit for monochromatic laser light emitters emitting in the near infrared ($\lambda = 905 \text{ nm}$), suitable for anti-inflammatory, regenerative and analgesic biostimulant treatment.

The **laserplus** device is intended for temporary use in outpatient settings, and is indicated with the monodiode handpiece (limited irradiation area), for the following therapeutic treatments:

- inflammatory pathologies affecting tendons and soft tissues (tendinitis, bursitis, enthesitis);
- insertional pathologies;
- superficial articular pain.

While with the tridiode handpiece (wider irradiation area), to perform the following therapeutic treatments:

- adjuvant in the treatment of pressure ulcers and sores;
- scarring forms and in the treatment of edema.

The effects of diffuse low-energy laser irradiation can be summarized below:

- Increased ATP production by stimulation of cellular mitochondria with consequent increase in cellular metabolism and RNA and DNA synthesis.
- Acceleration of protoplasmic cellular electrolyte turnover with consequent increase in metabolism and rates of cell mitosis.
- Stimulation of fibroblasts with consequent increase in the rate of formation of elastic fibers and collagen and therefore of connective tissue.
- Modification of intracapillary hydrostatic pressure with consequent greater absorption of interstitial fluids, reduction of edema and activation of tissue turnover.
- Neof ormation of vessels with consequent activation of vascularization.
- Stimulation of algotropic endings with consequent increase in the threshold of perception and production of β endorphin.

- Stimulation of immune systems with increased production of antibodies.

The laserplus device is fully managed and controlled by a microprocessor and has 1 output channel.

The device provides operation with parameter selection or with the use of programs stored inside. The internal programs automatically adapt according to the power of the connected handpiece.

The device has a safety control system (which can be activated) with automatic deactivation of laser irradiation in case of failure or incorrect positioning of the emitter handpiece with the part to be treated.

The device, supplied with a 30W monodiode handpiece, is able to drive different handpieces for power and number of laser emitters, for more information contact your dealer.

It is a device dispensing impulse bipolar square current, suitable for toning and rehabilitating muscle stimulation and for the treatment of pain syndromes acting as an analgesic stimulator.

The advantages of using bipolar impulse current, compared to monopolar currents, consist mainly in the elimination of unwanted electrolytic effects, thanks to the symmetrical bipolarity of the stimulating wave and therefore reduction of the possibility of burns; good tolerability.

The elestim 2 device is fully managed and controlled by a microprocessor and has 2 output channels.


The device provides operation with parameter selection or with the use of 12 programs (4 for beauty treatments – 5 for sports treatments – 3 for rehabilitation treatments) stored inside.

The device is designed for operation with adhesive electrodes that allow easier and faster arrangement on the points to be treated.

LEM MEANING

LEM is the brand owned by MORETTI SpA dedicated to the line of electro-medical devices for professional and home use in the healthcare environment.

BLOCK DIAGRAM

INTENDEDE USE

Active medical device and accessories for temporary and professional use, designed and manufactured for laser therapy applications in biostimulant, regenerative, analgesic, anti-inflammatory, and/or neuralgia and/or aesthetic blemishes therapeutic treatments.

INTENDED USER

Device for professional use: Medical degree with knowledge of technology applications and/or physiotherapy degree.

INTENDED PATIENT POPULATION

Adults- Male and Female (≥ 18 Years), with the exception of patients listed in the contraindications section.

LASER BIOSTIMULATION

The laser is an emitter of luminous electromagnetic radiation characterized by the following properties:

- Monochromicity = Emission of a single wavelength
- Phase coherence = Phase emission of radiation
- Directionality = Small angle of divergence of the emitted beam
- Brilliance = Emission with high energy density

In practice, the only characteristic that is common to all types of lasers, is the monochromicity that can be referred to one or more lines emitted.

The directionality characteristic is not present in diode lasers which have rather high divergence angles.

The laser emission wavelength is linked to the material used and the light outputs are always limited to a few percent of the energy supplied for the excitation of the material.

Some diode lasers can emit continuously with low power, others can emit only pulses but of high power because otherwise the energy dissipated by the system would result in the rapid destruction of the element. Continuous lasers (CW = Continuous Wave) can be operated in pulses with appropriate excitation control.

The emission wavelength is measured in Angstrom ($1\text{\AA} = 10^{-10}$) or in nanometres ($1\text{nm} = 10^{-9}$) 1 micron ($1\mu = 10^{-6}$ m).

In the medical therapeutic field the most used unit of measurement is the nanometres and therefore we will use this unit of measurement later.

History

The first observations of the biostimulant properties of laser radiation are due to TOMBERG who, already applying in 1961 a ruby laser ($\lambda = 694.3$ nm), observed on small irradiated white rats an increase in red blood cells, hemoglobin and platelets together with moderate leukocytosis with an increase in lymphocytes and a stimulus to the increase of bone marrow.

In 1966 KLEIN observed changes in haemagglutination in pigmented cells irradiated with ruby lasers. In 1969 MESTER presented a small statistic on the positive effects

of laser irradiation on the repetition of traumatic wounds and hair regrowth. At the same time, the results of experiments carried out in agronomy by KARLANDER and KRAUS were reported, which highlighted the difference in the effectiveness of coherent light and by RUBIN which highlighted the different effect of radiation from red laser radiation compared to the effect of radiation from ultraviolet laser radiation ($\lambda = 250 \text{ nm}$) which caused cell death.

In 1972 SCHUR reported positive results in the treatment of wounds with He-Ne laser.

Bioresonance Hypotheses

In 1923 the Russian Alexander GURVICH, working on cell cultivation, noticed that at the beginning of the cell cleavage in a vessel there was a signal that was transmitted to nearby cultures which in turn began the splitting, and that this signal was interrupted by a glass plate but not by a quartz plate. In 1973 MESTER reported its results in the study of increased collagen production in laser-irradiated wounds.

CHERKYROV reported positive results in the regeneration and formation of bone welding callus of fractures.

DURMANOV through irradiation with He-Ne documented a reduction of 40 to 50% of the healing time of burns.

SEMENOV described an experiment with 68 rats documenting hormonal changes in the adrenal glands with weight gain from 20 to 40% and reduction of lipid content in the adrenal cortex.

This change in internal adrenal secretion would justify the feeling of dizziness or drowsiness reported by some patients.

Profs. RACHISCHEN and TSOY observed the considerable acceleration in regeneration in the injured sciatic nerve with the improvement of electrophysiological functional parameters.

BENEDICENTI, through experiments with 904 nm diode laser, observed the increase of lymphatic drainage in the rat mesentery and the increase of

endocellular ATP in lymphocytes. Still Benedicenti reported the acceleration of the binary split in cultures of paramecium Aurelia.

JUDITH WALKER recorded evoked potentials in the human nervous system by irradiation of the ulnar nerve with He-Ne laser (1 mW) and detected marks near the shoulder with a latency time of 4 ms, demonstrating with this that laser can reduce pain through nerve stimulation similar to electrostimulation.

BENEDICENTI demonstrated increased circulating endocrine β in patients with trigeminal neuralgia undergoing diode laser treatment (904 nm).

PASSARIELLO and collaborators in 1983 confirmed the increase in circulating endocrine β in patients suffering from scapulohumeral peri-arthritis, rheumatic arthritis and lumbosciatica, subjected to laser treatment.

BALDONI and collaborators observed lively proliferated activity in the epithelial component of the mucosa in fragments of chicken embryonic intestine and better preservation than controls.

GIUBBOTTI deduced the increase in the transformation of prostaglandin PGC and PGH into prostacyclin PGI with consequent rebalancing of osmotic pressure, increase in blood circulation, platelet disaggregation, anti-edematous and anti-inflammatory action with pain reduction.

PASSARELLA, CATALANO and collaborators studying the effect of laser irradiation on rat liver mitochondria noted effects on the integrity and permeability of the mitochondrial membrane, on the metabolic processes involving carrier substrate transport, on the activity of mitochondrial enzymes and on the energy status of mitochondria.

GURVICH's observations were confirmed in 1940 by photomultiplier detections, and the Russian Prof. KAZNACHEV through his studies, concluded that contagious viral processes would be transmitted precisely by this message and not by the vermiculation of the virus.

The phenomenon of biological induction is an integral part of the biplane theory; according to the Russian INJUSHIN, the bioplasma obtains energy from the air that brings free electrons and how many and the conditions predisposing the concentration of diseases would be due to the decrease in bioplasmic energy.

According to the Polish Prof. SEDLAK the reason for the possibility of regulation of biological processes by laser irradiation lies in the influence of the same on electric and electromagnetic fields, and according to INIUSHIN the laser acts on the biological plasma by resonance effect restoring and strengthening the energy state. Prof. POPP affirms the generation by cells of electromagnetic waves similar to laser light and the existence of a complementary system to the cellular chemical component of light and acoustic waves that the cellular system needs for its internal and external regulation. The exchange of information would take place in the red and infrared spectrum and could be restored, when deficient, through laser irradiation.

According to Prof. POPP all plant and animal cells emit light radiation with a maximum of around 550 nm and a minimum of 350 nm. The latter wavelength would have a particular significance, and that would be the reason, why ultraviolet radiation of wavelengths of less than 350 nm is potentially carcinogenic.

PAIN CONTROL

Pain control theory using the Gate System by R. Melzack and P. Wall

The nerve impulses evoked by skin stimuli are transmitted through the fibers called A delta and C of small diameter and slow and through the vast myelin afferents Beta.

The impulses reach three systems in the spinal cord:

- cells of the gelatinous substance;
- the fibers of the dorsal column that project towards the brain;
- 1. the primary T transmission cells in the dorsal horn.

According to Gate's theory, the gelatinous substance functions as a gate control system that modulates the afferent beam in the dorsal spine system acts, at least in part, as a central trigger control that activates selective brain processes that affect the modulation properties of the gate control system; and T-cells activate neural mechanisms that include the system of action responsible for pain perception.

Painful phenomena are determined by the interaction of the above three systems. Gate's theory proposes that the gelatinous substance functions as a gate control system that modulates the synaptic transmission of nerve impulses from peripheral fibers to central cells.

The large diameter Beta fibers initially activate the T cells but later their effect is attenuated by a negative feedback system while the small diameter fibers A delta and C activate a positive reaction mechanism that elevates the effect of the impulses that arrive.

These reaction effects are radiated by the cells into the gelatinous substance.

Gate's theory indicates three characteristics as relevant to pain:

- the activity preceding the stimulus;
- the activity evoked by the stimulus;
- the ratio of the activity of large diameter fibers Beta to small diameter A delta and C.

In the absence of particular stimuli, nerve activity is carried by myelin and amyelin fibers of small diameter A delta C that adapt slowly and keep the gate open. When

a stimulus is applied, information is transmitted to the brain. Since in the absence of stimulus variations the Beta fibers are practically inactive, the stimulus produces greater increase in the beta fibers in relation to the A delta and C fibers, therefore the impulses of the Beta fibers in addition to triggering the T cells tend to close the presynaptic gate. If the stimulus is increased, the trigger frequency increases and therefore the positive and negative reactions compete resulting in a slow increase in the output of T cells. If the stimulus is prolonged the Beta fibers begin to adapt causing an increase in the effect of the A delta and C fibers, therefore the gate opens further and the exit from the Trigger cells increases. If at this point the activity of Beta fibers is artificially increased by counteracting their tendency to adapt, the exit from T cells is reduced. Therefore, the output from the T cells is a function of the number of fibers involved, the frequencies adopted and the relationship between the activity of the Beta type fibers to those of type A delta and C. Through electrical stimulation, or laser, at appropriate frequency, the activity of Beta fibers would be increased which would counteract, for the mechanisms indicated above with the activity of A delta and C fibers would produce a reduction in the output of T cells and therefore a lower sensation of CNS pain (see fig. A and B).


FIG. A


FIG. B

BIOSTIMULATION EFFECTS

The effects found with laser irradiation are to some extent similar to those found with means of physical stimulation of a different nature.

As for classical acupuncture and transcutaneous electrical stimulation, the monochromatic radiation emitted by the laser determines the remission of the algic crisis and the acceleration of the regenerative life cycle.

It has been found that both the stimulation produced by the manipulation of needles in acupuncture practices in the classical points dictated by Chinese texts,

as well as transcutaneous electrical stimulation by means of low-frequency electrical impulses induces, such as laser radiation, the secretion of enkephalins, including endocrine beta, which seems to preside over the control of the pain threshold. Also in common with electrical stimulation is the observation that high-frequency laser stimulation does not seem to cause an increase in circulating Beta endocrine. However, there is an analgesic action that is evidently established by different mechanisms and that persists for a shorter time.

The analgesic action seems to be a function of the frequency of the stimulation impulses both as regards the treatment times necessary to determine it, and the persistence times of algic remission.

It should also be noted that an increase in circulating Beta endorphins and consequent algic remission has been observed in some patients undergoing simulated treatment.

However, it is statistically established that the component of the placebo effect is not significant to such an extent that it can lead to the suspicion that acupuncture, transcutaneous electrostimulation or laser irradiation owe their effect mainly to placebo action. For the purposes of analgesic action there seems to be no effect or there is only a reduced effect by electrical stimulation with direct current, continuous laser irradiation or acupuncture without manipulation, while the regenerative effects seem to manifest themselves also with continuous stimulation and regardless of the frequency of stimulation used.

It is true that analgesic effects have been detected with continuous He-Ne laser irradiation with scanning equipment. It should be noted that, as a consequence of scanning, the continuous emission actually acts as pulsed emission as in successive moments the radiated area varies. Therefore, if a single area is considered, it will be irradiated with a frequency depending on the scanning speed and the area affected by the beam.

The same consideration applies to irradiation by high-frequency pulsed lasers with focused beam and destruction of the treated area.

However, it should be noted that to obtain analgesic effect with low frequencies of stimulation it is necessary to carry out the treatment for longer times.

Observations were also made through which it was possible to establish that the analgesic effect produced seems to be a function of the energy density delivered. In addition, it has been found that laser stimulation focused on classical acupuncture points, determines rapid onset of algic remission but the persistence of the state of analgesia is very reduced compared to the effect obtainable with stimulation spread directly in the algic areas or in the trigger points.

In summary, the decisive parameters to obtain analgesic action with laser emission seem to be:

1. Possibility of penetration into the tissue and therefore need for appropriate wavelengths of radiation and sufficient power.
2. Frequency of pulses emitted.
3. Choice of irradiated zone.

PARAMETERS OF LASER RADIATION

To have comparable data on the effect of the administration of laser radiation it is necessary to establish parameters that allow to relate the results achieved. The main parameters are as follows:

- (a) Wavelength ($\lambda = \text{nm}$)
- (b) Power density (radiation) ($D_p = \text{W}/\text{cm}^2$)
- (c) Pulse frequency ($F = \text{ips}$)

Wavelength

A fundamental difference in the method of administration of laser radiation to perform biostimulation, compared to photocoagulation or surgery, consists in the need to irradiate relatively large areas and carry out radiation in depth.

This entails first of all the need to use wavelengths for which there is poor absorption by hemoglobin and water; therefore, wavelengths between 600 nm and 1300 nm (Fig. C).


FIG. C

The emitters most commonly used for biostimulation are pulsed emission diode lasers with wavelengths of 904-905 nm (near infrared) with peak powers between 10 and 100 W and characterized by high emission divergence that allows to irradiate a sufficiently large area at the same time.

Power Density

From the power density it will be possible to have an evaluation of the penetration of the radiation into the tissue in sufficient quantity to obtain a biological stimulation effect.

The level of radiation penetration is difficult to establish because some of the radiant energy sent is reflected, and the energy penetrating the tissue is partly deflected due to the inhomogeneity of the tissue through.

To reduce the amount of reflected energy it is advisable, in addition to keeping the handpiece in contact, to clean the surface to be treated and direct the laser beam as much as possible in a direction orthogonal to the surface to be irradiated.


Taking into account the absorption curves of the most important components of human tissue: water and hemoglobin, we can calculate, at wavelength $\mu = 905$ nm, an absorption of about 15% for each millimeter of depth.

The power available in depth can be calculated with approximation by the relation:

$$P_d = P_o(1-A)^d$$

Where P_0 is the power density penetrated to the surface. A is the absorption and d the distance to the surface in millimeters.

To calculate power density (W/cm^2) we need to know the peak power emitted and the size of the area hit by the laser beam. By doing this simple calculation we can get an idea of the possibility of penetration of lasers with different emission powers. In fact, for the purposes of penetration obviously the power density on the radiated surface is important and not the power of the emitter.

For example, if we calculate the power density of a He-Ne laser with an emission power of 10 mW which invests with its beam, characterized by high directionality, an area with a diameter of one millimeter, we obtain an irradiated area equal to:

$$S = \pi D^2/4 = 0.785 \times 10^{-2} \text{ cm}^2$$

And therefore a power density:

$$Dp = P/S = 10^{-2}/0.785 \times 10^{-2} \times 1.27 \text{ W/cm}^2$$

If instead we use a diode laser with continuous emission, which is characterized by a high divergence (about 18) at the distance of about 13 millimeters of the emitter we have an irradiated area of about 1 cm^2 . In this case we would have a

power density equal to 10 mW/cm^2 and therefore a much lower chance of penetration.


If they still use a pulsed diode laser with peak power homogenized through optical fiber $\varnothing 4 \text{ mm}$ equal to 8 W and divergence characteristics equal to those of the previous diode laser, at a distance of 10 mm we will have a radiated surface of 1 cm^2 but the power density will be this time equal to 8 W/cm^2 and therefore much greater possibility of penetration.

At a wavelength of 905 nm , a power density of 5 W/cm^2 allows biological stimulation effect of up to 30 mm of soft tissue and about 2 mm of bone tissue.

Energy Density

While for emission lasers the value of the energy density expressed in Joule/cm^2 is numerically equal to the value of the power density, for pulsed lasers the pulse duration must be taken into account.

It has been established that biostimulation occurs with pulses of very short duration and some biological effects have been observed with during pulse less than nanosecond ($1 \text{ ns} = 10^{-9} \text{ sec}$).

Pulsed diode lasers normally have a pulse duration of 100 nsec and therefore the energy density for each pulse will be:

$$Dei = Dp \times \xi$$

Where Dp is the power density and ξ is the pulse duration time.

It has been established that biostimulation occurs with pulses of very short duration and some biological effects have been observed with during pulse less than nanosecond ($1ns=10^{-9}$ sec).

$$Dei = 8 \times 10^{-7} = 0,8 \times 10^{-6} \text{ J/cm}^2$$

Or $0,8 \mu\text{Joule/cm}^2$.

It seems appropriate to clarify that the energy of each photon is constant and independent of the emission power being related to the emission wavelength by Planck's constant and the speed of light.

$$E = C_{\lambda}h$$

Where $h=6,62 \times 10^{-34}$ Joule sec (Planck's constant) and $C=2,997 \times 10^8$ m/sec (Speed of light).

From studies carried out by E. Mester derives the opportunity not to exceed the dose of 3 J/cm^2 for each session.

Frequency of the impulses

The choice of pulse frequency depends, as mentioned above, on the main purpose for which the irradiation is carried out.

If the purpose is mainly analgesic, it is necessary to use rather low frequencies in order to obtain an action that persists sufficiently over time. If, on the other hand, the main purpose is to carry out a regenerative biostimulation, the frequency of the pulses may be higher.

In some cases, where both effects are desirable, it may be useful to radiate for a certain time at high frequencies and then progressively reduce the frequency of

the pulses. This technique seems to allow an algic remission that arises quickly and is maintained over time over time and at the same time with good regenerative stimulation.


Knowing the energy density per pulse and the pulsation frequency you can very simply calculate the energy delivered per unit area in one second.

$$E_s = E_{ei} \times f$$

In the case of a pulsed laser with energy of $1 \mu\text{J}/\text{cm}^2$ per impulse and an impulse frequency of 100 impulses per second, the following is obtained:

$$E_s = 100 \mu\text{J sec cm}^2$$

MEANINGS OF GRAPHIC SYMBOLS


Meaning of symbols:

The meaning of the graphic symbols printed on the equipment is as follows:

- | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ol style="list-style-type: none"> 1. Floating neutral electrode: patient circuit is isolated from earth at both high and low frequencies 2. The equipment is B class, protected against electric shock than that provided BY TYPE B APPLIED PARTS 3. Not Ionizing Radiation emitted. 4. Follow the instructions for use. 5. CE Mark (2017/745/EU) + Notified Number 0051 = IMQ Italy 6. The product mustn't be throwing in the containers for urban wastes, but it must be swallowed by a separate picking. 7. Manufacturer 8. Serial Number 9. Date of manufacture | <ol style="list-style-type: none"> 10. UDI Unique Device Identification 11. Medical Device 12. Distributor 13. Not to be serviced by users 14. REF Reference code number 15. Temperature Limits 16. Humidity Limits 17. Atmospheric pressure limits 18. This way up 19. FRAGILE – Handle with care 20. Keep away from sunlight 21. Keep Dry 22. Stacking limit by number 23. Weight 24. Dimensions 25. Number of pieces 26. Recycle 27. Laser STOP (to be used in case of emergency) 28. Explanatory laser label 29. Laser ready 30. Remote Interlock Connector 31. Model/Trade Name |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

BOX LABEL

With reference ISO15223-1 "Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied", and ISO780 "Packaging – Pictorial marking for handling of goods"
On label box of UNIT's carton box are present these indications:

CONTRAINDICATIONS

It is advisable to carry out the treatments after medical consultation.

The application of this device should be avoided in subjects:


- Wearers of pacemakers or other implanted electronic devices. These subjects should not be treated until the consent of the doctor who implanted the electronic device.
- Suffering from heart disorders.
- With severe imbalances in blood pressure.
- With the presence of suspected or overt carcinoma.
- With marked sensitivity to light.
- With skin lesions.
- With skin hypoesthesia.
- With neoplasms.
- With severe diseases of the nervous system.

- With bacterial infections in place.
- With severe kidney failure.
- Pregnant.
- Haemorrhagic.
- Epileptic.
- Wearers of metal prostheses.
- In a state of unconsciousness.
- Markedly asthenic.
- In infant-juvenile age.

Its application is also contraindicated:

- On tumors to prevent possible metastases.
- On areas of hemorrhage.

UNIT DESCRIPTION

1	Unit power supply socket
2	Power switch
3	Alphanumerical display
4	Control keyboard
4a	UP key to increase
4b	DWN key to decrease
4c	START/STOP key to start and interrupt
4d	ENTER key to confirm
4e	ESC key to go back
5	Warning/indication lights
5a	Output supplying warning light (yellow)
5b	Unit on indication light (green)
6	Output connector
6a	Laser danger plate

CONTROLS

POWER SOCKET UNIT (1)

On the back of the unit there is the power socket. Plug the supplied power cord into this outlet. The cable must be connected to a power socket with electrical characteristics compatible with the unit (see license plate data on the bottom of the unit).

POWER SWITCH (2)

On the left side of the unit is the power switch of the unit. Turning the switch to position I turns the unit on. The power on of the unit is indicated by the alphanumeric display and the power indicator light on the central part of the unit.

ALPHANUMERIC DISPLAY (3)

On the central part of the unit there is an alphanumeric display where all the indication parameters or variables in a given procedure are displayed.

CONTROL KEYBOARD (4)

The control keyboard consists of:

- **4a UP button** to increase the value, to select an option on the previous line of the display or to go to a possible previous page (indication on the display of the arrow at the top "↑").
- **4b DWN button** to decrease a value, to select an option on the next line of the display or to move to a possible next page (indication on the display of the arrow at the bottom "↓").
- **4c START/STOP button** to start/stop power supply.
- **4d ENTER button** to confirm a data on the display.
- **4e ESC button** to return to the previous selection page.

SPIES (5)

Under the keys there are three indicator lights that have the following meaning:

1	Dispensing output indicator lights (yellow).
2	
	Unit indicator light on (green color).

OUTPUT CONNECTOR (6)

In this area there is the connection point of the handpiece. Near the connector there is the laser danger plate (6a).

APPLICATION OF LASER INFORMATION PLATE

Together with the unit an information plate has been provided, which must be applied on the unit, below is an example of application, the user can choose a

position that he considers most suitable considering that the plate should be visible in the normal working position.

1. Remove the protection paper of the nameplate sticker.
2. Apply visibly.

TECHNICAL CHARACTERISTICS

Power Supply	230Vac / 50Hz – 10VA
Internal Power Fuse	T 500mA (5x20)
Wavelength Handpieces	905 nm \pm 10%
Maximum Power	Depends on the Handpiece
Output frequency	from 1 to 10000Hz
Pulse Duration	100 nsec
Output Channels	1
Treatment time	from 1 to 60 min
Output parameters Indicat.	Alphanumeric Display
Output indicator	Visual on the handpiece
Laser Classification	3B
Dimensions (mm)	270 L x 220 P x 110 H
Weight	1.3kg
Electrical Safety Class (EN60601-1)	II BF
Protection class (EN60529)	IP20
Classification MDR 2017/745/EU	IIB
Electromagnetic Environment of Utilization	Class A
Environmental Characteristics of Storage	
- Amb. Temperature	between -10°C and 50°C
- Relative Humidity	between 10 and 85%
- Atmospheric Pressure	between 50kPa and 106 kPa
Environmental Characteristics of Operation	
- Amb. Temperature	between 10 and 40 °C
- Relative Humidity	between 30 and 75%
- Atmospheric Pressure	between 60kPa and 106 kPa

Hardware Requirements

Microcontroller	16 bit
Clock Frequency	16 MHz
Rom	64 KB
Ram	2 KB
Peripherals	UART, I2C, Watch-dog timer
Visual	Display 7-Segments

STANDARD COMPOSITION

REF	Description	Quantity
00100.02	Power Cable	■/1
LEM271MN030	Monodiode handpiece	■/1
40100.99	One size laser glasses	■/2
LEM271PMDX1	Handpiece holder	■/1

LEM271MN030

00100.02

40100.99

LEM271PMDX1

PUTTING INTO SERVICE

- Unpack the appliance and accessories, check that they have not been damaged by transport. Any damage should be notified to the carrier immediately. In case of damage, do not use the appliance and / or accessories and have them checked by qualified personnel. In case of return of the appliance it is necessary to use the original packaging of the product or a packaging that guarantees equivalent safety for transport.
- Carefully study the documentation and operating instructions provided.
- Check that the available supply voltage complies with that required by the unit (see plate data on the bottom of the unit).
- Put the power switch on the left side to O (off).
- Plug the supplied power cord into the socket on the back.
- Insert the handpiece connector on the bottom front. Screw the ring nut on the handpiece connector. When not in use, it is recommended to store the handpiece in the case provided.
- Turn the unit on, setting the switch to position I (on). The successful ignition is indicated, in addition to the switching on of the alphanumeric display, by the lighting of the green power light. Now the device is ready to be used.
- Use the instructions in the "Applications" chapter.
- Start treatment according to the instructions in the chapter "Treatment".

APPLICATIONS

Any treatment protocols should be established by the attending physician. Below is a table with some treatments. For treatment points, please refer to the paragraph "POINTS OF APPLICATION".

The number of treatments needed depends on the pathology treated and varies from patient to patient. The interval between treatments for analgesic reasons, should be 24-48 hours. If these are carried out after a few hours, results are obtained more quickly. In general, treatments should be the most frequent initially and spaced subsequently according to the benefits obtained. In some pathologies, treatment cycles are necessary interspersed with a rest period after which treatment is resumed. In the event of the onset of a new algic crisis after a first cycle of treatment, it is normally sufficient to proceed with a reduced cycle of three / four sessions to fully restore the state of analgesia. Algic remission has variable trends from patient to patient.

APPLICATION TABLE

The P column shows the preset programs in the drive. For others, use the free program mode (see chapter "Treatment").

The duration of the treatments indicated refer to using handpieces with a maximum power of 30W, in the case of handpieces with lower or higher power the treatment time (of the preset programs) is automatically recalculated by the unit based on the maximum power of the handpiece.

<i>P</i>	<i>Treatment</i>	<i>FREQUENCY (Hz)</i>	<i>DURATION (min)</i>
X	ARTHRITIS	5000	30
X	CERVICAL ARTHROSIS	500	30
X	CERVICAL PAIN	500	30
X	COXO-ARTHROSIS	10000	20
X	DORSAL PAIN	500	20
X	EPICONDYLITIS	1000	20
X	EPITROCLEYTIS	1000	20
X	EXTENSION	1000	20
X	GONOARTRHOSIS	1000	35
X	HEMATOMA	1000	20
X	LUMBAGO	500	20
	MYALGIA TRAPEZIUS MUSCLE	8000	20
	OCCIPITAL NEURALGIA	10000	20
X	PAINFUL SHOULDER	500	20
	PARASCAPULAR PAIN	7000	30
X	SCIATICA	5000	20
	SUBDELTOID BURSITIS	10000	25
	TEMPOROMANDIBULAR PAIN	7000	20
X	TENDON PAIN	5000	20
X	TENDON VAGINITIS	8	30
X	TUNNEL CARPALE SYND	8	30

TIPS FOR BEST RESULTS

IMPORTANT: The handpiece must be used in contact with the area / point to be treated, it is advisable to use the contact sensor (see Treatment – Contact Control). This allows maximum power density to be transferred to the fabric, as well as minimizing the potential risk of accidental beam vision.

The contact method also allows to exert pressure on trigger or acupuncture points obtaining a synergy of effects.

When the contact method cannot be used (for example treatment on very painful parts or for pathologies where strict asepsis is required), the tip of the handpiece should not be kept at a distance greater than 0.5-1 cm.


Laser therapy should be considered adjuvant and not a substitute for the chemotherapy therapy or surgery indicated. The surface to be irradiated must be prepared, paying particular attention to cleaning and cleaning it from any residues of medicinal products. In case of inflammatory processes, associate laser treatment with antibiotic therapy. Any topical medications should be applied after laser therapy and not before. It is advisable to maintain the handpiece in such a way that the surface to be irradiated is normal to the radiation emergency plane so as to reduce the reflection of the beam.


In normal use of the handpieces supplied it is necessary to use glasses and it is absolutely necessary to avoid direct exposure to the laser beam in the eye area. The infrared laser beam ($\lambda=905$ nm) is invisible.


To reduce the amount of reflected energy it is advisable, in addition to keeping the handpiece in contact, to clean the surface to be treated and direct the laser beam as much as possible in a direction orthogonal to the surface to be irradiated.

POINTS OF APPLICATION


Below are some points of application.

TREATMENT

When switched on, the home screen appears on the display via the power switch on the side of the unit. The lower part displays the version of the installed software. Press the ENTER key.


The screen appears to select the program to use. The unit allows you to work with a free program, where the parameters to be used are entered from time to time, or with preset programs. To choose the program, press the UP (to go up) and DWN (to go down) keys, choose the desired program, press the ENTER key to move to the next screen.

IMPORTANT: In order to use a program (free or present) you must first have connected a handpiece, otherwise an error message is displayed.

WARNING: During use it is necessary to use laser screen glasses and direct exposure to the laser beam in the eye area is absolutely to be avoided. The infrared laser beam ($\lambda=905$ nm) is invisible.

Warning Mask

If the handpiece is not connected, the following screen is displayed:

To use the equipment you need to connect a handpiece.

Contact Control

The device has a safety control system (which can be activated) with automatic deactivation of laser irradiation in case of failure or incorrect positioning of the emitter handpiece with the part to be treated.


To enable or disable the control, follow these steps:

1. Turn on the unit, using the switch on the side, while holding down the ENTER button on the front.

After a few seconds, one of the following screens appears:


2. Release the ENTER button and using the UP and DWN buttons, select the second line ("Contact Control") and press ENTER again. One of the following screens appears:


3. Using the UP and DWN buttons, select the desired option and press ENTER again. You return to the home screen.

Free Program


Once the "Free Program" has been chosen, the screen is displayed to select the value of the working frequency. For the purpose of help, the unit stores all the latest used settings. To change the indicated value, press the UP (to increase the value) and Dwn (to decrease the value) keys, choose the desired value, press the ENTER key. Pressing the ESC key returns you to the previous screen.


After choosing the working frequency, the screen to select the working power is displayed. To change the indicated value, press the UP and DWN keys, choose the desired value, press the ENTER key. Pressing the ESC key returns you to the previous screen.


NOTE: The maximum usable power value depends on the type of handpiece connected.

After choosing the potency, the screen to select the duration of treatment is displayed. To change the value, press the UP and DWN keys. Press ENTER to confirm the value. Pressing the ESC key returns you to the previous screen.

Once the duration of the treatment has been chosen, the unit is ready to work, the word "PRESS START" flashes on the screen with an indication of the previously chosen settings below. To start the treatment press the Start/Stop button. Pressing the ESC key returns you to the previous screen of the selection of the duration of treatment.


In the Start state, the name of the program appears on the top, the arrow indicates the parameter that is changed by pressing the UP and DWN keys.

To switch between parameters, press the ENTER key.


The time of the duration of the treatment decreases only with the handpiece activated (red light on the handpiece on).

In case of use of the contact control system, bring the handpiece closer to the part to be treated, detected the contact the handpiece lights up and there is laser dispensing, moving the handpiece away the handpiece turns off and there is no dispensing, the time count also stops and the PAUSE indication appears on the display, approaching the handpiece again this turns on again and the time count begins again.

In case of non-use of the control system, the handpiece is stably lit. By pressing the Start/Stop button it is possible to stop the treatment.

On the screen appears the request to "Continue?", through the Up and Dwn buttons, select the desired option, in case of continuation of the treatment (SI), you return to the previous screen with the indication of the residual duration of the treatment. To resume treatment, press START again.

In case of NO answer or by pressing the ESC button, the treatment is finished.

At the natural end of the treatment time or after denial in a treatment interruption, the indication of the end of the treatment appears on the screen, pressing the ESC key, you return to the program selection screen.

Pre-Set Program

Through the Up and Dwn buttons you can select one of the programs in the unit.

PRESS START			
FREQUENCY	5000 Hz	FREQUENCY	5000 Hz
POWER	30 w	POWER	30 w
TIME	80:00 min	TIME	80:00 min

Once a program has been chosen, the unit is ready to work, the word "PRESS START" flashes on the screen with an indication of the parameters of the chosen program below.

NOTE: The duration of treatment may vary from the power of the connected handpiece.

To start the treatment press the Start/Stop button. Pressing the ESC key returns to the previous screen of the program selection.

FREQUENCY	5000 Hz
POWER	30 w
TIME	79:59 min

In the Start state, the name of the program appears at the top.


The time of the duration of the treatment decreases only with the handpiece activated (red light on the handpiece on).

In case of use of the contact control system, bring the handpiece closer to the part to be treated, detected the contact the handpiece lights up and there is laser dispensing, moving the handpiece away the handpiece turns off and there is no dispensing, the time count also stops and the PAUSE indication appears on the display, approaching the handpiece again this turns on again and the time count begins again.

In case of non-use of the control system, the handpiece is stably lit. By pressing the Start/Stop button it is possible to stop the treatment.

CONTINUE?
→YES
NO

On the screen appears the request to "Continue", through the Up and Dwn buttons, select the desired option, in case of continuation of the treatment (SI), you return to the previous screen with the indication of the residual duration of the treatment. To resume treatment, press START again.

CONTINUE?
YES
→ NO

In case of NO answer or by pressing the ESC button, the treatment is finished.

TREATMENT
ENDED

At the natural end of the treatment time or after denial in a treatment interruption, the indication of the end of the treatment appears on the screen, pressing the ESC key, you return to the program selection screen.

MAINTENANCE

PREVENTIVE MAINTENANCE

BEFORE PERFORMING ANY MAINTENANCE, DISCONNECT THE ELECTRICAL CONNECTIONS.

Preventive maintenance mainly consists of:

- checking the condition of equipment and accessories before each use to examine the integrity of all cables, insulation, casings, etc. designed to prevent access to live parts.
- periodic cleaning of the equipment to be carried out on a regular basis in order to keep the equipment in the best condition of functionality and appearance.

EQUIPMENT CLEANING

Clean the outer container and control panels with a damp cloth: it is recommended not to use abrasive products or solvents. Specifically, the control panel does not recommend the use of ethyl alcohol or detergents containing high percentages of ethyl alcohol. Do not immerse the equipment in liquids. In case of penetration of liquids do not use and contact professionally qualified personnel.

LASER HANDPIECES CLEANING AND DISINFECTION

For external cleaning of the handpieces use a soft cloth. Disinfect the parts in contact with the patient with an alcohol solution after each treatment.

Do not use aggressive detergents. Do not use abrasive products or solvents.

Gently clean the transparent plastic part, trying to avoid scratching it as much as possible. It is recommended, when not in use, to store the handpiece in the case provided.

Do not immerse the accessories in liquids, in case of penetration of liquids contact professionally qualified personnel.

NOTE: The appliance, periodically (at least once a year), should be checked by qualified personnel to verify the following electrical safety values:

- a measure of leakage currents.

CORRECTIVE MAINTENANCE

If a malfunction is found, it is suggested to first check that there is no error in the arrangement of the commands.

If the appliance(s) are subjected to external mechanical stress, e.g. after a serious fall, or if the apparatus or handpieces have been subjected to a dripping liquid, or if the apparatus or handpieces have been subjected to strong overheating (e.g. direct sunlight, fire), or if the functionality of the apparatus or handpieces appears to be altered or if parts of the housing of the apparatus or handpieces are broken, displaced or missing, or if any connector or cable shows signs of deterioration, the appliance and its accessories should be checked by professionally qualified personnel.

TROUBLESHOOTING GUIDE

The apparatus despite being with the switch on position I (on) does not Turn.	Verify that the power cord is inserted correctly. Make sure that your power outlet is actually working.
When turned on the unit does not seem to work properly.	Turn off the device, wait 30-45 seconds and try to turn it back on.
The device is starting but the handpiece does not turn on.	Check that the unit is not working with the detection system active (see Chapter Treatment - Contact Control). Bring the tip of the handpiece closer to the area to be treated. If the problem persists, try cleaning the clear plastic part (see Maintenance). If the transparent plastic part has scratches that need to be replaced.
With the contact system inserted, the light on the handpiece flashes and the unit sounds intermittently.	Bring the handpiece closer to the area to be treated, until the light is constantly switched on. If the problem persists, try cleaning the clear plastic part (see Maintenance). If the transparent plastic part has scratches, it is necessary Replace.

If, despite the interventions, you continue to have these problems, have the equipment and accessories checked by qualified personnel. Do not open the equipment for any reason, this in addition to being potentially dangerous, will drop any form of warranty.

WARRANTY

- The appliance is guaranteed for 24 months, any electric accumulators are guaranteed for 3 months, the accessory parts of use are excluded from the warranty.
- In the case of warranty service, a copy of the purchase document proving the date of purchase must be sent together with the device. If the user is not able to show the documents proving the date of purchase, the repair will still be carried out under warranty if on the date of intervention, of our personnel in charge, no more than 24 months have passed from the date of manufacture of the device.
- During the warranty period, all parts that at the sole discretion of the manufacturer are defective in material or workmanship will be replaced or repaired free of charge. The service is carried out at the manufacturer's service centre and the device is delivered to the user.
- Mobile or removable plastic parts, lamps, glass parts, external pipes and any accessories are not covered by the warranty.

The warranty does not apply and the service will be performed entirely to

Payment in the following cases:

- For imperfect or incomplete installation.
- For damage due to transport.
- For damage due to accidents (lightning, earthquakes, fires, etc.).
- If the serial number has been removed, deleted or altered.

- If the owner of the appliance is not the first buyer.
- For negligent use, carelessness or poor maintenance.
- For insufficient flow rate or abnormality of electrical systems.
- For repairs, modifications or tampering carried out by personnel not authorized by the manufacturer and in any case for reasons not dependent on the manufacturer.
- These conditions may not be modified by any other oral or written agreement.
- The warranty excludes any compensation for the period of inefficiency of the appliance.

EMC DECLARATIONS

EMISSIONS

<i>Guide and manufacturer's declaration – electromagnetic emissions</i>		
laserplus is intended to operate in the electromagnetic environment specified below. The Customer or user of the device should ensure that it is used in this environment.		
<i>Emission Test</i>	<i>Conformity</i>	<i>Electromagnetic environment</i>
RF CISPR 11 Emissions	Group 1	laserplus uses RF energy only for its internal functioning. As a result, its RF emissions are very low and probably does not cause any interference in nearby electronic devices.
RF CISPR 11 Emissions	Class A	laserplus is suitable for use in all environments, except domestic and those connected directly to a low-voltage public mains power supply that powers buildings used for domestic purposes.
Emission of harmonics on the electricity network IEC 61000-3-2	Class A	
Voltage fluctuations and flicker generation IEC 61000-3-3	Conform	

IMMUNITY

Electromagnetic immunity tests are intended to simulate disturbances typical of an electromagnetic environment. laserplus has been tested for immunity to disorders and their levels typical of a domestic, hospital and commercial environment.

Guidance and manufacturer's declaration - electromagnetic immunity		
laserplus is intended to operate in the electromagnetic environment specified below. The Customer or the user of the device should ensure that it is used in this environment.		
Immunity Test	Conformity	Electromagnetic Environment and Precautions
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Floors should be made by wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Transient/sequence of fast electrical pulses IEC 61000-4-4	±2 kV for power supply lines	The mains power supply quality and the magnetic fields at mains frequency should be the typical one used in domestic, commercial or hospital environment.
Impulse IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	
Voltage dips, short interruptions, and voltage variations on power input lines IEC 61000-4-11	<5 % of the nominal voltage (U_T) (Voltage dip >95 %) for half a cycle 70 % U_T (Voltage dip 30 %) for 25 cycles <5 % U_T (Voltage dip >95 %) for 5 sec	If the use of the device requires it to operate continuously during power outages, it is recommended to power it from an uninterruptible power supply (UPS) or batteries.

<p>Magnetic fields at mains frequency (50/60 Hz) IEC 61000-4-8</p>	<p>30 A/m</p>			<p>The magnetic fields at mains frequency should have characteristic levels of a typical location in a commercial or hospital environment.</p>
<p>Conducted RF Fields IEC 61000-4-6</p>	<p>3 Vrms from 150 kHz to 80 MHz</p>			<p>Electromagnetic interference may occur near the devices bearing the symbol below.</p> 
<p>Radiated RF Fields IEC61000-4-3</p>	<p>Field (V/m)</p>	<p>Frequency</p>	<p>Modulation</p>	
	<p>3</p>	<p>80MHz÷2700MHz</p>	<p>1kHz AM 80%</p>	
	<p>27</p>	<p>380MHz÷390MHz</p>	<p>18Hz PM 50%</p>	
	<p>28</p>	<p>430MHz÷470MHz</p>	<p>18Hz PM 50%</p>	
	<p>9</p>	<p>704MHz÷787MHz</p>	<p>217Hz PM 50%</p>	
	<p>28</p>	<p>800MHz÷960MHz</p>	<p>18Hz PM 50%</p>	
	<p>28</p>	<p>1700MHz÷1990MHz</p>	<p>217Hz PM 50%</p>	
	<p>28</p>	<p>2400MHz÷2570MHz</p>	<p>217Hz PM 50%</p>	
<p>9</p>	<p>5100MHz÷5800MHz</p>	<p>217Hz PM 50%</p>		

Information about elimination of this product (Applicable in the European Union and other European countries with separate collection systems)

On the end of the life, the present product mustn't be eliminated as urban refusal, but it must be eliminated in a separated collection.

If the product is eliminated in unsuitable way, it is possible that some parts of the product (for example some accumulators) could be negative for the environment and for the human health.

The symbol on the side (barred dustbin on wheel) denotes that the products mustn't throw into urban refuses container but it must be eliminated with separate collection.

In case of abusive elimination of this product, could be foreseen sanctions.

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