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CE I Class Medical Device
REGULATION (EU) 2017/745 OF THE EUROPEAN
PARLIAMENT AND OF THE COUNCIL
of 5 April 2017 concerning medical devices

CODICES
ST422/R RED TOURNIQUET
ST422/B BLUE TOURNIQUET

INTRODUCTION

Thank you for purchasing a tourniquet EASYRED by Moretti S.p.A. Its design and quality guarantees comfort, safety and reliability to the user. Easyred by Moretti S.p.A. tourniquets have been designed to satisfy all your demands. This instruction manual contains advices for your safety and for a correct use of the device. In case of doubts please contact your retailer.

INTENDED USE

The tourniquet is a device that blocks the blood flow. It can be used during blood sampling or for the treatment of external hemorrhage in first aid.

DECLARATION OF CONFORMITY

LMoretti S.p.A. declares under its sole responsibility that the product made and traded by Moretti S.p.A. and belonging to the group of TOURNIQUET complies with the provisions of the regulation 2017/745 on MEDICAL DEVICES of 5 April 2017. For this purpose, Moretti S.p.A. guarantees and declares under its sole responsibility what follows:

1. The devices satisfy the requirements of general safety and performance requested by the Annex I of regulation 2017/745 as laid down by the Annex IV of the above mentioned regulation.
2. The devices ARE NOT MEASURING INSTRUMENTS.
3. The devices ARE NOT MADE FOR CLINICAL TESTS.
4. The devices are packed in NON-STERILE BOX.
5. The devices belong to class I in accordance with the provisions of Annex VIII of the above mentioned regulation
6. Moretti S.p.A. provides to the Competent Authorities the technical documentation to prove the conformity to the 2017/745 regulation, for at least 10 years from the last lot production.

NOTE: Complete product codes, the manufacturer registration code (SRN), the UDI-DI code and any references to used regulations are included in the EU declaration of conformity that Moretti S.p.A. releases and makes available through its channels.

SYMBOLS

- | | |
|--|------------------------------|
| | Product code |
| | Unique Device Identification |
| | CE mark |
| | Manufacturer |
| | Batch Lot |
| | Read the instruction manual |
| | Medical Device |
| | Conditions of disposal |

HOW TO USE THE TOURNIQUETS

- 
1. Push the button on the buckle and open it.
 2. Put the tourniquet around the arm and close the buckle;
 3. Pull the elastic band outward until the arm is blocked;
 4. Push again the button on the buckle to open the tourniquet.

WARNING!

Read carefully this instruction manual for a correct use. Pay attention to the presence of children.

CLEANING

Use warm water and neutral soap to clean the product.

GENERAL WARNINGS

Do not use the product for a purpose not indicated in this manual. Moretti S.p.A. declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration of the product. The manufacturer can modify the device or the instruction manual without prior notification.

The user and/or the patient will have to report any serious accident that have occurred related to the device to the manufacturer and appropriate authority of the State which the user and/or patient belongs to.

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WARRANTY

The product has 2-year warranty from purchase date which is certified by the receipt or invoice. The distributor reserves the right either to repair or replace necessary parts to correct defects in the materials or workmanship without any charge. The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alteration or disassembly by unauthorized individuals. In case you need to get the product repaired and you fulfil the conditions of warranty, send it to the nearest retailer together with the purchase receipt.

WARRANTY CERTIFICATE

Product _____
Purchasing date _____
Authorized dealer _____
Street _____ Place _____
Sold to _____
Street _____ Place _____

Moretti S.p.a.
INTERNET: www.morettispa.com
E-MAIL: info@morettispa.com
MADE in P.R.C.

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