




























## Italiano

- | Simbolo   | Definizione                                  | Simbolo   | Definizione   |
|---|--|---|---|
|  | Parte applicata di tipo BF                   |  | Attenzione  |
|  | Seguire le istruzioni per l'uso              | <b>%SpO<sub>2</sub></b>   | Saturazione di ossigeno                                     |
| <b>PR bpm</b>   | Frequenza del polso (BPM)                    |  | Indicatore di bassa potenza                                 |
|  | Nessun allarme SpO <sub>2</sub>              | <b>SN</b>   | Numero di serie.  |
|  | Temperatura di stoccaggio e umidità relativa | <b>IP22</b>   | Il grado di protezione contro l'ingresso di polvere e acqua |
|  | Data di produzione                           |  | Dati del fabbricante  |
|  | Conformità alla Direttiva RAEE               | <b>?</b>  | Indica che il segnale non è stabile                         |
|  | Informazioni sull'importatore                | <b>MD</b>   | Dispositivo medico  |
| <b>PI%</b>  | Indice di perfusione                         | <b>RR</b>   | Frequenza respiratoria                                      |
|  | Omologazione Unione Europea                  | <b>EC REP</b>   | Rappresentante autorizzato nella comunità europea           |

## English

**Warnings!**

- Keep the oximeter away from young children. Small items such as the battery door, battery, and lanyard are choking hazards.
- Do not hang the lanyard from the device's electrical wire.

- | Symbol  | Definition                                | Symbol  | Definition  |
|---|---|---|---|
|  | Type BF applied part.                     |  | Storage temperature and relative humidity           |
|  | Attention                                 |  | Manufacturer's information                          |
|  | Date of Manufacture                       | <b>IP22</b>   | Protected against dripping water                    |
| <b>%SpO<sub>2</sub></b>   | Oxygen saturation                         |  | European union approval                             |
| <b>PR bpm</b>   | Pulse rate (BPM)                          |  | Authorized representative in the European community |
|  | Low power indication                      |  | consult accompanying documents.                     |
|  | No SpO <sub>2</sub> Alarm                 |  | Serial No.  |
|  | Waste electrical and electronic equipment | <b>?</b>  | Indicate the signal is not stable                   |
|  | Importer's information                    |  | Medical Device                                      |
| <b>PI%</b>  | Perfusion index                           | <b>RR</b>   | Respiration Rate                                    |



